


ENT Instruments – Micro Hartman Alligator Forceps

Hartmann Forceps Micro Black is used for the insertion or removal of Grommets / T Tubes and Foreign Bodies. The Forceps is black to reduce the glare of examination lights.



TWS 3690-01/04B	Micro Hartman Alligator Forceps, Grip Handle, Fine Tip, 90mm Length, 4mm Serrated Jaw, Straight, Black Coated
TWS 3690-01/05B	Micro Hartman Alligator Forceps, Grip Handle, Fine Tip, 90mm Length, 5mm Serrated Jaw, Straight, Black Coated
TWS 3690-01/06B	Micro Hartman Alligator Forceps, Grip Handle, Fine Tip, 90mm Length, 6mm Serrated Jaw, Straight, Black Coated
TWS 3690-01/07B	Micro Hartman Alligator Forceps, Grip Handle, Fine Tip, 90mm Length, 7mm Serrated Jaw, Straight, Black Coated

Crocodile Fine Jaw Micro Forceps is used for wax removal, foreign body removal and grommet placement



TWS 3691-01/03B	Micro Crocodile Forceps, 80mm Length, 3mm Serrated Jaw, Straight, Black Coated
TWS 3691-01/06B	Micro Crocodile Forceps, 80mm Length, 6mm Serrated Jaw, Straight, Black Coated