



**TW**  
Surgical

# SILICONE SHEETS AND BLOCKS



**UK  
CA** **CE**



TW Surgical

Treasure Wright (U.K.) Ltd.



[www.twrightmedical.com](http://www.twrightmedical.com)




[info@twrightmedical.com](mailto:info@twrightmedical.com)

# TW Surgical

## Medical Grade Silicone Sheets and Blocks

Firm, Plain and Reinforce



Shapes and Sizes



Treasure Wright (U.K.) Ltd.

128 -160 City Road, London, EC1V 2NX, United Kingdom  
[info@twrightmedical.com](mailto:info@twrightmedical.com)



[www.twrightmedical.com](http://www.twrightmedical.com)



## TWS Medical Grade Silicone Sheeting

TW Surgical Medical Grade Silicone Sheeting is a flexible silicone elastomer material designed for various medical applications.

Our Reinforced Sheeting is wire mesh embedded in silicone. The wire reinforcement minimizes stretching of the material, and allows for reliable attachment of sutures.




- Selection: Available in a variety of thicknesses in different Size and Shape
- Flexible: For ease of application
- Bio-compatibility: Our silicone has a high degree of bio-compatibility
- Adhesion: Reduces adhesion formation after the surgery and easily removed
- Twisted: Easily twisted and bent
- Supporting: To prevent a prolapse structural tissue
- Sterilization: Be cleaned and sterilized by steam autoclave
- Size and Shape Section: Be cut to the required size and shape
- Position Fixation: be sutured in the required wound position
- Packaging: Silicone sheet is enclosed between clean-room grade polyethylene sheets

Treasure Wright (U.K.) Ltd.  
128 -160 City Road, London, EC1V 2NX, United Kingdom  
[info@twrightmedical.com](mailto:info@twrightmedical.com)

## TWS Medical Grade Silicone Sheeting

- Easy to cut to required size and shape using scissors or knife
- Reinforcement minimizes stretching and allows for suture attachment
- Easy removal with minimal trauma
- Convenient storage of unused portion in a poly bag
- TW Surgical Medical Grade Silicone Sheeting is a flexible silicone elastomer material designed for various medical applications.
- Our Reinforced Sheeting is wire mesh embedded in silicone. The wire reinforcement minimizes stretching of the material, and allows for reliable attachment of sutures.



### APPLICATIONS:


- Any planned medical procedure deemed necessary by a medical professional
- Surgical repair of trauma injuries by medical physicians
- Cell growth and scar healing
- No control over the conditions of use, patient selection, surgical procedure
- The usage will be according to our instruction manual



## TWS Medical Grade Silicone Sheeting

TW Surgical Medical Grade Silicone Sheeting meet the safety test is produced as follows:

- |                             |                       |
|-----------------------------|-----------------------|
| 1. Hemolysis                | < 1 %                 |
| 2. Pyrogenicity             | Non pyrogenic         |
| 3. Intracutaneous Injection | Non Irritating        |
| 4. System Injection         | Non toxic             |
| 5. Intramuscular Implant    |                       |
| 3 days                      | Trauma free           |
| 10 days                     | Trauma free           |
| 90 days                     | Trauma free           |
| 6. Cell Tissue Culture      | No Cytopathic Effects |



Every production lot of elastomer is tested for levels of trace metals and for absence of cytopathic effects using a tissue culture test


The silicone Sheet is provided NON STERILE and must be sterilized before use.



TW Surgical

## TWS Medical Grade Silicone Sheeting – Square 0.5mm Thickness

### Transparent with Wire Pattern



TWS 3-33030/05T	TWS Medical Grade Reinforced Silicone Sheet, Square 3cm x 3cm x 0.5mm, Silicone and Dacron
TWS 3-35050/05T	TWS Medical Grade Reinforced Silicone Sheet, Square 5cm x 5cm x 0.5mm, Silicone and Dacron
TWS 3-37070/05T	TWS Medical Grade Reinforced Silicone Sheet, Square 7cm x 7cm x 0.5mm, Silicone and Dacron
TWS 3-37676/05T	TWS Medical Grade Reinforced Silicone Sheet, Square 7.6cm x 7.6cm x 0.5mm, Silicone and Dacron
TWS 3-3100100/05T	TWS Medical Grade Reinforced Silicone Sheet, Square 10cm x 10cm x 0.5mm, Silicone and Dacron
TWS 3-3150150/05T	TWS Medical Grade Reinforced Silicone Sheet, Square 15cm x 15cm x 0.5mm, Silicone and Dacron

**Different type of sizes and shapes to be consulted and produced**

**Treasure Wright (U.K.) Ltd.**


128 -160 City Road, London, EC1V 2NX, United Kingdom

[info@twrightmedical.com](mailto:info@twrightmedical.com)



## TWS Medical Grade Silicone Sheeting–Square 1mm Thickness

### Transparent with Wire Pattern



TWS 3-33030/10T	TWS Medical Grade Reinforced Silicone Sheet, Square 3cm x 3cm x 1mm, Silicone and Dacron
TWS 3-35050/10T	TWS Medical Grade Reinforced Silicone Sheet, Square 5cm x 5cm x 1mm, Silicone and Dacron
TWS 3-37070/10T	TWS Medical Grade Reinforced Silicone Sheet, Square 7cm x 7cm x 1mm, Silicone and Dacron
TWS 3-37676/10T	TWS Medical Grade Reinforced Silicone Sheet, Square 7.6cm x 7.6cm x 1mm, Silicone and Dacron
TWS 3-3100100/10T	TWS Medical Grade Reinforced Silicone Sheet, Square 10cm x 10cm x 1mm, Silicone and Dacron
TWS 3-3150150/10T	TWS Medical Grade Reinforced Silicone Sheet, Square 15cm x 15cm x 1mm, Silicone and Dacron


**Different type of sizes and shapes to be consulted and produced**

Treasure Wright (U.K.) Ltd.  
128 -160 City Road, London, EC1V 2NX, United Kingdom  
info@twrightmedical.com



## TWS Medical Grade Silicone Sheeting—Rectangle 0.5mm Thickness

### Transparent with Wire Pattern



TWS 3-32030/05T	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 2cm x 3cm x 0.5mm, Silicone and Dacron
TWS 3-32050/05T	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 2cm x 5cm x 0.5mm, Silicone and Dacron
TWS 3-32070/05T	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 2cm x 7cm x 0.5mm, Silicone and Dacron
TWS 3-320100/05T	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 2cm x 10cm x 0.5mm, Silicone and Dacron
TWS 3-320150/05T	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 2cm x 15cm x 0.5mm, Silicone and Dacron


**Different type of sizes and shapes to be consulted and produced**

Treasure Wright (U.K.) Ltd.  
128 -160 City Road, London, EC1V 2NX, United Kingdom  
[info@twrightmedical.com](mailto:info@twrightmedical.com)



## TWS Medical Grade Silicone Sheeting–Rectangle 0.5mm Thickness

### Transparent with Wire Pattern



TWS 3-33030/05T	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 3cm x 3cm x 0.5mm, Silicone and Dacron
TWS 3-33050/05T	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 3cm x 5cm x 0.5mm, Silicone and Dacron
TWS 3-33070/05T	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 3cm x 7cm x 0.5mm, Silicone and Dacron
TWS 3-330100/05T	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 3cm x 10cm x 0.5mm, Silicone and Dacron,
TWS 3-330150/05T	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 3cm x 15cm x 0.5mm, Silicone and Dacron

**Different type of sizes and shapes to be consulted and produced**

**Treasure Wright (U.K.) Ltd.**

128 -160 City Road, London, EC1V 2NX, United Kingdom

[info@twrightmedical.com](mailto:info@twrightmedical.com)



## TWS Medical Grade Silicone Sheeting–Rectangle 0.5mm Thickness

### Transparent with Wire Pattern



TWS 3-34030/05T	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 4cm x 3cm x 0.5mm, Silicone and Dacron
TWS 3-34050/05T	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 4cm x 5cm x 0.5mm, Silicone and Dacron
TWS 3-34070/05T	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 4cm x 7cm x 0.5mm, Silicone and Dacron
TWS 3-340100/05T	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 4cm x 10cm x 0.5mm, Silicone and Dacron
TWS 3-340150/05T	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 4cm x 15cm x 0.5mm, Silicone and Dacron


**Different type of sizes and shapes to be consulted and produced**

Treasure Wright (U.K.) Ltd.  
128 -160 City Road, London, EC1V 2NX, United Kingdom  
[info@twrightmedical.com](mailto:info@twrightmedical.com)



## TWS Medical Grade Silicone Sheeting–Rectangle 0.5mm Thickness

### Transparent with Wire Pattern



TWS 3-35030/05T	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5cm x 3cm x 0.5mm, Silicone and Dacron
TWS 3-35050/05T	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5cm x 5cm x 0.5mm, Silicone and Dacron
TWS 3-35070/05T	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5cm x 7cm x 0.5mm, Silicone and Dacron
TWS 3-350100/05T	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5cm x 10cm x 0.5mm, Silicone and Dacron
TWS 3-350150/05T	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5cm x 15cm x 0.5mm, Silicone and Dacron

**Different type of sizes and shapes to be consulted and produced**

Treasure Wright (U.K.) Ltd.  
128 -160 City Road, London, EC1V 2NX, United Kingdom  
info@twrightmedical.com



## TWS Medical Grade Silicone Sheeting–Rectangle 1mm Thickness

### Transparent with Wire Pattern



TWS 3-32030/10T	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 2cm x 3cm x 1mm, Silicone and Dacron
TWS 3-32050/10T	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 2cm x 5cm x 1mm, Silicone and Dacron
TWS 3-32070/10T	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 2cm x 7cm x 1mm, Silicone and Dacron
TWS 3-320100/10T	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 2cm x 10cm x 1mm, Silicone and Dacron
TWS 3-320150/10T	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 2cm x 15cm x 1mm, Silicone and Dacron

**Different type of sizes and shapes to be consulted and produced**

**Treasure Wright (U.K.) Ltd.**


128 -160 City Road, London, EC1V 2NX, United Kingdom  
info@twrightmedical.com



TW Surgical

## TWS Medical Grade Silicone Sheeting – Rectangle 1mm Thickness

### Transparent with Wire Pattern



TWS 3-33030/10T	TWS Medical Grade Reinforced Silicone Sheet 3cm x 3cm x 1mm, Silicone and Dacron
TWS 3-33050/10T	TWS Medical Grade Reinforced Silicone Sheet 3cm x 5cm x 1mm, Silicone and Dacron
TWS 3-33070/10T	TWS Medical Grade Reinforced Silicone Sheet 3cm x 7cm x 1mm, Silicone and Dacron
TWS 3-330100/10T	TWS Medical Grade Reinforced Silicone Sheet 3cm x 10cm x 1mm, Silicone and Dacron
TWS 3-330150/10T	TWS Medical Grade Reinforced Silicone Sheet 3cm x 15cm x 1mm, Silicone and Dacron

**Different type of sizes and shapes to be consulted and produced**

Treasure Wright (U.K.) Ltd.  
128 -160 City Road, London, EC1V 2NX, United Kingdom  
[info@twrightmedical.com](mailto:info@twrightmedical.com)



TW Surgical

## TWS Medical Grade Silicone Sheeting – Rectangle 1mm Thickness

### Transparent with Wire Pattern



TWS 3-34030/10T	TWS Medical Grade Reinforced Silicone Sheet 4cm x 3cm x 1mm, Silicone and Dacron
TWS 3-34050/10T	TWS Medical Grade Reinforced Silicone Sheet 4cm x 5cm x 1mm, Silicone and Dacron
TWS 3-34070/10T	TWS Medical Grade Reinforced Silicone Sheet 4cm x 7cm x 1mm, Silicone and Dacron
TWS 3-340100/10T	TWS Medical Grade Reinforced Silicone Sheet 4cm x 10cm x 1mm, Silicone and Dacron
TWS 3-340150/10T	TWS Medical Grade Reinforced Silicone Sheet 4cm x 15cm x 1mm, Silicone and Dacron


**Different type of sizes and shapes to be consulted and produced**

Treasure Wright (U.K.) Ltd.  
128 -160 City Road, London, EC1V 2NX, United Kingdom  
info@twrightmedical.com



## TWS Medical Grade Silicone Sheeting – Rectangle 1mm Thickness

### Transparent with Wire Pattern



TWS 3-35030/10T	TWS Medical Grade Reinforced Silicone Sheet 5cm x 3cm x 1mm, Silicone and Dacron
TWS 3-35050/10T	TWS Medical Grade Reinforced Silicone Sheet 5cm x 5cm x 1mm, Silicone and Dacron
TWS 3-35070/10T	TWS Medical Grade Reinforced Silicone Sheet 5cm x 7cm x 1mm, Silicone and Dacron
TWS 3-350100/10T	TWS Medical Grade Reinforced Silicone Sheet 5cm x 10cm x 1mm, Silicone and Dacron
TWS 3-350150/10T	TWS Medical Grade Reinforced Silicone Sheet 5cm x 15cm x 1mm, Silicone and Dacron

**Different type of sizes and shapes to be consulted and produced**

Treasure Wright (U.K.) Ltd.  
128 -160 City Road, London, EC1V 2NX, United Kingdom  
info@twrightmedical.com



TW Surgical

## TWS Medical Grade Silicone Sheeting – Circle 0.5mm Thickness

### Transparent with Wire Pattern



TWS 3-42020C/05T	TWS Medical Grade Reinforced Silicone Sheet, Circle 2cm Diameter x 0.5mm Thickness, Silicone and Dacron
TWS 3-43030C/05T	TWS Medical Grade Reinforced Silicone Sheet 3cm Diameter x 0.5mm Thickness, Silicone and Dacron
TWS 3-44040C/05T	TWS Medical Grade Reinforced Silicone Sheet 4cm Diameter x 0.5mm Thickness, Silicone and Dacron
TWS 3-45050C/05T	TWS Medical Grade Reinforced Silicone Sheet 5cm Diameter x 0.5mm Thickness, Silicone and Dacron
TWS 3-47070C/05T	TWS Medical Grade Reinforced Silicone Sheet 7cm Diameter x 0.5mm Thickness, Silicone and Dacron
TWS 3-4100100C/05T	TWS Medical Grade Reinforced Silicone Sheet 10cm Diameter x 0.5mm Thickness, Silicone and Dacron

**Different type of sizes and shapes to be consulted and produced**

Treasure Wright (U.K.) Ltd.  
128 -160 City Road, London, EC1V 2NX, United Kingdom  
[info@twrightmedical.com](mailto:info@twrightmedical.com)



TW Surgical

## TWS Medical Grade Silicone Sheeting – Circle 1mm Thickness

### Transparent with Wire Pattern



TWS 3-42020C/10T	TWS Medical Grade Reinforced Silicone Sheet, Circle 2cm Diameter x 1mm Thickness, Silicone and Dacron
TWS 3-43030C/10T	TWS Medical Grade Reinforced Silicone Sheet 3cm Diameter x 1mm Thickness, Silicone and Dacron
TWS 3-44040C/10T	TWS Medical Grade Reinforced Silicone Sheet 4cm Diameter x 1mm Thickness, Silicone and Dacron
TWS 3-45050C/10T	TWS Medical Grade Reinforced Silicone Sheet 5cm Diameter x 1mm Thickness, Silicone and Dacron
TWS 3-47070C/10T	TWS Medical Grade Reinforced Silicone Sheet 7cm Diameter x 1mm Thickness, Silicone and Dacron
TWS 3-4100100C/10T	TWS Medical Grade Reinforced Silicone Sheet 10cm Diameter x 1mm Thickness, Silicone and Dacron

**Different type of sizes and shapes to be consulted and produced**

Treasure Wright (U.K.) Ltd.  
128 -160 City Road, London, EC1V 2NX, United Kingdom  
[info@twrightmedical.com](mailto:info@twrightmedical.com)



## **TWS Medical Grade Silicone Sheeting – Rectangle 0.5mm Thickness**

### **Transparent with Wire Pattern**

### **A Rectangular Portion cut off from Rectangle Silicone Sheets**



TWS 3-520100/05T -0560	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5mm x 6cm off from 2cm x 10cm x 0.5mm Thickness, Silicone and Dacron
TWS 3-520150/05T -05110	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5mm x 11cm off from 2cm x 15cm x 0.5mm Thickness, Silicone and Dacron
TWS 3-520200/05T -05150	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5mm x 15cm off from 2cm x 20cm x 0.5mm Thickness, Silicone and Dacron
TWS 3-530100/05T -0560	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5mm x 6cm off from 3cm x 10cm x 0.5mm Thickness, Silicone and Dacron
TWS 3-530150/05T -05110	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5mm x 11cm off from 3cm x 15cm x 0.5mm Thickness, Silicone and Dacron
TWS 3-530200/05T -05150	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5mm x 15cm off from 3cm x 20cm x 0.5mm Thickness, Silicone and Dacron

**Different type of sizes and shapes to be consulted and produced**

**Treasure Wright (U.K.) Ltd.**  
128 -160 City Road, London, EC1V 2NX, United Kingdom  
[info@twrightmedical.com](mailto:info@twrightmedical.com)



TW Surgical

## TWS Medical Grade Silicone Sheeting – Rectangle 0.5mm Thickness

### Transparent with Wire Pattern

### A Rectangular Portion cut off from Rectangle Silicone Sheets



TWS 3-540100/05T -0560	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5mm x 6cm off from 4cm x 10cm x 0.5mm Thickness, Silicone and Dacron
TWS 3-540150/05T -05110	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5mm x 11cm off from 4cm x 15cm x 0.5mm Thickness, Silicone and Dacron
TWS 3-540200/05T -05150	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5mm x 15cm off from 4cm x 20cm x 0.5mm Thickness, Silicone and Dacron
TWS 3-550100/05T -0560	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5mm x 6cm off from 5cm x 10cm x 0.5mm Thickness, Silicone and Dacron
TWS 3-550150/05T -05110	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5mm x 11cm off from 5cm x 15cm x 0.5mm Thickness, Silicone and Dacron
TWS 3-550200/05T -05150	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5mm x 15cm off from 5cm x 20cm x 0.5mm Thickness, Silicone and Dacron

**Different type of sizes and shapes to be consulted and produced**

Treasure Wright (U.K.) Ltd.  
128 -160 City Road, London, EC1V 2NX, United Kingdom  
info@twrightmedical.com



## TWS Medical Grade Silicone Sheeting – Rectangle 1mm Thickness

### Transparent with Wire Pattern

### A Rectangular Portion cut off from Rectangle Silicone Sheets



TWS 3-520100/10T -0560	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5mm x 6cm off from 2cm x 10cm x 1mm Thickness, Silicone and Dacron
TWS 3-520150/10T -05110	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5mm x 11cm off from 2cm x 15cm x 1mm Thickness, Silicone and Dacron
TWS 3-520200/10T -05150	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5mm x 15cm off from 2cm x 20cm x 1mm Thickness, Silicone and Dacron
TWS 3-530100/10T -0560	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5mm x 6cm off from 3cm x 10cm x 1mm Thickness, Silicone and Dacron
TWS 3-530150/10T -05110	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5mm x 11cm off from 3cm x 15cm x 1mm Thickness, Silicone and Dacron
TWS 3-530200/10T -05150	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5mm x 15cm off from 3cm x 20cm x 1mm Thickness, Silicone and Dacron

**Different type of sizes and shapes to be consulted and produced**

Treasure Wright (U.K.) Ltd.  
128 -160 City Road, London, EC1V 2NX, United Kingdom  
info@twrightmedical.com



## TWS Medical Grade Silicone Sheeting – Rectangle 1mm Thickness

### Transparent with Wire Pattern

### A Rectangular Portion cut off from Rectangle Silicone Sheets



TWS 3-540100/10T -0560	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5mm x 6cm off from 4cm x 10cm x 1mm Thickness, Silicone and Dacron
TWS 3-540150/10T -05110	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5mm x 11cm off from 4cm x 15cm x 1mm Thickness, Silicone and Dacron
TWS 3-540200/10T -05150	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5mm x 15cm off from 4cm x 20cm x 1mm Thickness, Silicone and Dacron
TWS 3-550100/10T -0560	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5mm x 6cm off from 5cm x 10cm x 1mm Thickness, Silicone and Dacron
TWS 3-550150/10T -05110	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5mm x 11cm off from 5cm x 15cm x 1mm Thickness, Silicone and Dacron
TWS 3-550200/10T -05150	TWS Medical Grade Reinforced Silicone Sheet, Rectangle 5mm x 15cm off from 5cm x 20cm x 1mm Thickness, Silicone and Dacron

**Different type of sizes and shapes to be consulted and produced**

Treasure Wright (U.K.) Ltd.  
128 -160 City Road, London, EC1V 2NX, United Kingdom  
info@twrightmedical.com



## TWS Medical Grade Silicone Sheeting – Circle 0.5mm Thickness

### Transparent with Wire Pattern

### A Circular Portion cut off from Circular silicone Sheets



TWS 3-52020C/05T -1010	TWS Medical Grade Reinforced Silicone Sheet, Circle 1 cm off from 2cm Diameter x 0.5mm Thickness, Silicone and Dacron
TWS 3-53030C/05T -1010	TWS Medical Grade Reinforced Silicone Sheet, Circle 1 cm off from 3cm Diameter x 0.5mm Thickness, Silicone and Dacron
TWS 3-54040C/05T -1010	TWS Medical Grade Reinforced Silicone Sheet, Circle 1 cm off from 4cm Diameter x 0.5mm Thickness, Silicone and Dacron
TWS 3-55050C/05T -1010	TWS Medical Grade Reinforced Silicone Sheet, Circle 1 cm off from 5cm Diameter x 0.5mm Thickness, Silicone and Dacron
TWS 3-56060C/05T -1010	TWS Medical Grade Reinforced Silicone Sheet, Circle 1 cm off from 6cm Diameter x 0.5mm Thickness, Silicone and Dacron
TWS 3-57070C/05T -1010	TWS Medical Grade Reinforced Silicone Sheet, Circle 1 cm off from 7cm Diameter x 0.5mm Thickness, Silicone and Dacron

**Different type of sizes and shapes to be consulted and produced**

Treasure Wright (U.K.) Ltd.  
128 -160 City Road, London, EC1V 2NX, United Kingdom  
info@twrightmedical.com



TW Surgical

## TWS Medical Grade Silicone Sheeting – Circle 1mm Thickness

### Transparent with Wire Pattern

### A Circular Portion cut off from Circular silicone Sheets



TWS 3-52020C/10T -1010	TWS Medical Grade Reinforced Silicone Sheet, Circle 1 cm off from 2cm Diameter x 1mm Thickness, Silicone and Dacron
TWS 3-53030C/10T -1010	TWS Medical Grade Reinforced Silicone Sheet, Circle 1 cm off from 3cm Diameter x 1mm Thickness, Silicone and Dacron
TWS 3-54040C/10T -1010	TWS Medical Grade Reinforced Silicone Sheet, Circle 1 cm off from 4cm Diameter x 1mm Thickness, Silicone and Dacron
TWS 3-55050C/10T -1010	TWS Medical Grade Reinforced Silicone Sheet, Circle 1 cm off from 5cm Diameter x 1mm Thickness, Silicone and Dacron
TWS 3-56060C/10T -1010	TWS Medical Grade Reinforced Silicone Sheet, Circle 1 cm off from 6cm Diameter x 1mm Thickness, Silicone and Dacron
TWS 3-57070C/10T -1010	TWS Medical Grade Reinforced Silicone Sheet, Circle 1 cm off from 7cm Diameter x 1mm Thickness, Silicone and Dacron

**Different type of sizes and shapes to be consulted and produced**

Treasure Wright (U.K.) Ltd.  
128 -160 City Road, London, EC1V 2NX, United Kingdom  
info@twrightmedical.com

# TW Surgical Medical Grade Silicone Sheets and Blocks

---

Treasure Wright (U.K.) Ltd.  
128 -160 City Road, London, EC1V  
2NX, United Kingdom

