




TW
Surgical

ELECTRO SURGICAL



**UK
CA CE**



TW Surgical

Treasure Wright (U.K.) Ltd.





www.twrightmedical.com



info@twrightmedical.com

BIPOLAR FORCEPS



McPherson
(Straight)
TWS 84-2600/01



McPherson
(Angled)
TWS 84-2600/02



Coaptation
Nadler
TWS 84-2600/03



Coaptation
TWS 84-2600/04



Jeweler #5
TWS 84-2600/05

Cat #	Description	Tip (mm)	Size
TWS 84-2600/01	McPherson (Str.)	0.5	3½" (9.0 cm)
TWS 84-2600/02	McPherson (Ang.)	0.5	3½" (9.0 cm)
TWS 84-2600/03	Coaptation Nadler	0.5	4" (10.2 cm)
TWS 84-2600/04	Coaptation	0.5	4¼" (10.8 cm)
TWS 84-2600/05	Jeweler #5	Micro, 0.5	4½" (11.5 cm)

For your convenience add the following tip extension to the product code.
Needle Tip (-00), 0.5mm(-05), 1.0mm(-10) & 2.0mm (-20)

BIPOLAR FORCEPS



Jeweler #4



Tenzel



Jeweler #7



Iris



Iris

Cat #	Description	Tip (mm)	Size
TWS 84-2600/11	Jeweler #4	Micro, 0.5	4½" (11.5 cm)
TWS 84-2600/12	Tenzel	Micro, 0.5	4" (10.2 cm)
TWS 84-2600/13	Jeweler #7	Micro, 0.5	4½" (11.5 cm)
TWS 84-2600/14	Iris (Str.)	1.0	4¼" (10.8 cm)
TWS 84-2600/15	Iris (Ang.)	1.0	4¼" (10.8 cm)

For your convenience add the following tip extension to the product code.
 Needle Tip (-00), 0.5mm(-05), 1.0mm(-10) & 2.0mm (-20)

BIPOLAR FORCEPS



Adson
TWS 84-2608/64



Adson
TWS 84-2608/65



Adson
(Angled)
TWS 84-2608/66



Semkin
TWS 84-2608/67



Cushing
TWS 84-2608/68

Cat #	Description	Tip (mm)	Size
TWS 84-2608/64	Adson	1.0	4¼" (12.1 cm)
TWS 84-2608/65	Adson	1.0	6" (15.2 cm)
TWS 84-2608/66	Adson (Ang.)	1.0	6" (15.2 cm)
TWS 84-2608/67	Semkin	0.5	5½" (14.0 cm)
TWS 84-2608/68	Cushing	1.0, 1.5	6" (15.2 cm)
TWS 84-2600/16	Cushing	1.0, 1.5	6½" (16.5 cm)

For your convenience add the following tip extension to the product code.
Needle Tip (-00), 0.5mm(-05), 1.0mm(-10) & 2.0mm (-20)

BIPOLAR FORCEPS



Cushing



Cushing



Gerald



Gerald




Gerald Bayonet

Cat #	Description	Tip (mm)	Size
TWS 84-2600/22	Cushing	1.0, 2.0	7" (17.8 cm)
TWS 84-2600/23	Cushing (Gursh)	2.0	7" (17.8 cm)
TWS 84-2600/24	Gerald (Str.)	1.0, 1.5	7" (17.8 cm)
TWS 84-2600/25	Gerald (Ang.)	1.0, 1.5	7" (17.8 cm)
TWS 84-2600/26	Gerald Bayonet	1.0, 1.5	7¼" (18.4 cm)

For your convenience add the following tip extension to the product code.
 Needle Tip (-00), 0.5mm(-05), 1.0mm(-10) & 2.0mm (-20)

BIPOLAR FORCEPS



Jansen Bayonet
TWS 84-2608/74

Jansen Bayonet
TWS 84-2608/75

Jansen Bayonet
TWS 84-2608/76

Malis
TWS 84-2608/77

Malis
TWS 84-2608/78

Cat #	Description	Tip (mm)	Size
TWS 84-2608/74	Jansen Bayonet	0.5, 1.0, 1.5	5¼" (13.3 cm)
TWS 84-2608/75	Jansen Bayonet	0.5, 1.0, 1.5	6" (15.2 cm)
TWS 84-2608/76	Jansen Bayonet	0.5, 1.0, 1.5	7" (17.8 cm)
TWS 84-2600/30	Jansen Bayonet	0.5, 1.0, 1.5	7½" (19.1 cm)
TWS 84-2600/31	Jansen Bayonet	1.0, 1.5, 2.0	8" (20.4 cm)
TWS 84-2608/77	Malis	0.5, 1.0, 2.0	7" (17.8 cm)
TWS 84-2608/78	Malis	0.5, 1.0, 2.0	7" (17.8 cm)
TWS 84-2600/34	Malis	0.5, 1.0, 2.0	8" (20.4 cm)
TWS 84-2600/35	Malis	0.5, 1.0, 2.0	9" (22.8 cm)

For your convenience add the following tip extension to the product code.
 Needle Tip (-00), 0.5mm(-05), 1.0mm(-10) & 2.0mm (-20)

BIPOLAR FORCEPS



Yasargil

Yasargil

Scoville-Greenwood

Hardy Bayonet

Hardy Bayonet (Angled)

Cat #	Description	Tip (mm)	Size
TWS 84-2600/36	Scoville-Greenwood	1.0, 1.5	7½" (19.1 cm)
TWS 84-2600/37	Yasargil	1.0, 1.5	7¾" (20.0 cm)
TWS 84-2600/38	Yasargil	1.0, 1.5	6" (15.2 cm)
TWS 84-2600/39	Hardy Bayonet	1.0	8½" (21.6 cm)
TWS 84-2600/40	Hardy Bayonet (Ang.)	1.0	8½" (21.6 cm)

For your convenience add the following tip extension to the product code.
 Needle Tip (-00), 0.5mm(-05), 1.0mm(-10) & 2.0mm (-20)

BIPOLAR FORCEPS (NON-STICK)



McPherson
(Straight)



McPherson
(Angled)



Jeweler #5



Jeweler #4

Cat #	Description	Tip (mm)	Size
TWS 84-2600/41	McPherson (Str.)	0.5	3½" (9 cm)
TWS 84-2600/42	McPherson (Ang.)	0.5	3½" (9 cm)
TWS 84-2600/43	Jeweler #5	0.5, 0.7	4½" (11.5 cm)
TWS 84-2600/44	Jeweler #4	0.5, 0.7	4½" (11.5 cm)

Convert any regular bipolar forceps listed on pages 04 to 09 into non-stick bipolar forceps, use "NS" extension for conversion. Available Tip Size: 0.7mm, 1.0mm, 1.5mm & 2.0mm

BIPOLAR FORCEPS (NON-STICK)



Adson
TWS 84-2608/88



Semkin
TWS 84-2608/89



Cushing
TWS 84-2608/90



Cushing
TWS 84-2608/91

Cat #	Description	Tip (mm)	Size
TWS 84-2608/88	Adson	1.0	4¼" (12.1 cm)
TWS 84-2608/89	Semkin	0.5, 1.0	5½" (14.0 cm)
TWS 84-2608/90	Cushing	1.0, 1.5	6½" (16.5 cm)
TWS 84-2608/91	Cushing	1.0, 2.0	7" (17.8 cm)

Convert any regular bipolar forceps listed on pages 04 to 09 into non-stick bipolar forceps, use "NS" extension for conversion. Available Tip Size: 0.7mm, 1.0mm, 1.5mm & 2.0mm

BIPOLAR FORCEPS (NON-STICK)



Gerald Bayonet
TWS 84-2608/92



Jansen Bayonet
TWS 84-2608/93



Scoville-Greenwood
TWS 84-2608/94



Yasargil
TWS 84-2608/95



Hardy Bayonet
TWS 84-2608/96

Cat #	Description	Tip (mm)	Size
TWS 84-2608/92	Gerald Bayonet	1.0, 1.5	7¼" (18.4 cm)
TWS 84-2608/93	Jansen Bayonet	0.5, 1.0, 1.5	7" (17.8 cm)
TWS 84-2608/94	Scoville-Greenwood	1.0, 1.5	7½" (19.1 cm)
TWS 84-2608/95	Yasargil	1.0, 1.5	7¾" (20.0 cm)
TWS 84-2608/96	Hardy Bayonet	1.0, 1.5, 2.0	8½" (21.6 cm)

Convert any regular bipolar forceps listed on pages 04 to 09 into non-stick bipolar forceps, use "NS" extension for conversion.
Available Tip Size:
0.7mm, 1.0mm, 1.5mm & 2.0mm

BIPOLAR FORCEPS (ULTRA NON-STICK)



McPherson (Str.)



McPherson (Ang.)



Jeweler #5



Jeweler #4

Cat #	Description	Tip (mm)	Size
TWS 84-2600/54	McPherson (Str.)	0.5	3½" (9 cm)
TWS 84-2600/55	McPherson (Ang.)	0.5	3½" (9 cm)
TWS 84-2600/56	Jeweler #5	0.5, 0.7	4½" (11.5 cm)
TWS 84-2600/57	Jeweler #4	0.5, 0.7	4½" (11.5 cm)

Convert any regular bipolar forceps listed on pages 04 to 09 into ultra non-stick bipolar forceps, use "UNS" extension for conversion. Available Tip Size: 0.7mm, 1.0mm, 1.5mm & 2.0mm



TW Surgical

BIPOLAR FORCEPS (ULTRA NON-STICK)



Adson
TWS 84-2609/01



Semkin
TWS 84-2609/02



Cushing
TWS 84-2609/03



Cushing
TWS 84-2609/04

Cat #	Description	Tip (mm)	Size
TWS 84-2609/01	Adson	1.0	4¼" (12.1 cm)
TWS 84-2609/02	Semkin	0.5, 1.0	5½" (14.0 cm)
TWS 84-2609/03	Cushing	1.0, 1.5	6½" (16.5 cm)
TWS 84-2609/04	Cushing	1.0, 2.0	7" (17.8 cm)

Convert any regular bipolar forceps listed on pages 04 to 09 into ultra non-stick bipolar forceps, use "UNS" extension for conversion.
Available Tip Size:
0.7mm, 1.0mm, 1.5mm & 2.0mm

BIPOLAR FORCEPS (ULTRA NON-STICK)



Gerald Bayonet
TWS 84-2609/05



Jansen Bayonet
TWS 84-2609/06



Scoville-Greenwood
TWS 84-2609/07



Yasargil
TWS 84-2609/08



Hardy Bayonet
TWS 84-2609/09

Cat #	Description	Tip (mm)	Size
TWS 84-2609/05	Gerald Bayonet	1.0, 1.5	7¼" (18.4 cm)
TWS 84-2609/06	Jansen Bayonet	0.5, 1.0, 1.5	7" (17.8 cm)
TWS 84-2609/07	Scoville-Greenwood	1.0, 1.5	7½" (19.1 cm)
TWS 84-2609/08	Yasargil	1.0, 1.5	7¾" (20.0 cm)
TWS 84-2609/09	Hardy Bayonet	1.0, 1.5, 2.0	8½" (21.6 cm)

Convert any regular bipolar forceps listed on pages 04 to 09 into ultra non-stick bipolar forceps, use "UNS" extension for conversion.
Available Tip Size:
0.7mm, 1.0mm, 1.5mm & 2.0mm



TW Surgical

BIPOLAR FORCEPS (ULTRA NON-STICK)



Cat #	Description	Tip (mm)	Size
TWS 84-2600/67	Jeweler #5	0.5, 0.7	4 1/2" (11.5 cm)
TWS 84-2600/68	Jeweler #5 (Ang.)	0.5, 0.7	4 1/2" (11.5 cm)
TWS 84-2600/69	Straight	0.5, 0.7, 1.0, 1.5, 2.0	7" (17.8 cm)
TWS 84-2600/70	Straight	0.5, 0.7, 1.0, 1.5, 2.0	7 3/4" (20.0 cm)
TWS 84-2600/71	Straight	0.5, 0.7, 1.0, 1.5, 2.0	8 3/4" (22.0 cm)
TWS 84-2600/72	Straight	0.5, 0.7, 1.0, 1.5, 2.0	9 1/2" (24.0 cm)

Convert any regular bipolar forceps listed on pages 04 to 09 into ultra non-stick bipolar forceps, use "UNS" extension for conversion. Available Tip Size: 0.7mm, 1.0mm, 1.5mm & 2.0mm

BIPOLAR FORCEPS (ULTRA NON-STICK)



Straight
(Angled)
TWS 84-2609/15

Straight
(Angled)
TWS 84-2609/16

Straight
(Angled)
TWS 84-2609/17

Bayonet
TWS 84-2609/18

Bayonet
TWS 84-2609/19

Bayonet
TWS 84-2609/20

Cat #	Description	Tip (mm)	Size
TWS 84-2609/15	Straight (Ang.)	0.5, 0.7, 1.0, 1.5, 2.0	7" (17.8 cm)
TWS 84-2609/16	Straight (Ang.)	0.5, 0.7, 1.0, 1.5, 2.0	7 3/4" (20.0 cm)
TWS 84-2609/17	Straight (Ang.)	0.5, 0.7, 1.0, 1.5, 2.0	8 3/4" (22.0 cm)
TWS 84-2600/76	Straight (Ang.)	0.5, 0.7, 1.0, 1.5, 2.0	9 1/2" (24.0 cm)
TWS 84-2609/18	Bayonet	0.5, 0.7, 1.0, 1.5, 2.0	6 1/2" (16.5 cm)
TWS 84-2609/19	Bayonet	0.5, 0.7, 1.0, 1.5, 2.0	7" (17.8 cm)
TWS 84-2609/20	Bayonet	0.5, 0.7, 1.0, 1.5, 2.0	7 3/4" (20.0 cm)
TWS 84-2600/80	Bayonet	0.5, 0.7, 1.0, 1.5, 2.0	8 3/4" (22.0 cm)
TWS 84-2600/81	Bayonet	0.5, 0.7, 1.0, 1.5, 2.0	9 1/2" (24.0 cm)

Convert any regular bipolar forceps listed on pages 04 to 09 into ultra non-stick bipolar forceps, use "UNS" extension for conversion.
Available Tip Size:
0.7mm, 1.0mm, 1.5mm & 2.0mm



TW Surgical

BIPOLAR FORCEPS (ULTRA NON-STICK)



Bayonet (Angled Up)
TWS 84-2609/21

Bayonet (Angled Up)
TWS 84-2609/22

Bayonet (Angled Up)
TWS 84-2609/23

Bayonet (Angled Down)
TWS 84-2609/24


Bayonet (Angled Down)
TWS 84-2609/25

Bayonet (Angled Down)
TWS 84-2609/26

Cat #	Description	Tip (mm)	Size
TWS 84-2609/21	Bayonet (Ang. Up)	0.5, 0.7, 1.0, 1.5, 2.0	6 1/2" (16.5 cm)
TWS 84-2609/22	Bayonet (Ang. Up)	0.5, 0.7, 1.0, 1.5, 2.0	7" (17.8 cm)
TWS 84-2609/23	Bayonet (Ang. Up)	0.5, 0.7, 1.0, 1.5, 2.0	7 3/4" (20.0 cm)
TWS 84-2600/85	Bayonet (Ang. Up)	0.5, 0.7, 1.0, 1.5, 2.0	8 3/4" (22.0 cm)
TWS 84-2600/86	Bayonet (Ang. Up)	0.5, 0.7, 1.0, 1.5, 2.0	9 1/2" (24.0 cm)
TWS 84-2609/24	Bayonet (Ang. Down)	0.5, 0.7, 1.0, 1.5, 2.0	6 1/2" (16.5 cm)
TWS 84-2609/25	Bayonet (Ang. Down)	0.5, 0.7, 1.0, 1.5, 2.0	7" (17.8 cm)
TWS 84-2609/26	Bayonet (Ang. Down)	0.5, 0.7, 1.0, 1.5, 2.0	7 3/4" (20.0 cm)
TWS 84-2600/90	Bayonet (Ang. Down)	0.5, 0.7, 1.0, 1.5, 2.0	8 3/4" (22.0 cm)
TWS 84-2600/91	Bayonet (Ang. Down)	0.5, 0.7, 1.0, 1.5, 2.0	9 1/2" (24.0 cm)

Convert any regular bipolar forceps listed on pages 04 to 09 into ultra non-stick bipolar forceps, use "UNS" extension for conversion. Available Tip Size: 0.7mm, 1.0mm, 1.5mm & 2.0mm

BIPOLAR FORCEPS (HAND-SWITCH)



3 pin connector



Jeweler #4
TWS 84-2609/27



Semkin
TWS 84-2609/28



Jansen Bayonet
TWS 84-2609/29



Scoville-Greenwood
TWS 84-2609/30


Cat #	Description	Tip (mm)	Size
TWS 84-2609/27	Jeweler #4	0.5	4½" (11.5 cm)
TWS 84-2609/28	Semkin	0.75	5½" (14.0 cm)
TWS 84-2609/29	Jansen Bayonet	0.75	8" (20.4 cm)
TWS 84-2609/30	Scoville-Greenwood	1.5	7½" (19.1 cm)

More bipolar forceps hand-switch are available. Use "HS" extension with any regular bipolar forceps for conversion into bipolar forceps hand-switch.



TW Surgical

BIPOLAR FORCEPS (IRRIGATION)



Cushing
TWS 84-2609/31




Gerald Bayonet
TWS 84-2609/32



Jansen Bayonet
TWS 84-2609/33



Scoville-Greenwood
TWS 84-2609/34



Hardy Bayonet
TWS 84-2609/35

Cat #	Description	Tip (mm)	Size
TWS 84-2609/31	Cushing	1.0	6½" (16.5 cm)
TWS 84-2609/32	Gerald Bayonet	1.0	7¼" (18.4 cm)
TWS 84-2609/33	Jansen Bayonet	1.5	8" (20.4 cm)
TWS 84-2609/34	Scoville-Greenwood	1.5	7½" (19.1 cm)
TWS 84-2609/35	Hardy Bayonet	1.5	8½" (21.6 cm)

More bipolar forceps irrigation designs are available. Use "IR" extension with any regular bipolar forceps for conversion into bipolar forceps irrigation.

BIPOLAR FORCEPS (NON-INSULATED)



Jeweler #5
TWS 84-2609/36

Adson
TWS 84-2609/37

Malis
(Straight)
TWS 84-2609/38

Semkin
TWS 84-2609/39

Malis
TWS 84-2609/40

Hardy Bayonet
TWS 84-2609/41

Cat #	Description	Tip (mm)	Size
TWS 84-2609/36	Jeweler #5	Micro, 0.5	4½" (11.5 cm)
TWS 84-2609/38	Malis (Str.)	Micro, 0.5	6" (15.2 cm)
TWS 84-2609/37	Adson	1.0	4¾" (12.1 cm)
TWS 84-2609/39	Semkin	0.5, 0.7	5½" (14 cm)
TWS 84-2609/40	Malis	0.5, 1.0, 2.0	7" (17.8 cm)
TWS 84-2609/41	Hardy Bayonet	0.5, 1.0, 2.0	8½" (21.6 cm)

Convert any regular bipolar forceps listed on pages 04 to 09 into non-insulated bipolar forceps, use 120-series for conversion into non-insulated bipolar forceps.



TW Surgical

MONOPOLAR FORCEPS



Jeweler #4
TWS 84-2609/42



Jeweler #7
TWS 84-2609/43



Adson
TWS 84-2609/44



Potts-Smith
TWS 84-2609/45



Potts-Smith
TWS 84-2609/46



Malis
TWS 84-2609/47



Jansen Bayonet
TWS 84-2609/48



Jansen Bayonet
(Angled Up)
TWS 84-2601/18


Cat #	Description	Tip (mm)	Size
TWS 84-2609/42	Jeweler #4	Micro	4½" (11.5 cm)
TWS 84-2609/43	Jeweler #7	Micro	4½" (11.5 cm)
TWS 84-2609/44	Adson	1.0	4¾" (12.1 cm)
TWS 84-2609/45	Potts-Smith	2.0	7" (17.8 cm)
TWS 84-2609/46	Potts-Smith	2.0	8¼" (21.0 cm)
TWS 84-2601/12	Potts-Smith	2.0	9½" (24.1 cm)
TWS 84-2609/47	Malis	1.0, 1.5	7" (17.8 cm)

Cat #	Description	Tip (mm)	Size
TWS 84-2609/48	Jansen Bayonet	1.5	5¼" (13.3 cm)
TWS 84-2609/49	Jansen Bayonet	1.5	6¼" (15.9 cm)
TWS 84-2601/16	Jansen Bayonet	2.0	7½" (19.1 cm)
TWS 84-2601/17	Jansen Bayonet	2.0	8¼" (21.0 cm)
TWS 84-2601/18	Jansen Bayonet (Ang. Up)	2.0	7½" (19.1 cm)
TWS 84-2601/19	Jansen Bayonet (Ang. Up)	2.0	8¼" (21.0 cm)



TW Surgical

MONOPOLAR FORCEPS (HAND-SWITCH)



Jansen Bayonet
TWS 84-2609/50



Jansen Bayonet
TWS 84-2609/51



Potts-Smith
TWS 84-2609/52




Potts-Smith
TWS 84-2609/53

Cat #	Description	Tip (mm)	Size
TWS 84-2609/50	Jansen Bayonet	2.0	5½" (14.0 cm)
TWS 84-2609/51	Jansen Bayonet	2.0	7½" (19.0 cm)
TWS 84-2601/22	Jansen Bayonet	2.0	8¼" (21.0 cm)
TWS 84-2609/52	Potts-Smith	2.0	7" (17.8 cm)
TWS 84-2609/53	Potts-Smith	2.0	8¼" (21.0 cm)
TWS 84-2601/25	Potts-Smith	2.0	9½" (24.0 cm)



TW Surgical

MONOPOLAR FORCEPS (HAND-SWITCH NON-STICK)



Jansen Bayonet
TWS 84-2609/54



Jansen Bayonet
TWS 84-2609/55



Potts-Smith
TWS 84-2609/56




Potts-Smith
TWS 84-2609/57

Cat #	Description	Tip (mm)	Size
TWS 84-2609/54	Jansen Bayonet	2.0	5½" (14.0 cm)
TWS 84-2609/55	Jansen Bayonet	2.0	7½" (19.0 cm)
TWS 84-2601/28	Jansen Bayonet	2.0	8¼" (21.0 cm)
TWS 84-2609/56	Potts-Smith	2.0	7" (17.8 cm)
TWS 84-2609/57	Potts-Smith	2.0	8¼" (21.0 cm)
TWS 84-2601/31	Potts-Smith	2.0	9½" (24.0 cm)




TW Surgical

EUROPEAN BIPOLAR FORCEPS



Available Tip Widths (Always Specify Tips Width)



McPherson
TWS 84-2609/58



McPherson
(Angled)
TWS 84-2609/59



Jeweler #5
TWS 84-2609/60



Jeweler #5
(Angled)
TWS 84-2609/61



Jeweler #7
TWS 84-2609/62



Adson
TWS 84-2609/63

Cat #	Description	Tip (mm)	Size
TWS 84-2609/58	McPherson	0.5	4" (10.2 cm)
TWS 84-2609/59	McPherson (Ang.)	0.5	4" (10.2 cm)
TWS 84-2609/60	Jeweler #5	0.5, 0.7	4½" (11.5 cm)
TWS 84-2609/61	Jeweler #5 (Ang.)	0.5, 0.7	4½" (11.5 cm)
TWS 84-2609/62	Jeweler #7	0.5, 0.7	4½" (11.5 cm)
TWS 84-2609/63	Adson	1.0	4¾" (12.1 cm)
TWS 84-2601/38	Adson	1.0	6" (15.0 cm)

For your convenience add the following tip extension to the product code.
Needle Tip (-00), 0.5mm(-05), 1.0mm(-10) & 2.0mm (-20)



TW Surgical

EUROPEAN BIPOLAR FORCEPS



Straight
TWS 84-2609/64



Straight
TWS 84-2609/65



Straight
TWS 84-2609/66



Straight
TWS 84-2609/67

Cat #	Description	Tip (mm)	Size
TWS 84-2609/64	Straight	0.5, 1.0, 2.0	6" (15.0 cm)
TWS 84-2609/65	Straight	0.5, 1.0, 2.0	7" (17.8 cm)
TWS 84-2609/66	Straight	0.5, 1.0, 2.0	7¾" (20.0 cm)
TWS 84-2609/67	Straight	0.5, 1.0, 2.0	8¾" (22.0 cm)
TWS 84-2601/43	Straight	0.5, 1.0, 2.0	9½" (24.0 cm)

For your convenience add the following tip extension to the product code.
Needle Tip (-00), 0.5mm(-05), 1.0mm(-10) & 2.0mm (-20)



TW Surgical

EUROPEAN BIPOLAR FORCEPS



Straight (Angled)
TWS 84-2609/68



Straight (Angled)
TWS 84-2609/69



Straight (Angled)
TWS 84-2609/70



Straight (Angled)
TWS 84-2609/71

Cat #	Description	Tip (mm)	Size
TWS 84-2609/68	Straight (Ang.)	0.5, 1.0, 2.0	6" (15.0 cm)
TWS 84-2609/69	Straight (Ang.)	0.5, 1.0, 2.0	7" (17.8 cm)
TWS 84-2609/70	Straight (Ang.)	0.5, 1.0, 2.0	7¾" (20.0 cm)
TWS 84-2609/71	Straight (Ang.)	0.5, 1.0, 2.0	8¾" (22.0 cm)
TWS 84-2601/48	Straight (Ang.)	0.5, 1.0, 2.0	9½" (24.0 cm)

For your convenience add the following tip extension to the product code.
Needle Tip (-00), 0.5mm(-05), 1.0mm(-10) & 2.0mm (-20)




Bayonet
TWS 84-2609/72



Bayonet
TWS 84-2609/73



Bayonet
TWS 84-2609/74



Bayonet
TWS 84-2609/75

Cat #	Description	Tip (mm)	Size
TWS 84-2609/72	Bayonet	0.5, 1.0, 2.0	6½" (16.5 cm)
TWS 84-2609/73	Bayonet	0.5, 1.0, 2.0	7" (17.8 cm)
TWS 84-2609/74	Bayonet	0.5, 1.0, 2.0	7¾" (20.0 cm)
TWS 84-2609/75	Bayonet	0.5, 1.0, 2.0	8¾" (22.0 cm)
TWS 84-2601/53	Bayonet	0.5, 1.0, 2.0	9½" (24.0 cm)

For your convenience add the following tip extension to the product code.
Needle Tip (-00), 0.5mm(-05), 1.0mm(-10) & 2.0mm (-20)



TW Surgical

EUROPEAN BIPOLAR FORCEPS




Bayonet
(Angled Up)
TWS 84-2609/76



Bayonet
(Angled Up)
TWS 84-2609/77



Bayonet
(Angled Up)
TWS 84-2609/78



Bayonet
(Angled Up)
TWS 84-2609/79

Cat #	Description	Tip (mm)	Size
TWS 84-2609/76	Bayonet (Ang. Up)	0.5, 1.0, 2.0	6½" (16.5 cm)
TWS 84-2609/77	Bayonet (Ang. Up)	0.5, 1.0, 2.0	7" (17.8 cm)
TWS 84-2609/78	Bayonet (Ang. Up)	0.5, 1.0, 2.0	7¾" (20.0 cm)
TWS 84-2609/79	Bayonet (Ang. Up)	0.5, 1.0, 2.0	8¾" (22.0 cm)
TWS 84-2601/58	Bayonet (Ang. Up)	0.5, 1.0, 2.0	9½" (24.0 cm)

For your convenience add the following tip extension to the product code.
Needle Tip (-00), 0.5mm(-05), 1.0mm(-10) & 2.0mm (-20)




Bayonet
(Angled Down)
TWS 84-2609/80



Bayonet
(Angled Down)
TWS 84-2609/81



Bayonet
(Angled Down)
TWS 84-2609/82



Bayonet
(Angled Down)
TWS 84-2609/83

Cat #	Description	Tip (mm)	Size
TWS 84-2609/80	Bayonet (Ang. Down)	0.5, 1.0, 2.0	6½" (16.5 cm)
TWS 84-2609/81	Bayonet (Ang. Down)	0.5, 1.0, 2.0	7" (17.8 cm)
TWS 84-2609/82	Bayonet (Ang. Down)	0.5, 1.0, 2.0	7¾" (20.0 cm)
TWS 84-2609/83	Bayonet (Ang. Down)	0.5, 1.0, 2.0	8¾" (22.0 cm)
TWS 84-2601/83	Bayonet (Ang. Down)	0.5, 1.0, 2.0	9½" (24.0 cm)

For your convenience add the following tip extension to the product code.
 Needle Tip (-00), 0.5mm(-05), 1.0mm(-10) & 2.0mm (-20)



TW Surgical

EUROPEAN BIPOLAR FORCEPS (NON-STICK)



McPherson



McPherson (Angled)



Jeweler #5



Jeweler #5 (Angled)



Adson

Cat #	Description	Tip (mm)	Size
TWS 84-2601/64	McPherson	0.5	4" (10.2 cm)
TWS 84-2601/65	McPherson (Ang.)	0.5	4" (10.2 cm)
TWS 84-2601/66	Jeweler #5	0.5, 0.7	4½" (11.5 cm)
TWS 84-2601/67	Jeweler #5 (Ang.)	0.5, 0.7	4½" (11.5 cm)
TWS 84-2601/68	Adson	1.0	4¾" (12.1 cm)

Convert any European bipolar forceps listed on pages 25 to 30 into non-stick European bipolar forceps, use "NS" extension for conversion.
Available Tip Size:
0.7mm, 1.0mm, 1.5mm & 2.0mm

EUROPEAN BIPOLAR FORCEPS (NON-STICK)



Straight
TWS 84-2609/89



Straight
TWS 84-2609/90



Straight
TWS 84-2609/91



Straight
TWS 84-2609/92

Cat #	Description	Tip (mm)	Size
TWS 84-2609/89	Straight	0.5, 1.0, 2.0	6" (15.0 cm)
TWS 84-2609/90	Straight	0.5, 1.0, 2.0	7" (17.8 cm)
TWS 84-2609/91	Straight	0.5, 1.0, 2.0	7 3/4" (20.0 cm)
TWS 84-2609/92	Straight	0.5, 1.0, 2.0	8 3/4" (22.0 cm)
TWS 84-2601/73	Straight	0.5, 1.0, 2.0	9 1/2" (24.0 cm)

Convert any European bipolar forceps listed on pages 25 to 30 into non-stick European bipolar forceps, use "NS" extension for conversion.
Available Tip Size:
0.7mm, 1.0mm, 1.5mm & 2.0mm



TW Surgical

EUROPEAN BIPOLAR FORCEPS (NON-STICK)



Bayonet
TWS 84-2609/93



Bayonet
TWS 84-2609/94



Bayonet
TWS 84-2609/95



Bayonet
TWS 84-2609/96

Cat #	Description	Tip (mm)	Size
TWS 84-2609/93	Bayonet	0.5, 1.0, 2.0	6½" (16.5 cm)
TWS 84-2609/94	Bayonet	0.5, 1.0, 2.0	7" (17.8 cm)
TWS 84-2609/95	Bayonet	0.5, 1.0, 2.0	7¾" (20.0 cm)
TWS 84-2609/96	Bayonet	0.5, 1.0, 2.0	8¾" (22.0 cm)
TWS 84-2601/78	Bayonet	0.5, 1.0, 2.0	9½" (24.0 cm)

Convert any European bipolar forceps listed on pages 25 to 30 into non-stick European bipolar forceps, use "NS" extension for conversion.
Available Tip Size:
0.7mm, 1.0mm, 1.5mm & 2.0mm



McPherson



McPherson (Angled)



Jeweler #5



Jeweler #5 (Angled)



Adson

Cat #	Description	Tip (mm)	Size
TWS 84-2601/79	McPherson	0.5	4" (10.2 cm)
TWS 84-2601/80	McPherson (Ang.)	0.5	4" (10.2 cm)
TWS 84-2601/81	Jeweler #5	0.5, 0.7	4½" (11.5 cm)
TWS 84-2601/82	Jeweler #5 (Ang.)	0.5, 0.7	4½" (11.5 cm)
TWS 84-2601/83	Adson	1.0	4¾" (12.1 cm)

Convert any European bipolar forceps listed on pages 25 to 30 into ultra non-stick European bipolar forceps, use "UNS" extension for conversion.
Available Tip Size:
0.7mm, 1.0mm, 1.5mm & 2.0mm

EUROPEAN BIPOLAR FORCEPS (ULTRA NON-STICK)



Straight
TWS 84-2610/01



Straight
TWS 84-2610/02



Straight
TWS 84-2610/03




Straight
TWS 84-2610/04

Cat #	Description	Tip (mm)	Size
TWS 84-2610/01	Straight	0.5, 1.0, 2.0	6" (15.0 cm)
TWS 84-2610/02	Straight	0.5, 1.0, 2.0	7" (17.8 cm)
TWS 84-2610/03	Straight	0.5, 1.0, 2.0	7¾" (20.0 cm)
TWS 84-2610/04	Straight	0.5, 1.0, 2.0	8½" (22.0 cm)
TWS 84-2601/88	Straight	0.5, 1.0, 2.0	9½" (24.0 cm)

Convert any European bipolar forceps listed on pages 25 to 30 into ultra non-stick European bipolar forceps, use "UNS" extension for conversion.
Available Tip Size:
0.7mm, 1.0mm, 1.5mm & 2.0mm


EUROPEAN BIPOLAR FORCEPS (ULTRA NON-STICK)




Bayonet
TWS 84-2610/05



Bayonet
TWS 84-2610/06



Bayonet
TWS 84-2610/07



Bayonet
TWS 84-2610/08

Cat #	Description	Tip (mm)	Size
TWS 84-2610/05	Bayonet	0.5, 1.0, 2.0	6½" (16.5 cm)
TWS 84-2610/06	Bayonet	0.5, 1.0, 2.0	7" (17.8 cm)
TWS 84-2610/07	Bayonet	0.5, 1.0, 2.0	7¾" (20.0 cm)
TWS 84-2610/08	Bayonet	0.5, 1.0, 2.0	8¾" (22.0 cm)
TWS 84-2601/93	Bayonet	0.5, 1.0, 2.0	9½" (24.0 cm)

Convert any European bipolar forceps listed on pages 25 to 30 into ultra non-stick European bipolar forceps, use "UNS" extension for conversion.
Available Tip Size:
0.7mm, 1.0mm, 1.5mm & 2.0mm



TW Surgical

EUROPEAN BIPOLAR FORCEPS (ULTRA NON-STICK)



Jeweler #5



Jeweler #5 (Angled)



Straight



Straight



Straight

Cat #	Description	Tip (mm)	Size
TWS 84-2601/94	Jeweler #5	0.5, 0.7	4 1/2" (11.5 cm)
TWS 84-2601/95	Jeweler #5 (Ang.)	0.5, 0.7	4 1/2" (11.5 cm)
TWS 84-2601/96	Straight	0.5, 0.7, 1.0, 1.5, 2.0	7" (17.8 cm)
TWS 84-2601/97	Straight	0.5, 0.7, 1.0, 1.5, 2.0	7 3/4" (20.0 cm)
TWS 84-2601/98	Straight	0.5, 0.7, 1.0, 1.5, 2.0	8 3/4" (22.0 cm)
TWS 84-2601/99	Straight	0.5, 0.7, 1.0, 1.5, 2.0	9 1/2" (24.0 cm)

For your convenience add the following tip extension to the product code.
Needle Tip (-00), 0.5mm(-05), 1.0mm(-10) & 2.0mm (-20)



TW Surgical

EUROPEAN BIPOLAR FORCEPS (ULTRA NON-STICK)



Straight (Angled)
TWS 84-2610/14

Straight (Angled)
TWS 84-2610/15

Straight (Angled)
TWS 84-2610/16

Bayonet
TWS 84-2610/17

Bayonet
TWS 84-2610/18

Bayonet
TWS 84-2610/19

Cat #	Description	Tip (mm)	Size
TWS 84-2610/14	Straight (Ang.)	0.5, 0.7, 1.0, 1.5, 2.0	7" (17.8 cm)
TWS 84-2610/15	Straight (Ang.)	0.5, 0.7, 1.0, 1.5, 2.0	7 3/4" (20.0 cm)
TWS 84-2610/16	Straight (Ang.)	0.5, 0.7, 1.0, 1.5, 2.0	8 3/4" (22.0 cm)
TWS 84-2602/03	Straight (Ang.)	0.5, 0.7, 1.0, 1.5, 2.0	9 1/2" (24.0 cm)
TWS 84-2610/17	Bayonet	0.5, 0.7, 1.0, 1.5, 2.0	6 1/2" (16.5 cm)
TWS 84-2610/18	Bayonet	0.5, 0.7, 1.0, 1.5, 2.0	7" (17.8 cm)
TWS 84-2610/19	Bayonet	0.5, 0.7, 1.0, 1.5, 2.0	7 3/4" (20.0 cm)
TWS 84-2602/07	Bayonet	0.5, 0.7, 1.0, 1.5, 2.0	8 3/4" (22.0 cm)
TWS 84-2602/08	Bayonet	0.5, 0.7, 1.0, 1.5, 2.0	9 1/2" (24.0 cm)

For your convenience add the following tip extension to the product code.
Needle Tip (-00), 0.5mm(-05), 1.0mm(-10) & 2.0mm (-20)

EUROPEAN BIPOLAR FORCEPS (ULTRA NON-STICK)



Bayonet
(Angled Up)
TWS 84-2610/20

Bayonet
(Angled Up)
TWS 84-2610/21

Bayonet
(Angled Up)
TWS 84-2610/22

Bayonet
(Angled Down)
TWS 84-2610/23

Bayonet
(Angled Down)
TWS 84-2610/24

Bayonet
(Angled Down)
TWS 84-2610/25

Cat #	Description	Tip (mm)	Size
TWS 84-2610/20	Bayonet (Ang. Up)	0.5, 0.7, 1.0, 1.5, 2.0	6 ½" (16.5 cm)
TWS 84-2610/21	Bayonet (Ang. Up)	0.5, 0.7, 1.0, 1.5, 2.0	7" (17.8 cm)
TWS 84-2610/22	Bayonet (Ang. Up)	0.5, 0.7, 1.0, 1.5, 2.0	7 ¾" (20.0 cm)
TWS 84-2602/12	Bayonet (Ang. Up)	0.5, 0.7, 1.0, 1.5, 2.0	8 ¾" (22.0 cm)
TWS 84-2602/13	Bayonet (Ang. Up)	0.5, 0.7, 1.0, 1.5, 2.0	9 ½" (24.0 cm)
TWS 84-2610/23	Bayonet (Ang. Down)	0.5, 0.7, 1.0, 1.5, 2.0	6 ½" (16.5 cm)
TWS 84-2610/24	Bayonet (Ang. Down)	0.5, 0.7, 1.0, 1.5, 2.0	7" (17.8 cm)
TWS 84-2610/25	Bayonet (Ang. Down)	0.5, 0.7, 1.0, 1.5, 2.0	7 ¾" (20.0 cm)
TWS 84-2602/17	Bayonet (Ang. Down)	0.5, 0.7, 1.0, 1.5, 2.0	8 ¾" (22.0 cm)
TWS 84-2602/18	Bayonet (Ang. Down)	0.5, 0.7, 1.0, 1.5, 2.0	9 ½" (24.0 cm)

For your convenience add the following tip extension to the product code.
Needle Tip (-00), 0.5mm(-05), 1.0mm(-10) & 2.0mm (-20)



TW Surgical

EUROPEAN BIPOLAR FORCEPS (IRRIGATION)



Straight



Straight



Bayonet



Bayonet

Cat #	Description	Tip (mm)	Size
TWS 84-2620/1	Straight	0.5, 1.0, 2.0	6" (15.0 cm)
TWS 84-2620/2	Straight	0.5, 1.0, 2.0	8 1/2" (22.0 cm)
TWS 84-2620/3	Bayonet	0.5, 1.0, 2.0	6 1/2" (16.5 cm)
TWS 84-2620/4	Bayonet	0.5, 1.0, 2.0	7 3/4" (20.0 cm)

Convert any European bipolar forceps listed on pages 25 to 30 into European bipolar forceps irrigation. Use "IR" extension for conversion.



TW Surgical

EUROPEAN BIPOLAR FORCEPS (SUCTION)



Forceps (Straight)
TWS 84-2610/30



Forceps (Angled 30°)
TWS 84-2610/31



Forceps (Angled 45°)
TWS 84-2610/32



Forceps (Angled 60°)
TWS 84-2610/33

Cat #	Description	Tip (mm)	Size
TWS 84-2610/30	Forceps (Str.)	3.0	8¼" (21.0 cm)
TWS 84-2610/31	Forceps (Ang. 30°)	3.0	8¼" (21.0 cm)
TWS 84-2610/32	Forceps (Ang. 45°)	3.0	8¼" (21.0 cm)
TWS 84-2610/33	Forceps (Ang. 60°)	3.0	8¼" (21.0 cm)



TW Surgical

EUROPEAN BIPOLAR FORCEPS (TUNGSTEN CARBIDE)



Straight



Straight
(Angled)



Bayonet



Bayonet
(Angled Up)



Bayonet
(Angled Down)

Cat #	Description	Tip (mm)	Size
TWS 84-2602/23	Straight	2.0	7¼" (20.0 cm)
TWS 84-2602/24	Straight (Ang.)	2.0	7¼" (20.0 cm)
TWS 84-2602/25	Bayonet	2.0	7¼" (20.0 cm)
TWS 84-2602/26	Bayonet (Ang. Up)	2.0	7¼" (20.0 cm)
TWS 84-2602/27	Bayonet (Ang. Down)	2.0	7¼" (20.0 cm)


Treasure Wright (U.K.) Ltd.

128-160 City Road, London, ECLV 2NX, United Kingdom
info@twrightmedical.com



TW Surgical

EUROPEAN BIPOLAR FORCEPS



Straight
TWS 84-2610/39



Straight
TWS 84-2610/40



Straight
(Angled)
TWS 84-2610/41



Straight
(Angled)
TWS 84-2610/42

Cat #	Description	Tip (mm)	Size
TWS 84-2610/39	Straight	0.5, 1.0, 2.0	6" (15.0 cm)
TWS 84-2610/40	Straight	0.5, 1.0, 2.0	7" (17.8 cm)
TWS 84-2602/30	Straight	0.5, 1.0, 2.0	7½" (20.0 cm)
TWS 84-2602/31	Straight	0.5, 1.0, 2.0	8¼" (22.0 cm)
TWS 84-2602/32	Straight	0.5, 1.0, 2.0	9½" (24.0 cm)

Cat #	Description	Tip (mm)	Size
TWS 84-2610/41	Straight (Ang.)	0.5, 1.0, 2.0	6" (15.0 cm)
TWS 84-2610/42	Straight (Ang.)	0.5, 1.0, 2.0	7" (17.8 cm)
TWS 84-2602/35	Straight (Ang.)	0.5, 1.0, 2.0	7¾" (20.0 cm)
TWS 84-2602/36	Straight (Ang.)	0.5, 1.0, 2.0	8¾" (22.0 cm)
TWS 84-2602/37	Straight (Ang.)	0.5, 1.0, 2.0	9½" (24.0 cm)

For your convenience add the following tip extension to the product code.
Needle Tip (-00), 0.5mm(-05),
1.0mm(-10) & 2.0mm (-20)



TW Surgical


EUROPEAN BIPOLAR FORCEPS




Bayonet
TWS 84-2610/43



Bayonet
TWS 84-2610/44



Bayonet
(Angled Up)
TWS 84-2610/45



Bayonet
(Angled Down)
TWS 84-2610/46

Cat #	Description	Tip (mm)	Size
TWS 84-2610/43	Bayonet	0.5, 1.0, 2.0	5½" (14.0 cm)
TWS 84-2610/44	Bayonet	0.5, 1.0, 2.0	7" (17.8 cm)
TWS 84-2602/40	Bayonet	0.5, 1.0, 2.0	7¾" (20.0 cm)
TWS 84-2602/41	Bayonet	0.5, 1.0, 2.0	8¾" (22.0 cm)
TWS 84-2602/42	Bayonet	0.5, 1.0, 2.0	9½" (24.0 cm)

Cat #	Description	Tip (mm)	Size
TWS 84-2610/45	Bayonet (Ang. Up)	0.5, 1.0, 2.0	7¾" (20.0 cm)
TWS 84-2602/44	Bayonet (Ang. Up)	0.5, 1.0, 2.0	8¾" (22.0 cm)
TWS 84-2610/46	Bayonet (Ang. Down)	0.5, 1.0, 2.0	7¾" (20.0 cm)
TWS 84-2602/46	Bayonet (Ang. Down)	0.5, 1.0, 2.0	8¾" (22.0 cm)

For your convenience add the following tip extension to the product code. Needle Tip (-00), 0.5mm(-05), 1.0mm(-10) & 2.0mm (-20)

BIPOLAR FORCEPS SINGLE-USE (MOLDED)



Coaptation
TWS 84-2610/47



Jeweler #5
TWS 84-2610/48



Jeweler #5
(Angled)
TWS 84-2610/49



Cushing
TWS 84-2610/50




Jansen Bayonet
TWS 84-2610/51

Cat #	Description	Tip (mm)	Size
TWS 84-2610/47	Coaptation	0.5	4" (10.2 cm)
TWS 84-2610/48	Jeweler #5	0.5	4" (10.2 cm)
TWS 84-2610/49	Jeweler #5 (Ang.)	0.5	4" (10.2 cm)
TWS 84-2610/50	Cushing	1.0	8" (20.4 cm)
TWS 84-2610/51	Jansen Bayonet	1.0	8" (20.4 cm)



TW Surgical

BIPOLAR FORCEPS SINGLE-USE



Jeweler #5
TWS 84-2610/52



Adson
TWS 84-2610/53



Semkin
(Straight)
TWS 84-2610/54



Cushing
(Straight)
TWS 84-2610/55



Gerald
(Straight)
TWS 84-2610/56

Cat #	Description	Tip (mm)	Size
TWS 84-2610/52	Jeweler #5	Micro, 0.5	4½" (11.5 cm)
TWS 84-2610/53	Adson	1.0	4¾" (12.1 cm)
TWS 84-2610/54	Semkin (Str.)	0.5	5½" (14.0 cm)
TWS 84-2610/55	Cushing (Str.)	1.0, 1.5	6½" (16.5 cm)
TWS 84-2610/56	Gerald (Str.)	1.0, 1.5	7" (17.8 cm)

Forceps are supplied with 3m cable.
Available in Sterile/Non-Sterile packing.
Sterilization by Ethylene Oxide (EO).



TW Surgical

BIPOLAR FORCEPS SINGLE-USE



Gerald Bayonet
TWS 84-2610/57



Jansen Bayonet (Irrigation)
TWS 84-2610/58



Jansen Bayonet
TWS 84-2610/59



Malis
TWS 84-2610/60



Hardy Bayonet (Irrigation)
TWS 84-2610/61


Cat #	Description	Tip (mm)	Size
TWS 84-2610/57	Gerald Bayonet	1.0, 1.5	7¼" (18.4 cm)
TWS 84-2610/58	Jansen Bayonet (Irrigation)	0.5, 0.7, 1.0, 1.5, 2.0	7" (17.8 cm)
TWS 84-2610/59	Jansen Bayonet	1.0, 1.5, 2.0	8" (20.4 cm)
TWS 84-2610/60	Malis	1.0, 1.5, 2.0	7" (17.8 cm)
TWS 84-2610/61	Hardy Bayonet (Irrigation)	1.0	8½" (21.6 cm)

Forceps are supplied with 3m cable.
Available in Sterile/Non-Sterile packing.
Sterilization by Ethylene Oxide (EO).



TW Surgical

BIPOLAR FORCEPS SINGLE-USE (NON-STICK)



Jeweler #5
TWS 84-2610/62



Adson
TWS 84-2610/63



Semkin
(Straight)
TWS 84-2610/64



Cushing
(Straight)
TWS 84-2610/65



Gerald
(Straight)
TWS 84-2610/66

Cat #	Description	Tip (mm)	Size
TWS 84-2610/62	Jeweler #5	Micro, 0.5	4½" (11.5 cm)
TWS 84-2610/63	Adson	1.0	4¾" (12.1 cm)
TWS 84-2610/64	Semkin (Str.)	0.5	5½" (14.0 cm)
TWS 84-2610/65	Cushing (Str.)	1.0, 1.5	6½" (16.5 cm)
TWS 84-2610/66	Gerald (Str.)	1.0, 1.5	7" (17.8 cm)

Forceps are supplied with 3m cable.
Available in Sterile/Non-Sterile packing.
Sterilization by Ethylene Oxide (EO).

BIPOLAR FORCEPS SINGLE-USE (NON-STICK)



Gerald Bayonet
TWS 84-2610/67



Jansen Bayonet (Irrigation)
TWS 84-2610/68



Jansen Bayonet
TWS 84-2610/69




Malis
TWS 84-2610/70



Hardy Bayonet (Irrigation)
TWS 84-2610/71

Cat #	Description	Tip (mm)	Size
TWS 84-2610/67	Gerald Bayonet	1.0, 1.5	7¼" (18.4 cm)
TWS 84-2610/68	Jansen Bayonet (Irrigation)	0.5, 0.7, 1.0, 1.5, 2.0	7" (17.8 cm)
TWS 84-2610/69	Jansen Bayonet	1.0, 1.5, 2.0	8" (20.4 cm)
TWS 84-2610/70	Malis	1.0, 1.5, 2.0	7" (17.8 cm)
TWS 84-2610/71	Hardy Bayonet (Irrigation)	1.0	8½" (21.6 cm)

Forceps are supplied with 3m cable.
Available in Sterile/Non-Sterile packing.
Sterilization by Ethylene Oxide (EO).



Jeweler #5
TWS 84-2610/72



Adson
TWS 84-2610/73



Semkin
(Straight)
TWS 84-2610/74



Cushing
(Straight)
TWS 84-2610/75



Gerald
(Straight)
TWS 84-2610/76

Cat #	Description	Tip (mm)	Size
TWS 84-2610/72	Jeweler #5	Micro, 0.5	4½" (11.5 cm)
TWS 84-2610/73	Adson	1.0	4¾" (12.1 cm)
TWS 84-2610/74	Semkin (Str.)	0.5	5½" (14.0 cm)
TWS 84-2610/75	Cushing (Str.)	1.0, 1.5	6½" (16.5 cm)
TWS 84-2610/76	Gerald (Str.)	1.0, 1.5	7" (17.8 cm)

Forceps are supplied with 3m cable.
Available in Sterile/Non-Sterile packing.
Sterilization by Ethylene Oxide (EO).

BIPOLAR FORCEPS SINGLE-USE (ULTRA NON-STICK)



Gerald Bayonet
TWS 84-2610/77



Jansen Bayonet
(Irrigation)
TWS 84-2610/78



Jansen Bayonet
TWS 84-2610/79




Malis
TWS 84-2610/80



Hardy Bayonet
(Irrigation)
TWS 84-2610/81

Cat #	Description	Tip (mm)	Size
TWS 84-2610/77	Gerald Bayonet	1.0, 1.5	7¼" (18.4 cm)
TWS 84-2610/78	Jansen Bayonet (Irrigation)	0.5, 0.7, 1.0, 1.5, 2.0	7" (17.8 cm)
TWS 84-2610/79	Jansen Bayonet	1.0, 1.5, 2.0	8" (20.4 cm)
TWS 84-2610/80	Malis	1.0, 1.5, 2.0	7" (17.8 cm)
TWS 84-2610/81	Hardy Bayonet (Irrigation)	1.0	8½" (21.6 cm)

Forceps are supplied with 3m cable.
Available in Sterile/Non-Sterile packing.
Sterilization by Ethylene Oxide (EO).



McPherson
(Straight)
TWS 84-2610/82

Jeweler #5
TWS 84-2610/83

Adson
TWS 84-2610/84

Semkin
TWS 84-2610/85

Jansen Bayonet
TWS 84-2610/86

Malis
TWS 84-2610/87


Cat #	Description	Tip (mm)	Size
TWS 84-2610/82	McPherson (Str.)	0.5	3½" (9.0 cm)
TWS 84-2610/83	Jeweler #5	Needle, 0.5	4½" (11.5 cm)
TWS 84-2610/84	Adson	0.5, 1.0	4¾" (12.1 cm)
TWS 84-2610/85	Semkin	0.5, 0.7, 1.0	5½" (14.0 cm)
TWS 84-2610/86	Jansen Bayonet	0.5, 1.0, 1.5, 2.0	7" (17.8 cm)
TWS 84-2610/87	Malis	0.5, 1.0, 1.5	7" (17.8 cm)

EUROPEAN BIPOLAR FORCEPS SINGLE-USE

Series supplied with
Banana plug



Series supplied with
Valleylab plug



McPherson
(Angled)
TWS 84-2610/88

Jeweler #5
TWS 84-2610/89

Jeweler #5
(Angled)
TWS 84-2610/90

Straight
TWS 84-2610/91

Bayonet
TWS 84-2610/92

Bayonet
TWS 84-2610/93


Cat #	Description	Tip (mm)	Size
TWS 84-2610/88	McPherson (Ang.)	0.5	4" (10.2 cm)
TWS 84-2610/89	Jeweler #5	0.5, 0.7	4½" (11.5 cm)
TWS 84-2610/90	Jeweler #5 (Ang.)	0.5, 0.7	4½" (11.5 cm)
TWS 84-2610/91	Straight	0.5, 1.0, 1.5	7¾" (20.0 cm)
TWS 84-2610/92	Bayonet	0.5, 1.0, 1.5	7" (17.8 cm)
TWS 84-2610/93	Bayonet	0.5, 1.0, 1.5	7¾" (20.0 cm)

Forceps are supplied with 3m cable.
Available in Sterile/Non-Sterile packing.
Sterilization by Ethylene Oxide (EO).


EUROPEAN BIPOLAR FORCEPS

SINGLE-USE (NON-STICK)

Series supplied with
Banana plug



Series supplied with
Valleylab plug



McPherson
(Angled)
TWS 84-2610/94

Jeweler #5
TWS 84-2610/95

Jeweler #5
(Angled)
TWS 84-2610/96

Straight
TWS 84-2610/97

Bayonet
TWS 84-2610/98

Bayonet
TWS 84-2610/99


Cat #	Description	Tip (mm)	Size
TWS 84-2610/94	McPherson (Ang.)	0.5	4" (10.2 cm)
TWS 84-2610/95	Jeweler #5	0.5, 0.7	4½" (11.5 cm)
TWS 84-2610/96	Jeweler #5 (Ang.)	0.5, 0.7	4½" (11.5 cm)
TWS 84-2610/97	Straight	0.5, 1.0, 1.5	7¾" (20.0 cm)
TWS 84-2610/98	Bayonet	0.5, 1.0, 1.5	7" (17.8 cm)
TWS 84-2610/99	Bayonet	0.5, 1.0, 1.5	7¾" (20.0 cm)

Forceps are supplied with 3m cable.
Available in Sterile/Non-Sterile packing.
Sterilization by Ethylene Oxide (EO).


EUROPEAN BIPOLAR FORCEPS SINGLE-USE (ULTRA NON-STICK)



Series supplied with
Banana plug



Series supplied with
Valleylab plug



McPherson
(Angled)
TWS 84-2611/00

Jeweler #5
TWS 84-2611/01

Jeweler #5
(Angled)
TWS 84-2611/02

Straight
TWS 84-2611/03

Bayonet
TWS 84-2611/04

Bayonet
TWS 84-2611/05

Cat #	Description	Tip (mm)	Size
TWS 84-2611/00	McPherson (Ang.)	0.5	4" (10.2 cm)
TWS 84-2611/01	Jeweler #5	0.5, 0.7	4½" (11.5 cm)
TWS 84-2611/02	Jeweler #5 (Ang.)	0.5, 0.7	4½" (11.5 cm)
TWS 84-2611/03	Straight	0.5, 1.0, 1.5	7¾" (20.0 cm)
TWS 84-2611/04	Bayonet	0.5, 1.0, 1.5	7" (17.8 cm)
TWS 84-2611/05	Bayonet	0.5, 1.0, 1.5	7¾" (20.0 cm)

Forceps are supplied with 3m cable.
Available in Sterile/Non-Sterile packing.
Sterilization by Ethylene Oxide (EO).



Jeweler #5
TWS 84-2611/06



Jeweler #5
(Angled)
TWS 84-2611/07



Adson
TWS 84-2611/08



Straight
TWS 84-2611/09



Straight
(Angled)
TWS 84-2611/10

Cat #	Description	Tip (mm)	Size
TWS 84-2611/06	Jeweler #5	0.5, 0.7	4½" (11.5 cm)
TWS 84-2611/07	Jeweler #5 (Ang.)	0.5, 0.7	4½" (11.5 cm)
TWS 84-2611/08	Adson	1.0	4¾" (12.1 cm)
TWS 84-2611/09	Straight	1.0, 1.5, 2.0	7¾" (20.0 cm)
TWS 84-2611/10	Straight (Ang.)	1.0, 1.5, 2.0	7¾" (20.0 cm)

Available in Sterile/Non-Sterile packing.
Sterilization by Ethylene Oxide (EO).

More designs are available in European bipolar forceps (single-use).



EUROPEAN BIPOLAR FORCEPS SINGLE-USE (NON-STICK)



Jeweler #5
TWS 84-2611/11



Jeweler #5
(Angled)
TWS 84-2611/12



Adson
TWS 84-2611/13



Straight
TWS 84-2611/14



Straight
(Angled)
TWS 84-2611/15

Cat #	Description	Tip (mm)	Size
TWS 84-2611/11	Jeweler #5	0.5, 0.7	4½" (11.5 cm)
TWS 84-2611/12	Jeweler #5 (Ang.)	0.5, 0.7	4½" (11.5 cm)
TWS 84-2611/13	Adson	1.0	4¾" (12.1 cm)
TWS 84-2611/14	Straight	1.0, 1.5, 2.0	7¾" (20.0 cm)
TWS 84-2611/15	Straight (Ang.)	1.0, 1.5, 2.0	7¾" (20.0 cm)

Available in Sterile/Non-Sterile packing.
Sterilization by Ethylene Oxide (EO).

More designs are available in non-stick
European bipolar forceps (single-use).

EUROPEAN BIPOLAR FORCEPS

SINGLE-USE (ULTRA NON-STICK)



Jeweler #5
TWS 84-2611/16



Jeweler #5
(Angled)
TWS 84-2611/17



Adson
TWS 84-2611/18



Straight
TWS 84-2611/19




Straight
(Angled)
TWS 84-2611/20

Cat #	Description	Tip (mm)	Size
TWS 84-2611/16	Jeweler #5	0.5, 0.7	4½" (11.5 cm)
TWS 84-2611/17	Jeweler #5 (Ang.)	0.5, 0.7	4½" (11.5 cm)
TWS 84-2611/18	Adson	1.0	4¾" (12.1 cm)
TWS 84-2611/19	Straight	1.0, 1.5, 2.0	7¾" (20.0 cm)
TWS 84-2611/20	Straight (Ang.)	1.0, 1.5, 2.0	7¾" (20.0 cm)

Available in Sterile/Non-Sterile packing.
Sterilization by Ethylene Oxide (EO).

More designs are available in ultra non-stick
European bipolar forceps (single-use).

Silicone Bipolar Cable
10' (3.0 m)



Artery Sealer
TWS 84-2611/21



Artery Sealer
(Non-Stick)
TWS 84-2611/22





Artery Sealer
TWS 84-2611/23



Artery Sealer
(Non-Stick)
TWS 84-2611/24

Cat #	Description	Size
TWS 84-2611/21	Bipolar Artery Sealer	10½" (26.0 cm)
TWS 84-2611/22	Bipolar Artery Sealer (Non-Stick)	10½" (26.0 cm)
TWS 84-2611/23	Bipolar Artery Sealer	9½" (24.0 cm)
TWS 84-2611/24	Bipolar Artery Sealer (Non-Stick)	9½" (24.0 cm)


Silicone Bipolar Cable
10' (3.0 m)



Artery Sealer
TWS 84-2611/26



Artery Sealer
(Non-Stick)
TWS 84-2611/27




Artery Sealer
TWS 84-2611/28



Artery Sealer
(Non-Stick)
TWS 84-2611/29

Cat #	Description	Size
TWS 84-2611/26	Bipolar Artery Sealer	6½" (16.0 cm)
TWS 84-2611/27	Bipolar Artery Sealer (Non-Stick)	6½" (16.0 cm)
TWS 84-2611/28	Bipolar Artery Sealer	7½" (19.0 cm)
TWS 84-2611/29	Bipolar Artery Sealer (Non-Stick)	7½" (19.0 cm)

Silicone Bipolar Cable
10' (3.0 m)



Artery Sealer
TWS 84-2611/30



Artery Sealer
(Non-Stick)
TWS 84-2611/31



Artery Sealer
TWS 84-2611/33



Artery Sealer
(Non-Stick)
TWS 84-2611/32

Cat #	Description	Size
TWS 84-2611/30	Bipolar Artery Sealer	8½" (22.0 cm)
TWS 84-2611/31	Bipolar Artery Sealer (Non-Stick)	8½" (22.0 cm)
TWS 84-2611/33	Bipolar Artery Sealer	9½" (24.0 cm)
TWS 84-2611/32	Bipolar Artery Sealer (Non-Stick)	9½" (24.0 cm)
TWS 84-2603/34	Bipolar Artery Sealer	10½" (26.0 cm)
TWS 84-2603/36	Bipolar Artery Sealer (Non-Stick)	10½" (26.0 cm)


BIPOLAR SCISSORS



TW Surgical



Bipolar Scissors
(Curved)
TWS 84-2611/34



Bipolar Scissors
(Curved)
TWS 84-2611/35




Bipolar Scissors
(Curved)
TWS 84-2611/36

Cat #	Description	Size
TWS 84-2611/34	Bipolar Scissors (Crv.)	7" (18.0 cm)
TWS 84-2611/35	Bipolar Scissors (Crv.)	8 1/4" (21.0 cm)
TWS 84-2611/36	Bipolar Scissors (Crv.)	9" (23.0 cm)



Ethicon Bipolar Cable
10' (3.0 m)



Bipolar Clamp Scissors
TWS 84-2611/38




Bipolar Clamp Scissors
TWS 84-2611/39

Cat #	Description	Size
TWS 84-2611/38	Bipolar Clamp Scissors	6" (15.0 cm)
TWS 84-2611/39	Bipolar Clamp Scissors	7½" (19.0 cm)

TONSILLECTOMY BIPOLAR CLAMP



TW Surgical



Silicone Bipolar Cable
10' (3.0 m)



Tonsillectomy Bipolar Clamp
TWS 84-2611/41




Tonsillectomy Bipolar Clamp (Single-Use)
TWS 84-2611/42

Cat #	Description	Size
TWS 84-2611/41	Tonsillectomy Bipolar Clamp	7¼" (18.5 cm)
TWS 84-2611/42	Tonsillectomy Bipolar Clamp (Single-Use)	7¼" (18.5 cm)

ENDOSCOPIC BIPOLAR FORCEPS

Series supplied with
Banana plug



Endoscopic Bipolar
Forceps (Straight)
(Reusable)
26.0 cm
TWS 84-2611/43

Endoscopic Bipolar
Forceps (Angled)
(Reusable)
26.0 cm
TWS 84-2611/44

Endoscopic Bipolar
Forceps (Straight)
(Single-Use)
27.0 cm
TWS 84-2611/45

Endoscopic Bipolar
Forceps (Angled)
(Single-Use)
27.0 cm
TWS 84-2611/46


Endoscopic Bipolar
Forceps (Straight)
(Single-Use)
27.0 cm
TWS 84-2611/47

Endoscopic Bipolar
Forceps (Angled)
(Single-Use)
27.0 cm
TWS 84-2611/47

POLYPECTOMY SNARE



TW Surgical



Polypectomy Snare
TWS 84-2611/49




Polypectomy Snare (Single-Use)
TWS 84-2611/50



TW Surgical

DIATHERMY INSTRUMENTS



Tuner-Warwick
(Straight)




Wilson Hey Forceps
(Straight)




Tuner-Warwick
(Curved)



Wilson Hey Forceps
(Curved)



Surgical Scissors
(Straight)




Artery Forceps
(Straight)



Surgical Scissors
(Curved)



Artery Forceps
(Curved)



Richie Forceps
(Straight)



Abbey Needle

English Monopolar
(Male 4.8 mm connector)

European Monopolar
(Female 4.0 mm connector)

Cat #	Cat #	Description	Size
TWS 84-2603/84	TWS 84-2603/64	Tuner-Warwick Scissors (Str.)	7" (17.8 cm)
TWS 84-2603/85	TWS 84-2603/65	Tuner-Warwick Scissors (Crv.)	7" (17.8 cm)
TWS 84-2603/86	TWS 84-2603/66	Tuner-Warwick Scissors (Str.)	9" (23.0 cm)
TWS 84-2603/87	TWS 84-2603/67	Tuner-Warwick Scissors (Crv.)	9" (23.0 cm)
TWS 84-2603/88	TWS 84-2603/68	Surgical Scissors (Str.)	7" (17.8 cm)
TWS 84-2603/89	TWS 84-2603/69	Surgical Scissors (Crv.)	7" (17.8 cm)
TWS 84-2603/90	TWS 84-2603/70	Surgical Scissors (Str.)	9" (23.0 cm)
TWS 84-2603/91	TWS 84-2603/71	Surgical Scissors (Crv.)	9" (23.0 cm)
TWS 84-2603/92	TWS 84-2603/72	Wilson Hey Forceps (Str.)	7" (17.8 cm)
TWS 84-2603/93	TWS 84-2603/73	Wilson Hey Forceps (Crv.)	7" (17.8 cm)
TWS 84-2603/94	TWS 84-2603/74	Wilson Hey Forceps (Str.)	9" (23.0 cm)
TWS 84-2603/95	TWS 84-2603/75	Wilson Hey Forceps (Crv.)	9" (23.0 cm)
TWS 84-2603/96	TWS 84-2603/76	Artery Forceps (Str.)	7" (17.8 cm)
TWS 84-2603/97	TWS 84-2603/77	Artery Forceps (Crv.)	7" (17.8 cm)
TWS 84-2603/98	TWS 84-2603/78	Artery Forceps (Str.)	9" (23.0 cm)
TWS 84-2603/99	TWS 84-2603/79	Artery Forceps (Crv.)	9" (23.0 cm)
TWS 84-2604/00	TWS 84-2603/80	Richie Forceps (Str.)	7" (17.8 cm)
TWS 84-2604/01	TWS 84-2603/81	Richie Forceps (Str.)	8" (20.0 cm)
TWS 84-2604/02	TWS 84-2603/82	Richie Forceps (Str.)	9" (23.0 cm)
TWS 84-2604/03	TWS 84-2603/83	Abbey Needle	6¾" (17.5 cm)

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com

DIATHERMY INSTRUMENTS



TW Surgical



Adson



McIndoe



Gillies (1x2)



Dissecting Forceps



Gerald (1x2)



Adson Bayonet



Remington Hobb



DeBakey Tissue Forceps

English Monopolar
(Male 4.8 mm connector)

European Monopolar
(Female 4.0 mm connector)

European Monopolar
(Male 4.0 mm connector)

Monopolar
(Without connector)

Cat #	Cat #	Cat #	Cat #	Description	Tip (mm)	Size
TWS 84-2604/53	TWS 84-2604/04	TWS 84-2604/13	TWS 84-2604/32	Adson	1.0	4 3/4" (12.0 cm)
TWS 84-2604/54	TWS 84-2604/05	TWS 84-2604/14	TWS 84-2604/33	McIndoe	1.0, 2.0	6" (15.0 cm)
TWS 84-2604/55	TWS 84-2604/06	TWS 84-2604/15	TWS 84-2604/35	Gillies (1x2)	-	6" (15.0 cm)
TWS 84-2604/56	TWS 84-2604/07	TWS 84-2604/16	TWS 84-2604/36	Gillies (1x2)	-	7" (17.8 cm)
TWS 84-2604/57	TWS 84-2604/08	TWS 84-2604/17	TWS 84-2604/37	Gillies (1x2)	-	8" (20.0 cm)
TWS 84-2604/58	TWS 84-2604/09	TWS 84-2604/18	TWS 84-2604/38	Gillies (1x2)	-	9" (23.0 cm)
TWS 84-2604/59	TWS 84-2604/10	TWS 84-2604/19	TWS 84-2604/39	Dissecting Forceps	1.0, 2.0	6" (15.0 cm)
TWS 84-2604/60	TWS 84-2604/11	TWS 84-2604/20	TWS 84-2604/40	Dissecting Forceps	1.0, 2.0	7" (17.8 cm)
TWS 84-2604/61	TWS 84-2604/12	TWS 84-2604/21	TWS 84-2604/41	Dissecting Forceps	1.0, 2.0	8 1/2" (21.0 cm)
TWS 84-2604/62	TWS 84-2604/13	TWS 84-2604/22	TWS 84-2604/42	Dissecting Forceps	1.0, 2.0	9 1/2" (24.1 cm)
TWS 84-2604/63	TWS 84-2604/14	TWS 84-2604/23	TWS 84-2604/43	Gerald (1x2)	-	7" (17.8 cm)
TWS 84-2604/64	TWS 84-2604/15	TWS 84-2604/24	TWS 84-2604/44	Adson Bayonet	1.0, 2.0	6 1/2" (16.0 cm)
TWS 84-2604/65	TWS 84-2604/16	TWS 84-2604/25	TWS 84-2604/45	Adson Bayonet	1.0, 2.0	7 1/2" (19.1 cm)
TWS 84-2604/66	TWS 84-2604/17	TWS 84-2604/26	TWS 84-2604/46	Adson Bayonet (Ang. Up)	1.0, 2.0	6 1/2" (16.0 cm)
TWS 84-2604/67	TWS 84-2604/18	TWS 84-2604/27	TWS 84-2604/47	Adson Bayonet (Ang. Up)	1.0, 2.0	7 1/2" (19.1 cm)
TWS 84-2604/68	TWS 84-2604/19	TWS 84-2604/28	TWS 84-2604/48	Adson Bayonet (Ang. Up)	1.0, 2.0	8 1/2" (21.0 cm)
TWS 84-2604/69	TWS 84-2604/10	TWS 84-2604/29	TWS 84-2604/49	Remington Hob	1.0, 2.0	9 1/2" (24.1 cm)
TWS 84-2604/70	TWS 84-2604/11	TWS 84-2604/30	TWS 84-2604/50	DeBakey Tissue Forceps	1.5, 2.0	7 3/4" (19.5 cm)
TWS 84-2604/71	TWS 84-2604/12	TWS 84-2604/31	TWS 84-2604/51	DeBakey Tissue Forceps	1.5, 2.0	9 1/2" (24.1 cm)



TW Surgical

BIPOLAR ELECTRODES



Straight
TWS 84-2611/68



Angled
(Upward Tip)
TWS 84-2611/69



Angled
TWS 84-2611/70



Angled 90°
TWS 84-2611/71



Angled 45°
TWS 84-2611/72

Cat #	Description	Size
TWS 84-2611/68	Straight	18.0 cm
TWS 84-2604/73	Straight	20.0 cm
TWS 84-2604/74	Straight	22.0 cm
TWS 84-2604/75	Straight	27.0 cm
TWS 84-2611/70	Angled	18.0 cm
TWS 84-2604/77	Angled	20.0 cm
TWS 84-2604/78	Angled	22.0 cm
TWS 84-2611/69	Angled (Upward Tip)	16.5 cm
TWS 84-2604/80	Angled (Upward Tip)	20.0 cm
TWS 84-2611/72	Angled 45°	19.0 cm
TWS 84-2611/71	Angled 90°	16.0 cm

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com

BIPOLAR ELECTRODES SINGLE-USE



TW Surgical



Cat #	Description	Size
TWS 84-2611/73	Straight	18.0 cm
TWS 84-2604/84	Straight	20.0 cm
TWS 84-2604/85	Straight	22.0 cm
TWS 84-2604/86	Straight	27.0 cm
TWS 84-2611/75	Angled	18.0 cm
TWS 84-2604/88	Angled	20.0 cm
TWS 84-2604/89	Angled	22.0 cm
TWS 84-2611/74	Angled (Upward Tip)	16.5 cm
TWS 84-2604/91	Angled (Upward Tip)	20.0 cm
TWS 84-2611/77	Angled 45°	19.0 cm
TWS 84-2611/76	Angled 90°	16.0 cm



TW Surgical

BIPOLAR ELECTRODES



Teflon Coated



Bipolar Eraser
TWS 84-2611/78



Bipolar Eraser
(Single-Use)
TWS 84-2611/79



Bi-Turbinate Probe
TWS 84-2611/80



Bi-Turbinate Probe
(Single-Use)
TWS 84-2611/81



Bi-Turbinate Probe
TWS 84-2611/82



Bi-Turbinate Probe
(Single-Use)
TWS 84-2611/83


Cat #	Description	Size
TWS 84-2611/78	Bipolar Eraser (18 Ga) Flat Tip	11.5 cm
TWS 84-2611/79	Bipolar Eraser (18 Ga) Flat Tip (Single-Use)	11.5 cm
TWS 84-2604/96	Bipolar Eraser (18 Ga) Tapered Tip	11.5 cm
TWS 84-2604/97	Bipolar Eraser (18 Ga) Tapered Tip (Single-Use)	11.5 cm

Cat #	Description	Size
TWS 84-2611/80	Bi-Turbinate Probe	17.5 cm
TWS 84-2604/99	Bi-Turbinate Probe Teflon Coated	17.5 cm
TWS 84-2611/81	Bi-Turbinate Probe (Single-Use)	17.5 cm
TWS 84-2605/01	Bi-Turbinate Probe Teflon Coated (Single-Use)	17.5 cm
TWS 84-2611/82	Bi-Turbinate Probe	19.5 cm
TWS 84-2605/03	Bi-Turbinate Probe Teflon Coated	19.5 cm
TWS 84-2611/83	Bi-Turbinate Probe (Single-Use)	19.5 cm
TWS 84-2605/05	Bi-Turbinate Probe Teflon Coated (Single-Use)	19.5 cm



TW Surgical


BIPOLAR CABLES



Silicone Bipolar Cable

Cat #	Size
TWS 84-2620/30	10' (3.0 m)
TWS 84-2620/50	16' (5.0 m)


Instrument Connector: USA 2-Pin Plug
 Generator Connector: 4 mm Banana Plug
 Compatibility: Valleylab, Conmed, Codman, Olympus & most Bipolar Generator



Silicone Bipolar Cable

Cat #	Size
TWS 84-2605/06	10' (3.0 m)
TWS 84-2605/07	16' (5.0 m)


Instrument Connector: USA 2-Pin Plug
 Generator Connector: 4 mm Banana Plug
 Compatibility: Valleylab, Conmed, Codman, Olympus & most Bipolar Generator



Silicone European Bipolar Cable

Cat #	Size
TWS 84-2605/08	10' (3.0 m)
TWS 84-2605/09	16' (5.0 m)

Instrument Connector: European Flat Plug
 Generator Connector: 4 mm Banana Plug
 Compatibility: Valleylab, Conmed, Codman, Olympus & most Bipolar Generator



Silicone European Bipolar Cable

Cat #	Size
TWS 84-2605/10	10' (3.0 m)
TWS 84-2605/11	16' (5.0 m)


Instrument Connector: European Flat Plug
 Generator Connector: Martin Type
 Compatibility: Aesculap, Martin, Berchtold, Wolf



Silicone European Bipolar Cable

Cat #	Size
TWS 84-2605/12	10' (3.0 m)
TWS 84-2605/13	16' (5.0 m)


Instrument Connector: European Flat Plug
 Generator Connector: Martin Type
 Compatibility: Aesculap, Martin, Berchtold, Wolf



Silicone Bipolar Cable

Cat #	Size
TWS 84-2605/14	10' (3.0 m)
TWS 84-2605/15	16' (5.0 m)

Instrument Connector: USA 2-Pin Plug
 Generator Connector: Martin Type
 Compatibility: Aesculap, Martin, Berchtold, Wolf



Silicone European Bipolar Cable

Cat #	Size
TWS 84-2605/16	10' (3.0 m)
TWS 84-2605/17	16' (5.0 m)

Instrument Connector: Aesculap
 Generator Connector: 4 mm Banana Plug
 Compatibility: Aesculap, Valleylab, Conmed, Codman & most Bipolar Generator



TW Surgical

BIPOLAR CABLES



Silicone European Bipolar Cable

Cat #	Size
TWS 84-2605/18 10'	(3.0 m)
TWS 84-2605/19 16'	(5.0 m)

Instrument Connector: European Flat Plug
Generator Connector: Erbe Type

Compatibility: Erbe, Downs, Siemens




Silicone European Bipolar Cable

Cat #	Size
TWS 84-2605/20 10'	(3.0 m)
TWS 84-2605/21 16'	(5.0 m)

Instrument Connector: European Flat Plug
Generator Connector: Erbe Type

Compatibility: Erbe, Downs, Siemens




Silicone Bipolar Cable

Cat #	Size
TWS 84-2605/22 10'	(3.0 m)
TWS 84-2605/23 16'	(5.0 m)

Instrument Connector: USA 2-Pin Plug
Generator Connector: Erbe Type

Compatibility: Erbe, Downs, Siemens




Silicone European Bipolar Cable

Cat #	Size
TWS 84-2605/24 10'	(3.0 m)
TWS 84-2605/25 16'	(5.0 m)

Instrument Connector: European Flat Plug
Generator Connector: Pin Distance 29 mm

Compatibility: Valleylab, EMC, Dolley




Silicone Bipolar Cable

Cat #	Size
TWS 84-2605/26 10'	(3.0 m)
TWS 84-2605/27 16'	(5.0 m)

Instrument Connector: USA 2-Pin Plug
Generator Connector: Pin Distance 29 mm

Compatibility: Valleylab, EMC, Dolley

4 mm




Silicone Bipolar Cable

Cat #	Size
TWS 84-2605/28 10'	(3.0 m)
TWS 84-2605/29 16'	(5.0 m)

Instrument Connector: 4 mm Banana Plug
Generator Connector: 4 mm Banana Plug

Compatibility: Valleylab, Conmed, Codman, Olympus & most Bipolar Generator

2.5 mm




Silicone Bipolar Cable

Cat #	Size
TWS 84-2605/30 10'	(3.0 m)
TWS 84-2605/31 16'	(5.0 m)

Instrument Connector: Ethicon / 2.5 mm
Generator Connector: 4 mm Banana Plug

Compatibility: Valleylab, Conmed, Codman, Olympus & most Bipolar Generator

BIPOLAR CABLES



Silicone European Bipolar Cable

Cat #	Size
TWS 84-2605/32	10' (3.0 m)
TWS 84-2605/33	16' (5.0 m)

Instrument Connector: Aesculap
Generator Connector: Martin Type

Compatibility: Aesculap, Martin, Berchtold, Wolf




Silicone European Bipolar Cable

Cat #	Size
TWS 84-2605/34	10' (3.0 m)
TWS 84-2605/35	16' (5.0 m)

Instrument Connector: Aesculap
Generator Connector: Erbe Type

Compatibility: Erbe, Downs, Siemens




Silicone European Bipolar Cable

Cat #	Size
TWS 84-2605/36	10' (3.0 m)
TWS 84-2605/37	16' (5.0 m)

Instrument Connector: Aesculap
Generator Connector: Pin Distance 29 mm

Compatibility: Valleylab, EMC, Dolley




Silicone Bipolar Cable

Cat #	Size
TWS 84-2605/39	10' (3.0 m)
TWS 84-2605/38	16' (5.0 m)

Instrument Connector: USA 2-Pin Plug
Generator Connector: Universal 18 - 28 mm

Compatibility: Valleylab, Conmed, Codman, Olympus & most Bipolar Generator




Silicone European Bipolar Cable

Cat #	Size
TWS 84-2605/41	10' (3.0 m)
TWS 84-2605/40	16' (5.0 m)

Instrument Connector: European Flat Plug

Generator Connector: Universal 18 - 28 mm

Compatibility: Valleylab, Conmed, Codman, Olympus & most Bipolar Generator



Silicone European Bipolar Cable

Cat #	Size
TWS 84-2605/42	10' (3.0 m)
TWS 84-2605/43	16' (5.0 m)

Instrument Connector: Aesculap


Generator Connector: Universal 18 - 28 mm

Compatibility: Valleylab, Conmed, Codman, Olympus & most Bipolar Generator



TW Surgical

BIPOLAR CABLES




Silicone Bipolar Cable

Cat #	Size
TWS 84-2605/44	10' (3.0 m)
TWS 84-2605/45	16' (5.0 m)

Instrument Connector: USA 2-Pin Plug

Generator Connector: Ellman Plug

Compatibility: Ellman




Silicone Bipolar Cable

Cat #	Size
TWS 84-2605/47	10' (3.0 m)
TWS 84-2605/48	16' (5.0 m)

Instrument Connector: USA 2-Pin Plug

Generator Connector: Ellman Plug

Compatibility: Ellman



Silicone European Bipolar Cable

Cat #	Size
TWS 84-2605/49	10' (3.0 m)
TWS 84-2605/50	16' (5.0 m)


Instrument Connector: European Flat Plug

Generator Connector: Ellman Plug

Compatibility: Ellman



4 mm



Silicone Bipolar Cable

Cat #	Size
TWS 84-2605/51	10' (3.0 m)
TWS 84-2605/52	16' (5.0 m)


Instrument Connector: 4 mm Banana Plug

Generator Connector: Ellman Plug

Compatibility: Ellman



2.5 mm



Silicone Bipolar Cable


Cat #	Size
TWS 84-2605/53	10' (3.0 m)
TWS 84-2605/54	16' (5.0 m)

Instrument Connector: Ethicon / 2.5 mm

Generator Connector: Ellman Plug

Compatibility: Ellman


MONOPOLAR CABLES



Silicone Monopolar Cable

Cat #	Size
TWS 84-2605/5510'	(3.0 m)
TWS 84-2605/5616'	(5.0 m)


Instrument Connector: 4 mm Male Fitting (European Diathermy)
 Generator Connector: Standard Bovie® 8 mm Plug
 Compatibility: Valleylab, Conmed, Bovie®



Silicone Monopolar Cable

Cat #	Size
TWS 84-2605/5710'	(3.0 m)
TWS 84-2605/5816'	(5.0 m)

Instrument Connector: Screw Type
 Generator Connector: Standard Bovie® 8 mm Plug
 Compatibility: Valleylab, Conmed, Bovie®



Silicone Monopolar Endoscopic Cable

Cat #	Size
TWS 84-2605/5910'	(3.0 m)
TWS 84-2605/6016'	(5.0 m)

Instrument Connector: 4 mm Female Fitting
 Generator Connector: Standard Bovie® 8 mm Plug
 Compatibility: Valleylab, Conmed, Bovie®



Silicone Monopolar Endoscopic Cable

Cat #	Size
TWS 84-2605/6110'	(3.0 m)
TWS 84-2605/6216'	(5.0 m)

Instrument Connector: 4 mm Female Fitting
 Generator Connector: Standard 4 mm Banana Plug
 Compatibility: Most ES Generators



Silicone Diathermy Cable


Cat #	Size
TWS 84-2605/6310'	(3.0 m)
TWS 84-2605/6416'	(5.0 m)

Instrument Connector: 4 mm Male Fitting (European Diathermy)
 Generator Connector: Standard 4 mm Banana Plug
 Compatibility: Most ES Generators



TW Surgical

MONOPOLAR CABLES



Silicone Diathermy Cable

Cat #	Size
TWS 84-2605/6510'	(3.0 m)
TWS 84-2605/6616'	(5.0 m)


Instrument Connector: 4.8 mm Female Fitting (English Diathermy)
 Generator Connector: Standard Bovie® 8 mm Plug
 Compatibility: Valleylab, Conmed, Bovie®



Silicone Diathermy Cable

Cat #	Size
TWS 84-2605/6710'	(3.0 m)
TWS 84-2605/6816'	(5.0 m)

Instrument Connector: 4.8 mm Female Fitting (English Diathermy)
 Generator Connector: Standard 4 mm Banana Plug
 Compatibility: Most ES Generators



Silicone Urology (Double Stem) Cable

Cat #	Size
TWS 84-2605/6910'	(3.0 m)
TWS 84-2605/7016'	(5.0 m)

Instruments Connector: Resectoscope
 Generator Connector: Standard 4 mm Banana Plug
 Compatibility: Most ES Generators




Silicone Urology (Single Stem) Cable

Cat #	Size
TWS 84-2605/7110'	(3.0 m)
TWS 84-2605/7216'	(5.0 m)

Instruments Connector: Resectoscope
 Generator Connector: Standard 4 mm Banana Plug
 Compatibility: Most ES Generators

BIPOLAR CABLES SINGLE-USE




Bipolar Cable

Cat #	Size
-------	------

TWS 84-2605/73 10' (3.0 m)
TWS 84-2605/74 16' (5.0 m)

Instrument Connector: USA 2-Pin Plug
Generator Connector: 4 mm Banana Plug

Compatibility: Valleylab, Conmed, Codman, Olympus & most Bipolar Generator




European Bipolar Cable

Cat #	Size
-------	------

TWS 84-2605/75 10' (3.0 m)
TWS 84-2605/76 16' (5.0 m)

Instrument Connector: European Flat Plug
Generator Connector: 4 mm Banana Plug

Compatibility: Valleylab, Conmed, Codman, Olympus & most Bipolar Generator




European Bipolar Cable

Cat #	Size
-------	------

TWS 84-2605/77 10' (3.0 m)
TWS 84-2605/78 16' (5.0 m)

Instrument Connector: European Flat Plug
Generator Connector: Martin Type

Compatibility: Aesculap, Martin, Berchtold, Wolf



Bipolar Cable

Cat #	Size
-------	------

TWS 84-2605/79 10' (3.0 m)
TWS 84-2605/80 16' (5.0 m)

Instrument Connector: USA 2-Pin Plug
Generator Connector: Martin Type

Compatibility: Aesculap, Martin, Berchtold, Wolf




European Bipolar Cable

Cat #	Size
-------	------

TWS 84-2605/81 10' (3.0 m)
TWS 84-2605/82 16' (5.0 m)

Instrument Connector: European Flat Plug
Generator Connector: Erbe Type

Compatibility: Erbe, Downs, Siemens



Bipolar Cable

Cat #	Size
-------	------

TWS 84-2605/83 10' (3.0 m)
TWS 84-2605/84 16' (5.0 m)

Instrument Connector: USA 2-Pin Plug
Generator Connector: Erbe Type

Compatibility: Erbe, Downs, Siemens



TW Surgical

BIPOLAR CABLES SINGLE-USE



European Bipolar Cable

Cat #	Size
TWS 84-2605/85	10' (3.0 m)
TWS 84-2605/86	16' (5.0 m)

Instrument Connector: European Flat Plug
Generator Connector: Pin Distance 29 mm

Compatibility: Valleylab, EMC, Dolley




Bipolar Cable

Cat #	Size
TWS 84-2605/87	10' (3.0 m)
TWS 84-2605/88	16' (5.0 m)

Instrument Connector: USA 2-Pin Plug
Generator Connector: Pin Distance 29 mm

Compatibility: Valleylab, EMC, Dolley




European Bipolar Cable

Cat #	Size
TWS 84-2605/89	10' (3.0 m)
TWS 84-2605/90	16' (5.0 m)

Instrument Connector: Aesculap
Generator Connector: 4 mm Banana Plug

Compatibility: Aesculap, Valleylab, Conmed, Codman & most Bipolar Generator




European Bipolar Cable

Cat #	Size
TWS 84-2605/91	10' (3.0 m)
TWS 84-2605/92	16' (5.0 m)

Instrument Connector: Aesculap
Generator Connector: Erbe Type

Compatibility: Erbe, Downs, Siemens



European Bipolar Cable

Cat #	Size
TWS 84-2605/93	10' (3.0 m)
TWS 84-2605/94	16' (5.0 m)

Instrument Connector: Aesculap
Generator Connector: Pin Distance 29 mm

Compatibility: Valleylab, EMC, Dolley



Bipolar Cable

Cat #	Size
TWS 84-2605/95	10' (3.0 m)
TWS 84-2605/96	16' (5.0 m)

Instrument Connector: USA 2-Pin Plug
Generator Connector: Ellman Plug

Compatibility: Ellman



European Bipolar Cable

Cat #	Size
TWS 84-2605/97	10' (3.0 m)
TWS 84-2605/98	16' (5.0 m)

Instrument Connector: European Flat Plug
Generator Connector: Ellman Plug

Compatibility: Ellman

MONOPOLAR CABLES SINGLE-USE



Monopolar Cable

Cat #	Size
TWS 84-2605/99	10' (3.0 m)
TWS 84-2606/00	16' (5.0 m)

Instrument Connector: 4 mm Male Fitting (European Diathermy)

Generator Connector: Standard Bovie® 8 mm Plug

Compatibility: Valleylab, Conmed, Bovie®



Monopolar Cable

Cat #	Size
TWS 84-2606/01	10' (3.0 m)
TWS 84-2606/02	16' (5.0 m)

Instrument Connector: Screw Type

Generator Connector: Standard Bovie® 8 mm Plug

Compatibility: Valleylab, Conmed, Bovie®



Monopolar Cable

Cat #	Size
TWS 84-2606/03	10' (3.0 m)
TWS 84-2606/04	16' (5.0 m)

Instrument Connector: 4 mm Female Fitting

Generator Connector: Standard Bovie® 8 mm Plug

Compatibility: Valleylab, Conmed, Bovie®




Monopolar Cable

Cat #	Size
TWS 84-2606/05	10' (3.0 m)
TWS 84-2606/06	16' (5.0 m)

Instrument Connector: 4 mm Female Fitting

Generator Connector: Standard 4 mm Banana Plug

Compatibility: Most ES Generators



Monopolar Cable

Cat #	Size
TWS 84-2606/07	10' (3.0 m)
TWS 84-2606/08	16' (5.0 m)

Instrument Connector: 4 mm Male Fitting (European Diathermy)

Generator Connector: Standard 4 mm Banana Plug

Compatibility: Most ES Generators



Monopolar Cable

Cat #	Size
TWS 84-2606/09	10' (3.0 m)
TWS 84-2606/10	16' (5.0 m)

Instrument Connector: 4.8 mm Female Fitting (English Diathermy)

Generator Connector: Standard Bovie® 8 mm Plug

Compatibility: Valleylab, Conmed, Bovie®



Monopolar Cable

Cat #	Size
TWS 84-2606/11	10' (3.0 m)
TWS 84-2606/12	16' (5.0 m)

Instrument Connector: 4.8 mm Female Fitting (English Diathermy)

Generator Connector: Standard 4 mm Banana Plug

Compatibility: Most ES Generators



TW Surgical

ELECTROSURGICAL PENCILS



Finger Switch Pencil (Single-Use)

Description	Regular	With Tip Cleaner	With Holster	With Holster & With Tip Cleaner
Finger Switch Pencil (3m) Cable	TWS 84-2606/13	TWS 84-2606/15	TWS 84-2606/16	TWS 84-2606/17
Finger Switch Pencil (5m) Cable	TWS 84-2606/14	TWS 84-2606/18	TWS 84-2606/19	TWS 84-2606/20



Finger Switch Pencil - Non-Stick (Single-Use)

Description	Regular	With Tip Cleaner	With Holster	With Holster & With Tip Cleaner
Finger Switch Pencil (3m) Cable	TWS 84-2606/21	TWS 84-2606/23	TWS 84-2606/25	TWS 84-2606/27
Finger Switch Pencil (5m) Cable	TWS 84-2606/22	TWS 84-2606/24	TWS 84-2606/26	TWS 84-2606/28



Packed Finger Switch Pencil



Tip Cleaner

TWS 84-2611/85



Holster

TWS 84-2611/86

Products are available in Sterile/Non-Sterile packaging. Sterilization by Ethylene Oxide (EO).


ELECTROSURGICAL PENCILS



Rocker-Switch Pencil (Single-Use)
TWS 84-2611/87 Cable)



Foot Control Pencil (Single-Use)
TWS 84-2611/88 Cable)
TWS 84-2611/89 Cable)



Finger Switch Pencil - Reusable
TWS 84-2611/90 Silicone Cable)



Electrode Handle - Reusable
(3m Silicone Cable)

Cat #	Size
TWS 84-2606/29	2.4 mm
TWS 84-2606/30	1.6 mm

Products are available in Sterile/Non-Sterile packaging.
Sterilization by Ethylene Oxide (EO).



TW Surgical

ELECTRODES (2.4 mm)



Blade Electrode
TWS 84-2611/91 (Single-Use)
TWS 84-2611/92 (Reusable)



Blade Electrode
TWS 84-2611/93 (Single-Use)
TWS 84-2611/94 (Reusable)



Blade Electrode (Curved)
TWS 84-2611/95 (Single-Use)
TWS 84-2611/96 (Reusable)



Blade Electrode (Curved)
TWS 84-2611/97 (Single-Use)
TWS 84-2611/98 (Reusable)



Knife Electrode
TWS 84-2611/99 (Single-Use)
TWS 84-2612/00 (Reusable)



Knife Electrode (Curved)
TWS 84-2612/01 (Single-Use)
TWS 84-2612/02 (Reusable)



Needle Electrode
TWS 84-2612/03 (Single-Use)
TWS 84-2612/04 (Reusable)



Needle Electrode (Curved)
TWS 84-2612/05 (Single-Use)
TWS 84-2612/06 (Reusable)



Ball Electrode (2 mm)
TWS 84-2612/07 (Single-Use)
TWS 84-2612/08 (Reusable)



Ball Electrode (Curved) (2 mm)
TWS 84-2612/09 (Single-Use)
TWS 84-2612/10 (Reusable)



Ball Electrode (3 mm)
TWS 84-2612/11 (Single-Use)
TWS 84-2612/12 (Reusable)



Ball Electrode (Curved) (3 mm)
TWS 84-2612/13 (Single-Use)
TWS 84-2612/14 (Reusable)



Ball Electrode (4 mm)
TWS 84-2612/15 (Single-Use)
TWS 84-2612/16 (Reusable)



Ball Electrode (Curved) (4 mm)
TWS 84-2612/17 (Single-Use)
TWS 84-2612/18 (Reusable)



Ball Electrode (5 mm)
TWS 84-2612/19 (Single-Use)
TWS 84-2612/20 (Reusable)



Ball Electrode (Curved) (5 mm)
TWS 84-2612/21 (Single-Use)
TWS 84-2612/22 (Reusable)



Blade Electrode (Shaft 13 cm)
TWS 84-2612/23 (Single-Use)
TWS 84-2612/24 (Reusable)



Needle Electrode (Shaft 13 cm)
TWS 84-2612/25 (Single-Use)
TWS 84-2612/26 (Reusable)



Ball Electrode 3 mm (Shaft 13 cm)
TWS 84-2612/27 (Single-Use)
TWS 84-2612/28 (Reusable)



Ball Electrode 5 mm (Shaft 13 cm)
TWS 84-2612/29 (Single-Use)
TWS 84-2612/30 (Reusable)



Circle Loop Electrode (4 mm)
TWS 84-2612/31 (Reusable)



Circle Loop Electrode (6 mm)
TWS 84-2612/32 (Reusable)



Circle Loop Electrode (8 mm)
TWS 84-2612/33 (Reusable)



Diamond Electrode
TWS 84-2612/34 (Reusable)



Fine Needle Electrode
TWS 84-2612/35 (Reusable)

ELECTRODES (2.4 mm) NON-STICK



Blade Electrode (6 cm)
TWS 84-2612/38 (Single-Use)



Blade Electrode (6 cm)
TWS 84-2612/36 (Single-Use)



Blade Electrode (7 cm)
TWS 84-2612/39 (Single-Use)



Blade Electrode (7 cm)
TWS 84-2612/37 (Single-Use)



Blade Electrode (10 cm)
TWS 84-2612/40 (Single-Use)



Blade Electrode (10 cm)
TWS 84-2612/41 (Single-Use)



Blade Electrode (16 cm)
TWS 84-2612/42 (Single-Use)



Blade Electrode (16 cm)
TWS 84-2612/43 (Single-Use)



Needle Electrode (7 cm)
TWS 84-2612/44 (Single-Use)



Needle Electrode (7 cm)
TWS 84-2612/45 (Single-Use)



Needle Electrode (10 cm)
TWS 84-2612/46 (Single-Use)



Needle Electrode (10 cm)
TWS 84-2612/47 (Single-Use)



Needle Electrode (15 cm)
TWS 84-2612/48 (Single-Use)



Needle Electrode (15 cm)
TWS 84-2612/49 (Single-Use)



Ball Electrode (5 cm)
TWS 84-2612/50 (Single-Use)



Ball Electrode (13 cm)
TWS 84-2612/51 (Single-Use)



TUNGSTEN NEEDLE ELECTRODES (2.4 mm)

TW Surgical



Micro Needle Electrode (5 cm)
TWS 84-2612/52 (Single-Use)
TWS 84-2612/53 (Reusable)



Micro Needle Electrode Angled (5 cm)
TWS 84-2612/54 (Single-Use)
TWS 84-2612/55 (Reusable)



Micro Needle Electrode (6 cm)
TWS 84-2612/56 (Single-Use)
TWS 84-2612/57 (Reusable)



Micro Needle Electrode Angled (6 cm)
TWS 84-2612/58 (Single-Use)
TWS 84-2612/59 (Reusable)



Micro Needle Electrode (7 cm)
TWS 84-2612/60 (Single-Use)
TWS 84-2612/61 (Reusable)



Micro Needle Electrode Angled (7 cm)
TWS 84-2612/62 (Single-Use)
TWS 84-2612/63 (Reusable)



Micro Needle Electrode (9 cm)
TWS 84-2612/64 (Single-Use)
TWS 84-2612/65 (Reusable)



Micro Needle Electrode Angled (9 cm)
TWS 84-2612/66 (Single-Use)
TWS 84-2612/67 (Reusable)



Micro Needle Electrode (12 cm)
TWS 84-2612/68 (Single-Use)
TWS 84-2612/69 (Reusable)



Micro Needle Electrode Angled (12 cm)
TWS 84-2612/70 (Single-Use)
TWS 84-2612/71 (Reusable)



Micro Needle Electrode (15 cm)
TWS 84-2612/72 (Single-Use)
TWS 84-2612/73 (Reusable)




Micro Needle Electrode Angled (15 cm)
TWS 84-2612/74 (Single-Use)
TWS 84-2612/75 (Reusable)

Tungsten Needles are widely used in a variety of surgical procedures like;

- Plastic Surgery
- ENT Surgery
- Head & Neck Surgery
- Oral & Maxillo Surgery

Tungsten has extremely high melting point (>3400°C) which helps maintain a sharpness compared to stainless steel tips that dull rapidly. All our tungsten needles are compatible with Colorado needles.



Micro Needle Electrode (Reusable)

LOOP ELECTRODES (2.4 mm)



TW Surgical

All loop electrodes are available in 8.5 cm shaft length also, use:



Size	1.0 x 0.4 cm	1.0 x 0.5 cm	1.0 x 0.7 cm	1.0 x 0.8 cm	1.0 x 1.0 cm
Single-Use	TWS 84-2606/31	TWS 84-2606/33	TWS 84-2606/35	TWS 84-2606/37	TWS 84-2606/39
Reusable	TWS 84-2606/32	TWS 84-2606/34	TWS 84-2606/36	TWS 84-2606/38	TWS 84-2606/40



Size	1.0 x 0.5 cm	1.0 x 0.7 cm	1.0 x 0.8 cm	1.0 x 1.0 cm	1.0 x 1.5 cm
Single-Use	TWS 84-2606/47	150-175-07	TWS 84-2606/45	TWS 84-2606/43	TWS 84-2606/41
Reusable	TWS 84-2606/48	152-175-07	TWS 84-2606/46	TWS 84-2606/44	TWS 84-2606/42



Size	1.5 x 0.5 cm	1.5 x 0.7 cm	1.5 x 0.8 cm	1.5 x 1.0 cm	1.5 x 1.5 cm	1.5 x 2.0 cm
Single-Use	TWS 84-2606/49	TWS 84-2606/51	TWS 84-2606/53	TWS 84-2606/55	TWS 84-2606/57	TWS 84-2606/59
Reusable	TWS 84-2606/50	TWS 84-2606/52	TWS 84-2606/54	TWS 84-2606/56	TWS 84-2606/58	TWS 84-2606/60



Size	2.0 x 0.8 cm	2.0 x 1.0 cm	2.0 x 1.2 cm	2.0 x 1.3 cm	2.0 x 1.5 cm	2.0 x 2.0 cm
Single-Use	TWS 84-2606/61	TWS 84-2606/62	TWS 84-2606/63	TWS 84-2606/64	TWS 84-2606/65	TWS 84-2606/66
Reusable	TWS 84-2606/67	TWS 84-2606/68	TWS 84-2606/69	TWS 84-2606/70	TWS 84-2606/65	TWS 84-2606/71



Size	2.5 x 1.0 cm	2.5 x 1.2 cm	2.5 x 1.3 cm	2.5 x 1.5 cm	2.5 x 2.0 cm	2.5 x 2.5 cm
Single-Use	TWS 84-2606/72	TWS 84-2606/73	TWS 84-2606/74	TWS 84-2606/75	TWS 84-2606/76	TWS 84-2606/77
Reusable	TWS 84-2606/78	TWS 84-2606/79	TWS 84-2606/74	TWS 84-2606/80	TWS 84-2606/81	TWS 84-2606/83



TW Surgical

LOOP ELECTRODES (2.4 mm)



Size	1.0 x 0.4 cm	1.0 x 0.5 cm	1.0 x 0.7 cm	1.0 x 0.8 cm	1.0 x 1.0 cm
Single-Use	TWS 84-2606/84	TWS 84-2606/85	TWS 84-2606/86	TWS 84-2606/87	TWS 84-2606/88
Reusable	TWS 84-2606/89	TWS 84-2606/90	TWS 84-2606/91	TWS 84-2606/92	TWS 84-2606/93



Size	1.0 x 0.5 cm	1.0 x 0.7 cm	1.0 x 0.8 cm	1.0 x 1.0 cm
Single-Use	TWS 84-2606/97	TWS 84-2606/96	TWS 84-2606/95	TWS 84-2606/94
Reusable	TWS 84-2606/98	TWS 84-2606/99	TWS 84-2607/00	TWS 84-2607/01



Size	1.5 x 0.5 cm	1.5 x 0.7 cm	1.5 x 0.8 cm	1.5 x 1.0 cm	1.5 x 1.2 cm	1.5 x 1.5 cm
Single-Use	TWS 84-2607/02	TWS 84-2607/03	TWS 84-2607/04	TWS 84-2607/05	TWS 84-2607/06	TWS 84-2607/07
Reusable	TWS 84-2607/08	TWS 84-2607/09	TWS 84-2607/10	TWS 84-2607/11	TWS 84-2607/12	TWS 84-2607/13



Size	2.0 x 0.8 cm	2.0 x 1.0 cm	2.0 x 1.2 cm	2.0 x 1.3 cm	2.0 x 1.5 cm	2.0 x 2.0 cm
Single-Use	TWS 84-2607/14	TWS 84-2607/15	TWS 84-2607/16	TWS 84-2607/17	TWS 84-2607/18	TWS 84-2607/19
Reusable	TWS 84-2607/20	TWS 84-2607/21	TWS 84-2607/22	TWS 84-2607/23	TWS 84-2607/24	TWS 84-2607/25



TW Surgical

CONIZATION ELECTRODES (2.4 mm)



Conization Electrode (Small)
TWS 84-2612/76



Conization Electrode (Medium)
TWS 84-2612/77



Conization Electrode (Large)
TWS 84-2612/78



Conization Electrode (Medium Extended)
TWS 84-2612/79



Conization Electrode (Large Shallow)
TWS 84-2612/80



Electrode Extender (7 cm)
TWS 84-2612/81 (Reusable)



Electrode Extender (13 cm)
TWS 84-2612/82 (Reusable)



Electrode Extender (18 cm)
TWS 84-2612/83 (Reusable)

Cat #	Description	Cone Size
TWS 84-2607/26	Conization Electrode	8 x 15 mm
TWS 84-2607/27	Conization Electrode	11 x 15 mm
TWS 84-2607/28	Conization Electrode	13 x 15 mm
TWS 84-2607/29	Conization Electrode	11 x 21 mm
TWS 84-2607/30	Conization Electrode	14 x 8 mm

Cat #	Description	Size
TWS 84-2607/31	Electrode Extender	7.0 cm
TWS 84-2607/32	Electrode Extender	13.0 cm
TWS 84-2607/33	Electrode Extender	15.0 cm
TWS 84-2607/34	Electrode Extender	18.0 cm



TW Surgical

ELECTRODES (4.0 mm)



Blade Electrode
TWS 84-2612/84 (Reusable)



Blade Electrode (Shaft 13 cm)
TWS 84-2612/96 (Reusable)



Blade Electrode (Curved)
TWS 84-2612/85(Reusable)



Needle Electrode (Shaft 13 cm)
TWS 84-2612/97(Reusable)



Needle Electrode
TWS 84-2612/86(Reusable)



Ball Electrode 3 mm (Shaft 13 cm)
TWS 84-2612/98 (Reusable)



Needle Electrode (Curved)
TWS 84-2612/87 (Reusable)



Ball Electrode 4 mm (Shaft 13 cm)
TWS 84-2612/99 (Reusable)



Ball Electrode (2 mm)
TWS 84-2612/88 (Reusable)



Ball Electrode 5 mm (Shaft 13 cm)
TWS 84-2613/00(Reusable)



Ball Electrode (2 mm) Curved
TWS 84-2612/89(Reusable)



Circle Loop Electrode 7 cm (4 mm)
TWS 84-2613/01 (Reusable)



Ball Electrode (3 mm)
TWS 84-2612/90 (Reusable)



Circle Loop Electrode 7 cm (6 mm)
TWS 84-2613/02 (Reusable)



Ball Electrode (3 mm) Curved
TWS 84-2612/91(Reusable)



Circle Loop Electrode 7 cm (8 mm)
TWS 84-2613/03 (Reusable)



Ball Electrode (4 mm)
TWS 84-2612/92 (Reusable)



Circle Loop Electrode 7 cm (10 mm)
TWS 84-2613/04(Reusable)



Ball Electrode (4 mm) Curved
TWS 84-2612/93(Reusable)



Plate Electrode (Small)
TWS 84-2613/05(Reusable)



Ball Electrode (5 mm)
TWS 84-2612/94 (Reusable)



Plate Electrode (Large)
TWS 84-2613/06(Reusable)



Ball Electrode (5 mm) Curved
TWS 84-2612/95 (Reusable)

LOOP ELECTRODES (4.0 mm)



TW Surgical



Size	1.0 x 0.4 cm	1.0 x 0.5 cm	1.0 x 0.7 cm	1.0 x 0.8 cm	1.0 x 1.0 cm
Reusable	TWS 84-2607/35	TWS 84-2607/36	TWS 84-2607/37	TWS 84-2607/38	TWS 84-2607/39



Size	1.0 x 0.5 cm	1.0 x 0.7 cm	1.0 x 0.8 cm	1.0 x 1.0 cm	1.0 x 1.5 cm
Reusable	TWS 84-2607/40	TWS 84-2607/41	TWS 84-2607/42	TWS 84-2607/43	TWS 84-2607/44



Size	1.5 x 0.5 cm	1.5 x 0.7 cm	1.5 x 0.8 cm	1.5 x 1.0 cm	1.5 x 1.2 cm	1.5 x 1.5 cm
Reusable	TWS 84-2607/45	TWS 84-2607/46	TWS 84-2607/47	TWS 84-2607/48	TWS 84-2607/49	TWS 84-2607/50



Size	2.0 x 0.8 cm	2.0 x 1.0 cm	2.0 x 1.2 cm	2.0 x 1.3 cm	2.0 x 1.5 cm	2.0 x 2.0 cm
Reusable	TWS 84-2607/51	TWS 84-2607/52	TWS 84-2607/53	TWS 84-2607/54	TWS 84-2607/55	TWS 84-2607/56



TW Surgical

ELECTRODES (1.6 mm) SINGLE-USE



Blade Electrode
TWS 84-2613/07



Blade Electrode (Crv.)
TWS 84-2613/09



Ball Electrode (2 mm)
TWS 84-2613/11



Ball Electrode (Crv.) (2 mm)
TWS 84-2613/13



Ball Electrode (3 mm)
TWS 84-2613/14



Needle Electrode

Cat #	Description	Size
TWS 84-2607/57	Needle Electrode	5.0 cm
TWS 84-2607/58	Needle Electrode	6.0 cm
TWS 84-2607/59	Needle Electrode	7.0 cm
TWS 84-2607/60	Needle Electrode	11.0 cm



Needle Electrode (Crv.)

Cat #	Description	Size
TWS 84-2607/61	Needle Electrode (Crv.)	5.0 cm
TWS 84-2607/62	Needle Electrode (Crv.)	6.0 cm
TWS 84-2607/63	Needle Electrode (Crv.)	7.0 cm
TWS 84-2607/64	Needle Electrode (Crv.)	11.0 cm




Tungsten Needle Electrode

Cat #	Description	Size
TWS 84-2607/65	Tungsten Needle Electrode	5.0 cm
TWS 84-2607/66	Tungsten Needle Electrode	6.0 cm
TWS 84-2607/67	Tungsten Needle Electrode	7.0 cm
TWS 84-2607/68	Tungsten Needle Electrode	9.0 cm




Tungsten Needle Electrode (Crv.)


Cat #	Description	Size
TWS 84-2607/69	Tungsten Needle Electrode (Crv.)	5.0 cm
TWS 84-2607/70	Tungsten Needle Electrode (Crv.)	6.0 cm
TWS 84-2607/71	Tungsten Needle Electrode (Crv.)	7.0 cm
TWS 84-2607/72	Tungsten Needle Electrode (Crv.)	9.0 cm



Ball Electrode (Crv.) (3 mm)
TWS 84-2613/08




Ball Electrode (4 mm)
TWS 84-2613/10



Ball Electrode (Crv.) (4 mm)
TWS 84-2613/12



Ball Electrode (5 mm)
TWS 84-2613/15



Ball Electrode (Crv.) (5 mm)
TWS 84-2613/16



Circle Loop Electrode

Cat #	Description	Size
TWS 84-2607/73	Circle Loop Electrode	4.0 mm
TWS 84-2607/74	Circle Loop Electrode	6.0 mm
TWS 84-2607/75	Circle Loop Electrode	8.0 mm
TWS 84-2607/76	Circle Loop Electrode	10.0 mm



Diamond Loop Electrode

Cat #	Description	Size
TWS 84-2607/77	Diamond Loop Electrode	4.0 mm
TWS 84-2607/78	Diamond Loop Electrode	6.0 mm
TWS 84-2607/79	Diamond Loop Electrode	8.0 mm



Triangle Loop Electrode

Cat #	Description	Size
TWS 84-2607/80	Triangle Loop Electrode	4.0 mm
TWS 84-2607/81	Triangle Loop Electrode	6.0 mm
TWS 84-2607/82	Triangle Loop Electrode	8.0 mm



Oval Loop Electrode


Cat #	Description	Size
TWS 84-2607/83	Oval Loop Electrode	4.0 mm
TWS 84-2607/84	Oval Loop Electrode	6.0 mm
TWS 84-2607/85	Oval Loop Electrode	8.0 mm

PATIENT RETURN PLATES



TW Surgical

Size: 21 x 17 cm



Patient Return Plate (Reusable)
TWS 84-2613/17 (5m cable)

Size: 21 x 10 cm



Bipolar Patient Plate (Single-Use)
TWS 84-2613/19 (Adult Size)
TWS 84-2613/20 (Child Size)

Size: 21 x 10 cm



Monopolar Patient Plate (Single-Use)
TWS 84-2613/21 (Adult Size)
TWS 84-2613/22 (Child Size)

Size: 16 x 12 cm



S.S. Patient Return Plate
TWS 84-2613/18



Patient Return Cable
TWS 84-2613/23 (5m cable)
TWS 84-2613/24 (3m cable)



Patient Return Cable
TWS 84-2613/25 (5m cable)
TWS 84-2613/26 (3m cable)




Patient Return Cable (Reusable)
TWS 84-2613/27 (5m cable)
TWS 84-2613/28 (3m cable)



Patient Return Cable (Reusable)
TWS 84-2613/29 (5m cable)
TWS 84-2613/30 (3m cable)

Size: 13 x 9 cm



Patient Return Plate
(3m cable)
TWS 84-2613/31



TW Surgical


PATIENT RETURN PLATES SILICONE

Size: 20 x 15 cm




Monopolar
TWS 84-2613/32


Size: 20 x 15 cm



Bipolar
TWS 84-2613/33




Erbe
TWS 84-2613/34



Valleylab
TWS 84-2613/35

Size: 25 x 15 cm




Silicone Return Plate
Erbe Connector
TWS 84-2613/36 (5m cable)

Size: 25 x 15 cm



Erbe Connector
TWS 84-2613/37 (5m cable)



Silicone Patient Return
Valleylab Connector
TWS 84-2613/38 (5m cable)

LAPAROSCOPIC ELECTRODES



TW Surgical

Without Suction



L-Hook Electrode (36 cm)
TWS 84-2613/39 (Reusable)



L-Hook Electrode (36 cm)
TWS 84-2613/40 (Single-Use)



J-Hook Electrode (36 cm)
TWS 84-2613/45 (Reusable)



J-Hook Electrode (36 cm)
TWS 84-2613/41 (Single-Use)



Needle Electrode (36 cm)
TWS 84-2613/46 (Reusable)



Needle Electrode (36 cm)
TWS 84-2613/42 (Single-Use)



Spatula Electrode (36 cm)
TWS 84-2613/47 (Reusable)



Spatula Electrode (36 cm)
TWS 84-2613/43 (Single-Use)



Ball Electrode (36 cm)
TWS 84-2613/48 (Reusable)



Ball Electrode (36 cm)
TWS 84-2613/44 (Single-Use)

With Suction



L-Hook Electrode (36 cm)
TWS 84-2613/49 (Reusable)



J-Hook Electrode (36 cm)
TWS 84-2613/50 (Reusable)



Needle Electrode (36 cm)
TWS 84-2613/51 (Reusable)



Spatula Electrode (36 cm)
TWS 84-2613/52 (Reusable)



Ball Electrode (36 cm)
TWS 84-2613/53 (Reusable)





TW Surgical

ARTHROSCOPIC ELECTRODES



Dagger (12 cm)
TWS 84-2613/54



Dagger (18 cm)
TWS 84-2613/55



Hook (12 cm)
TWS 84-2613/56




Hook (18 cm)
TWS 84-2613/57



Hockey (12 cm)
TWS 84-2613/58



Hockey (18 cm)
TWS 84-2613/59



Nasal Coagulation Cannula
TWS 84-2613/60 (10 cm) (Reusable)
TWS 84-2613/61 (26 cm) (Reusable)

Nasal Coagulation Cannula
TWS 84-2613/62 (10 cm) (Single-Use)
TWS 84-2613/63 (26 cm) (Single-Use)



Frazier Suction Tube 11 French (5 1/4", 13.5 cm)
TWS 84-2613/64



Frazier Suction Tube 11 French (8", 20.5 cm)
TWS 84-2613/65




Frazier Suction Tube 11 French (16", 40.6 cm)
TWS 84-2613/66



Suction Coagulator (Single-Use)
TWS 84-2613/67 (3m cable)

ELECTROSURGICAL INSTRUMENTS FOR GYNECOLOGY



Disposable Tube



Graves Vaginal Speculum
With Fixed or Disposable Smoke Tube

Fixed Tube Cat #	Disposable Tube Cat #	Size	Length	Width
TWS 84-2607/86	TWS 84-2607/87	X-Large	150 mm	40 mm
TWS 84-2607/89	TWS 84-2607/88	Large	115 mm	35 mm
TWS 84-2607/90	TWS 84-2607/91	Medium	95 mm	35 mm
TWS 84-2607/92	TWS 84-2607/93	Small	75 mm	20 mm

VU-More (Graves) Speculum
With Fixed or Disposable Smoke Tube

Fixed Tube Cat #	Disposable Tube Cat #	Size	Length	Width
TWS 84-2607/95	TWS 84-2607/94	Large	115 mm	35 mm
TWS 84-2607/97	TWS 84-2607/96	Medium	100 mm	35 mm
TWS 84-2607/99	TWS 84-2607/98	Small	80 mm	30 mm



VU-Max Graves Speculum
With Fixed or Disposable Smoke Tube

Fixed Tube Cat #	Disposable Tube Cat #	Size	Length	Width
TWS 84-2608/00	TWS 84-2608/01	X-Large	135 mm	37 mm
TWS 84-2608/02	TWS 84-2608/04	Large	110 mm	35 mm
TWS 84-2608/03	TWS 84-2608/05	Medium	100 mm	35 mm



Pederson Speculum
With Fixed or Disposable Smoke Tube

Fixed Tube Cat #	Disposable Tube Cat #	Size	Length	Width
TWS 84-2608/06	TWS 84-2608/10	X-Large	165 mm	25 mm
TWS 84-2608/07	TWS 84-2608/11	Large	115 mm	25 mm
TWS 84-2608/08	TWS 84-2608/12	Medium	100 mm	22 mm
TWS 84-2608/09	TWS 84-2608/13	Small	75 mm	14 mm



VU-Max Pederson Speculum
With Fixed or Disposable Smoke Tube

Fixed Tube Cat #	Disposable Tube Cat #	Size	Length	Width
TWS 84-2608/14	TWS 84-2608/16	Large	115 mm	25 mm
TWS 84-2608/15	TWS 84-2608/17	Medium	100 mm	22 mm



Chinese Pattern Speculum
With Smoke Tube

Cat #	Size	Length	Width
TWS 84-2608/18	Medium	95 mm	28 mm
TWS 84-2608/19	Small	70 mm	35 mm



TW Surgical

ELECTROSURGICAL INSTRUMENTS FOR GYNECOLOGY



Cusco Speculum
With Smoke Tube

Cat #	Size	Length	Width
TWS 84-2608/20	Large	95 mm	36 mm
TWS 84-2608/21	Medium	85 mm	35 mm
TWS 84-2608/22	Small	75 mm	32 mm



Collin Speculum
With Smoke Tube

Cat #	Size	Length	Width
TWS 84-2608/23	Large	110 mm	40 mm
TWS 84-2608/24	Medium	100 mm	35 mm
TWS 84-2608/25	Small	95 mm	30 mm



Weisman Speculum (Right Open)
With Smoke Tube

Cat #	Size	Length	Width
TWS 84-2608/26	Large	115 mm	35 mm
TWS 84-2608/27	Medium	95 mm	35 mm



Weisman Speculum (Left Open)
With Smoke Tube

Cat #	Size	Length	Width
TWS 84-2608/28	Large	115 mm	35 mm
TWS 84-2608/29	Medium	95 mm	35 mm



Weisman Pederson Speculum 90° (Right Open)
With Smoke Tube

Cat #	Size	Length	Width
TWS 84-2608/30	Large	115 mm	25 mm
TWS 84-2608/31	Medium	100 mm	25 mm



Weisman Pederson Speculum 90° (Left Open)
With Smoke Tube

Cat #	Size	Length	Width
TWS 84-2608/32	Large	115 mm	25 mm
TWS 84-2608/33	Medium	100 mm	25 mm

ELECTROSURGICAL INSTRUMENTS FOR GYNECOLOGY



Weisman Speculum 90° (Right Open)
With Smoke Tube

Cat #	Size	Length	Width
TWS 84-2608/34	X-Large	135 mm	40 mm
TWS 84-2608/35	Large	115 mm	35 mm
TWS 84-2608/36	Medium	95 mm	35 mm



Weisman Speculum 90° (Left Open)
With Smoke Tube

Cat #	Size	Length	Width
TWS 84-2608/37	X-Large	135 mm	40 mm
TWS 84-2608/38	Large	115 mm	35 mm
TWS 84-2608/39	Medium	95 mm	35 mm



CUSCO Rocket Speculum
With Smoke Tube

Cat #	Size	Length	Width
TWS 84-2608/40	Virgin	86 mm	16 mm
TWS 84-2608/41	Small	95 mm	27 mm
TWS 84-2608/42	Small Long	107 mm	27 mm
TWS 84-2608/43	Small Extra Long	125 mm	27 mm
TWS 84-2608/44	Medium	95 mm	30 mm
TWS 84-2608/45	Medium Long	110 mm	30 mm
TWS 84-2608/46	Medium Extra Long	125 mm	30 mm
TWS 84-2608/47	Large	110 mm	35 mm
TWS 84-2608/48	Large Long	120 mm	35 mm
TWS 84-2608/49	Extra Large Extra Long	130 mm	36 mm



Winterton Speculum
With Smoke Tube

Cat #	Size	Length	Width
TWS 84-2608/50	Large	125 mm	28 mm



TW Surgical

ELECTROSURGICAL INSTRUMENTS FOR GYNECOLOGY

Dressing Forceps

TWS 84-2613/68(21 cm)

TWS 84-2613/69(25 cm)



Dressing Forceps (Fully Coated)

TWS 84-2613/70 (21 cm)

TWS 84-2613/71 (25 cm)



Tissue Forceps (Exposed Tip)

TWS 84-2613/72 (21 cm)

TWS 84-2613/73 (25 cm)



Tissue Forceps

TWS 84-2613/74 (21 cm)

TWS 84-2613/75 (25 cm)



Iris Hook

TWS 84-2613/76 (23 cm)



Emmet Hook

TWS 84-2613/77 (23 cm)



Straight Hook

TWS 84-2613/78 (25 cm)



Two-Prong Hook

TWS 84-2613/79 (25 cm)



Three-Prong Hook

TWS 84-2613/80 (25 cm)



ELECTROSURGICAL INSTRUMENTS FOR GYNECOLOGY



TW Surgical



Lateral Vaginal Retractor
TWS 84-2613/81



Lateral Vaginal Retractor
Wide Open
TWS 84-2613/82



Lateral Vaginal Retractor
1/2" Open
TWS 84-2613/83



Lateral Vaginal Retractor
1" Open
TWS 84-2613/84



Gelpi Retractor
TWS 84-2613/88



Rigby Retractor
TWS 84-2613/89



Kogan Endocervical Retractor



3mm TWS 84-2613/85



4mm TWS 84-2613/86



5mm TWS 84-2613/87


Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical


ELECTROSURGICAL INSTRUMENTS FOR GYNECOLOGY




Schroeder Tenaculum (24 cm)
TWS 84-2613/90




Schroeder Tenaculum (24 cm)
TWS 84-2613/91



Single Tooth Tenaculum (24 cm)
TWS 84-2613/92



Ring Forceps (22 cm)
TWS 84-2613/93



Campion Forceps (22 cm)
TWS 84-2613/94

ELECTROSURGICAL INSTRUMENTS FOR GYNECOLOGY



Sims Speculum (Large)
TWS 84-2613/95




Sims Speculum (Medium)
TWS 84-2613/96




Sims Speculum (Small)
TWS 84-2613/97




Auvar Weighted Speculum
TWS 84-2613/98



US Army Retractor (Set of 2) 8 1/2" (21.5 cm)
TWS 84-2613/99



Farabeuf Retractor (Set of 2) 6" (15 cm)
TWS 84-2614/00





TW Surgical

RETRACTORS AND SKIN HOOKS



Jager Lid Plate 4½" (11.0 cm)
TWS 84-2614/01



Desmarres Lid Retractor 5½" (14.0 cm)

Cat #	Size
TWS 84-2608/51	10 mm
TWS 84-2608/52	12 mm
TWS 84-2608/53	14 mm
TWS 84-2608/54	16 mm




Fomon Retractor - Double Ball End 6¼" (16.0 cm)
TWS 84-2614/02



Senn-Muller Retractor (Sharp) 6¼" (16.0 cm)
TWS 84-2614/03



Senn-Muller Retractor (Blunt) 6¼" (16.0 cm)
TWS 84-2614/04



Desmarres Lid Retractor 8/10 mm 5½" (14.0 cm)
TWS 84-2614/05




Double Eye Lid Retractor 5½" (14.0 cm)
TWS 84-2614/06




Ragnell-Davis Retractor 6" (15.0 cm)
TWS 84-2614/07

STERILIZATION TRAYS




Sterilization Tray (Stainless Steel)
TWS 84-2614/08




Sterilization Tray (Stainless Steel)
TWS 84-2614/09

Cat #	Description	Length	Width	Height
TWS 84-2614/08	Sterilization Tray (3 pcs.) Stainless Steel	30.5 cm	8 cm	4.5 cm
TWS 84-2614/09	Sterilization Tray (6 pcs.) Stainless Steel	30.5 cm	15 cm	4.5 cm



Sterilization Tray (Aluminum)
TWS 84-2614/10



Sterilization Tray (Aluminum)
TWS 84-2614/11

Cat #	Description	Length	Width	Height
TWS 84-2614/10	Sterilization Tray (2 pcs.) Aluminum	25 cm	6.5 cm	5 cm
TWS 84-2614/11	Sterilization Tray (5 pcs.) Aluminum	30 cm	14 cm	5 cm



TW
Surgical

Treasure Wright (U.K.) Ltd.

**Address: 128-160 , City Road,
London EC1V 2NX. United Kingdom**

www.twrightmedical.com

info@twrightmedical.com



TW Surgical



**UK
CA CE**