

TW Surgical Esmarch Tourniquet Bandages - Latex-free

TW Surgical Esmarch Tourniquet Bandages are soft and very comfortable to use, hence preferred by medical professionals for their patients.

A popular feature of these bandages is that they can be sterilized by Steam, Gamma or ETO.



Utility: TW Surgical Esmarch Tourniquet Bandages are used to control bleeding by applying it around a limb in such a way that blood is expelled from it. TW Surgical Esmarch Tourniquet Bandages is ideal for orthopedic blood-less surgery.



Specifications of TW Surgical Esmarch Tourniquet Bandages Standard Sizes:

- Width: 60 mm, 75 mm, 100 mm, 120 mm, 150 mm (2.5"/ 3"/ 4"/ 5"/ 6").
- Length: 6 Feet, 8 Feet, 10 Feet, 12 Feet



Surface:

- The latex-free esmarch bandages are available with a Smooth or Textured surface.



TW Surgical Esmarch Tourniquet Bandages Latex-free

Material: Rubber and Silicone based

		
Conforms to Biocompatibility Standards	Excellent Stretch Recovery	Plain or Textured Surface
		
Silicone Treated Surface	Range of Wide and Length Available	Can be Sterilized –Steam, Gamma, ETO

The **TW Surgical Esmarch Tourniquet Bandages** can also be manufactured in the different sizes, surface & supply form of the latex-free tourniquets.

Order Information:

TWS23060-06	Tourniquet Bandages 60mm (2.5”) wide x 6 Feet Long with Smooth surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
TWS23060-08	Tourniquet Bandages 60mm (2.5”) wide x 8 Feet Long with Smooth surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
TWS23060-10	Tourniquet Bandages 60mm (2.5”) wide x 10 Feet Long with Smooth surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
TWS23060-12	Tourniquet Bandages 60mm (2.5”) wide x 12 Feet Long with Smooth surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
TWS23075-06	Tourniquet Bandages 75mm (3”) wide x 6 Feet Long with Smooth surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
TWS23075-08	Tourniquet Bandages 75mm (3”) wide x 8 Feet Long with Smooth surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
TWS23075-10	Tourniquet Bandages 75mm (3”) wide x 10 Feet Long with Smooth surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
TWS23075-12	Tourniquet Bandages 75mm (3”) wide x 12 Feet Long with Smooth surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma




TW Surgical Esmarch Tourniquet Bandages Latex-free

Order Information:

23075-20	Esmarch Bandages 75mm (3") wide x 2M Long with Smooth surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
23075-20T	Esmarch Bandages 75mm (3") wide x 2M Long with Textured surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
23075-30	Esmarch Bandages 75mm (3") wide x 3M Long with Smooth surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
23075-30T	Esmarch Bandages 75mm (3") wide x 3M Long with Textured surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
23075-40	Esmarch Bandages 75mm (3") wide x 4M Long with Smooth surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
23075-40T	Esmarch Bandages 75mm (3") wide x 4M Long with Textured surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
23075-50	Esmarch Bandages 75mm (3") wide x 5M Long with Smooth surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
23075-50T	Esmarch Bandages 75mm (3") wide x 5M Long with Textured surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma


23100-20	Esmarch Bandages 100mm (4") wide x 2M Long with Smooth surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
23100-20T	Esmarch Bandages 100mm (4") wide x 2M Long with Textured surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
23100-30	Esmarch Bandages 100mm (4") wide x 3M Long with Smooth surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
23100-30T	Esmarch Bandages 100mm (4") wide x 3M Long with Textured surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
23100-40	Esmarch Bandages 100mm (4") wide x 4M Long with Smooth surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
23100-40T	Esmarch Bandages 100mm (4") wide x 4M Long with Textured surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
23100-50	Esmarch Bandages 100mm (4") wide x 5M Long with Smooth surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
23100-50T	Esmarch Bandages 100mm (4") wide x 5M Long with Textured surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma



TW Surgical Esmarch Tourniquet Bandages Latex-free

Order Information:

TWS23120-06	Tourniquet Bandages 120mm (5") wide x 6 Feet Long with Smooth surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
TWS23120-08	Tourniquet Bandages 120mm (5") wide x 8 Feet Long with Smooth surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
TWS23120-10	Tourniquet Bandages 120mm (5") wide x 10 Feet Long with Smooth surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
TWS23120-12	Tourniquet Bandages 120mm (5") wide x 12 Feet Long with Smooth surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
TWS23150-06	Tourniquet Bandages 150mm (6") wide x 6 Feet Long with Smooth surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
TWS23150-08	Tourniquet Bandages 150mm (6") wide x 8 Feet Long with Smooth surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
TWS23150-10	Tourniquet Bandages 150mm (6") wide x 10 Feet Long with Smooth surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma
TWS23150-12	Tourniquet Bandages 150mm (6") wide x 12 Feet Long with Smooth surface, Latex-free Washable and Can be sterilized by Steam, ETO, or Gamma

 Conforms to Biocompatibility Standards	 Excellent Stretch Recovery	 Plain or Textured Surface
 Silicone Treated Surface	 Range of Wide and Length Available	 Can be Sterilized –Steam, Gamma, ETO

