


## Hand and Foot Instruments –Kleinert-Ragnell Retractor



TWS 35-6140/0310	<p>Kleinert Ragnell Double End Retractor          140 mm Long, Oval Handle, Smooth Curved Blades          One End: 3 mm Wide and 10 mm Reach, Round Blades          One End: 2 mm Wide blade and 5 mm Reach, Round Blades          Stainless Steel</p>
TWS 35-6140/0515	<p>Kleinert Ragnell Double End Retractor          140 mm Long, Oval Handle, Smooth Curved Blades          One End: 5 mm Wide and 15 mm Reach, Round Blades          One End: 3 mm Wide blade and 8 mm Reach, Round Blades          Stainless Steel</p>
TWS 35-6140/0820	<p>Kleinert Ragnell Double End Retractor          140 mm Long, Oval Handle, Smooth Curved Blades          One End: 8 mm Wide and 20 mm Reach, Round Blades          One End: 5 mm Wide blade and 11 mm Reach, Round Blades          Stainless Steel</p>