

# General Instrument – Louisville Elevators

Louisville Elevators will assist in lifting the tissue flaps.



TWS 52-1700/2434	Louisville Double End Elevators, Hexagonal Handle, Slightly Curved Tips, 2mm x 4mm and 3mm x 4mm Tips 170mm Long, Material: Stainless Steel
TWS 52-1700/3545	Louisville Double End Elevators, Hexagonal Handle, Slightly Curved Tips, 3mm x 5mm and 4mm x 5mm Tips 170mm Long, Material: Stainless Steel
TWS 52-1700/5869	Louisville Double End Elevators, Hexagonal Handle, Slightly Curved Tips, 5mm x 8mm and 6mm x 9mm Tips 170mm Long, Material: Stainless Steel



TWS 52-1701/3537	Louisville Double End Elevators, Hexagonal Handle, Slightly Curved Tip and Pointed Tip 3mm x 5mm Curved Tip and 3mm Pointed Tip, 170mm Long Material: Stainless Steel
TWS 52-1700/5759	Louisville Double End Elevators, Hexagonal Handle, Slightly Curved Tip and Pointed Tip 5mm x 7mm Curved Tip and 5mm Pointed Tip, 170mm Long Material: Stainless Steel

**Treasure Wright (U.K.) Ltd.**

128-160 City Road, London, EC1V 2NX, UK  
info@twrightmedical.com