

Orthopedic Instrument – TWS Spatula

TWS Spatula is useful in exploration and handling of delicate nerves. This instrument features a slim handle with an elongated spatula blade portion that widens at that end and features a slightly curvature. These details make this instrument well suited for use in delicate tissue and prevent damage or trauma to surrounding areas.

Prevents Tearing: The wide, flat shape of a spatula distributes pressure evenly.

Protects Blood Supply: Peeling the tissue off the bone in one clean sheet keeps those tiny blood vessels intact.

Bone is a Safe Shield: Surgeon can push firmly against it with the blunt spatula without any risk of cutting into deeper, dangerous areas.



TWS 15-2300/08TM	TWS Spatula Curved Up 8mm Width Delicate Round Tip (Smooth and Blunt), 23cm Long Material: Stainless Steel
TWS 15-2301/08TM	TWS Spatula Curved Up 8mm Width Delicate Round Tip (Flap and Sharp), 23cm Long Material: Stainless Steel