

Orthopaedic General Instrument –Stille Bone Cutting Forceps

Stille Bone cutter is a double-action instrument that provides greater pressure on the bones without as much effort by the Surgeon. The jaws of this instrument have sharp blades that square when closed. Stille Bone cutter specifically most useful in orthopedic surgeries that require cutting or clipping of smaller, lower density bones. The Bone Cutter handle curves slightly at the joints, and then again more strongly and in the opposite direction at the jaws. This product is re-usable and are designed to perform specific functions such as cutting, clamping, dissecting, and other similar procedures.



TWS 50-1530/27	Stille Horsley Bone Cutting Forceps, Double Action 27cm Length, Angled Jaws, Stainless Steel
TWS 50-1530/32	Stille Horsley Bone Cutting Forceps, Double Action 32cm Length, Angled Jaws, Stainless Steel
TWS 50-1530/36	Stille Horsley Bone Cutting Forceps, Double Action 36cm Length, Angled Jaws, Stainless Steel
TWS 50-1590/27	Stille Horsley Bone Cutting Forceps, Double Action 27cm Length, Angled 90° Jaws, Stainless Steel
TWS 50-1590/32	Stille Horsley Bone Cutting Forceps, Double Action 32cm Length, Angled 90° Jaws, Stainless Steel
TWS 50-1590/36	Stille Horsley Bone Cutting Forceps, Double Action 36cm Length, Angled 90° Jaws, Stainless Steel