



TW
Surgical

SUCTION TUBE



TW Surgical

Treasure Wright (U.K.) Ltd.



www.twrightmedical.com



info@twrightmedical.com

UK
CA

CE



SUCTION TUBE

TW Surgical supplies suction instrument for neurosurgery to very fine suction instruments for narrow approaches, offers a portfolio of suction cannulas to meet your needs.

- Extremely fine suction control due to tear drop shaped suction interrupter
- Individual adaption of the suction tube
- Conical suction tube prevents plugging
- Atraumatic distal tip
- Conical suction tube
- Tailor-made size



TWS200-4712/33	Round Hollow Suction Tube For Suction With Suction Hose Length: 20cm, Wall Thickness: 1mm, Diameter: 3.3mm Oval Shape Head: 4.7mm Height x 1.5mm Width x 20mm Long Material: Stainless Steel
TWS200-4720/33	Round Hollow Suction Tube For Suction With Suction Hose Length: 20cm, Wall Thickness: 1mm, Diameter: 3.3mm Oval Shape Head: 4.7mm Height x 2.0mm Width x 20mm Long Material: Stainless Steel
TWS200-4725/33	Round Hollow Suction Tube For Suction With Suction Hose Length: 20cm, Wall Thickness: 1mm, Diameter: 3.3mm Oval Shape Head: 4.7mm Height x 2.5mm Width x 20mm Long Material: Stainless Steel
TWS200-4730/33	Round Hollow Suction Tube For Suction With Suction Hose Length: 20cm, Wall Thickness: 5mm, Diameter: 3.3mm Oval Shape Head: 7mm Height x 3.0mm Width x 20mm Long Material: Stainless Steel
TWS200-4740/33	Round Hollow Suction Tube For Suction With Suction Hose Length: 20cm, Wall Thickness: 5mm, Diameter: 3.3mm Oval Shape Head: 7mm Height x 4.0mm Width x 20mm Long Material: Stainless Steel



TW
Surgical

Treasure Wright (U.K.) Ltd.

**Address: 128-160 , City Road,
London EC1V 2NX. United Kingdom**

www.twrightmedical.com

info@twrightmedical.com