



TW
Surgical

ASPIRATION NEEDLES



**UK
CA** **CE**



TW Surgical

Treasure Wright (U.K.) Ltd.



www.twrightmedical.com



info@twrightmedical.com

Reusable Long Aspiration Needles



Lancet point

- Nickel plated big hub with super smooth lapped finished cannula
- Available from G 10 to G 24 in any length combinations
- Pack of 10 pieces

- **Intended use** : For Laboratory use.



TW Surgical

CODE	SIZE
TWS 12381	10G x 2" (50mm x 3.3mm)
TWS 12382	10G x 2.1/2" (60mm x 3.3mm)
TWS 12383	10G x 3" (75mm x 3.3mm)
TWS 12384	10G x 4" (100mm x 3.3mm)
TWS 12385	10G x 5" (130mm x 3.3mm)
TWS 12386	10G x 6" (150mm x 3.3mm)
TWS 12387	10G x 8" (200mm x 3.3mm)
TWS 12388	11G x 2" (50mm x 3.00mm)
TWS 12389	11G x 2.1/2" (60mm x 3.00mm)
TWS 12390	11G x 3" (75mm x 3.00mm)
TWS 12391	11G x 4" (100mm x 3.00mm)
TWS 12392	11G x 5" (130mm x 3.00mm)
TWS 12393	11G x 6" (150mm x 3.00mm)
TWS 12394	12G x 2" (50mm x 2.8mm)
TWS 12395	12G x 2.1/2" (60mm x 2.8mm)
TWS 12397	12G x 3" (75mm x 2.8mm)
TWS 12396	12G x 4" (100mm x 2.8mm)
TWS 12398	12G x 5" (130mm x 2.8mm)
TWS 12399	12G x 6" (150mm x 2.8mm)
TWS 12400	12G x 8" (200mm x 2.8mm)
TWS 12401	12G x 10" (250mm x 2.8mm)
TWS 12402	13G x 2" (50mm x 2.5mm)
TWS 12403	13G x 2.1/2" (60mm x 2.5mm)
TWS 12404	13G x 3" (75mm x 2.5mm)
TWS 12405	13G x 4" (100mm x 2.5mm)
TWS 12406	13G x 5" (130mm x 2.5mm)
TWS 12407	13G x 6" (150mm x 2.5mm)
TWS 12408	13G x 8" (200mm x 2.5mm)
TWS 12409	13G x 10" (250mm x 2.5mm)
TWS 12410	14G x 2" (50mm x 2.00mm)
TWS 12411	14G x 2.1/2" (60mm x 2.00mm)
TWS 12412	14G x 3" (75mm x 2.00mm)
TWS 12413	14G x 4" (100mm x 2.00mm)
TWS 12414	14G x 5" (130mm x 2.00mm)
TWS 12415	14G x 6" (150mm x 2.00mm)

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical

CODE	SIZE
TWS 12416	14G x 8" (200mm x 2.00mm)
TWS 12417	14G x 10" (250mm x 2.00mm)
TWS 12418	15G x 2" (50mm x 1.8mm)
TWS 12419	15G x 2.1/2" (60mm x 1.8mm)
TWS 12420	15G x 3" (75mm x 1.8mm)
TWS 12421	15G x 4" (100mm x 1.8mm)
TWS 12422	15G x 5" (130mm x 1.8mm)
TWS 12423	15G x 6" (150mm x 1.8mm)
TWS 12424	15G x 8" (200mm x 1.8mm)
TWS 12425	15G x 10" (250mm x 1.8mm)
TWS 12426	16G x 2" (50mm x 1.6mm)
TWS 12427	16G x 2.1/2" (60mm x 1.6mm)
TWS 12428	16G x 3" (75mm x 1.6mm)
TWS 12429	16G x 4" (100mm x 1.6mm)
TWS 12430	16G x 5" (130mm x 1.6mm)
TWS 12431	16G x 6" (150mm x 1.6mm)
TWS 12432	16G x 8" (200mm x 1.6mm)
TWS 12433	16G x 10" (250mm x 1.6mm)
TWS 12434	16G x 12" (300mm x 1.6mm)
TWS 12435	17G x 2" (50mm x 1.5mm)
TWS 12436	17G x 2.1/2" (60mm x 1.5mm)
TWS 12437	17G x 3" (75mm x 1.5mm)
TWS 12438	17G x 4" (100mm x 1.5mm)
TWS 12439	17G x 5" (130mm x 1.5mm)
TWS 12440	17G x 6" (150mm x 1.5mm)
TWS 12441	17G x 8" (200mm x 1.5mm)
TWS 12442	17G x 10" (250mm x 1.5mm)
TWS 12443	17G x 12" (250mm x 1.5mm)
TWS 12444	18G x 2" (50mm x 1.3mm)
TWS 12445	18G x 2.1/2" (60mm x 1.3mm)
TWS 12446	18G x 3" (75mm x 1.3mm)

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical

CODE	SIZE
TWS 12671	18G x 4" (100mm x 1.3mm)
TWS 12670	18G x 5" (130mm x 1.3mm)
TWS 12672	18G x 6" (150mm x 1.3mm)
TWS 12673	18G x 8" (200mm x 1.3mm)
TWS 12674	18G x 10" (250mm x 1.3mm)
TWS 12675	18G x 12" (300mm x 1.3mm)
TWS 12676	19G x 2" (50mm x 1.00mm)
TWS 12677	19G x 2.1/2" (60mm x 1.00mm)
TWS 12678	19G x 3" (75mm x 1.00mm)
TWS 12679	19G x 4" (100mm x 1.00mm)
TWS 12680	19G x 5" (130mm x 1.00mm)
TWS 12681	19G x 6" (150mm x 1.00mm)
TWS 12682	19G x 8" (200mm x 1.00mm)
TWS 12683	19G x 10" (250mm x 1.00mm)
TWS 12684	19G x 12" (300mm x 1.00mm)
TWS 12685	20G x 2" (50mm x 0.9mm)
TWS 12686	20G x 2.1/2" (60mm x 0.9mm)
TWS 12687	20G x 3" (75mm x 0.9mm)
TWS 12688	20G x 4" (100mm x 0.9mm)
TWS 12689	20G x 5" (130mm x 0.9mm)
TWS 12690	20G x 6" (150mm x 0.9mm)
TWS 12691	20G x 8" (200mm x 0.9mm)
TWS 12692	20G x 10" (250mm x 0.9mm)
TWS 12693	20G x 12" (300mm x 0.9mm)
TWS 12694	21G x 2" (50mm x 0.8mm)
TWS 12695	21G x 2.1/2" (60mm x 0.8mm)
TWS 12696	21G x 3" (75mm x 0.8mm)
TWS 12697	21G x 4" (100mm x 0.8mm)
TWS 12698	21G x 5" (130mm x 0.8mm)
TWS 12699	21G x 6" (150mm x 0.8mm)
TWS 12700	21G x 8" (200mm x 0.8mm)
TWS 12701	21G x 10" (250mm x 0.8mm)
TWS 12702	22G x 2" (50mm x 0.7mm)
TWS 12703	22G x 2.1/2" (60mm x 0.7mm)

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical

CODE	SIZE
TWS 12704	22G x 3" (75mm x 0.7mm)
TWS 12705	22G x 4" (100mm x 0.7mm)
TWS 12706	22G x 5" (130mm x 0.7mm)
TWS 12707	22G x 6" (150mm x 0.7mm)
TWS 12708	22G x 8" (200mm x 0.7mm)
TWS 12709	22G x 10" (250mm x 0.7mm)
TWS 12710	23G x 2" (50mm x 0.6mm)
TWS 12711	23G x 2.1/2"(60mm x 0.6mm)
TWS 12712	23G x 3" (75mm x 0.6mm)
TWS 12713	23G x 4" (100mm x 0.6mm)
TWS 12714	23G x 5" (130mm x 0.6mm)
TWS 12715	23G x 6" (150mm x 0.6mm)
TWS 12716	23G x 8" (200mm x 0.6mm)
TWS 12717	23G x 10" (250mm x 0.6mm)
TWS 12718	24G x 2" (50mm x 0.57mm)
TWS 12719	24G x 2.1/2"(60mm x 0.57mm)
TWS 12720	24G x 3" (75mm x 0.57mm)
TWS 12721	24G x 4" (100mm x 0.57mm)
TWS 12722	24G x 5" (130mm x 0.57mm)
TWS 12723	24G x 6"(150mm x 0.57mm)
TWS 12724	24G x 8"(200mm x 0.57mm)
TWS 12725	24G x 10"(250mm x 0.57mm)

NOTE: Customized sizes are also available on request.

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com

Bone Marrow Biopsy- Jamshidi type



Pyramidal Trocar



Taper point

Bone marrow biopsy needle designed for puncture and extraction of histological tissue from the iliac crest.

- Sturdy metal device for better Insertion and Withdrawal.
- Three part Device: Needle, Stylet, Obturator.
- Ergonomic handle gripto distribute pressure in an optimal way applied by the operator.
- Device is offered in two tip styles
 1. Classical tapered tip with e)tended chisel stylet.
 2. Tapered diamond cut tip with pyramidal trocar tip.
- Marked obturator to measure and retrieve sample.
- Luer lock connection for syringe.

Bone Marrow Biopsy- Jamshidi type - Sizes

Code	Diameter	Length	Compatible Gauge
TWS 12727	4.0 mm	100/150 mm	8G
TWS 12726	3.5 mm	100/150 mm	10G
TWS 12728	3.0 mm	100 mm	11G
TWS 20000	2.4 mm	100 mm	13G

**Custom sizes available on request*

Spinal Needle



Whitacre Pencil point




Quincke Bevel

- Professional standard quality hub
- Accurately matched stylet wire for smooth gliding
- Made from AISI SS 304
- Available in Quincke Bevel/Pencil Point.
- Supplied unsterile in individual peel pouch pack.

Spinal Needle- Sizes

Product Code	Size
TWS 12729	14G x 3½"
TWS 12730	16G x 3½"
TWS 12731	18G x 3½"
TWS 12732	19G x 3½"
TWS 12733	20G x 3½"
TWS 12734	21G x 3½"
TWS 12735	22G x 3½"
TWS 12736	23G x 3½"
TWS 12737	24G x 3½"
TWS 12738	25G x 3½"
TWS 12739	26G x 3½"
TWS 12740	27G x 3½"

PELD Needle



Needle for Percutaneous Endoscopic Lumbar Discectomy

- Cannula made from AISI 304 Stainless Steel
- Perfect matching between cannula and stylet offers a safe penetration.
- Tip shape is designed to facilitate needle insertion smoothly.
- Lumen of the needle accepts a standard guide wire.
- Centimeter markings for measuring depth of insertion.
- Luer Lock connection for syringe.


PELD Needle- Sizes

Code	Diameter	Length	Compatible Gauge
TWS 12741	1.25 mm	200 mm	18G
TWS 12742	1.25 mm	250 mm	18G

All 18G needles are compatible with \varnothing 0.95 mm guide wire.

**Custom sizes available on request*

Extra Strong Guide Wire



Extra Strong Straightened Guide Wire for Introduction of cannulated instruments

- Blunt on both ends
- Different working lengths available
- Made from Surgical grade Stainless Steel 300 series
- Package of 5
- Custom length and diameter options available

Extra Strong Guide Wire- Sizes

Code	Diameter	Length	Compatible Gauge
TWS 12743	1 mm / 1.2mm	400mm	16G
TWS 12744	0.8 mm	300mm	18G
TWS 12745	0.8 mm	410mm	18G
TWS 12746	0.8mm	450mm	18G
TWS 12747	0.95mm	410mm	18G
TWS 12748	0.95mm	450mm	18G

**Custom sizes available on request*


Pudendal Nerve Block Needle & Guide (Paracervical)



- Two Needle Instrument with sharp point for better insertion.
- Packed in a single unit paper pouch.
- Supplied Un sterile and has to be pre cleaned and sterile before use.
- In case of re-use, inspect-clean & sterile

CODE	SIZE
TWS 12749	G20 x 55 mm
TWS 12750	G25 x 85 mm

Chiba Needle



Matched Bevel

- Fine Aspiration Chiba Needle.
- Packed in a single unit paper pouch.
- Supplied Un sterile and has to be pre cleaned and sterile before use.
- In case of re-use, inspect-clean & sterile

Chiba Needle- Sizes

Code	Size
TWS 12751	23G x 8"
TWS 12752	22G x 8"
TWS 12753	21G x 8"
TWS 12754	19G x 8"
TWS 12755	18G x 8"
TWS 12756	17G x 8"
TWS 12757	16G x 8"
TWS 12758	G8 to G20, Biopsy Needles


Franklin-Silverman's Biopsy Needle



- Silverman Biopsy for soft tissue biopsy with Franklin's modification.
- Made from medical grade SS AISI 304.
- Supplied in single piece packing and Unsterile.

Code	Size	Suitability
TWS 12759	G14	Adult use
TWS 12760	G16	Pediatric use

Bone Marrow Biopsy- Osgood



- For Bone Marrow Biopsy or transfusion.
- Has Round Luer lock hub.
- Sharp pointed cannula.
- Fitted stylet with matched bevel point.
- Supplied in paper pouch pack & Unsterile.

- **Intended use** : For treatment use.

Bone Marrow Biopsy- Osgood- Sizes

CODE	SIZE
TWS 12761	16G x 35mm
TWS 12762	16G x 50mm
TWS 12763	18G x 35mm
TWS 12764	18G x 50mm

NOTE: Customized sizes are also available on request

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com

Menghini Liver Biopsy Needle- Sizes

Code/Length	Length	Compatible Gauge
TWS 12765/L	70/100/120 mm	19G
TWS 12766/L	70/100/120 mm	18G
TWS 12767/L	70/100/120 mm	17G
TWS 12768/L	70/100/120 mm	16G

Vertebroplasty Needle Set



- Sturdy metal device with ergonomic handle for better insertion & withdrawal
- Comprises of 1. Beveled point insertion needle 2. Blunt point cement filling needle 3. Cement pusher rod.
- Centimetre graduation markings to guide penetration depth
- Luer lock connection for syringe
- Reusable and supplied Unsterile
- **Intended use** : For Vertebroplasty procedure.

Vertebroplasty Needle Set- Sizes

Product Code	Size
TWS 12769	10G x 150mm / 100mm
TWS 12770	11G x 150mm / 100mm
TWS 12771	13G x 150mm / 100mm

Bone Marrow Aspiration Needles



- Large spool Luer Hub for secure grip
- Sharp pointed cannula made from AISI S 304
- Adjustable stopper for advanced depth control
- Supplied Unsterile
- **Intended use** : For medical procedures.

Code/Length	Length	Compatible Gauge
TWS 12772	50 mm	14G
TWS 12773	50 mm	16G
TWS 12774	50 mm	18G

Abram's Lung Biopsy



- A three part fine biopsy instrument, with polished trocar point for easy insertion.
- Made of surgical grade stainless steel 304, having side port opening for cutting tissue.
- A spiral grooved hub with stillate adapter and blunt obturator is provided.
- Ergonomic grip is provided by two hexagonal fittings.
- Suitable for Luer tip syringes.
- Reusable, supplied unsterile and can be steam sterilized and reused.

Code	Compatible Gauge
TWS 12775	8G
TWS 12776	10G
TWS 12777	11G

Epidural Tuohy Needle



- Epidural needle with Tuohy's point.
- Available with or without stylet wire, fitted with adapter.
- Matched point bevel grinding.
- Reusable – All metal construction.
- Can be reused by autoclaving / sterilization.
- Available in single piece pre-sterile, ready to use pack.

Epidural Tuohy Needle- Sizes

Product Code	Size
TWS 12778E	14G x 3" (76mm)
TWS 12779E	14G x 3.1/2" (89mm)
TWS 12780E	16G x 3" (76mm)
TWS 12781E	16G x 3.1/2" (89mm)
TWS 12782E	16G x 4.1/2" (115mm)
TWS 12783E	16G x 5" (130mm)
TWS 12784E	16G x 7" (180mm)
TWS 12785E	17G x 3" (76mm)
TWS 12786E	17G x 3.1/2" (89mm)
TWS 12787E	17G x 4" (100mm)
TWS 12788E	17G x 7" (180mm)
TWS 12789E	18G x 3" (76mm)
TWS 12790E	18G x 3.1/2" (89mm)
TWS 12791E	19G x 3" (76mm)
TWS 12792E	20G x 3" (76mm)
TWS 12793E	21G x 3" (76mm)
TWS 12794E	22G x 3" (76mm)
TWS 12795E	23G x 3" (76mm)
TWS 12796E	24G x 3" (76mm)


Cope's Pleural Biopsy Needle



- Sturdy three part needle set to perform pleural biopsy
- Made from medical grade SS 304
- Can be steamed and sterilized
- Connection for Luer slip syringes
- Outer needle 11G size

Code	Compatible Gauge
TWS 12797	11G x 75mm

Menghini Liver Biopsy Needle



- A five part fine biopsy instrument, with pencil trocar point for easy insertion.
- Made of surgical grade stainless steel tubing
- Three Stillates Provided :
 - 1.Pencil Point
 - 2.Rounded Head
 - 3.Biopsy Stop
- Suitable for luer tip / luer lock syringes.
- Reusable, supplied unsterile.
- Can be steam sterilized and reused.

Ophthalmic Needle



- Blunt point needle with Stainless Steel AISI 304 cannula
 - Standard and professional quality Luer lock Hub
 - Precision tip for controlled flow
 - Nickel plated brass Hub
 - Sold in a pack of 12 pieces
-
- **Intended use** : For Ophthalmic use.

Ophthalmic Needle

Product Code	Size
TWS 12798	14G
TWS 12799	16G
TWS 12800	18G
TWS 12801	19G
TWS 12802	20G
TWS 12804	21G
TWS 12805	22G
TWS 12806	23G
TWS 12807	24G
TWS 12808	25G
TWS 12809	26G
TWS 12810	27G
TWS 12811	28G
TWS 12812	30G


NOTE: Needle length will be provided as per customer's request.

Dental Irrigation Needle



- Luer lock fitting, Cannula made of medical grade stainless steel.
- Available in 26G x 40 mm
- For irrigation of root canal
- Autoclavable
- Double Vent of 0.25 mm each
- Closed Rounded Tip
- 12pcs in Plastic Case

Liposuction Cannula



- Used for laser assisted lipo suction.
- Cannula made of medical grade stainless steel - 304
- Handle made of anodized aluminum
- Chemical and corrosion resistant.
- Protective wooden box packing.
- Best suitable for laser fibers.

Hand Pcs set for Surgical Lasers



- Cannula made of medical grade stainless steel - 304
- Handle made of anodized aluminum
- Chemical and corrosion resistant.
- Protective wooden box packing.

ENT Suction Cannula



- The cannula is made from SS 304.
- Styles available :
 - Rounded blunt tip / bent / bent distance from hub is 20mm for bigger gauges & 10mm for thin gauges..
- Packed in a single unit paper pouch.
- Supplied Un sterile.

Code	Size	Length
TWS 12813E	G10, G12, G14	150 mm
TWS 12814E	G14, G17, G19, G22, G26	70 mm

Tru-cut Biopsy



Designed for manual capture of high quality tissue samples with minimal trauma to the patient.

- All metal part biopsy, Aluminum anodized, to make it resistance to all chemicals, except caustic soda.
- Water and whether proof.
- A light weight instrument with good handle grip.
- A scalpel sharp surgical steel cutting edge, retrieves clearly cut high quality specimens.
- The tru-cut also offers a 15mm specimen notch.
- Can be steam sterilized and reused.

Code	Size	Suitability
TWS 12815	G14 x 150 mm	Adult use
TWS 12816	G18 x 150 mm	Pediatric use

Treasure Wright (U.K.) Ltd.

Reusable Hypodermic Needles



Back-cut point



- Professional standard quality hub
- Conical section, 6% Luer Lock
- Cannula made from AISI SS 304 and hub made from Nickle plated . Brass
- Back Cut cannula



TW Surgical

CODE	SIZE
TWS 12817	12G x ¾" (20mm x 2.8mm)
TWS 12818	12G x 1" (25mm x 2.8mm)
TWS 12819	12G x 1.1/4" (30mm x 2.8mm)
TWS 12820	12G x 1.1/2" (40mm x 2.8mm)
TWS 12821	13G x 1" (25mm x 2.5mm)
TWS 12822	13G x 1.1/2" (40mm x 2.5mm)
TWS 12823	14G x ½" (13mm x 2.00mm)
TWS 12824	14G x 5/8" (15mm x 2.00mm)
TWS 12825	14G x ¾" (20mm x 2.00mm)
TWS 12826	14G x 1" (25mm x 2.00mm)
TWS 12827	14G x 1.1/4" (30mm x 2.00mm)
TWS 12828	14G x 1.1/2" (40mm x 2.00mm)
TWS 12829	14G x 3/8" (9.5 mm x 2.00mm)
TWS 12830	15G x ½" (13mm x 1.8mm)
TWS 12831	15G x 5/8" (15mm x 1.8mm)
TWS 12832	15G x ¾" (20mm x 1.8mm)
TWS 12833	15G x 1" (25mm x 1.8mm)
TWS 12834	15G x 1.1/4" (30mm x 1.8mm)
TWS 12835	15G x 1.1/2" (40mm x 1.8mm)
TWS 12836	15G x 3/8" (9.5 mm x 1.8mm)
TWS 12837	16G x ½" (13mm x 1.6mm)
TWS 12838	16G x 5/8" (15mm x 1.6mm)
TWS 12839	16G x ¾" (20mm x 1.6mm)
TWS 12840	16G x 1" (25mm x 1.6 mm)
TWS 12841	16G x 1.1/4" (30mm x 1.6mm)
TWS 12842	16G x 1.1/2" (40mm x 1.6mm)
TWS 12843	16G x 3/8" (9.5 mm x 1.6mm)
TWS 12844	17G x ½" (13mm x 1.5mm)
TWS 12845	17G x 5/8" (15mm x 1.5mm)
TWS 12846	17G x ¾" (20mm x 1.5mm)
TWS 12847	17G x 1" (25mm x 1.5mm)

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com

CODE	SIZE
TWS 12848	17G x 1.1/4" (30mm x 1.5mm)
TWS 12849	17G x 1.1/2" (40mm x 1.5mm)
TWS 12850	17G x 3/8" (9.5 mm x 1.5mm)
TWS 12851	18G x 1/2" (13mm x 1.3mm)
TWS 12852	18G x 5/8" (15mm x 1.3mm)
TWS 12853	18G x 3/4" (20mm x 1.3mm)
TWS 12854	18G x 1" (25mm x 1.3mm)
TWS 12855	18G x 1.1/4" (30mm x 1.3mm)
TWS 12856	18G x 1.1/2" (40mm x 1.3mm)
TWS 12857	18G x 3/8" (9.5 mm x 1.3mm)
TWS 12858	18G x 1/4" (6mm x 1.3mm)
TWS 12859	19G x 1/2" (13mm x 1.00mm)
TWS 12860	19G x 5/8" (15mm x 1.00mm)
TWS 12861	19G x 3/4" (20mm x 1.00mm)
TWS 12862	19G x 1" (25mm x 1.00mm)
TWS 12863	19G x 1.1/4" (30mm x 1.00mm)
TWS 12864	19G x 1.1/2" (40mm x 1.00mm)
TWS 12865	19G x 3/8" (9.5 mm x 1.00mm)
TWS 12866	20G x 1/2" (13mm x 0.9mm)
TWS 12867	20G x 3/4"(20mm x 0.9mm)
TWS 12868	20G x 1" (25mm x 0.9mm)
TWS 12869	20G x 1.1/4" (30mm x 0.9mm)
TWS 12870	20G x 1.1/2"(40mm x 0.9mm)
TWS 12871	20G x 5/8" (15 mm x 0.9mm)
TWS 12872	21G x 1/2" (13mm x 0.8mm)
TWS 12873	21G x 3/4" (20mm x 0.8mm)
TWS 12874	21G x 1" (25mm x 0.8mm)
TWS 12875	21G x 1.1/4" (30mm x 0.8mm)
TWS 12876	21G x 1.1/2" (40mm x 0.8mm)
TWS 12877	21G x 5/8" (15 mm x 0.8mm)



CODE	SIZE
TWS 12878	22G x ½" (13mm x 0.7mm)
TWS 12879	22G x ¾" (20mm x 0.7mm)
TWS 12880	22G x 1" (25mm x 0.7mm)
TWS 12881	22G x 1.1/4" (30mm x 0.7mm)
TWS 12882	22G x 1.1/2" (40mm x 0.7mm)
TWS 12883	22G x 5/8" (15 mm x 0.7mm)
TWS 12884	22G x 1/4" (6 mm x 0.7mm)
TWS 12885	23G x 1" (25mm x 0.6mm)
TWS 12886	23G x 1.1/4" (32mm x 0.6mm)
TWS 12887	23G x 1.1/2" (40mm x 0.6mm)
TWS 12888	23G x 1/2" (13mm x 0.6mm)
TWS 12889	23G x 1/4" (6 mm x 0.6mm)
TWS 12890	23G x 5/8" (15mm x 0.6mm)
TWS 12891	24G x 1" (25mm x 0.5mm)
TWS 12892	24G x 1.1/4" (32mm x 0.5mm)
TWS 12893	24G x 1.1/2" (40mm x 0.5mm)
TWS 12894	24G x 1/2" (13mm x 0.5mm)
TWS 12895	25G x 1" (25mm x 0.5mm)
TWS 12896	25G x 1.1/4" (32mm x 0.5mm)
TWS 12897	25G x 1.1/2" (40mm x 0.5mm)
TWS 12898	26G x ½" (13mm x 0.45mm)
TWS 12899	26G x 1" (25mm x 0.45mm)
TWS 12900	26G x 1.5/8"(32mm x 0.45mm)
TWS 12901	26G x 1.1/2" (40mm x 0.45mm)
TWS 12902	27G x ½" (13mm x 0.4mm)
TWS 12903	27G x 1" (25mm x 0.4mm)
TWS 12904	27G x 1.1/4" (32mm x 0.4mm)
TWS 12905	27G x 1.1/2" (40mm x 0.4mm)

NOTE: Customized sizes are also available on request.

Animal Feeding Needles



Curved



Ball point



Oval point



Straight

- Allows fast, accurate feeding/medication
- Superior & export quality stainless steel AISI 304
- Round ball/oval tip
- Luer lock & Leak proof
- Set of 6 pieces
- Available in Straight & Curved

Reusable Animal Feeding Needles- Sizes

CODE	SIZE
TWS 12906	8G x 105mm
TWS 12907	10G x 105mm
TWS 12908	12G x 105mm
TWS 12909	14G x 95mm
TWS 12910	16G x 95mm
TWS 12911	18G x 55mm
TWS 12912	20G x 95mm
TWS 12913	20G x 55mm

Reusable Luer Lock Syringes



- Metal Luer Lock
- Non-slip grip and premium quality handle
- Plain metal plunger for higher strength
- Silicon washer ring for higher accuracy and durability (1 extra piece included)
- Boilable and Reusable

Reusable Luer Lock Syringes- Sizes

CODE	SIZE
TWS 12914	5 ml
TWS 12915	10 ml
TWS 12916	20 ml
TWS 12917	30 ml
TWS 12918	50 ml
TWS 12919	100 ml

Rumen Puncture Needles



- Triple beveled sharpened point for easy penetration
- Made with SS 304
- With removable stylet for ease

CODE	SIZE
TWS 12920	12G X 100mm
TWS 12921	10G X 100mm
TWS 12922	8G X 160mm (without adapter)

**Custom sizes available on request*

Teat Cannula set



Milk
Siphon

Teat tumor
extractor

Teat catheter

- Used for treatment of Mastities
- Made from AISI SS 304
- Two lateral eye opening in Milk Siphon & Teat Catheter

Teat Cannula set- Sizes

CODE	SIZE
TWS 12923	13 G x 100mm (Total length) (2.5mm x 100mm)
TWS 12924	13 G x 120mm (Total length) (2.5mm x 120mm)
TWS 12925	13G, With S.S. handle and Cleaning wire, Without Hub

Teat Siphon- Heavy



- Used for treatment of Mastities
- Made from AISI SS 304
- Two lateral eye openings.
- Comes with S.S Mandrill / cleaning wire.

CODE	SIZE
TWS 12926	13G, With S.S. handle and Cleaning wire, Without Hub

Tumor Extractor




- Heavy walled tubes for high strength
- Radius Tip to prevent injury
- One loop ring for easy grip
- Used for removing / extracting tumors

Teat Catheter



- Heavy walled tubes for high strength
- Radius Tip to prevent injury
- Two lateral eye openings for easy flow of medicine
- Supplied with Stylet (S.S 304 Cleaning wire) and curved handle grip

Artificial Insemination Gun Auto lock type




- Automatic lock- Spring type locking mechanism for cut or un-cut sheaths
- Suitable for French medium and mini type straws (0.25cc & 0.5cc)
- Universal type made of SS 304
- Cone shape head with smart design and easy to use

Artificial Insemination Gun Auto lock type- Sizes


CODE	SIZE
TWS 12927	0.50cc. (Suitable for French - Medium Straws)
TWS 12928	0.25cc. (Suitable for French - Mini Straws)
TWS 12929	Universal Type (Suitable for Both Medium and Mini Straws)
TWS 12930	Universal Type (Suitable for Both Medium and Mini Straws)

Artificial Insemination Gun Ring lock type



- Ring type locking mechanism for pre-cut sheaths
- Suitable for French medium and mini type straws (0.25cc & 0.5cc)
- Universal type made of SS 304
- Cone shape head with smart design and easy to use

Intra Uterine Catheters



- Made from AISI SS 304
- Luer slip hub
- Can be used with Sheath with ring lock
-

Feeding Nozzle



- Feeding nozzle for drenching
- 7/16 threaded hub or Luer lock hub
- Ball tip with leak proof joint
- Made from Chrome plated brass
- Available with 30ml Plastic steel Luer lock syringe

CODE	SIZE
TWS 12931	Threaded 6.5mm x 4mm x 12 cm, Curved
TWS 12933	Luer Lock 6.5mm x 4mm x 12 cm, Curved

**Custom sizes available on request*

Reusable Hypodermic Needle, Record



- Professional standard quality hub
- Hex section, Record taper
- Made from AISI SS 304
- Back Cut/ Lancet point cannula
- Available in dozen pc plastic box packing and outer box of one gross needles.

Reusable Hypodermic Needle, Record- Sizes

CODE	SIZE
TWS 12933	Record Needle 17G X 30mm
TWS 12934	Record Needle 18G x 30mm
TWS 12935	Record Needle 19G x 30mm
TWS 12936	Record Needle 14G X 35mm
TWS 12937	Record Needle 14G x 70mm

NOTE: Customized sizes are also available on request.

Ear tag applicator needle



- Applicator pin needle with nut, length: 1.26 inch
- Wide suitable range fitted for many pliers
- High quality accessory with high strength
- Easy to use
- Made from steel
- High tensile strength and it never bends

Reusable Long Aspiration Needles



Lancet point

- Nickel plated big hub with super smooth lapped finished cannula
- Available from G 10 to G 24 in any length combinations
- Pack of 10 pieces



TW Surgical

CODE	SIZE
TWS 12938	10G x 2" (50mm x 3.3mm)
TWS 12939	10G x 2.1/2" (60mm x 3.3mm)
TWS 12940	10G x 3" (75mm x 3.3mm)
TWS 12941	10G x 4" (100mm x 3.3mm)
TWS 12942	10G x 5" (130mm x 3.3mm)
TWS 12943	10G x 6" (150mm x 3.3mm)
TWS 12944	10G x 8" (200mm x 3.3mm)
TWS 12945	11G x 2" (50mm x 3.00mm)
TWS 12946	11G x 2.1/2" (60mm x 3.00mm)
TWS 12947	11G x 3" (75mm x 3.00mm)
TWS 12948	11G x 4" (100mm x 3.00mm)
TWS 12949	11G x 5" (130mm x 3.00mm)
TWS 12950	11G x 6" (150mm x 3.00mm)
TWS 12951	12G x 2" (50mm x 2.8mm)
TWS 12952	12G x 2.1/2" (60mm x 2.8mm)
TWS 12953	12G x 3" (75mm x 2.8mm)
TWS 12954	12G x 4" (100mm x 2.8mm)
TWS 12955	12G x 5" (130mm x 2.8mm)
TWS 12956	12G x 6" (150mm x 2.8mm)
TWS 12957	12G x 8" (200mm x 2.8mm)
TWS 12958	12G x 10" (250mm x 2.8mm)
TWS 12959	13G x 2" (50mm x 2.5mm)
TWS 12960	13G x 2.1/2" (60mm x 2.5mm)
TWS 12961	13G x 3" (75mm x 2.5mm)
TWS 12962	13G x 4" (100mm x 2.5mm)
TWS 12963	13G x 5" (130mm x 2.5mm)
TWS 12964	13G x 6" (150mm x 2.5mm)
TWS 12965	13G x 8" (200mm x 2.5mm)
TWS 12966	13G x 10" (250mm x 2.5mm)
TWS 12967	14G x 2" (50mm x 2.00mm)
TWS 12968	14G x 2.1/2" (60mm x 2.00mm)
TWS 12969	14G x 3" (75mm x 2.00mm)
TWS 12970	14G x 4" (100mm x 2.00mm)
TWS 12971	14G x 5" (130mm x 2.00mm)
TWS 12972	14G x 6" (150mm x 2.00mm)

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical

CODE	SIZE
TWS 12974	14G x 8" (200mm x 2.00mm)
TWS 12975	14G x 10" (250mm x 2.00mm)
TWS 12976	15G x 2" (50mm x 1.8mm)
TWS 12977	15G x 2.1/2" (60mm x 1.8mm)
TWS 12978	15G x 3" (75mm x 1.8mm)
TWS 12979	15G x 4" (100mm x 1.8mm)
TWS 12980	15G x 5" (130mm x 1.8mm)
TWS 12981	15G x 6" (150mm x 1.8mm)
TWS 12982	15G x 8" (200mm x 1.8mm)
TWS 12983	15G x 10" (250mm x 1.8mm)
TWS 12984	16G x 2" (50mm x 1.6mm)
TWS 12985	16G x 2.1/2" (60mm x 1.6mm)
TWS 12986	16G x 3" (75mm x 1.6mm)
TWS 12987	16G x 4" (100mm x 1.6mm)
TWS 12988	16G x 5" (130mm x 1.6mm)
TWS 12989	16G x 6" (150mm x 1.6mm)
TWS 12990	16G x 8" (200mm x 1.6mm)
TWS 12991	16G x 10" (250mm x 1.6mm)
TWS 12992	16G x 12" (300mm x 1.6mm)
TWS 12993	17G x 2" (50mm x 1.5mm)
TWS 12994	17G x 2.1/2" (60mm x 1.5mm)
TWS 12995	17G x 3" (75mm x 1.5mm)
TWS 12996	17G x 4" (100mm x 1.5mm)
TWS 12997	17G x 5" (130mm x 1.5mm)
TWS 12998	17G x 6" (150mm x 1.5mm)
TWS 12999	17G x 8" (200mm x 1.5mm)
TWS 13000	17G x 10" (250mm x 1.5mm)
TWS 13001	17G x 12" (250mm x 1.5mm)
TWS 13002	18G x 2" (50mm x 1.3mm)
TWS 13003	18G x 2.1/2" (60mm x 1.3mm)
TWS 13004	18G x 3" (75mm x 1.3mm)

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical

CODE	SIZE
TWS 13005	18G x 4" (100mm x 1.3mm)
TWS 13006	18G x 5" (130mm x 1.3mm)
TWS 13007	18G x 6" (150mm x 1.3mm)
TWS 13008	18G x 8" (200mm x 1.3mm)
TWS 13009	18G x 10" (250mm x 1.3mm)
TWS 13010	18G x 12" (300mm x 1.3mm)
TWS 13011	19G x 2" (50mm x 1.00mm)
TWS 13012	19G x 2.1/2" (60mm x 1.00mm)
TWS 13013	19G x 3" (75mm x 1.00mm)
TWS 13014	19G x 4" (100mm x 1.00mm)
TWS 13015	19G x 5" (130mm x 1.00mm)
TWS 13016	19G x 6" (150mm x 1.00mm)
TWS 13017	19G x 8" (200mm x 1.00mm)
TWS 13018	19G x 10" (250mm x 1.00mm)
TWS 13019	19G x 12" (300mm x 1.00mm)
TWS 13020	20G x 2" (50mm x 0.9mm)
TWS 13022	20G x 2.1/2" (60mm x 0.9mm)
TWS 13023	20G x 3" (75mm x 0.9mm)
TWS 13024	20G x 4" (100mm x 0.9mm)
TWS 13025	20G x 5" (130mm x 0.9mm)
TWS 13026	20G x 6" (150mm x 0.9mm)
TWS 13027	20G x 8" (200mm x 0.9mm)
TWS 13028	20G x 10" (250mm x 0.9mm)
TWS 13029	20G x 12" (300mm x 0.9mm)
TWS 13030	21G x 2" (50mm x 0.8mm)
TWS 13032	21G x 2.1/2" (60mm x 0.8mm)
TWS 13033	21G x 3" (75mm x 0.8mm)
TWS 13034	21G x 4" (100mm x 0.8mm)
TWS 13035	21G x 5" (130mm x 0.8mm)
TWS 13036	21G x 6" (150mm x 0.8mm)
TWS 13037	21G x 8" (200mm x 0.8mm)
TWS 13038	21G x 10" (250mm x 0.8mm)
TWS 13039	22G x 2" (50mm x 0.7mm)
TWS 13040	22G x 2.1/2" (60mm x 0.7mm)

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical

CODE	SIZE
TWS 13041	22G x 3" (75mm x 0.7mm)
TWS 13042	22G x 4" (100mm x 0.7mm)
TWS 13043	22G x 5" (130mm x 0.7mm)
TWS 13044	22G x 6" (150mm x 0.7mm)
TWS 13045	22G x 8" (200mm x 0.7mm)
TWS 13046	22G x 10" (250mm x 0.7mm)
TWS 13047	23G x 2" (50mm x 0.6mm)
TWS 13048	23G x 2.1/2"(60mm x 0.6mm)
TWS 13049	23G x 3" (75mm x 0.6mm)
TWS 13050	23G x 4" (100mm x 0.6mm)
TWS 13051	23G x 5" (130mm x 0.6mm)
TWS 13052	23G x 6" (150mm x 0.6mm)
TWS 13053	23G x 8" (200mm x 0.6mm)
TWS 13054	23G x 10" (250mm x 0.6mm)
TWS 13055	24G x 2" (50mm x 0.57mm)
TWS 13056	24G x 2.1/2"(60mm x 0.57mm)
TWS 13057	24G x 3" (75mm x 0.57mm)
TWS 13058	24G x 4" (100mm x 0.57mm)
TWS 13059	24G x 5" (130mm x 0.57mm)
TWS 13060	24G x 6"(150mm x 0.57mm)
TWS 13061	24G x 8"(200mm x 0.57mm)
TWS 13062	24G x 10"(250mm x 0.57mm)

NOTE: Customized sizes are also available on request.

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com

GAUGE SIZE (INCH)	GAUGE SIZE (mm)
18G x 36"	1.3mm x 900mm
16G x 36"	1.6mm x 900mm
13G x 48"	2.5mm x 1200mm
14G x 36"	2.00mm x 900mm
14G x 48"	2.00mm x 1200mm
18G x 24"	1.3mm x 600mm
16G x 24"	1.6mm x 600mm
16G x 36"	1.6mm x 900mm
14G x 24"	2.00mm x 600mm
16G x 20"	1.6mm x 500mm
18G x 16"	1.3mm x 400mm
16G x 20"	1.6mm x 500mm
18G x 12"	1.3mm x 300mm
20G x 12"	0.9mm x 300mm
22G x 12"	0.7mm x 300mm
15G x 12"	1.8mm x 300mm
16G x 12"	1.6mm x 300mm
14G x 12"	2.00mm x 300mm
19G x 12"	1.00mm x 300mm

GAUGE SIZE (INCH)	GAUGE SIZE (mm)
23G x 20"	0.6mm x 500mm
22G x 20"	0.7mm x 500mm
19G x 20"	1.00mm x 500mm
22G x 20"	0.7mm x 500mm
17G x 20"	1.5mm x 500mm
14G x 20"	2.00mm x 500mm
17G x 20"	1.5mm x 500mm
18G x 16"	1.3mm x 400mm
23G x 12"	0.6mm x 300mm
17G x 16"	1.5mm x 400mm
18G x 20"	1.3mm x 500mm
17G x 12"	1.5mm x 300mm
21G x 12"	0.8mm x 300mm
14G x 24"	2.00mm x 600mm
20G x 12"	0.9mm x 300mm
18G x 20"	1.3mm x 500mm
20G x 16"	0.9mm x 400mm
20G x 20"	0.9mm x 500mm
19G x 20"	1.00mm x 500mm
21G x 20"	0.8mm x 500mm

**Custom sizes available on request*

Reusable S.S. Transfer Needles- Double Beveled



- Stainless Steel Transfer Needles
- Both side sharp bevel.
- Suitable for transfer between bottles.
- Can be steam sterilized and reused.

Industrial Dispensing Needle



Blunt tip



- Blunt point with Stainless Steel AISI 304 cannula
- Standard and professional quality Luer lock Hub
- Precision tip for controlled flow
- Sold in a pack of 12 pieces

Industrial Dispensing Needle- Sizes

GAUGE	O.D. mm	I.D. mm
5	5.62	4.84
6	5.25	4.47
7	4.49	3.65
8	4.10	3.19
9	n/a	n/a
10	3.53	2.92
11	2.97	2.45
12	2.72	2.00
13	2.40	2.00
14	2.00	1.50
15	1.80	1.40
16	1.60	1.20
17	1.40	1.00
18	1.20	0.80
19	1.00	0.65
20	0.90	0.60
21	0.80	0.50
22	0.70	0.40
23	0.60	0.30
24	0.55	0.28
25	0.50	0.25
26	0.45	0.25
27	0.40	0.20
28	0.36	0.14
30	0.30	0.10

Specialty Needles



- Special needle and cannula
- Available curved, bent and straight.
- Sharp pointed or blunt
- Any gauge and length combination
- As per drawing or samples
- Material : Nickle plated brass or stainless steel 304.
- Specific OD and ID available on request.



TW
Surgical

Treasure Wright (U.K.) Ltd.

**Address: 128-160 , City Road,
London EC1V 2NX. United Kingdom**

www.twrightmedical.com

info@twrightmedical.com



TW Surgical



**UK
CA CE**