



TW
Surgical

CAPSULORHEXIS FORCEP



**UK
CA** **CE**



TW Surgical

Treasure Wright (U.K.) Ltd.



www.twrightmedical.com




info@twrightmedical.com

CAPSULORHEXIS FORCEPS



TW Surgical


TWS 85-1201/47
INAMURA CAPSULORHEXIS
FORCEPS CURVED, TITANIUM



Overall length 120 mm, 4.7 inches

- Delicate shafts 11 mm long
- TC sharp tips
- Round handle


TWS 85-1202/47
INAMURA CAPSULORHEXIS
FORCEPS CURVED,
TITANIUM



Overall length 120 mm, 4.7 inches

- Delicate shafts 11 mm long
- TC grasping tips
- Round handle


TWS 85-1203/47
INAMURA CAPSULORHEXIS
FORCEPS CURVED, TITANIUM



Overall length 120 mm, 4.7 inches

- Delicate shafts 11 mm long
- TC 1 blunt and 1 sharp tips
- Round handle


TWS 85-1204/47
INAMURA CAPSULORHEXIS
FORCEPS ANGLED, TITANIUM



Overall length 120 mm, 4.7 inches

- Delicate shafts 11 mm long
- TC sharp tips
- Round handle


TWS 85-1205/47
INAMURA CAPSULORHEXIS
FORCEPS ANGLED, TITANIUM



Overall length 120 mm, 4.7 inches

- Delicate shafts 11 mm long
- TC grasping tips
- Round handle


TWS 85-1206/47
INAMURA CAPSULORHEXIS
FORCEPS ANGLED, TITANIUM



Overall length 120 mm, 4.7 inches

- Delicate shafts 11 mm long
- TC 1 blunt and 1 sharp tips
- Round handle

TWS 85-1207/47
INAMURA CAPSULORHEXIS
FORCEPS ANGLED, TITANIUM



Overall length 119 mm, 4.7 inches


- Delicate shafts 11 mm long
- Sharp tips
- Round handle

CAPSULORHEXIS FORCEPS



TW Surgical

TWS 85-2120/47
INAMURA CAPSULORHEXIS FORCEPS
CURVED, STAINLESS




- Delicate shafts 11 mm long
- Sharp tips

Overall length 120 mm, 4.7 inches

- Round handle

TWS 85-2121/47
INAMURA CAPSULORHEXIS FORCEPS
CURVED, STAINLESS




- Delicate shafts 11 mm long
- Grasping tips

Overall length 120 mm, 4.7 inches

- Round handle

TWS 85-2122/47
INAMURA CAPSULORHEXIS FORCEPS
CURVED, STAINLESS




- Delicate shafts 11 mm long
- 1 blunt and 1 sharp tips

Overall length 120 mm, 4.7 inches

- Round handle

TWS 85-2123/47
INAMURA CAPSULORHEXIS FORCEPS
ANGLED, STAINLESS




- Delicate shafts 11 mm long

Overall length 120 mm, 4.7 inches

- Round handle

TWS 85-2124/47
INAMURA CAPSULORHEXIS
FORCEPS ANGLED, STAINLESS




- Delicate shafts 11 mm long
- Grasping tips

Overall length 120 mm, 4.7 inches

- Round handle

TWS 85-2125/47
CAPSULORHEXIS FORCEPS CURVED,
INAMURA TYPE, STAINLESS




- Delicate shafts 11 mm long
- 1 blunt and 1 sharp tips

Overall length 120 mm, 4.7 inches

- Round handle

TWS 85-2125/47
CAPSULORHEXIS FORCEPS CURVED,
INAMURA TYPE, STAINLESS



- Delicate shafts 12.5 mm long
- Sharp tips

Overall length 115 mm, 4.5 inches


- Round handle

CAPSULORHEXIS FORCEPS



TW Surgical



TWS 85-3101/39
INAMURA CAPSULORHEXIS
FORCEPS CURVED, TITANIUM



Overall length 101 mm, 3.9 inches

- Delicate shafts 10.5 mm long
- TC grasping tips
- Flat handle


TWS 85-3100/39
INAMURA CAPSULORHEXIS
FORCEPS CURVED, TITANIUM



Overall length 100 mm, 3.9 inches


- Delicate shafts 12 mm long with 6 notches
- 1 mm distance between notches
- TC grasping tips
- Flat handle

TWS 85-3102/39
INAMURA CAPSULORHEXIS
FORCEPS CURVED, TITANIUM



- Delicate shafts 10 mm long
- TC sharp tips
- Flat handle


TWS 85-3103/39
INAMURA CAPSULORHEXIS
FORCEPS CURVED, ROUND HANDLE



Overall length 116 mm, 4.5 inches

- Delicate shafts 10.8 mm long, stainless
- Ultra-thin profile and limiter for capsulorhexis even through 1.8 mm incision
- TC sharp tips
- Flat handle, titanium


TWS 85-3104/39
INAMURA CAPSULORHEXIS FORCEPS,
CURVED, TITANIUM/ STAINLESS



Overall length 119 mm, 4.7 inches

- Delicate shafts 11 mm long, stainless
- TC grasping tips
- Flat handle, titanium


TWS 85-3104/39
INAMURA CAPSULORHEXIS FORCEPS, CURVED,
TITANIUM/ STAINLESS



Overall length 118 mm, 4.6 inches

- Delicate shafts 11 mm long
- TC sharp tips
- Flat handle

TWS 85-3104/39
INAMURA CAPSULORHEXIS
FORCEPS ANGLED, TITANIUM



Overall length 118 mm, 4.6 inches

- Delicate shafts 11 mm long
- Angled 45 degree
- TC sharp tips
- Flat handle


CAPSULORHEXIS FORCEPS



TW Surgical

TWS 85-0481/46

INAMURA CAPSULORHEXIS FORCEPS CURVED, TITANIUM




Overall length 118 mm, 4.6 inches

- Delicate shafts 11 mm long
- TC grasping tips
- Flat handle

TWS 85-0432/48

C ALLADINE-INAMURA CAPSULORHEXIS FORCEPS ANGLED, TITANIUM





Overall length 123 mm, 4.8 inches

- Delicate shafts 14.7 mm long
- 11 mm long from tip to screw
- Pointed serrated interlocking tips
- Maximum opening 3.7 mm wide
- Round handle

TWS 85-0433/48

C ALLADINE-INAMURA CAPSULORHEXIS FORCEPS, SLIGHTLY CURVED, TITANIUM




Overall length 123 mm, 4.8 inches

- Delicate shafts 14.7 mm long
- 11 mm long from tip to screw
- Pointed serrated interlocking tips
- Maximum opening 3.7 mm wide
- Round handle

TWS 85-4102/04

UTRATA CAPSULORHEXIS FORCEPS CURVED, TITANIUM



Overall length 102 mm, 4 inches

- Delicate shafts 11 mm long
- TC grasping tips
- Flat handle

TWS 85-4103/04

UTRATA CAPSULORHEXIS FORCEPS CURVED, TITANIUM




Overall length 102 mm, 4 inches

- Delicate shafts 11 mm long
- TC sharp tips
- Flat handle

TWS 85-4104/04

UTRATA CAPSULORHEXIS FORCEPS CURVED, TITANIUM




Overall length 102 mm, 4 inches

- Delicate shafts 11 mm long
- TC grasping tips with groove
- Flat handle

TWS 85-4105/04

DR. DAVID ISRAELI'S UTRATA CAPSULORHEXIS FORCEPS CURVED, TITANIUM




Overall length 102 mm, 4 inches

- Delicate shafts 11 mm long
- TC grasping tips
- Reversed tips
- Flat handle

Treasure Wright (U.K.) Ltd.
128 City Road, London, EC1V 2NX, UK
info@twrightmedical.com

TWS 85-5102/04

UTRATA CAPSULORHEXIS FORCEPS ANGLED, TITANIUM



- Delicate shafts 11 mm long
- TC grasping tips


Overall length 102 mm, 4 inches

- Flat handle



TWS 85-5103/04

UTRATA CAPSULORHEXIS FORCEPS ANGLED, TITANIUM



- Delicate shafts 11 mm long
- TC sharp tips


Overall length 102 mm, 4 inches

- Flat handle



TWS 85-5104/04

UTRATA CAPSULORHEXIS FORCEPS ANGLED, TITANIUM



- Delicate shafts 11 mm long
- TC grasping tips with groove


Overall length 102 mm, 4 inches

- Flat handle



TWS 85-5105/04

JAFFE UTRATA CAPSULORHEXIS FORCEPS ANGLED, TITANIUM



- Delicate shafts 11 mm long
- TC round grasping tips


Overall length 102 mm, 4 inches

- Flat handle



TWS 85-5106/04

DR. DAVID ISRAELI'S UTRATA CAPSULORHEXIS FORCEPS ANGLED, TITANIUM



- Delicate shafts 11 mm long
- TC grasping/reversed tips


Overall length 102 mm, 4 inches

- Flat handle



TWS 85-5100/41

UTRATA CAPSULORHEXIS FORCEPS ANGLED, TITANIUM



- Delicate shafts 11.5 mm long
- TC grasping tips with groove


Overall length 105 mm, 4.1 inches

- Round handle



TWS 85-1050/41

UTRATA CAPSULORHEXIS FORCEPS ANGLED, TITANIUM



- Delicate shafts 11.5 mm long
- TC tips

Overall length 105 mm, 4.1 inches

- Round handle



TWS 85-1051/41

UTRATA CAPSULORHEXIS FORCEPS ANGLED, TITANIUM



Overall length 105 mm, 4.1 inches
6 engravings at 1, 2, 3, 4, 5, 6 mm

- Delicate shafts 11.5 mm long
- TC grasping tips
- Round handle

TWS 85-1052/41

UTRATA CAPSULORHEXIS FORCEPS ANGLED, TITANIUM




Overall length 105 mm, 4.1 inches

- Delicate shafts 11.5 mm long
- TC sharp tips
- Round handle

TWS 85-1053/41

UTRATA CAPSULORHEXIS FORCEPS ANGLED, TITANIUM




Overall length 105 mm, 4.1 inches

- Delicate shafts 11.5 mm long
- TC grasping tips with groove
- Round handle

TWS 85-1054/41

UTRATA CAPSULORHEXIS FORCEPS CURVED, TITANIUM



Overall length 105 mm, 4.1 inches
6 engravings at 1, 2, 3, 4, 5, 6 mm

- Delicate shafts 11.5 mm long
- TC grasping tips
- Round handle

TWS 85-1055/41

UTRATA CAPSULORHEXIS FORCEPS CURVED, TITANIUM




Overall length 105 mm, 4.1 inches

- Delicate shafts 11.5 mm long
- TC sharp tips
- Round handle

TWS 85-1056/41

UTRATA CAPSULORHEXIS FORCEPS ANGLED, TITANIUM




Overall length 105 mm, 4.1 inches

- Delicate shafts 11.5 mm long
- Round handle

TWS 85-1057/41

CAPSULORHEXIS FORCEPS ANGLED, TITANIUM



Overall length 105 mm, 4.1 inches

- Delicate shafts 11.5 mm long
- TC 1 blunt and 1 sharp tips
- Round handle

CAPSULORHEXIS FORCEPS



TW Surgical

TWS 85-7105/41

CAPSULORHEXIS FORCEPS CURVED, TITANIUM



Overall length 105 mm, 4.1 inches

- Delicate shafts 11.5 mm long
- TC 1 blunt and 1 sharp tips
- Round handle

TWS 85-7104/41

UTRATA CAPSULORHEXIS FORCEPS CURVED, TITANIUM OF 162.01 ANGLED



Overall length 104 mm, 4 inches

- Delicate shafts 11.5 mm long
- TC sharp tips
- Light round handle

TWS 85-7115/45

UTRATA CAPSULORHEXIS FORCEPS CURVED, TITANIUM




Overall length 115 mm, 4.5 inches

- Delicate shafts 11.5 mm long
- TC sharp tips
- Light round handle

TWS 85-0383/32

UTRATA CAPSULORHEXIS FORCEPS ANGLED, TITANIUM / STAINLESS




Overall length 83 mm, 3.2 inches

- Delicate shafts 11.5 mm long, stainless
- Grasping tips
- Flat handle, titanium

TWS 85-0384/32

UTRATA CAPSULORHEXIS FORCEPS ANGLED, TITANIUM / STAINLESS




Overall length 83 mm, 3.2 inches

- Delicate shafts 11.5 mm long, stainless
- Very sharp tips
- Flat handle, titanium

TWS 85-0388/32

UTRATA CAPSULORHEXIS FORCEPS ANGLED, TITANIUM MS 300S STAINLESS




Overall length 83 mm, 3.2 inches

- Delicate shafts 12.3 mm long, stainless
- Very sharp tips
- Flat handle, titanium

TWS 85-0387/32

UTRATA CAPSULORHEXIS FORCEPS ANGLED, TITANIUM



Overall length 83 mm, 3.2 inches

- Very delicate shafts 12.5 mm long
- TC grasping tips
- Flat handle

Treasure Wright (U.K.) Ltd.
128 City Road, London, EC1V 2NX, UK
info@twrightmedical.com


CAPSULORHEXIS FORCEPS



TW Surgical

TWS 85-0883/32

UTRATA CAPSULORHEXIS FORCEPS ANGLED, STAINLESS



- Very delicate shafts 12.5 mm long
- Grasping tips


Overall length 83 mm, 3.2 inches

- Flat handle



TWS 85-0884/33

UTRATA CAPSULORHEXIS FORCEPS ANGLED, TITANIUM



- Delicate shafts 11.5 mm long
- TC sharp tips


Overall length 83 mm, 3.3 inches

- Flat handle



TWS 85-0885/33

UTRATA CAPSULORHEXIS FORCEPS, CURVED, TITANIUM
MS375 01 STAINLESS, NO TC



- Delicate shafts 12.3 mm long
- TC grasping tips

Overall length 83 mm, 3.3 inches

- Flat handle



TWS 85-8105/41

UTRATA CAPSULORHEXIS FORCEPS CURVED, STAINLESS



- Delicate shafts 11 mm long
- Sharp tips


Overall length 105 mm, 4.1 inches

- Flat handle



TWS 85-0870/27

UTRATA CAPSULORHEXIS FORCEPS ANGLED, STAINLESS



- Delicate shafts 11 mm long
- Sharp tips


Overall length 70 mm, 2.7 inches

- Flat handle




TWS 85-8104/41

MICROCOAXIAL CAPSULORHEXIS FORCEPS ANGLED,
STAINLESS




- Delicate shafts 10.5 mm long
- Ultra-thin profile and limiter for capsulorhexis even through 1.8 mm incision
- TC sharp tips
- Round handle

Overall length 104 mm, 4.1 inches




TWS 85-8106/41

MICROCOAXIAL CAPSULORHEXIS FORCEPS CURVED,
STAINLESS



- Delicate shafts 10.5 mm long
- Ultra-thin profile and limiter for capsulorhexis even through 1.8 mm incision
- TC sharp tips
- Round handle

Overall length 105 mm, 4.1 inches



Treasure Wright (U.K.) Ltd.

128 City Road, London, EC1V 2NX, UK

info@twrightmedical.com


CAPSULORHEXIS FORCEPS



TW Surgical

TWS 85-9104/41

MICROCOAXIAL CAPSULORHEXIS FORCEPS ANGLED,
TITANIUM



Overall length 104 mm, 4.1 inches

- Delicate shafts 10.5 mm long
- Ultra-thin profile and limiter for capsulorhexis even through 1.8 mm incision
- TC sharptips
- Round handle

TWS 85-9105/41

MICROCOAXIAL CAPSULORHEXIS FORCEPS CURVED,
TITANIUM



Overall length 105 mm, 4.1 inches

- Delicate shafts 10.5 mm long
- Ultra-thin profile and limiter for capsulorhexis even through 1.8 mm incision
- TC sharptips
- Round handle

TWS 85-9106/41

MICROCOAXIAL CAPSULORHEXIS FORCEPS CURVED,
STAINLESS MS 261.02 TITANIUM



Overall length 105 mm, 4.1 inches

- Delicate shafts 10.5 mm long
- Ultra-thin profile and limiter for capsulorhexis even through 1 mm incision
- TC sharptips
- Round handle

TWS 85-1050/41

MICROCOAXIAL CAPSULORHEXIS FORCEPS
CURVED, STAINLESS
MS 261.03





Overall length 105 mm, 4.1 inches

- Delicate shafts 10.5 mm long
- Ultra-thin profile and limiter for capsulorhexis even through 1 mm incision
- TC sharptips
- Round handle

TWS 85-1052/41

MICROCOAXIAL CAPSULORHEXIS FORCEPS
CURVED, STAINLESS



Overall length 105 mm, 4.1 inches

- Delicate shafts 10.5 mm long
- Ultra-thin profile and limiter for capsulorhexis even through 1.5 mm incision
- TC sharp tips
- Round handle

TWS 85-0999/39

CAPSULORHEXIS FORCEPS, CURVED, TITANIUM MS
385.01 STAINLESS, NO TC





Overall length 99 mm, 3.9 inches

- Delicate shafts 11 mm long
- Maximum opening 0.8 mm wide
- TC grasping tips
- Round handle

All forceps also available in stainless steel / TC – Tungsten Carbide

TWS 85-9118/46

MICROCOAXIAL CAPSULORHEXIS FORCEPS
CURVED, TITANIUM



Overall length 118 mm, 4.6 inches

- Delicate shafts 10.5 mm long
- Ultra-thin profile and limiter for capsulorhexis even through 1.5 mm incision
- TC sharptips
- Round handle

Treasure Wright (U.K.) Ltd.
128 City Road, London, EC1V 2NX, UK
info@twrightmedical.com


CAPSULORHEXIS FORCEPS



TW Surgical

TWS 85-1090/35

MICROCOAXIAL CAPSULORHEXIS FORCEPS
CURVED, TITANIUM





Overall length 90 mm, 3.5 inches

- Delicate shafts 10.5 mm long
- Ultra-thin profile and limiter for capsulorhexis even through 1.5 mm incision
- TC sharptips
- Round handle

TWS 85-8110/46

MICROCOAXIAL UTRATA CAPSULORHEXIS FORCEPS
CURVED, TITANIUM




Overall length 118 mm, 4.6 inches

- Delicate shafts 10 mm long
- Ultra-thin profile and limiter for capsulorhexis even through 1.5 mm incision
- TC sharptips
- Round handle

TWS 85-1097/38

DOUBLE CROSS ACTION CAPSULORHEXIS FORCEPS
FOR SMALL INSICIONS (MICS), TITANIUM



Overall length 97 mm, 3.8 inches

- Extremely delicate shafts 11.6 mm long
- Serrated grasping tips 0.4mm wide
- Maximum opening 0.5mm
- Flat handle

TWS 85-9710/38

DOUBLE CROSS ACTION CAPSULORHEXIS FORCEPS
FOR SMALL INSICIONS (MICS), STAINLESS




Overall length 97 mm, 3.8 inches

- Extremely delicate shafts 11.6 mm long
- Serrated grasping tips 0.4mm wide
- Maximum opening 0.5mm
- Flat handle

TWS 85-1099/39

DOUBLE CROSS ACTION CAPSULORHEXIS FORCEPS
FOR SMALL INSICIONS (MICS), TITANIUM




Overall length 99 mm, 3.9 inches

- Extremely delicate shafts 11.6 mm long
- Serrated grasping tips 0.4mm wide
- Maximum opening 0.5mm
- Flat handle

TWS 85-9910/39

DOUBLE CROSS ACTION CAPSULORHEXIS FORCEPS FOR
SMALL INSICIONS (MICS), STAINLESS



Overall length 99 mm, 3.9 inches

- Extremely delicate shafts 11.6 mm long
- Serrated grasping tips 0.4mm wide
- Maximum opening 0.5mm
- Flat handle


PRECHOPPERS



TW Surgical

TWS 85-1121/48

PRECHOPPER FOR HARD NUCLEUS STRAIGHT,
TITANIUM
MS 150.01 STAINLESS




- Maximum opening 3 mm
- Cross action

Overall length 121 mm, 4.8 inches

- Round handle

TWS 85-1122/48

PRECHOPPER FOR SOFT NUCLEUS STRAIGHT,
TITANIUM



- Maximum opening 3 mm
- Cross action

Overall length 121 mm, 4.8 inches

- Round handle

TWS 85-1217

PRECHOPPER FOR HARD NUCLEUS ANGLED,
TITANIUM
MS 150.01 STAINLESS




- Maximum opening 3 mm
- Cross action

Overall length 117 mm, 4.6 inches

- Round handle

TWS 85-1172/47

PRECHOPPER FOR SOFT NUCLEUS ANGLED,
TITANIUM





- Maximum opening 3 mm
- Cross action

Overall length 117 mm, 4.6 inches

- Round handle

TWS 85-1192/47

PRECHOPPER FOR HARD NUCLEUS CURVED,
TITANIUM
MS 150.01 STAINLESS




- Maximum opening 3 mm
- Cross action

Overall length 119 mm, 4.7 inches

- Round handle

TWS 85-9112/47

PRECHOPPER FOR SOFT NUCLEUS
CURVED, TITANIUM




- Maximum opening 3 mm
- Cross action

Overall length 119 mm, 4.7 inches

- Round handle

TWS 85-1123/48

AKAHOSHI COMBO PRECHOPPER STRAIGHT,
STAINLESS



- Maximum opening 2.6 mm
- Cross action

Overall length 121 mm, 4.8 inches

- Round handle

Treasure Wright (U.K.) Ltd.
128 City Road, London, EC1V 2NX, UK
info@twrightmedical.com


PRECHOPPERS



TW Surgical

TWS 85-1213/48

AKAHOSHI COMBO PRECHOPPER STRAIGHT, TITANIUM



- Maximum opening 2.6 mm
- Cross action


Overall length 121 mm, 4.8 inches

- TC tips
- Round handle



TWS 85-1214/48

AKAHOSHI COMBO PRECHOPPER STRAIGHT, STAINLES



- Maximum opening 2.6 mm
- Cross action


Overall length 121 mm, 4.8 inches

- Round handle



TWS 85-1200/47

COMBO PRECHOPPER FOR SUB-1.5 mm COAXIAL MICROPHACO STRAIGHT, STAINLESS



- Maximum opening 1.5 mm
- Specially designed to split nuclei through small 1.5 mm incision

Overall length 120 mm, 4.7 inches

- Cross action
- Round handle



TWS 85-1201/47

AKAHOSHI HYBRID COMBO PRECHOPPER STRAIGHT STAINLESS



- Maximum opening 2.9 mm
- Groove enables easy rotation of the nuclei during prechopping


Overall length 120 mm, 4.7 inches

- Cross action
- Round handle



TWS 85-1210/48

AKAHOSHI NUCLEUS SPLITTER/ COMPLETER/ SEPARATOR, STRAIGHT, TITANIUM



- Blunt, atraumatic tips ideal for completing separation of nucleus after initiation of femtosecond
- Used to divide the nucleus after grooving


Overall length 121 mm, 4.8 inches

- Maximum opening 3.2 mm
- Blunt ring tip 1.6 mm dia
- Cross action
- Round handle



TWS 85-1211/48

AKAHOSHI HYBRID COMBO PRECHOPPER, TITANIUM



- Used for prechopping grade 1 to 4 nuclei
- Rounded tips with sharpened angular edge

Overall length 121 mm, 4.8 inches

- Maximum opening 2.9 mm
- Cross action
- Round handle




SPLITTERS



TW Surgical


TWS 85-1375/29
NUCLEUS SPLITTER FORCEPS, TITANIUM



Overall length 75 mm, 2.9 inches

- TC external surface of paddle shape tips
- To crush nucleus small segments for phaco
- Cross action criss-cross jaws
- 1.1 × 2.5 mm tips
- Flat handle


TWS 85-7513/29
NUCLEUS SPLITTER FORCEPS, TITANIUM



Overall length 75 mm, 2.9 inches

- TC external surface of paddle shape tips
- To crush nucleus small segments for phaco
- Cross action criss-cross jaws
- 1.7 × 2.8 mm tips
- Flat handle


TWS 85-1307/42
NUCLEUS SPLITTER FORCEPS, TITANIUM



Overall length 107 mm, 4.2 inches

- TC external surface of paddle shape tips
- To crush nucleus small segments for phaco
- Maximum opening 7mm
- Cross action criss-cross jaws
- 1.1 × 2.5 mm tips
- Flat handle

TWS 85-7013/42
ERNEST NUCLEUS CRACKER, TITANIUM



Overall length 107 mm, 4.2 inches

- 2.5 × 1.3 mm tips
- Maximum opening 4mm
- Cross action criss-cross jaws
- Round handle


FIXATION FORCEPS



TW Surgical

TWS 85-1480/31

PIERSE TYPE FIXATION FORCEPS STRAIGHT, TITANIUM




Overall length 80 mm, 3.1 inches

- 1×1 Pierse tips
- Tips 0.2 mm wide
- 0.1 mm groove
- Flat handle

TWS 85-1481/31

PIERSE TYPE FIXATION FORCEPS STRAIGHT, TITANIUM




Overall length 80 mm, 3.1 inches

- 1×1 Pierse tips
- 0.1 mm groove
- Tips 0.2 mm wide
- TC tying platform 2.5 mm long
- Flat handle

TWS 85-1410/31

PIERSE TYPE FIXATION FORCEPS CURVED, TITANIUM




Overall length 110 mm, 4.3 inches

- 1×1 Pierse tips
- 0.1 mm groove
- Tips 0.2 mm wide
- TC tips
- Flat handle

TWS 85-1498/39

PIERSE TYPE FIXATION FORCEPS STRAIGHT, TITANIUM




Overall length 98 mm, 3.9 inches

- 1×1 Pierse tips
- 0.1 mm groove
- Tips 0.2 mm wide
- TC tips
- Round handle

TWS 85-1490/36

PIERSE TYPE FIXATION FORCEPS ANGLED, TITANIUM




Overall length 90 mm, 3.6 inches

- 1×1 Pierse tips
- 0.1 mm groove
- Tips 0.4 mm wide
- Flat handle

TWS 85-1400/39

HOFMAN - POLACK FIXATION FORCEPS ANGLED, TITANIUM




Overall length 100 mm, 3.9 inches

- 0.1mm Pierse type tips
- U-shaped 2 mm spread
- Angled 2 point fixation
- Round handle

TWS 85-1402/39

FORCEPS ANGLED, TITANIUM BORES CORNEAL FIXATION



Overall length 100 mm, 3.9 inches

- 1×2 teeth, 0.12 mm height
- Angled 2 point fixation
- 0.4×0.5mm tips
- U-shaped 2 mm spread
- TC teeth
- Round handle

Treasure Wright (U.K.) Ltd.
128 City Road, London, EC1V 2NX, UK
info@twrightmedical.com


SPLITTERS



TW Surgical

TWS 85-1500/39

BORES CORNEAL FORCEPS ANGLED, TITANIUM



- 1×2 teeth, 0.12 mm height
- Angled 2 point fixation
- 0.3×0.3mm tips




Overall length 100 mm, 3.9 inches

- U-shaped 2 mm spread
- TC teeth
- Round handle

TWS 85-1595/37

THORPE CORNEAL FIXATION FORCEPS STRAIGHT, TITANIUM



- 2×3 teeth, 0.2 mm height
- 90 degree teeth
- Tips 0.9 mm wide




Overall length 95 mm, 3.7 inches

- 3.5 mm tying platforms
- TC teeth and platforms
- Light round handle

TWS 85-1504/04

GRANT CORNEAL FIXATION FORCEPS CURVED WITH LOCK, TITANIUM



- 2×3 teeth, 0.2 mm height
- 90 degree teeth
- Tips 0.9 mm wide




Overall length 104 mm, 4 inches

- 3.5 mm tying platforms
- TC teeth and platforms
- Flat handle

TWS 85-1045/04

WRIGHT LOCKED MUSCLE FORCEPS STRAIGHT, TITANIUM



- 2×3 teeth, 0.5 mm height
- 45 degree teeth
- 6 mm tying platforms




Overall length 104 mm, 4 inches

- Tips 1.4 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-1046/04

LOCKED MUSCLE FORCEPS CURVED (MOODY), TITANIUM



- 2×3 teeth, 0.5 mm height
- 45 degree teeth
- Tips 1.4 mm wide




Overall length 104 mm, 4 inches

- TC teeth
- Flat handle

TWS 85-1047/04

LOCKED MUSCLE FORCEPS CURVED, TITANIUM



- 1×2 teeth, 0.5 mm height
- 45 degree teeth
- Tips 0.8 mm wide




Overall length 104 mm, 4 inches

- TC teeth
- Flat handle

TWS 85-1048/04

WRIGHT LOCKED MUSCLE FORCEPS CURVED, TITANIUM



- 1×2 teeth, 0.3 mm height
- 45 degree teeth
- Tips 0.5 mm wide



Overall length 104 mm, 4 inches

- TC teeth
- Flat handle


Treasure Wright (U.K.) Ltd.

128 City Road, London, EC1V 2NX, UK

info@twrightmedical.com

SPLITTERS

TWS 85-1604/04
WRIGHT LOCKED MUSCLE FORCEPS
CURVED, TITANIUM




- 1×2 teeth, 0.4 mm height
- 45 degree teeth
- Tips 0.6 mm wide
- TC teeth
- Flat handle



Overall length 104 mm, 4 inches

TWS 85-1699/39
JAMESON MUSCLE FORCEPS SMALL LEFT,
TITANIUM




- 4 pins
- 8 mm angled jaws
- Sliding lock
- Flat handle



Overall length 99 mm, 3.9 inches

TWS 85-9916/39
JAMESON MUSCLE FORCEPS SMALL RIGHT,
TITANIUM




- 4 pins
- 8 mm angled jaws
- Sliding lock
- Flat handle



Overall length 99 mm, 3.9 inches

TWS 85-1610/39
JAMESON MUSCLE FORCEPS LARGE LEFT,
TITANIUM




- 6 pins
- 10.6 mm angled jaws
- Sliding lock
- Flat handle



Overall length 100 mm, 3.9 inches

TWS 85-1016/39
JAMESON MUSCLE FORCEPS LARGE RIGHT,
TITANIUM




- 6 pins
- 10.6 mm angled jaws
- Sliding lock
- Flat handle



Overall length 100 mm, 3.9 inches

TWS 85-1695/37
PUTTERMAN CLAMP, TITANIUM




- 6 pins
- Serrated jaws
- Sliding lock
- Flat handle



Overall length 95 mm, 3.7 inches

FIXATION FORCEPS


TWS 85-1704/41
APT'S CLAMP, STAINLESS



Overall length 104 mm, 4.1 inches

- Opening 2, 4, 6 mm
- 11 mm serrated working part angled 77 degree
- One shaft with a groove
- Second shaft with a projection (a cog)
- Round handle


TWS 85-1779/31
APT'S CLAMP, STAINLESS



Overall length 79 mm, 3.1 inches

- Opening 2, 4, 6 mm
- 10 mm serrated working part angled 77 degree
- One shaft with a groove
- Second shaft with a projection (a cog)
- Round handle


TWS 85-1785/33
CILIA FORCEPS NARROW, TITANIUM



Overall length 85 mm, 3.3 inches

- Narrowed smooth jaws
- Tips 0.4 mm wide
- TC tips 6 mm long
- Flat handle


TWS 85-1791/36
ZIEGLER CILIA FORCEPS, TITANIUM



Overall length 91 mm, 3.6 inches

- 4.5 mm rounded platforms
- TC platforms
- Flat handle

TWS 85-1776/03
CILIA FORCEPS WIDE, TITANIUM




Overall length 76 mm, 3 inches

- Platform 1×3 mm wide
- TC platforms
- Flat handle

CORNEAL TYING FORCEPS

TWS 85-1875/29

CORNEAL FORCEPS COLIBRI TYPE, TITANIUM



- 1×2 teeth, 0.1 mm height
- 45 degree teeth
- 3 mm tying platforms




Overall length 75 mm, 2.9 inches

- Tips 0.3 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-1870/29

CORNEAL FORCEPS COLIBRI TYPE, TITANIUM



- 1×2 teeth, 0.1 mm height
- 90 degree teeth
- 3 mm tying platforms




Overall length 75 mm, 2.9 inches

- Tips 0.3 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-1871/29

TROUTMAN-BARRAQUER CORNEAL UTILITY FORCEPS (COLIBRI), TITANIUM



- 1×2 teeth, 0.12 mm height
- 45 degree teeth
- 3 mm tying platforms




Overall length 75 mm, 2.9 inches

- Tips 0.3 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-1872/29

TROUTMAN-BARRAQUER CORNEAL UTILITY FORCEPS (COLIBRI), TITANIUM



- 1×2 teeth, 0.12 mm height
- 90 degree teeth
- 3 mm tying platforms




Overall length 75 mm, 2.9 inches

- Tips 0.3 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-1873/29

TROUTMAN-BARRAQUER CORNEAL UTILITY FORCEPS (COLIBRI), TITANIUM



- 1×2 teeth, 0.2 mm height
- 45 degree teeth
- 4 mm tying platforms




Overall length 75 mm, 2.9 inches

- Tips 0.5 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-1874/29

TROUTMAN-BARRAQUER CORNEAL UTILITY FORCEPS (COLIBRI), TITANIUM



- 1×2 teeth, 0.2 mm height
- 90 degree teeth
- 4 mm tying platforms




Overall length 75 mm, 2.9 inches

- Tips 0.5 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-1876/29

TROUTMAN-BARRAQUER CORNEAL FORCEPS COLIBRI TYPE, TITANIUM



- 1×2 teeth, 0.4 mm height
- 45 degree teeth
- 5 mm tying platforms



Overall length 75 mm, 2.9 inches

- Tips 0.6 mm wide
- TC teeth and platforms
- Flat handle


CORNEAL TYING FORCEPS



TW Surgical

TWS 85-1975/29

TROUTMAN- BARRAQUER CORNEAL FORCEPS
COLIBRI TYPE, TITANIUM



Overall length 75 mm, 2.9 inches

- 1×2 teeth, 0.4 mm height
- 90 degree teeth
- 5 mm tying platforms
- Tips 0.6 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-1978/03

CORNEAL FORCEPS COLIBRI
TYPE, STAINLESS




Overall length 78 mm, 3 inches

- 1×2 teeth, 0.2 mm height
- 75 degree teeth
- 5 mm tying platforms
- Tips 0.3 mm wide
- Flat handle

TWS 85-1986/34

CORNEAL FORCEPS COLIBRI TYPE,
TITANIUM




Overall length 86 mm, 3.4 inches

- 1x2 teeth, 0.12 mm height
- 90 degree teeth
- 5.0mm tying platforms
- Tips 0.35 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-1904/41

CORNEAL FORCEPS COLIBRI TYPE,
TITANIUM



Overall length 104 mm, 4.1 inches

- 1×2 teeth, 0.1 mm height
- 45 degree teeth
- 5.5 mm tying platforms
- Tips 0.3 mm wide
- TC teeth and platforms
- Light round handle

TWS 85-1920/47

CORNEAL FORCEPS COLIBRI TYPE,
TITANIUM




Overall length 120 mm, 4.7 inches

- 1×2 teeth, 0.1 mm height
- 45 degree teeth
- 6 mm tying platforms
- Tips 0.3 mm wide
- TC teeth and platforms
- Light round handle

TWS 85-1921/40

CORNEAL FORCEPS PIERCE TYPE,
TITANIUM




Overall length 122 mm, 4.8 inches

- 1×2 teeth, 0.12 mm height
- Pierce tip 0.3 mm wide with groove 0.06 mm
- 4 mm tying platforms
- TC teeth and platforms
- Round handle

TWS 85-1906/42

CORNEAL FORCEPS COLIBRI TYPE,
TITANIUM



Overall length 106 mm, 4.2 inches

- 1×2 teeth, 0.12 mm height
- 45 degree teeth
- 6 mm tying platforms
- Tips 0.3 mm wide
- TC teeth and platforms
- Flat handle


Treasure Wright (U.K.) Ltd.

128 City Road, London, EC1V 2NX, UK

info@twrightmedical.com

CORNEAL TYING FORCEPS

TWS 85-2016/42
CORNEAL FORCEPS COLIBRI TYPE,
TITANIUM




- 1×2 teeth, 0.2 mm height
- 45 degree teeth
- 6 mm tying platforms

Overall length 106 mm, 4.2 inches

- Tips 0.3 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-2100/39
CORNEAL FORCEPS COLIBRI TYPE,
TITANIUM




- 1×2 teeth, 0.1 mm height
- 45 degree teeth
- Tips 0.3 mm wide

Overall length 100 mm, 3.9 inches

- TC teeth
- Round handle

TWS 85-2101/39
OGAWA CORNEAL FORCEPS
COLIBRI TYPE, TITANIUM




- 1×2 teeth, 0.12 mm height
- 45 degree teeth
- Tips 0.3 mm wide

Overall length 100 mm, 3.9 inches

- TC teeth
- Round handle

TWS 85-2102/39
CORNEAL FORCEPS COLIBRI
TYPE, TITANIUM





- 1×2 teeth, 0.1 mm height
- 90 degree teeth
- Tips 0.3 mm wide

Overall length 100 mm, 3.9 inches

- TC teeth
- Round handle

TWS 85-2103/39
OGAWA CORNEAL FORCEPS
COLIBRI TYPE, TITANIUM




- 1×2 teeth, 0.12 mm height
- 90 degree teeth
- Tips 0.3 mm wide

Overall length 100 mm, 3.9 inches

- TC teeth
- Round handle

TWS 85-2104/39
GIRARD CORNEAL FORCEPS
COLIBRI TYPE, TITANIUM




- 1×2 teeth, 0.12 mm height
- 90 degree teeth
- Tips 0.3 mm wide

Overall length 100 mm, 3.9 inches

- 6 mm tying platforms
- TC teeth and platforms
- Round handle

TWS 85-2105/35
CATALANO CORNEAL FORCEPS,
TITANIUM



- 1×2 teeth, 0.1 mm height
- 30 degree teeth
- Tips 0.3 mm wide

Overall length 90 mm, 3.5 inches


- TC teeth
- Round handle

CORNEAL TYING FORCEPS



TW Surgical

TWS 85-2184/33
CORNEAL FORCEPS, TITANIUM




- 1×2 teeth, 0.3 mm height
- 90 degree teeth
- Tips 0.5 mm wide

Overall length 84 mm, 3.3 inches

- TC teeth
- Round handle



TWS 85-2190/35
PIERSE TYPE CORNEAL
COLIBRI FORCEPS, TITANIUM




- 1×1 Pierse tip
- 0.1 mm groove

Overall length 90 mm, 3.5 inches

- Tips 0.4 mm wide
- Round handle



TWS 85-2175/29
HOSKIN CORNEAL FORCEPS,
TITANIUM




- 1×1 Pierse tip
- 0.1 mm groove

Overall length 75 mm, 2.9 inches

- Tips 0.4 mm wide
- Flat handle



TWS 85-2170/27
HOSKIN CORNEAL FORCEPS,
TITANIUM



- 1×1 Pierse tip
- 0.1 mm groove
- 6 mm tying platforms


Overall length 70 mm, 2.7 inches

- Tips 0.4 mm wide
- Flat handle



SCLERAL FORCEPS

TWS 85-2285/33
PAUFIQUE TYING SCLERAL FORCEPS
STRAIGHT, TITANIUM
MS074.100,12 mm




- 1×2 teeth, 0.3 mm height
- 45 degree teeth
- 6 mm tying platforms



Overall length 85 mm, 3.3 inches

- Tips 0.5 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-8522/33
PAUFIQUE TYING SCLERAL FORCEPS
STRAIGHT, TITANIUM
MS074.060,1 mm MS074.070,2 mm




- 1×2 teeth, 0.3 mm height
- 90 degree teeth
- 6 mm tying platforms



Overall length 85 mm, 3.3 inches

- TC teeth and platforms
- Tips 0.5 mm wide
- Flat handle

TWS 85-8584/33
TYING FORCEPS, STAINLESS




- 1x2 teeth, 0.1 mm height
- TC teeth and tip 6 mm long
- 90 degree teeth



Overall length 84 mm, 3.3 inches

- No platform
- Tips 0.5 x 0.35 mm wide
- Flat handle

TWS 85-2286/33
PAUFIQUE TYING SCLERAL FORCEPS
STRAIGHT, TITANIUM




- 1×2 teeth, 0.5 mm height
- 45 degree teeth
- 6 mm tying platforms



Overall length 85 mm, 3.3 inches

- Tips 0.8 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-2287/33
PAUFIQUE TYING SCLERAL FORCEPS
STRAIGHT, TITANIUM




- 1×2 teeth, 0.5 mm height
- 90 degree teeth
- 6 mm tying platforms



Overall length 85 mm, 3.3 inches

- Tips 0.8 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-2288/33
PAUFIQUE TYING SCLERAL FORCEPS
STRAIGHT, TITANIUM




- 1×2 teeth, 0.8 mm height
- 45 degree teeth
- 6 mm tying platforms



Overall length 85 mm, 3.3 inches

- Tips 1.2 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-2289/33
PAUFIQUE TYING SCLERAL FORCEPS
STRAIGHT, TITANIUM



- 1×2 teeth, 0.8 mm height
- 90 degree teeth
- 6 mm tying platforms




Overall length 85 mm, 3.3 inches

- Tips 1.2 mm wide
- TC teeth and platforms
- Flat handle

SCLERAL FORCEPS


TWS 85-2315/41
WRIGHT TYING SCLERAL
FORCEPS STRAIGHT, TITANIUM



Overall length 105 mm, 4.1 inches

- 1×2 teeth, 0.75 mm height
- 60 degree curved teeth
- 6 mm tying platforms
- Tips 1.2 mm wide
- TC teeth and platforms
- Flat handle


TWS 85-1523/41
TYING SCLERAL FORCEPS
STRAIGHT, TITANIUM



Overall length 105 mm, 4.1 inches

- 1×2 teeth, 0.65 mm height
- 60 degree curved teeth
- 6 mm tying platforms
- Tips 1.15 mm wide
- TC teeth and platforms
- Flat handle


TWS 85-2311/43
CASTROVIEJO TYING FORCEPS
STRAIGHT, TITANIUM



Overall length 110 mm, 4.3 inches

- 1×2 teeth, 0.9 mm height
- 6 mm tying platforms
- Tips 1 mm wide
- TC teeth and platforms
- Flat handle


TWS 85-1123/43
CASTROVIEJO TYING FORCEPS
STRAIGHT, TITANIUM



Overall length 110 mm, 4.3 inches

- 1×2 teeth, 1.25 mm height
- 6 mm tying platforms
- Tips 1.2 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-2317/42
SCLERAL PLUG FORCEPS
CROSS ACTION, TITANIUM




Overall length 107 mm, 4.2 inches

- Angled shafts with groove 0.3 mm
- Flat handle
- Designed to securely hold scleral plugs during insertion and removal. The cross action handle is designed to reduce hand fatigue while grasping plug

TYING FORCEPS

TWS 85-2487/35
CASTROVIEJO TYING FORCEPS
STRAIGHT, TITANIUM




- 1×2 teeth, 0.1 mm height
- 45 degree teeth
- 6 mm tying platforms

Overall length 87 mm, 3.5 inches

- Tips 0.3 mm wide
- TC teeth and platforms
- Flat handle



TWS 85-8724/35
CASTROVIEJO TYING FORCEPS
STRAIGHT, TITANIUM




- 1×2 teeth, 0.2 mm height
- 45 degree teeth
- 6 mm tying platforms

Overall length 87 mm, 3.5 inches

- Tips 0.4 mm wide
- TC teeth and platforms
- Flat handle



TWS 85-2484/33
CASTROVIEJO TYING FORCEPS
ANGLED, TITANIUM




- 1×2 teeth, 0.1 mm height
- 45 degree teeth
- 6 mm tying platforms

Overall length 84 mm, 3.3 inches

- Tips 0.3 mm wide
- TC teeth and platforms
- Flat handle



TWS 85-8424/33
CASTROVIEJO TYING FORCEPS
ANGLED, TITANIUM




- 1×2 teeth, 0.2 mm height
- 45 degree teeth
- 6 mm tying platforms

Overall length 84 mm, 3.3 inches

- Tips 0.4 mm wide
- TC teeth and platforms
- Flat handle



TWS 85-2478/35
CASTROVIEJO TYING FORCEPS
CURVED, TITANIUM




- 1×2 teeth, 0.1 mm height
- 45 degree teeth
- 6 mm tying platforms

Overall length 87 mm, 3.5 inches

- Tips 0.3 mm wide
- TC teeth and platforms
- Flat handle



TWS 85-8724/35
CASTROVIEJO TYING FORCEPS
CURVED, TITANIUM




- 1×2 teeth, 0.2 mm height
- 45 degree teeth
- 6 mm tying platforms

Overall length 87 mm, 3.5 inches

- Tips 0.4 mm wide
- TC teeth and platforms
- Flat handle



TWS 85-2493/36
TYING FORCEPS STRAIGHT, STAINLESS



- 1×2 teeth, 0.1 mm height
- 90 degree curved teeth
- 5 mm tying platforms

Overall length 93 mm, 3.6 inches


- Tips 0.3 mm wide
- TC teeth and platforms
- Flat handle



TYING FORCEPS

TWS 85-2593/36

TYING FORCEPS STRAIGHT, TITANIUM




Overall length 93 mm, 3.6 inches

- 1×2 teeth, 0.12 mm height
- 90 degree curved teeth
- 5 mm tying platforms
- Tips 0.3 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-2500/39

BONN TYING FORCEPS STRAIGHT, TITANIUM




Overall length 100 mm, 3.9 inches

- 1×2 teeth, 0.12 mm height
- 45 degree teeth
- 6 mm tying platforms
- Tips 0.3 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-1025/39

BONN TYING FORCEPS STRAIGHT, TITANIUM




Overall length 100 mm, 3.9 inches

- 1×2 teeth, 0.12 mm height
- 90 degree teeth
- 6 mm tying platforms
- Tips 0.4 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-9825/38

BONN TYING FORCEPS ANGLED, TITANIUM




Overall length 98 mm, 3.8 inches

- 1×2 teeth, 0.12 mm height
- 45 degree teeth
- 6 mm tying platforms
- Tips 0.3 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-9925/39

BONN TYING FORCEPS CURVED, TITANIUM




Overall length 99 mm, 3.9 inches

- 1×2 teeth, 0.12 mm height
- 45 degree teeth
- 6 mm tying platforms
- Tips 0.3 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-1052/41

CASTROVIEJO TYING FORCEPS STRAIGHT, TITANIUM




Overall length 105 mm, 4.1 inches

- 1×2 teeth, 0.12 mm height
- 45 degree teeth
- 6 mm tying platforms
- Tips 0.3 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-1205/41

CASTROVIEJO TYING FORCEPS STRAIGHT, TITANIUM



Overall length 105 mm, 4.1 inches

- 1×2 teeth, 0.3 mm height
- 45 degree teeth
- 6 mm tying platforms
- Tips 0.5 mm wide
- TC teeth and platforms
- Flat handle

Treasure Wright (U.K.) Ltd.


128 City Road, London, EC1V 2NX, UK

info@twrightmedical.com

TYING FORCEPS

TWS 85-2615/41

CASTROVIEJO TYING FORCEPS STRAIGHT, TITANIUM



- 1×2 teeth, 0.5 mm height
- 45 degree teeth
- 6 mm tying platforms




Overall length 105 mm, 4.1 inches

- Tips 0.7 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-1526/41

CASTROVIEJO TYING FORCEPS STRAIGHT, TITANIUM



- 1×2 teeth, 0.6 mm height
- 45 degree teeth
- 6 mm tying platforms




Overall length 105 mm, 4.1 inches

- Tips 0.8 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-2602/04

CASTROVIEJO TYING FORCEPS ANGLED, TITANIUM



- 1×2 teeth, 0.12 mm height
- 45 degree teeth
- 6 mm tying platforms




Overall length 102 mm, 4 inches

- Tips 0.3 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-1026/04

CASTROVIEJO TYING FORCEPS ANGLED, TITANIUM



- 1×2 teeth, 0.3 mm height
- 45 degree teeth
- 6 mm tying platforms




Overall length 102 mm, 4 inches

- Tips 0.5 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-2605/41

CASTROVIEJO TYING FORCEPS CURVE TITANIUM



- 1×2 teeth, 0.12 mm height
- 45 degree teeth
- 6 mm tying platforms




Overall length 105 mm, 4.1 inches

- Tips 0.3 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-1052/41

CASTROVIEJO TYING FORCEPS CURVED, TITANIUM



- 1×2 teeth, 0.3 mm height
- 45 degree teeth
- 6 mm tying platforms




Overall length 105 mm, 4.1 inches

- Tips 0.5 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-1053/41

CASTROVIEJO TYING FORCEPS CURVED, TITANIUM



- 1×2 teeth, 0.5 mm height
- 45 degree teeth
- 6 mm tying platforms




Overall length 105 mm, 4.1 inches

- Tips 0.7 mm wide
- TC teeth and platforms
- Flat handle

TYING FORCEPS


TWS 85-2715/41
CASTROVIEJO TYING FORCEPS
STRAIGHT, TITANIUM



Overall length 105 mm, 4.1 inches

- 1×2 teeth, 0.12 mm height
- 90 degree teeth
- 6 mm tying platforms
- Tips 0.3 mm wide
- TC teeth and platforms
- Flat handle


TWS 85-1527/41
CASTROVIEJO TYING FORCEPS
STRAIGHT, TITANIUM



Overall length 105 mm, 4.1 inches

- 1×2 teeth, 0.3 mm height
- 90 degree teeth
- 6 mm tying platforms
- Tips 0.5 mm wide
- TC teeth and platforms
- Flat handle


TWS 85-1027/41
CASTROVIEJO TYING FORCEPS
STRAIGHT, TITANIUM



Overall length 105 mm, 4.1 inches

- 1×2 teeth, 0.5 mm height
- 90 degree teeth
- 6 mm tying platforms
- Tips 0.7 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-1127/43
CASTROVIEJO TYING FORCEPS
STRAIGHT, TITANIUM



Overall length 110 mm, 4.3 inches

- 1×2 teeth, 0.12 mm height
- 90 degree curved teeth
- 6 mm tying platforms
- Tips 0.3 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-2711/43
CASTROVIEJO TYING FORCEPS
STRAIGHT, TITANIUM



Overall length 110 mm, 4.3 inches

- 1×2 teeth, 0.3 mm height
- 90 degree curved teeth
- 6 mm tying platforms
- Tips 0.5 mm wide
- TC teeth and platforms
- Flat handle


TWS 85-1270/43
CASTROVIEJO TYING FORCEPS
STRAIGHT, TITANIUM



Overall length 110 mm, 4.3 inches

- 1×2 teeth, 0.5 mm height
- 90 degree curved teeth
- 6 mm tying platforms
- Tips 0.55 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-1057/41
TYING FORCEPS STRAIGHT,
TITANIUM




Overall length 105 mm, 4.1 inches

- 1×2 teeth, 0.5 mm height
- 45 degree teeth
- 6 mm tying platforms
- Tips 0.8 mm wide
- TC teeth and platforms
- Flat handle

TYING FORCEPS

TWS 85-2880/31
TYING FORCEPS ANGLED, STAINLESS




- 1×2 teeth, 0.1 mm height
- 60 degree teeth
- Pointed platforms 11 mm long



- Overall length 80 mm, 3.1 inches**
- Tips 0.3 mm wide
 - Flat handle

TWS 85-8028/31
TYING FORCEPS ANGLED, STAINLESS




- 1×2 teeth, 0.1 mm height
- 45 degree teeth
- Pointed platforms 11 mm long



- Overall length 80 mm, 3.1 inches**
- Tips 0.3 mm wide
 - Flat handle

TWS 85-2808/31
TYING FORCEPS ANGLED, STAINLESS




- 1×2 teeth, 0.2 mm height
- 45 degree teeth
- Pointed platforms 10 mm long



- Overall length 80 mm, 3.1 inches**
- Tips 0.5 mm wide
 - Flat handle

TWS 85-0828/31
TYING FORCEPS ANGLED, STAINLESS




- 1×2 teeth, 0.2 mm height
- 40 degree teeth
- Pointed platforms 10 mm long



- Overall length 80 mm, 3.1 inches**
- Tips 0.5 mm wide
 - Flat handle

TWS 85-2814/41
TYING FORCEPS STRAIGHT, TITANIUM




- 1×2 teeth, 0.12 mm height
- 45 degree teeth
- 6 mm tying platforms



- Overall length 104 mm, 4.1 inches**
- Tips 0.3 mm wide
 - TC teeth and platforms
 - Light round handle

TWS 85-2815/41
TENNANT TYING FORCEPS STRAIGHT,
TITANIUM




- 1×2 teeth, 0.1 mm height
- 45 degree teeth
- 6 mm tying platforms



- Overall length 105 mm, 4.1 inches**
- Tips 0.3 mm wide
 - TC teeth and platforms
 - Round handle

TWS 85-1528/41
TENNANT TYING FORCEPS STRAIGHT,
TITANIUM




- 1×2 teeth, 0.12 mm height
- 45 degree teeth
- 6 mm tying platforms



- Overall length 105 mm, 4.1 inches**
- Tips 0.3 mm wide
 - TC teeth and platforms
 - Round handle


TYING FORCEPS

TWS 85-2999/38
TENNANT TYING FORCEPS
STRAIGHT, TITANIUM




- 1×2 teeth, 0.15 mm height
- 90 degree teeth
- 6.3 mm tying platforms
- Overall length 99 mm, 3.8 inches
- Tips 0.3 mm wide
- TC teeth and platforms
- Round handle

TWS 85-1292/40
TENNANT TYING FORCEPS
ANGLED, TITANIUM




- 1×2 teeth, 0.1 mm height
- 45 degree teeth
- 6 mm tying platforms
- Overall length 102 mm, 4.0 inches
- Tips 0.3 mm wide
- TC teeth and platforms
- Round handle

TWS 85-1294/41
TENNANT TYING FORCEPS
CURVED, TITANIUM




- 1×2 teeth, 0.1 mm height
- 30 degree teeth
- 6 mm tying platforms
- Overall length 104 mm, 4.1 inches
- Tips 0.3 mm wide
- TC teeth and platforms
- Round handle

TWS 85-2974/29
McPHERSON TYING FORCEPS
STRAIGHT, TITANIUM




- 4 mm tying platforms
- Tips 0.2 mm wide
- Overall length 74 mm, 2.9 inches
- TC platforms
- Flat handle

TWS 85-2972/28
BONN TYING FORCEPS
ANGLED, TITANIUM




- 4 mm tying platforms
- Tips 0.2 mm wide
- Overall length 72 mm, 2.8 inches
- TC platforms
- Flat handle

TWS 85-2973/28
OSHER TYING FORCEPS
CURVED, TITANIUM



- 4 mm tying platforms
- Tips 0.2 mm wide
- Overall length 73 mm, 2.8 inches
- TC platforms
- Flat handle

TWS 85-2980/31
CATALANO TYING FORCEPS
ANGLED, STAINLESS




- 10 mm V-groove tying platforms
- Tips 0.5 mm wide
- Overall length 80 mm, 3.1 inches
- TC platforms
- Flat handle

TYING FORCEPS

TWS 85-3087/34

McPHERSON TYING FORCEPS STRAIGHT, TITANIUM




- 4 mm tying platforms
- Tips 0.3 mm wide

Overall length 87 mm, 3.4 inches

- TC platforms
- Flat handle

TWS 85-8730/34

McPHERSON TYING FORCEPS STRAIGHT, TITANIUM




- 6 mm tying platforms
- Tips 0.3 mm wide

Overall length 87 mm, 3.4 inches

- TC platforms
- Flat handle

TWS 85-3807/34

McPHERSON TYING FORCEPS STRAIGHT, TITANIUM




- 8 mm tying platforms
- Tips 0.3 mm wide

Overall length 87 mm, 3.4 inches

- TC platforms
- Flat handle

TWS 85-0387/34

McPHERSON TYING FORCEPS ANGLED, TITANIUM




- 10 mm tying platforms
- Tips 0.3 mm wide

Overall length 87 mm, 3.4 inches

- TC platforms
- Flat handle

TWS 85-3084/33

McPHERSON TYING FORCEPS STRAIGHT, TITANIUM




- 4 mm tying platforms
- Tips 0.3 mm wide

Overall length 84 mm, 3.3 inches

- TC platforms
- Flat handle

TWS 85-8430/33

McPHERSON TYING FORCEPS ANGLED, TITANIUM




- 6 mm tying platforms
- Tips 0.3 mm wide

Overall length 84 mm, 3.3 inches

- TC platforms
- Flat handle

TWS 85-3084/33

McPHERSON TYING FORCEPS ANGLED, TITANIUM



- 8 mm tying platforms
- Tips 0.3 mm wide


Overall length 84 mm, 3.3 inches

- TC platforms
- Flat handle

TYING FORCEPS

TWS 85-3184/33

McPHERSON TYING FORCEPS
ANGLED, TITANIUM



- 10 mm tying platforms
- Tips 0.3 mm wide


Overall length 84 mm, 3.3 inches

- TC platforms
- Flat handle



TWS 85-3187/24

McPHERSON TYING FORCEPS CURVED,
TITANIUM



- 4 mm tying platforms
- Tips 0.3 mm wide


Overall length 87 mm, 3.4 inches

- TC platforms
- Flat handle



TWS 85-8731/24

McPHERSON TYING FORCEPS CURVED,
TITANIUM



- 6 mm tying platforms
- Tips 0.3 mm wide


Overall length 87 mm, 3.4 inches

- TC platforms
- Flat handle



TWS 85-3187/34

McPHERSON TYING FORCEPS CURVED,
TITANIUM



- 8 mm tying platforms
- Tips 0.3 mm wide


Overall length 87 mm, 3.4 inches

- TC platforms
- Flat handle



TWS 85-3184/33

McPHERSON TYING FORCEPS STRAIGHT,
TITANIUM



- 4.5 mm tying platforms
- Tips 0.4 mm wide


Overall length 84 mm, 3.3 inches

- TC platforms
- Flat handle



TWS 85-8431/33

McPHERSON TYING FORCEPS ANGLED,
TITANIUM



- 3.5 mm tying platforms
- Tips 0.45 mm wide


Overall length 84 mm, 3.3 inches

- TC platforms
- Flat handle



TWS 85-3485/33

McPHERSON TYING FORCEPS ANGLED,
TITANIUM



- 5 mm tying platforms
- Tips 0.3 mm wide

Overall length 85 mm, 3.3 inches


- TC platforms
- Flat handle



TYING FORCEPS

TWS 85-3286/34

McPHERSON TYING FORCEPS ANGLED, TITANIUM




Overall length 86 mm, 3.4 inches

- 8 mm tying platforms
- Tips 0.5 mm wide
- TC platforms
- Flat handle

TWS 85-3284/33

McPHERSON TYING FORCEPS ANGLED, TITANIUM




Overall length 84 mm, 3.3 inches

- 8.5 mm tying platforms
- Tips 0.3 mm wide
- TC platforms
- Flat handle

TWS 85-3200/39

McPHERSON TYING FORCEPS STRAIGHT, TITANIUM




Overall length 100 mm, 3.9 inches

- 4 mm tying platforms
- Tips 0.3 mm wide
- TC platforms
- Flat handle

TWS 85-0032/39

McPHERSON TYING FORCEPS STRAIGHT, TITANIUM




Overall length 100 mm, 3.9 inches

- 6 mm tying platforms
- Tips 0.3 mm wide
- TC platforms
- Flat handle

TWS 85-3298/39

McPHERSON TYING FORCEPS ANGLED, TITANIUM




Overall length 98 mm, 3.9 inches

- 4 mm tying platforms
- Tips 0.3 mm wide
- TC platforms
- Flat handle

TWS 85-9832/39

McPHERSON TYING FORCEPS ANGLED, TITANIUM




Overall length 98 mm, 3.9 inches

- 6 mm tying platforms
- Tips 0.3 mm wide
- TC platforms
- Flat handle

TWS 85-3299/39

McPHERSON TYING FORCEPS CURVED, TITANIUM




Overall length 99 mm, 3.9 inches

- 4 mm tying platforms
- Tips 0.3 mm wide
- TC platforms
- Flat handle

TYING FORCEPS

TWS 85-3399/39

McPHERSON TYING FORCEPS
CURVED, TITANIUM




Overall length 99 mm, 3.9 inches

- TC platforms
- Flat handle

TWS 85-3300/39

CATALANO TYING FORCEPS
STRAIGHT, TITANIUM




Overall length 100 mm, 3.9 inches

- 6 mm V- groove tying platforms
- TC platforms
- Flat handle
- Tips 0.5 mm wide

TWS 85-0033/39

TYING FORCEPS STRAIGHT,
TITANIUM MS 367,02 STAINLESS




Overall length 100 mm, 3.9 inches

- 5 mm tying platforms
- TC platforms
- Tips 0.4 x 0.25 mm wide
- Flat handle with holes

TWS 85-0303/39

TYING FORCEPS ANGLED,
TITANIUM MS 367,03
STAINLESS



Overall length 100 mm, 3.9 inches

- 5 mm tying platforms
- TC platforms
- Tips 0.3 x 0.25 mm wide
- Flat handle with holes

TWS 85-3315/41

McPHERSON TYING FORCEPS
CURVED, TITANIUM




CURVED, TITANIUM

Overall length 100 mm, 3.9 inches

- 6 mm V- groove tying platforms
- TC platforms
- Flat handle
- Tips 0.5 mm wide

TWS 85-3232/39

CATALANO TYING FORCEPS
STRAIGHT, TITANIUM




Overall length 105 mm, 4.1 inches

- 4 mm tying platforms
- TC platforms
- Flat handle
- Tips 0.3 mm wide

TWS 85-1533/41

McPHERSON TYING FORCEPS
STRAIGHT, TITANIUM




Overall length 105 mm, 4.1 inches

- 6 mm tying platforms
- TC platforms
- Flat handle
- Tips 0.3 mm wide

TYING FORCEPS

TWS 85-3405/41
McPHERSON TYING FORCEPS
STRAIGHT, TITANIUM




- 8 mm tying platforms
- Tips 0.3 mm wide

Overall length 105 mm, 4.1 inches

- TC platforms
- Flat handle

TWS 85-3412/04
McPHERSON TYING FORCEPS ANGLED,
TITANIUM




- 4 mm tying platforms
- Tips 0.3 mm wide

Overall length 102 mm, 4 inches

- TC platforms
- Flat handle

TWS 85-1034/04
McPHERSON TYING FORCEPS ANGLED,
TITANIUM




- 6 mm tying platforms
- Tips 0.3 mm wide

Overall length 102 mm, 4 inches

- TC platforms
- Flat handle

TWS 85-3410/04
McPHERSON TYING FORCEPS ANGLED,
TITANIUM




- 8 mm tying platforms
- Tips 0.3 mm wide

Overall length 102 mm, 4 inches

- TC platforms
- Flat handle

TWS 85-3415/41
McPHERSON TYING FORCEPS CURVED,
TITANIUM




- 4 mm tying platforms
- Tips 0.3 mm wide

Overall length 105 mm, 4.1 inches

- TC platforms
- Flat handle

TWS 85-1534/41
McPHERSON TYING FORCEPS CURVED,
TITANIUM




- 6 mm tying platforms
- Tips 0.3 mm wide

Overall length 105 mm, 4.1 inches

- TC platforms
- Flat handle

TWS 85-3451/41
McPHERSON TYING FORCEPS CURVED,
TITANIUM




- 8 mm tying platforms
- Tips 0.3 mm wide

Overall length 105 mm, 4.1 inches

- TC platforms
- Flat handle

TYING FORCEPS


TWS 85-3515/41
TYING FORCEPS STRAIGHT,
TITANIUM



Overall length 105 mm, 4.1 inches

- 6 mm V-groove tying platforms
- Tips 0.5 mm wide
- TC platforms
- Flat handle


TWS 85-3512/04
TYING FORCEPS ANGLED,
TITANIUM



Overall length 102 mm, 4 inches

- 6 mm V-groove tying platforms
- Tips 0.5 mm wide
- TC platforms
- Flat handle


TWS 85-1035/41
TYING FORCEPS CURVED,
TITANIUM



Overall length 105 mm, 4.1 inches

- 6 mm V-groove tying platforms
- Tips 0.5 mm wide
- TC platforms
- Flat handle


TWS 85-3510/41
CATALANO TYING FORCEPS
STRAIGHT, TITANIUM



Overall length 105 mm, 4.1 inches

- 2 mm V-groove tying platforms
- Tips 0.35 mm wide
- TC platforms
- Flat handle


TWS 85-3501/41
CATALANO TYING FORCEPS
ANGLED, TITANIUM



Overall length 105 mm, 4.1 inches

- 2 mm V-groove tying platforms
- Tips 0.35 mm wide
- TC platforms
- Flat handle


TWS 85-3517/42
McPHERSON TYING FORCEPS
ANGLED, TITANIUM



Overall length 107 mm, 4.2 inches

- 11 mm tying platforms
- Tips 0.3 mm wide
- TC platforms
- Flat handle

TWS 85-3511/43
TYING FORCEPS STRAIGHT,
TITANIUM MS368&STAINLESS,
NO TC OF 388.02




Overall length 110 mm, 4.3 inches

- 6.5 mm tying platforms
- Tips 0.5 mm wide
- TC platforms
- Flat handle with holes

TYING FORCEPS

TWS 85-3610/39

CATALANO TYING FORCEPS
CURVED, TITANIUM




Overall length 100 mm, 3.9 inches

- 6 mm V-groove tying platforms
- Tips 0.5 mm wide
- TC platforms
- Round handle

TWS 85-1036/39

CATALANO TYING FORCEPS
STRAIGHT, TITANIUM




Overall length 100 mm, 3.9 inches

- 6 mm V-groove tying platforms
- Tips 0.5 mm wide
- TC platforms
- Round handle

TWS 85-0136/39

CATALANO TYING FORCEPS
ANGLED, TITANIUM




Overall length 100 mm, 3.9 inches

- 6 mm V-groove tying platforms
- Tips 0.5 mm wide
- TC platforms
- Round handle

TWS 85-3611/39

TENNANT TYING FORCEPS
STRAIGHT, TITANIUM





Overall length 101 mm, 3.9 inches

- 10 mm tying platforms
- Tips 0.25 mm wide
- TC platforms
- Round handle

TWS 85-1136/39

TENNANT TYING FORCEPS
CURVED, TITANIUM




Overall length 101 mm, 3.9 inches

- 10 mm tying platforms
- Tips 0.25 mm wide
- TC platforms
- Round handle

TWS 85-3615/41

GIRARD TYING FORCEPS
STRAIGHT, TITANIUM




Overall length 105 mm, 4.1 inches

- 4 mm tying platforms
- Tips 0.25 mm wide
- TC platforms
- Round handle

TWS 85-1536/41

GIRARD TYING FORCEPS
STRAIGHT, TITANIUM




Overall length 105 mm, 4.1 inches

- 6 mm tying platforms
- Tips 0.3 mm wide
- TC platforms
- Round handle

TYING FORCEPS

TWS 85-3715/41

GIRARD TYING FORCEPS STRAIGHT, TITANIUM




Overall length 105 mm, 4.1 inches

- 8 mm tying platforms
- Tips 0.35 mm wide
- TC platforms
- Round handle

TWS 85-3712/04

GIRARD TYING FORCEPS ANGLED, TITANIUM




Overall length 102 mm, 4 inches

- 4 mm tying platforms
- Tips 0.25 mm wide
- TC platforms
- Round handle

TWS 85-1237/04

GIRARD TYING FORCEPS ANGLED, TITANIUM





Overall length 102 mm, 4 inches

- 6 mm tying platforms
- Tips 0.3 mm wide
- TC platforms
- Round handle

TWS 85-2037/04

GIRARD TYING FORCEPS ANGLED, TITANIUM




Overall length 102 mm, 4 inches

- 8 mm tying platforms
- Tips 0.35 mm wide
- TC platforms
- Round handle

TWS 85-3713/04

TYING FORCEPS STRAIGHT, TITANIUM





Overall length 103 mm, 4 inches

- 6.2 mm tying platforms
- Tips 0.27 mm wide
- TC platforms
- Round handle

TWS 85-3702/04

TYING FORCEPS CURVED, TITANIUM




Overall length 102 mm, 4 inches

- 6.2 mm tying platforms
- Tips 0.27 mm wide
- TC platforms
- Round handle

TWS 85-1537/41

GIRARD TYING FORCEPS CURVED, TITANIUM




Overall length 105 mm, 4.1 inches

- 4 mm tying platforms
- Tips 0.25 mm wide
- TC platforms
- Round handle


TYING FORCEPS

TWS 85-3815/41
GIRARD TYING FORCEPS CURVED,
TITANIUM




- Overall length 105 mm, 4.1 inches
- 6 mm tying platforms
- Tips 0.3 mm wide
- TC platforms
- Round handle

TWS 85-1538/41
GIRARD TYING FORCEPS CURVED,
TITANIUM




- Overall length 105 mm, 4.1 inches
- 8 mm tying platforms
- Tips 0.35 mm wide
- TC platforms
- Round handle

TWS 85-3813/04
CATALANO TYING FORCEPS
STRAIGHT, TITANIUM




- Overall length 103 mm, 4 inches
- 2 mm V-groove tying platforms
- Tips 0.26 mm wide
- TC platforms
- Round handle

TWS 85-1338/04
CATALANO TYING FORCEPS
ANGLED, TITANIUM





- Overall length 103 mm, 4 inches
- 2 mm V-groove tying platforms
- Tips 0.35 mm wide
- TC platforms
- Round handle

TWS 85-3810/04
CATALANO TYING FORCEPS
CURVED, TITANIUM




- Overall length 103 mm, 4 inches
- 2 mm V-groove tying platforms
- Tips 0.35 mm wide
- TC platforms
- Round handle

TWS 85-3811/41
CATALANO TYING FORCEPS
CURVED, TITANIUM



- Overall length 105 mm, 4.1 inches
- 4 mm V-groove tying platforms
- Tips 0.25 mm wide
- TC platforms
- Round handle


TWS 85-3800/39
TYING FORCEPS FOR FEMTO / SMILE,
TITANIUM



- Overall length 100 mm, 3.9 inches
- 8.0mm delicate serrated tips
- Tips 0.3 mm wide
- Light round handle


TYING FORCEPS

TWS 85-3904/41
TYING FORCEPS STRAIGHT, TITANIUM




- Overall length 104 mm, 4.1 inches**
- 6 mm V-groove tying platforms
 - Tips 0.5 mm wide
 - TC platforms
 - Light round handle

TWS 85-3902/04
TYING FORCEPS ANGLD, TITANIUM




- Overall length 102 mm, 4 inches**
- 6 mm V-groove tying platforms
 - Tips 0.5 mm wide
 - TC platforms
 - Light round handle

TWS 85-3903/04
TYING FORCEPS CURVED, TITANIUM




- Overall length 103 mm, 4 inches**
- 6 mm V-groove tying platforms
 - Tips 0.5 mm wide
 - TC platforms
 - Light round handle

TWS 85-3940/41
TYING FORCEPS CURVED, TITANIUM




- Overall length 104 mm, 4.1 inches**
- 2 mm V-groove tying platforms
 - Tips 0.35 mm wide
 - Light round handle

TWS 85-3915/45
TYING FORCEPS CURVED, TITANIUM




- Overall length 115 mm, 4.5 inches**
- No platforms
 - Tips 0.3 mm wide
 - TC tips 10 mm long
 - Round handle

TWS 85-3551/45
TYING FORCEPS STRAIGHT, TITANIUM



- Overall length 115 mm, 4.5 inches**
- No platforms
 - Tips 0.3 mm wide
 - TC tips 10 mm long
 - Round handle


TWS 85-3923/48
TYING FORCEPS ANGLED, TITANIUM



- Overall length 123 mm, 4.8 inches**
- No platforms
 - Tips 0.3 mm wide
 - TC tips 9 mm long
 - Round handle

IRIS FORCEPS

TWS 85-4074/29
BISHOP-HARMON IRIS FORCEPS
STRAIGHT, TITANIUM



Overall length 74 mm, 2.9 inches

- 4 mm delicate serrated tips
- Tips 0.3 mm wide
- Flat handle


TWS 85-4072/28
BISHOP-HARMON IRIS FORCEPS ANGLED,
TITANIUM



Overall length 72 mm, 2.8 inches

- 4 mm delicate serrated tips
- Tips 0.3 mm wide
- Flat handle


TWS 85-7440/28
BISHOP-HARMON IRIS FORCEPS CURVED,
TITANIUM



Overall length 74 mm, 2.9 inches

- 4 mm delicate serrated tips
- Tips 0.3 mm wide
- Flat handle


TWS 85-7240/28
DVALI IRIS FORCEPS ANGLED, TITANIUM



Overall length 72 mm, 2.8 inches

- 1×2 lateral teeth 0.063 mm height
- TC teeth
- Flat handle


TWS 85-4075/29
BOTVIN IRIS FORCEPS FOR PERIPHERAL
IRIDECTOMY, TITANIUM



Overall length 75mm, 2.9 inches

- Angled 45 degree very delicate tips
- 1×2 teeth, 0.1 mm height
- TC teeth
- Tips 0.3 mm wide
- Flat handle


TWS 85-7540/29
BOTVIN IRIS FORCEPS FOR PERIPHERAL
IRIDECTOMY, TITANIUM



Overall length 75mm, 2.9 inches

- Angled 45 degree very delicate serrated tips 2 mm long
- Tips 0.4 mm wide
- Flat handle

TWS 85-5774/29
BONN IRIS FORCEPS STRAIGHT,
TITANIUM




Overall length 74 mm, 2.9 inches

- 1×2 teeth, 0.063 mm height
- 45 degree teeth
- 4 mm tying platforms
- Tips 0.2 mm wide
- TC teeth and platforms
- Flat handle

IRIS FORCEPS

TWS 85-4174/29
BONN IRIS FORCEPS STRAIGHT,
TITANIUM




- 1×2 teeth, 0.063 mm height
- 90 degree teeth
- 4 mm tying platforms

Overall length 74 mm, 2.9 inches

- Tips 0.2 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-7441/29
BONN IRIS FORCEPS STRAIGHT,
TITANIUM




- 1×2 teeth, 0.12 mm height
- 45 degree teeth
- 4 mm tying platforms

Overall length 74 mm, 2.9 inches

- Tips 0.3 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-7414/29
BONN IRIS FORCEPS STRAIGHT,
TITANIUM




- 1×2 teeth, 0.12 mm height
- 90 degree teeth
- 4 mm tying platforms

Overall length 74 mm, 2.9 inches

- Tips 0.3 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-4172/28
BONN IRIS FORCEPS ANGLED,
TITANIUM




- 1×2 teeth, 0.063 mm height
- 45 degree teeth
- 4 mm tying platforms

Overall length 72 mm, 2.8 inches

- Tips 0.2 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-4173/28
BONN IRIS FORCEPS CURVED,
TITANIUM




- 1×2 teeth, 0.063 mm height
- 45 degree teeth
- 4 mm tying platforms

Overall length 73 mm, 2.8 inches

- Tips 0.2 mm wide
- TC teeth and platforms
- Flat handle

TWS 85-4110/39
BONN IRIS FORCEPS STRAIGHT,
TITANIUM




- 1×2 teeth, 0.063 mm height
- 45 degree teeth
- 4 mm tying platforms

Overall length 100 mm, 3.9 inches

- Tips 0.2 mm wide
- TC teeth and platforms
- Flat handle

JEWELERS FORCEPS


TWS 85-4210/43
JEWELERS FORCEPS STRAIGHT,
TITANIUM



Overall length 110 mm, 4.3 inches

- Fine shafts with 0.2 mm sharp pointed tips
- TC shafts 6 mm
- Flat handle


TWS 85-4220/47
JEWELERS FORCEPS STRAIGHT,
TITANIUM



Overall length 120 mm, 4.7 inches

- Fine shafts with 0.2 mm sharp pointed tips
- TC shafts 6 mm
- Flat handle


TWS 85-4201/43
JEWELERS FORCEPS STRAIGHT,
TITANIUM



Overall length 110 mm, 4.3 inches

- Fine shafts with 0.2 mm sharp pointed tips
- TC tying platforms 6 mm
- Flat handle

TWS 85-4202/47
JEWELERS FORCEPS STRAIGHT,
TITANIUM



Overall length 120 mm, 4.7 inches

- Fine shafts with 0.2 mm sharp pointed tips
- TC tying platforms 6 mm
- Flat handle


TWS 85-4203/47
JEWELERS FORCEPS, TITANIUM



Overall length 120 mm, 4.7 inches

- Fine narrowed shafts with 0.2 mm sharp pointed tips
- TC shafts 6 mm
- Flat handle


TWS 85-4104/43
JEWELERS FORCEPS CURVED,
TITANIUM



Overall length 110 mm, 4.3 inches

- Fine shafts with 0.2 mm sharp pointed tips
- TC shafts 6 mm
- Flat handle

TWS 85-4015/43
JEWELERS FORCEPS CURVED,
TITANIUM




Overall length 110 mm, 4.3 inches

- Fine shafts with 0.2 mm sharp pointed tips
- TC tying platforms 6 mm
- Flat handle

JEWELERS FORCEPS


TWS 85-4312/47
JEWELERS FORCEPS, CURVED,
TITANIUM



Overall length 120 mm, 4.7 inches

- Fine narrowed shafts with 0.2 mm sharp pointed tips
- TC tying platforms 6 mm
- Flat handle


TWS 85-4310/43
JEWELERS FORCEPS
STRAIGHT, TITANIUM



Overall length 110 mm, 4.3 inches

- Fine narrowed shafts with 0.2 mm sharp pointed tips
- TC shafts 6 mm
- Flat handle


TWS 85-4311/43
JEWELERS FORCEPS
STRAIGHT, TITANIUM



Overall length 110 mm, 4.3 inches

- Fine narrowed shafts with 0.2 mm sharp pointed tips
- TC tying platforms 6 mm
- Flat handle


TWS 85-1143/43
JEWELERS FORCEPS
STRAIGHT, TITANIUM



Overall length 110 mm, 4.3 inches

- Fine narrowed shafts with 0.3 mm sharp pointed tips
- TC shafts 6 mm
- Flat handle


TWS 85-1043/43
JEWELERS FORCEPS STRAIGHT,
TITANIUM



Overall length 110 mm, 4.3 inches

- Fine narrowed shafts with 0.3 mm sharp pointed tips
- TC tying platforms 6 mm
- Flat handle


TWS 85-1034/43
JEWELERS FORCEPS STRAIGHT,
TITANIUM



Overall length 110 mm, 4.3 inches

- Fine narrowed shafts with 0.15 mm very sharp pointed tips
- TC shafts 6 mm
- Flat handle

TWS 85-1104/43
JEWELERS FORCEPS STRAIGHT,
TITANIUM




Overall length 110 mm, 4.3 inches

- Fine narrowed shafts with 0.15 mm very sharp pointed tips
- TC shafts 6 mm
- Flat handle

IOL FORCEPS


TWS 85-4400/39
McDONALD IOL IMPLANTATION
FORCEPS, TITANIUM



Overall length 100 mm, 3.9 inches

- Angled curved jaws
- Flat handle


TWS 85-4110/39
HALLENS INSERTION FORCEPS ANGLED,
TITANIUM
MS 103L WITH LOCK



Overall length 100 mm, 3.9 inches

- Round handle


TWS 85-4132/52
ICL CARTRIDGE LOADING
FORCEPS, ANGLED, TITANIUM



Overall length 132 mm, 5.2 inches

- Tips 0.45 × 3 mm wide
- Shafts 17 mm long
- Designed to remove the ICL from the container and push it inside the cartridge
- Round handle


TWS85-4415/41
IOL IMPLANTATION FORCEPS
ANGLED, TITANIUM



Overall length 105 mm, 4.1 inches

- Shafts 6.6 mm long
- Flat handle


TWS 85-4412/04
FAULKNER IOL
IMPLANTATION FORCEPS
ANGLED, TITANIUM



Overall length 102 mm, 4 inches

- Tips 0.3 mm dia
- Round shafts 11 mm long
- Flat handle


TWS85-4414/41
GODFREY IOL
IMPLANTATION FORCEPS
ANGLED, TITANIUM



Overall length 104 mm, 4.1 inches

- Tips 0.8 × 0.3 mm wide
- Shafts 12.3 mm long
- Flat handle

TWS 85-4411/43
CARTRIDGE LOADING
FORCEPS STRAIGHT,
TITANIUM




Overall length 110 mm, 4.3 inches

- Tips 1.2 mm wide
- Flat handle
- Designed for inserting Acrylic IOL into A, B, C, D cartridges
- To be used with injectors OJ 004, OJ 005, OJ 008, OJ 009

IOL FORCEPS

TWS 85-4555/60

KOCH UNIVERSAL PREFOLDER WITH LOCK, TITANIUM




Overall length 155 mm, 6.0 inches

- Round handle

TWS 85-4510/43

SIEBEL IOL FOLDING FORCEPS ANGLED, TITANIUM




Overall length 110 mm, 4.3 inches

- Flat handle

TWS 85-4505/41

BURATTO IOL IMPLANTATION FORCEPS ANGLED, TITANIUM




Overall length 105 mm, 4.1 inches

- Round handle

TWS 85-4510/43

IOL FOLDING FORCEPS ANGLED, TITANIUM




Overall length 110 mm, 4.3 inches

- Round handle

TWS 85-1145/43

IOL IMPLANTATION FORCEPS ANGLED, TITANIUM




Overall length 110 mm, 4.3 inches

- Flat handle

TWS 85-1545/41

BURATTO IOL IMPLANTATION FORCEPS ANGLED, TITANIUM




Overall length 105 mm, 4.1 inches

- Flat handle

TWS 85-1145/45

FORCEPS, TITANIUM



Overall length 115 mm, 4.5 inches

- Tips 0.65 × 0.6 mm wide
- Maximum opening - 27 mm
- Flat handle
- Screw mechanism

CHALAZION FORCEPS

TWS 85-4610/39

FRANCIS CHALAZION FORCEPS, TITANIUM



Overall length 100 mm, 3.9 inches

- Lower plate 20 mm × 12 mm, external size
- Upper plate window 17 mm × 9 mm, inner size
- Flat handle

TWS 85-1046/39

DESMARRES CHALAZION FORCEPS, TITANIUM



Overall length 100 mm, 3.9 inches

- Lower plate 34 mm × 24 mm, external size
- Upper plate window 31 mm × 21 mm, inner size
- Flat handle

TWS 85-4611/39

HUNT CHALAZION FORCEPS, TITANIUM



Overall length 100 mm, 3.9 inches

- Lower plate 14 mm dia, external size
- Upper plate window 11.6 mm, inner size
- Flat handle

TWS 85-1406/39

DESMARRES CHALAZION FORCEPS, TITANIUM



Overall length 100 mm, 3.9 inches

- Lower plate 12 mm × 11 mm, external size
- Upper plate window 10 mm × 9 mm, inner size
- Flat handle

TWS 85-4698/38

CHALAZION FORCEPS, TITANIUM



Overall length 98 mm, 3.8 inches

- Lower plate 17 mm × 15 mm, external size
- Upper plate window 14 mm × 12 mm, inner size
- Flat handle

TWS 85-9846/38

FRANCIS CHALAZION FORCEPS, TITANIUM




Overall length 98 mm, 3.8 inches

- Lower plate 15.5 mm d, external size
- Upper plate window 12.5 mm dia, inner size
- Flat handle

CONJUNTIVAL FORCEPS


TWS 85-4719/43
FECHTNER CONJUNCTIVAL
FORCEPS STRAIGHT, TITANIUM



Overall length 109 mm, 4.3 inches

- Delicate ring jaws
- Tips 0.75 mm external dia, 0.4 mm inner dia
- 6 mm tying platforms
- TC platforms
- Flat handle


TWS 85-4791/43
FECHTNER CONJUNCTIVAL
FORCEPS CURVED, TITANIUM



Overall length 109 mm, 4.3 inches

- Delicate ring jaws
- Tips 0.75 mm external dia, 0.4 mm inner dia
- 6 mm tying platforms
- TC platforms
- Flat handle


TWS 85-1094/43
FECHTNER CONJUNCTIVAL
FORCEPS STRAIGHT, TITANIUM



Overall length 109 mm, 4.3 inches

- Delicate ring jaws
- Two holes 0.75 mm external dia, 0.4 mm inner dia
- 6.8 mm tying platforms
- TC platforms
- Flat handle


TWS 85-4709/43
SOLISH FORCEPS STRAIGHT,
TITANIUM



Overall length 109 mm, 4.3 inches

- Tips 2 mm wide with two holes 0.4 mm dia
- 6 mm tying platforms
- TC platforms
- Flat handle

TWS 85-4917/43
SOLISH FORCEPS STRAIGHT,
TITANIUM




Overall length 109 mm, 4.3 inches

- Tips 2 mm wide with two holes 0.4 mm dia
- 6 mm tying platforms
- TC platforms
- Flat handle


DRESSING, LASIK, DSAEK FORCEPS

TWS 85-4813/04
MENDEZ FORCEPS FOR LASIK ANGLED,
TITANIUM




- Overall length 103 mm, 4 inches
- Tips 0.25 mm wide
- Curved tips 10 mm long
- Round handle

TWS 85-4818/46
OGAWA FORCEPS FOR DSAEK, ANGLED,
TITANIUM




- Overall length 118 mm, 4.6 inches
- Shafts 11.2 mm long
- Tips 0.85 mm wide
- 0.9 mm tying platforms
- Round handle

TWS 85-4875/29
FORCEPS COLIBRI TYPE, STAINLESS
MS 192.01T TITANIUM TC platform



- Overall length 75 mm, 2.9 inches
- Boot-shaped sharp tips 2 mm wide
- 2.5 mm tying platforms
- Used to remove sutures
- Flat handle



TWS 85-4816/45
MICROSURGICAL FORCEPS STRAIGHT,
TITANIUM



- Overall length 116 mm, 4.5 inches
- 6 mm working parts
- Tips 0.9 mm wide
- TC surface
- Round handle




TWS 85-4815/45
MICROSURGICAL FORCEPS CURVED, TITANIUM




- Overall length 115 mm, 4.5 inches
- 6 mm working parts
- Tips 0.9 mm wide
- TC surface
- Round handle

TWS 85-4814/41
NUGENT UTILITY FORCEPS ANGLED, TITANIUM



- Overall length 104 mm, 4.1 inches
- 9 mm cross-wise platforms
- Tips 1.1 mm wide
- Flat handle


TWS 85-4810/41
BONACCOLTO UTILITY
FORCEPS STRAIGHT, TITANIUM



- Overall length 104 mm, 4.1 inches
- 16 mm tying platforms
- Tips 1.2 mm wide
- Toothed tips 2.5 mm long
- Flat handle

DRESSING, DIFFERENT FORCEPS

TWS 85-4900/39
MICROSURGICAL FORCEPS
STRAIGHT, TITANIUM




- Serrated tips 12 mm long
- Tips 0.4 mm wide

Overall length 100 mm, 3.9 inches

- Flat handle

TWS 85-1049/39
MICROSURGICAL FORCEPS
STRAIGHT, TITANIUM




- Serrated tips 11.5 mm long
- Tips 0.37 mm wide

Overall length 100 mm, 3.9 inches

- Flat handle

TWS 85-4901/39
DRESSING FORCEPS
CURVED, TITANIUM




- Serrated tips 9 mm long
- Tips 0.7 mm wide

Overall length 101 mm, 3.9 inches

- Flat handle

TWS 85-4980/31
DRESSING FORCEPS
STRAIGHT, TITANIUM




- Delicate serrated tips 4 mm long
- Tips 0.4 mm wide

Overall length 80 mm, 3.1 inches

- Flat handle

TWS 85-4990/35
DRESSING FORCEPS
STRAIGHT, TITANIUM




- Delicate serrated tips 10 mm long
- Tips 0.6 mm wide

Overall length 90 mm, 3.5 inches

- Round handle

TWS 85-4905/41
DRESSING FORCEPS
STRAIGHT, TITANIUM




- Delicate serrated tips 8 mm long
- Tips 0.6 mm wide

Overall length 105 mm, 4.1 inches

- Flat handle

TWS 85-4986/34
DRESSING FORCEPS
STRAIGHT, TITANIUM




- Cross serrated tips 4 mm long
- Tips 0.8 x 0.8 mm wide

Overall length 86 mm, 3.4 inches

- Flat handle with holes

DRESSING, DIFFERENT FORCEPS


TWS 85-5000/43
DRESSING FORCEPS ANGLED, TITANIUM



- Delicate serrated tips 8 mm long
- Tips 0.6 mm wide
- Flat handle



TWS 85-5011/43
TUBE INSERTION FORCEPS FOR AHMED VALVES, TITANIUM




- Tips with groove for secure holding valve tube
- Tips 1.0 mm wide
- Flat handle

Overall length 111 mm, 4.3 inches



TWS 85-5008/42
TUBE INSERTION FORCEPS FOR AHMED VALVES, TITANIUM




- Tips with groove for secure holding valve tube
- Tip 0.95 mm wide
- Flat handle

Overall length 108 mm, 4.2 inches



TWS 85-5023/48
BROWN-ADSON FORCEPS, TITANIUM




- Brown teeth 0.5 mm height
- Serrated tips 4 mm long
- Flat handle

Overall length 123 mm, 4.8 inches



TWS 85-5020/47
ADSON TISSUE FORCEPS, TITANIUM




- 1 x 2 teeth 0.7 mm height
- TC teeth
- Cross-serrated tying platforms 15 mm long
- Flat handle

Overall length 120 mm, 4.7 inches



TWS 85-5085/33
FORCEPS FOR MEIBOMIAN GLANDS, TITANIUM




- Oval working tips 1 x 1.8 mm wide
- Tips 0.3 mm height
- Flat handle

Overall length 85 mm, 3.3 inches



TWS 85-5009/43
COMPRESSING LID FORCEPS (FOR MEIBOMIAN GLANDS), TITANIUM




- Working tips 4 x 10 mm wide
- TC tips
- Flat handle

Overall length 109 mm, 4.3 inches



DRESSING, DIFFERENT FORCEPS

TWS 85-5184/33
TROCAR PORTS GRASPING
FORCEPS, ANGLED, TITANIUM



Overall length 84 mm, 3.3 inches

- Delicate shafts 11 mm long
- Angled 0.4×0.4mm tips
- Flat handle
- Used for easy removing of trocar port from sclera

NEEDLE HOLDERS



TW Surgical

TWS 85-1000/39
COHAN NEEDLE HOLDER CURVED,
TITANIUM
OH 007.01 STRAIGHT



- Fine jaws 7.5 mm long
- Tips 0.4 mm wide

Overall length 100 mm, 3.9 inches

- TC jaws
- Round handle

TWS 85-1001/39
BARRAQUER NEEDLE HOLDER
CURVED, TITANIUM
OH 008.01 STRAIGHT



- Fine jaws 8.5 mm long
- Tips 0.4 mm wide

Overall length 100 mm, 3.9 inches

- TC jaws
- Round handle

TWS 85-1002/39
BARRAQUER NEEDLE HOLDER
CURVED, TITANIUM
OH015.01 STRAIGHT



- Fine jaws 9.5 mm long
- Tips 0.4 mm wide

Overall length 100 mm, 3.9 inches

- TC jaws
- Round handle

TWS 85-1003/39
NEEDLE HOLDER CURVED,
TITANIUM
MS100 STRAIGHT



- Fine jaws 5.5 mm long
- Tips 0.3 mm wide

Overall length 100 mm, 3.9 inches

- TC jaws
- Round handle

TWS 85-1040/41
NEEDLE HOLDER CURVED,
TITANIUM
OH 047.02 U STRAIGHT



- Fine jaws 11 mm long
- Tips 0.3 mm wide

Overall length 104 mm, 4.1 inches

- TC jaws
- Light round handle

TWS 85-1041/41
NEEDLE HOLDER CURVED
WITH UNIVERSAL LOCK,
TITANIUM



- Fine jaws 11 mm long
- Tips 0.3 mm wide

Overall length 104 mm, 4.1 inches

- TC jaws
- Light round handle

NEEDLE HOLDERS

TWS 85-1005/39
BECHERT- SINSKEY NEEDLE
HOLDER STRAIGHT,
TITANIUM



- Overall length 100 mm, 3.9 inches**
- Fine jaws 7.5 mm long
 - Tips 0.4 mm wide
 - TC jaws
 - Round handle

TWS 85-1050/41
NEEDLE HOLDER STRAIGHT,
TITANIUM
MS 025.01 CURVED



- Overall length 105 mm, 4.1 inches**
- Standard jaws 9 mm long
 - Tips 0.4 mm wide
 - TC jaws
 - Round handle double action

TWS 85-1051/41
BARRAQUER NEEDLE
HOLDER CURVED,
TITANIUM



- Overall length 105 mm, 4.1 inches**
- Standard jaws 7.7 mm long
 - Tips 0.5 mm wide
 - TC jaws
 - Round handle

TWS85-1060/42
NEEDLE HOLDER CURVED,
TITANIUM



- Overall length 106 mm, 4.2 inches**
- Delicate jaws 7 mm long
 - Tips 0.3 mm wide
 - TC jaws
 - Round handle

TWS 85-1101/43
McPHERSON NEEDLE
HOLDER CURVED, TITANIUM
OH 023 STRAIGHT



- Overall length 110 mm, 4.3 inches**
- Standard jaws 10.5 mm long
 - Tips 0.4 mm wide
 - TC jaws
 - Round handle

TWS 85-1251/43
CASTROVIEJO NEEDLE
HOLDER CURVED, TITANIUM



- Overall length 125 mm, 4.9 inches**
- Fine jaws 7.5 mm long
 - Tips 0.4 mm wide
 - TC jaws
 - Round handle

NEEDLE HOLDERS



TW Surgical

TWS 85-1252/49

CASTROVIEJO NEEDLE HOLDER
STRAIGHT, TITANIUM



- Fine jaws 7.5 mm long
- Tips 0.4 mm wide

Overall length 125 mm, 4.9 inches

- TC jaws
- Round handle

TWS 85-1253/49

BARRAQUER NEEDLE HOLDER
CURVED, TITANIUM



- Fine jaws 9.5 mm long
- Tips 0.4 mm wide

Overall length 125 mm, 4.9 inches

- TC jaws
- Round handle

TWS 85-1254/49

BECHERT-SINSKEY NEEDLE
HOLDER STRAIGHT, TITANIUM



- Fine jaws 8 mm long
- Tips 0.4 mm wide

Overall length 125 mm, 4.9 inches

- TC jaws
- Round handle

TWS 85-1255/49

BECHERT-SINSKEY NEEDLE
HOLDER STRAIGHT, TITANIUM



- Fine jaws 8 mm long
- Tips 0.4 mm wide

Overall length 125 mm, 4.9 inches

- TC jaws
- Round handle

TWS 85-1256/49

McPHERSON NEEDLE
HOLDER CURVED, TITANIUM



- Fine jaws 10.5 mm long
- Tips 0.4 mm wide

Overall length 125 mm, 4.9 inches

- TC jaws
- Round handle

TWS 85-1257/49

BARRAQUER NEEDLE
HOLDER CURVED
TITANIUM MS020.01
STRAIGHT



- Fine jaws 9.5 mm long
- Tips 0.4 mm wide

Overall length 125 mm, 4.9 inches

- TC jaws
- Round handle

Treasure Wright (U.K.) Ltd.
128 City Road, London, EC1V 2NX, UK
info@twrightmedical.com

NEEDLE HOLDERS



TW Surgical

TWS 85-1258/49
CASTROVIE NEEDLE HOLDER
CURVED, TITANIUM
OH 026.01 STRAIGHT



- Strong jaws 9.5 mm long
- Tips 0.5 mm wide

Overall length 125 mm, 4.9 inches

- TC jaws
- Round handle

TWS 85-1259/49
NEEDLE HOLDER CURVED,
TITANIUM
MS 056.01 STRAIGHT



- Delicate jaws 6 mm long
- Tips 0.3 mm wide

Overall length 125 mm, 4.9 inches

- TC jaws
- Round handle

TWS 85-5210/49
BARRAQUER NEEDLE HOLDER
CURVED WITH UNIVERSAL LOCK,
TITANIUM



- Strong jaws 9 mm long
- Tips 0.7 mm wide

Overall length 125 mm, 4.9 inches

- TC jaws
- Round handle

TWS 85-7210/50
NEEDLE HOLDER CURVED,
TITANIUM



- Strong jaws 10.5 mm long
- Tips 1.2 mm wide

Overall length 127 mm, 5 inches

- TC jaws
- Flat handle

TWS 85-9210/51
CASTROVIEJO NEEDLE HOLDER
CURVED, TITANIUM
MS 010 STRAIGHT



- Standard jaws 11 mm long
- Tips 0.4 mm wide

Overall length 129 mm, 5.1 inches

- TC jaws
- Round handle

TWS 85-5310/53
TENNANT NEEDLE HOLDER
CURVED, TITANIUM



- Standard jaws 11 mm long
- Tips 0.4 mm wide

Overall length 135 mm, 5.3 inches

- TC jaws
- Round handle

NEEDLE HOLDERS

TWS 85-1350/53
 BARRAQUER NEEDLE HOLDER
 CURVED WITH UNIVERSAL LOCK,
 TITANIUM



- Overall length 135 mm, 5.3 inches**
- Standard jaws 9.5 mm long
 - Tips 0.4 mm wide
 - TC jaws
 - Round handle

TWS 85-1351/53
 BARRAQUER NEEDLE HOLDER
 CURVED, TITANIUM
 OH 021.01 STRAIGHT OH 021.02
 LIGHT CURVED



- Overall length 135 mm, 5.3 inches**
- Strong jaws 8.5 mm long
 - Tips 0.6 mm wide
 - TC jaws
 - Round handle

TWS 85-1352/53
 CASTROVIEJO NEEDLE HOLDER
 STRAIGHT, TITANIUM



- Overall length 135 mm, 5.3 inches**
- Standard jaws 11 mm long
 - Tips 0.4 mm wide
 - TC jaws
 - Round handle

TWS 85-1353/53
 CASTROVIEJO NEEDLE HOLDER
 CURVED, TITANIUM
 OH 027.01 STRAIGHT



- Overall length 135 mm, 5.3 inches**
- Strong jaws 10.5 mm long
 - Tips 0.6 mm wide
 - TC jaws
 - Round handle

TWS 85-1354/54
 CASTROVIEJO NEEDLE HOLDER
 CURVED, TITANIUM



- Overall length 135 mm, 5.3 inches**
- Standard jaws 10.5 mm long
 - Tips 0.4 mm wide
 - TC jaws
 - Round handle

TWS 85-1355/53
 NEEDLE HOLDER CURVED,
 TITANIUM
 OH 045.01 STRAIGHT



- Overall length 135 mm, 5.3 inches**
- Strong jaws 9 mm long
 - Tips 0.5 mm wide
 - TC jaws
 - Flat handle

NEEDLE HOLDERS



TW Surgical

TWS 85-1356/53

NEEDLE HOLDER CURVED WITH
UNIVERSAL LOCK, TITANIUM



Overall length 135 mm, 5.3 inches

- Strong jaws 10 mm long
- Tips 0.8 mm wide
- TC jaws
- Light round handle

TWS 85-1357/53

CASTROVIEJO NEEDLE HOLDER
CURVED, TITANIUM
OH 048.01 STRAIGHT



Overall length 135 mm, 5.3 inches

- Strong jaws 11 mm long
- Tips 0.95 mm wide
- TC jaws
- Flat handle

TWS 85-1358/53

NEEDLE HOLDER CURVED,
TITANIUM



Overall length 135 mm, 5.3 inches

- Fine jaws 7.5 mm long
- Tips 0.4 mm wide
- TC jaws
- Round handle

TWS 85-1359/53

NEEDLE HOLDER CURVED,
TITANIUM



Overall length 135 mm, 5.3 inches

- Fine jaws 11 mm long
- Tips 0.3 mm wide
- TC jaws
- Round handle

TWS 85-5311/53

BARRAQUER NEEDLE HOLDER
CURVED WITH UNIVERSAL LOCK,
TITANIUM



Overall length 135 mm, 5.3 inches

- Strong jaws 10 mm long
- Tips 0.9 mm wide
- TC jaws
- Round handle

TWS 85-5312/53

WRIGHT NEEDLE HOLDER CURVED,
TITANIUM
MS 050.01 STRAIGHT



Overall length 135 mm, 5.3 inches

- Strong jaws 11 mm long
- Tips 0.95 mm wide
- TC jaws
- Round handle

NEEDLE HOLDERS

TWS 85-5313/53
KALT NEEDLE HOLDER
STRAIGHT, TITANIUM
OH 005.01 CURVED



Overall length 135 mm, 5.3 inches

- Strong jaws 10.5 mm long
- Tips 0.8 mm wide
- TC jaws
- Flat locking handle

TWS 85-5314/53
SNEAD NEEDLE HOLDER -
SCISSORS, STAINLESS



Overall length 135 mm, 5.3 inches

- Fine 2.5 mm long needle holder jaws
- Curved 6.5 mm long scissors blades
- TC jaws
- Round handle

TWS 85-1400/55
SNEAD NEEDLE HOLDER -
SCISSORS, STAINLESS



Overall length 140 mm, 5.5 inches

- Fine 3.5 mm long needle holder jaws
- Curved 7.5 mm long scissors blades
- TC jaws
- Round handle

TWS 85-1401/55
SNEAD NEEDLE HOLDER
-SCISSORS, TITANIUM
OH 029.01 CURVED



Overall length 140 mm, 5.5 inches

- Fine 3.5 mm long needle holder jaws
- Straight 7.5 mm long scissors blades
- TC jaws
- Round handle

TWS 85-1402/47
BARRAQUER NEEDLE HOLDER,
CURVED WITH UNIVERSAL LOCK,
TITANIUM



Overall length 120 mm, 4.7 inches

- Fine jaws 9.5 mm long
- Tips 0.4 mm wide
- TC jaws
- Light round handle

TWS 85-1403/46
DR. ZLOTY'S NEEDLE HOLDER -
FORCEPS CURVED TITANIUM
MS 041.01 STRAIGHT



Overall length 118 mm, 4.6 inches

- Needle Holder tips:
- Standard jaws 9 mm long
- Tips 0.8 mm wide
- TC jaws
- Forceps tips:
- Curved tips 0.4 mm wide
- TC V-grooved platform 6 mm
- Round handle swing joint

NEEDLE HOLDERS



TW Surgical

TWS 85-4210/49
APT'S NEEDLE HOLDER CURVED
WITH LOCK, TITANIUM
MS 058 L STRAIGHT



Overall length 124 mm, 4.9 inches

- Strong jaws 17 mm long
- Tips 0.75 mm wide
- TC jaws
- Round handle

TWS 85-9210/50
APT'S NEEDLE HOLDER CURVED
WITH LOCK, TITANIUM



Overall length 129 mm, 5.0 inches

- Fine jaws 23.3 mm long
- Tips 0.7 mm wide
- TC jaws
- Round handle

TWS 85-0031/51
APT'S NEEDLE HOLDER CURVED
WITH LOCK, TITANIUM



Overall length 130 mm, 5.1 inches

- Fine jaws 10.2 mm long
- Tips 0.6 mm wide
- TC jaws
- Round handle

TWS 85-0051/59
NEEDLE HOLDER CURVED,
TITANIUM



Overall length 150 mm, 5.9 inches

- Extra delicate jaws 12.3 mm long
- Tips 0.3 mm wide
- TC jaws
- Round handle

TWS 85-0050/59
NEEDLE HOLDER STRAIGHT,
TITANIUM



Overall length 150 mm, 5.9 inches

- Extra delicate jaws 12.3 mm long
- Tips 0.3 mm wide
- TC jaws
- Round handle

WESTCOTT SCISSORS



TW Surgical

TWS 85-0210/47

WESTCOTT TENOTOMY SCISSORS
CURVED, STAINLESS

MS 00 01 Ov. length 117mm, 4.6
inches, 13 mm blades



- 16.5 mm blades
- Blunt tips

Overall length 120 mm, 4.7 inches

- Mat finish
- Flat handle

TWS 85-0211/47

WESTCOTT TENOTOMY SCISSORS
CURVED, STAINLESS
OS 039.01 FLAT HANDLE



- 15 mm blades
- Blunt tips

Overall length 120 mm, 4.7 inches

- Mat finish
- Round handle

TWS 85-0212/47

WESTCOTT TENOTOMY SCISSORS
CURVED, STAINLESS-TITANIUM
MS 051.01 SHARP POINTED TIPS



- 15.5 mm blades, stainless
- Blunt tips

Overall length 120 mm, 4.7 inches

- Mat finish
- Flat handle, titanium

TWS 85-0213/47

WESTCOTT TENOTOMY SCISSORS
CURVED, STAINLESS-TITANIUM



- 16 mm blades, stainless
- Blunt tips

Overall length 120 mm, 4.7 inches

- Mat finish
- Flat handle, titanium

TWS 85-0215/47

WESTCOTT TENOTOMY SCISSORS
STRAIGHT, STAINLESS-TITANIUM
OS 054 CURVED

OS 054.02 STRAIGHT, sharp tips



- 17 mm blades, stainless
- Blunt tips

Overall length 120 mm, 4.7 inches

- Mat finish
- Round handle, titanium

TWS 85-0217/47

WESTCOTT STITCH SCISSORS
STRAIGHT, STAINLESS-TITANIUM



- 17 mm blades, stainless
- Sharp pointed tips

Overall length 120 mm, 4.7 inches

- Mat finish
- Flat handle, titanium

WESTCOTT SCISSORS



TW Surgical

TWS 85-0218/47
WESTCOTT STITCH SCISSORS
CURVED, STAINLESS-TITANIUM




- 16.5 mm blades, stainless
- Sharp pointed tips

Overall length 120 mm, 4.7 inches

- Mat finish
- Flat handle, titanium

TWS 85-1170/46
WESTCOTT STITCH
SCISSORS CURVED,
STAINLESS



- 12.5 mm blades
- Sharp pointed tips

Overall length 117 mm, 4.6 inches

- Mat finish
- Flat handle

TWS 85-0140/55
WESTCOTT STITCH SCISSORS
CURVED, STAINLESS



- 14.5 mm blades
- Sharp pointed tips

Overall length 140 mm, 5.5 inches

- Mat finish
- Flat handle

TWS 85-2102/47
WESTCOTT TENOTOMY
SCISSORS CURVED, TITANIUM



- 16.5 mm blades
- Blunt tips

Overall length 120 mm, 4.7 inches

- Mat finish
- Round handle

TWS 85-0101/43
WESTCOTT TENOTOMY
SCISSORS CURVED, STAINLESS



- 20 mm blades
- Blunt tips

Overall length 110 mm, 4.3 inches

- Mat finish
- Flat handle

TWS 85-0102/43
WESTCOTT STITCH SCISSORS
CURVED, STAINLESS-TITANIUM



- 11 mm blades, stainless
- Sharp pointed tips

Overall length 110 mm, 4.3 inches

- Mat finish
- Flat handle, titanium

WESTCOTT SCISSORS



TW Surgical

TWS 85-4101/45
WESTCOTT STITCH
SCISSORS CURVED,
STAINLESS



Overall length 114 mm, 4.5 inches

- 20.5 mm blades
- Sharp pointed tips
- Mat finish
- Flat handle


TWS 85-0120/47
WESTCOTT STITCH SCISSORS
CURVED, STAINLESS



Overall length 120 mm, 4.7 inches

- 16 mm blades
- Sharp pointed tips
- Mat finish
- Flat handle

TWS 85-0012/47
WESTCOTT TENOTOMY
SCISSORS CURVED, STAINLESS



Overall length 120 mm, 4.7 inches

- 16 mm blades
- Blunt tips
- Mat finish
- Flat handle

TWS 85-0011/43
WESTCOTT SCISSORS,
CURVED, STAINLESS



Overall length 110 mm, 4.3 inches

- 14 mm blades
- Mat finish
- Blunt tips
- Flat handle

TWS 85-1021/47
WESTCOTT STITCH
SCISSORS CURVED,
STAINLESS



Overall length 120 mm, 4.7 inches

- 17 mm blades
- Sharp pointed tips
- Mat finish
- Flat handle

TWS 85-1710/46
WESTCOTT TENOTOMY
SCISSORS CURVED,
STAINLESS



Overall length 117 mm, 4.6 inches

- 8.5 mm blades
- Blunt tips
- Mat finish
- Flat handle

VANNAS SCISSORS



TW Surgical

TWS 85-8400/33
VANNAS SCISSORS
STRAIGHT, STAINLESS



- 9.5 mm blades
- Sharp pointed tips

Overall length 84 mm, 3.3 inches

- Mat finish
- Flat handle

TWS 85-8000/31
VANNAS SCISSORS
STRAIGHT, STAINLESS



- 5.5 mm blades
- Sharp pointed tips

Overall length 80 mm, 3.1 inches

- Mat finish
- Flat handle

TWS 85-8001/31
VANNAS SCISSORS,
STRAIGHT, STAINLESS-
TITANIUM



- Blades 6.5 mm, stainless
- Mat finish

Overall length 80 mm, 3.1 inches

- Sharp pointed tips
- Flat handle, titanium

TWS 85-8401/33
VANNAS SCISSORS
CURVED, STAINLESS



- 9.5 mm blades
- Sharp pointed tips

Overall length 84 mm, 3.3 inches

- Mat finish
- Flat handle

TWS 85-8002/31
VANNAS SCISSORS
CURVED, STAINLESS



- 5.5 mm blades
- Sharp pointed tips

Overall length 80 mm, 3.1 inches

- Mat finish
- Flat handle

TWS 85-8100/32
VANNAS SCISSORS
ANGLED, STAINLESS



- 9.5 mm blades
- Sharp pointed tips

Overall length 81 mm, 3.2 inches

- Mat finish
- Flat handle

VANNAS SCISSORS

TWS 85-8101/32
VANNAS SCISSORS ANGLED,
STAINLESS



- 5.5 mm blades
- Sharp pointed tips

Overall length 81 mm, 3.2 inches

- Mat finish
- Flat handle

TWS 85-8102/32
VANNAS SCISSORS ANGLED
TO SIDE, STAINLESS



- 9.5 mm blades
- Sharp pointed tips

Overall length 81 mm, 3.2 inches

- Mat finish
- Flat handle

TWS 85-8003/31
VANNAS SCISSORS ANGLED TO
SIDE, STAINLESS



- 5.5 mm blades
- Sharp pointed tips

Overall length 80 mm, 3.1 inches

- Mat finish
- Flat handle

TWS 85-9000/35
VANNAS SCISSORS ANGLED TO
SIDE, STAINLESS
MTITANIUM



- 10 mm blades
- Sharp pointed tips

Overall length 90 mm, 3.5 inches

- Mat finish
- Flat handle

TWS 85-8400/33
VANNAS SCISSORS CURVED,
STAINLESS




- 8.5 mm blades
- Sharp pointed tips

Overall length 84 mm, 3.3 inches

- Mat finish
- Round handle

TWS 85-8103/31
VANNAS SCISSORS CURVED,
STAINLESS
S 016.02 STRAIGHT



- 6 mm blades
- Sharp pointed tips

Overall length 81 mm, 3.1 inches

- Mat finish
- Round handle

VANNAS SCISSORS

TWS 85-9001/35
VANNAS SCISSORS CURVED,
STAINLESS



- 8.5 mm blades
- Sharp pointed tips

Overall length 90 mm, 3.5 inches

- Mat finish
- Flat handle

TWS 85-1150/46
VANNAS SCISSORS STRAIGHT,
STAINLESS



- 9 mm blades
- Sharp pointed tips

Overall length 115 mm, 4.6 inches

- Mat finish
- Flat handle

TWS 85-1151/46
VANNAS SCISSORS CURVED,
STAINLESS



- 9 mm blades
- Sharp pointed tips

Overall length 115 mm, 4.6 inches

- Mat finish
- Flat handle

TWS 85-1152/46
VANNAS SCISSORS ANGLED,
STAINLESS



- 9 mm blades
- Sharp pointed tips

Overall length 115 mm, 4.6 inches

- Mat finish
- Flat handle

TWS 85-8401/33
VANNAS SCISSORS
CURVED, TITANIUM



- 8.5 mm blades
- Sharp pointed tips

Overall length 84 mm, 3.3 inches

- Mat finish
- Round handle

TWS 85-8402/33
VANNAS SCISSORS CURVED,
TITANIUM



- 9 mm blades
- Sharp pointed tips

Overall length 84 mm, 3.3 inches

- Mat finish
- Flat handle

VANNAS SCISSORS

TWS 85-8102/31
VANNAS SCISSORS,
STRAIGHT, TITANIUM



- 6mm blades
- Blunt tips

Overall length 81mm, 3.1 inches

- Mat finish
- Flat handle

TWS 85-4100/41
VANNAS SCISSORS CURVED,
STAINLESS-TITANIUM



- 8.5 mm blades, stainless
- Sharp pointed tips

Overall length 104 mm, 4.1 inches

- Mat finish
- Round handle, titanium

TWS 85-8004/31
VANNAS SCISSORS CURVED,
STAINLESS-TITANIUM



- 5.2 mm blades, stainless
- Sharp pointed tips

Overall length 80 mm, 3.1 inches

- Mat finish
- Flat handle, titanium

TWS 85-8005/31
VANNAS SCISSORS CURVED,
STAINLESS-TITANIUM



- 5.2 mm blades, stainless
- Blunt tips

Overall length 80 mm, 3.1 inches

- Mat finish
- Flat handle, titanium

TWS 85-8006/31
VANNAS SCISSORS STRAIGHT,
STAINLESS- TITANIUM



- 5.2 mm blades, stainless
- Sharp pointed tips

Overall length 80 mm, 3.1 inches

- Mat finish
- Flat handle, titanium

TWS 85-6011/42
VANNAS SCISSORS STRAIGHT,
STAINLESS-TITANIUM



- 6 mm blades, stainless
- Sharp pointed tips

Overall length 106 mm, 4.2 inches

- Mat finish
- Round double action handle, titanium

VANNAS SCISSORS

TWS 85-1060/42
VANNAS SCISSORS CURVED,
STAINLESS-TITANIUM
OS 027.02 ANGLED, blades 5.5 mm



- Overall length 106 mm, 4.2 inches**
- 6 mm blades, stainless
 - Sharp pointed tips
 - Mat finish
 - Round double-action handle, titanium

TWS 85-1107/43
VANNAS SCISSORS STRAIGHT,
STAINLESS- TITANIUM



- Overall length 110 mm, 4.3 inches**
- 7.2 mm blades, stainless
 - Sharp pointed tips
 - Mat finish
 - Flat handle , titanium

TWS 85-8008/31
VANNAS SCISSORS STRAIGHT,
STAINLESS
OS 068.01 SLIGHTLY CURVED



- Overall length 80 mm, 3.1 inches**
- 6.5 mm blades
 - Very sharp pointed tips
 - Mat finish
 - Flat handle

TWS 85-9601/38
VANNAS SCISSORS, CURVED,
TITANIUM/STAINLESS
MS 076 (TITANIUM)



- Overall length 96 mm, 3.8 inches**
- Blades 8 mm, stainless
 - Mat finish
 - Sharp tips
 - Flat handle with hole, titanium

TWS 85-8405/33
VANNAS CAPSULOTOMY
SCISSORS STRAIGHT, STAINLESS



- Overall length 84 mm, 3.3 inches**
- 6 mm blades
 - Sharp pointed tips
 - Mat finish
 - Flat handle

TWS 85-8500/33
VANNAS CAPSULOTOMY SCISSORS
STRAIGHT. STAINLESS
MS 322.01 CURVED



- Overall length 85 mm, 3.3 inches**
- 6 mm blades
 - Sharp pointed tips
 - Mat finish
 - Flat handle

NOYES SCISSORS

TWS 85-1051/41
NOYES SCISSORS STRAIGHT,
STAINLESS



- 12 mm blades
- Sharp pointed tips
- Mat finish
- Flat handle

Overall length 105 mm, 4.1 inches

TWS 85-1052/41
NOYES COMBINED SCISSORS
STRAIGHT, STAINLESS



- 12 mm blades
- 1 sharp and 1 blunt blades
- Mat finish
- Flat handle

Overall length 105 mm, 4.1 inches

TWS 85-1053/41
NOYES TENOTOMY SCISSORS
STRAIGHT, STAINLESS



- 12 mm blades
- Blunt tips
- Mat finish
- Flat handle

Overall length 105 mm, 4.1 inches

TWS 85-1054/41
TENOTOMY SCISSORS
STRAIGHT, STAINLESS



- 12 mm blades
- Blunt tips
- Mat finish
- Flat handle

Overall length 105 mm, 4.1 inches

TWS 85-9900/39
NOYES SCISSORS SLIGHTLY
CURVED, STAINLESS



- 10 mm blades
- Sharp pointed tips
- Mat finish
- Flat handle

Overall length 99 mm, 3.9 inches

CORNEAL SCISSORS

TWS 85-1059/41
 CASTROVIEJO MICRO SCISSORS
 CURVED, STAINLESS



Overall length 105 mm, 4.1 inches

- Strongly curved 9 mm blades
- Sharp pointed tips
- Mat finish
- Round handle

TWS 85-0001/39
 CASTROVIEJO MINI
 SCISSORS CURVED, STAINLESS



Overall length 100 mm, 3.9 inches

- 8 mm blades
- Sharp pointed tips
- Mat finish
- Flat handle

TWS 85-0002/39
 CASTROVIEJO MINI SCISSORS
 CURVED, TITANIUM



Overall length 100 mm, 3.9 inches

- 8 mm blades
- Sharp pointed tips
- Mat finish
- Flat handle

TWS 85-9405/37
 CASTROVIEJO MINI SCISSORS
 CURVED, STAINLESS-TITANIUM OS
 045.01 BLUNT TIPS



Overall length 94 mm, 3.7 inches

- 7 mm blades, stainless
- Sharp pointed tips
- Mat finish
- Flat handle, titanium

TWS 85-1057/41
 CASTROVIEJO UNIVERSAL CORNEAL
 SCISSORS CURVED, STAINLESS



Overall length 105 mm, 4.1 inches

- 14 mm blades
- Blunt tips
- Mat finish
- Flat handle

TWS 85-0011/43
 CASTROVIEJO UNIVERSAL
 CORNEAL SCISSORS CURVED,
 STAINLESS
 IMS 018.02 Flat handle
 MS 018.03 11.5mm blades



Overall length 110 mm, 4.3 inches

- 13.5 mm blades
- Blunt tips
- Mat finish
- Round handle

CORNEAL SCISSORS

TWS 85-1107/43
 CASTROVIEJO UNIVERSAL
 CORNEAL SCISSORS
 CURVED, TITANIUM



- 13.5 mm blades
- Blunt tips

Overall length 110 mm, 4.3 inches

- Mat finish
- Round handle

TWS 85-1106/43
 CASTROVIEJO UNIVERSAL
 CORNEAL SCISSORS CURVED,
 STAINLESS-TITANIUM



- 13.5 mm blades, stainless
- Blunt tips

Overall length 110 mm, 4.3 inches

- Mat finish
- Round handle, titanium

TWS 85-1057/41
 TROUTMAN-CASTROVIEJO
 CORNEAL SCISSORS RIGHT,
 STAINLESS



- 9.5 mm blades
- Blunt tips

Overall length 105 mm, 4.1 inches

- Mat finish
- Flat handle

TWS 85-1056/41
 TROUTMAN-CASTROVIEJO CORNEAL
 SCISSORS
 LEFT, STAINLESS



- 9.5 mm blades
- Blunt tips

Overall length 105 mm, 4.1 inches

- Mat finish
- Flat handle

TWS 85-0015/39
 CASTROVIEJO CORNEAL
 SCISSORS RIGHT, STAINLESS



- 11mm blades
- Blunt tips

Overall length 100 mm, 3.9 inches

- Mat finish
- Flat handle

TWS 85-0016/39
 CASTROVIEJO CORNEAL SCISSORS
 LEFT, STAINLESS



- 11mm blades
- Blunt tips

Overall length 100 mm, 3.9 inches

- Mat finish
- Flat handle

CORNEAL SCISSORS

TWS 85-0014/39

CASTROVIEJO CORNEAL SCISSORS
RIGHT, TITANIUM



- 11 mm blades
- Blunt tips

Overall length 100 mm, 3.9 inches

- Mat finish
- Flat handle

TWS 85-0015/39

CASTROVIEJO CORNEAL SCISSORS
LEFT, TITANIUM



- 11 mm blades
- Blunt tips

Overall length 100 mm, 3.9 inches

- Mat finish
- Flat handle

TWS 85-1057/41

CASTROVIEJO CORNEAL SCISSORS
RIGHT, STAINLESS



- 8.5 mm blades
- Blunt tips

Overall length 105 mm, 4.1 inches

- Mat finish
- Flat handle

TWS 85-1055/41

CASTROVIEJO CORNEAL SCISSORS
LEFT, STAINLESS



- 8.5mm blades
- Blunt tips

Overall length 105 mm, 4.1 inches

- Mat finish
- Flat handle

TWS 85-1042/41

TROUTMAN-CASTROVIEJO
CORNEAL SCISSORS RIGHT,
STAINLESS-TITANIUM



- 9.5 mm blades, stainless
- Blunt tips

Overall length 104 mm, 4.1 inches

- Mat finish
- Flat handle, titanium

TWS 85-1043/41

TROUTMAN-CASTROVIEJO
CORNEAL SCISSORS LEFT,
STAINLESS-TITANIUM



- 9.5 mm blades, stainless
- Blunt tips

Overall length 104 mm, 4.1 inches

- Mat finish
- Flat handle, titanium

CORNEAL SCISSORS

TWS 85-1020/40
 TROUTMAN-CASTROVIEJO CORNEAL
 SCISSORS RIGHT, STAINLESS-TITANIUM



- 6.5 mm blades, stainless
- Mat finish

Overall length 102 mm, 4 inches

- Blunt tips
- Flat handle, titanium

TWS 85-1021/40
 TROUTMAN-CASTROVIEJO CORNEAL
 SCISSORS LEFT, STAINLESS-TITANIUM



- 6.5mm blades, stainless
- Mat finish

Overall length 102 mm, 4 inches

- Blunt tips
- Flat handle, titanium

TWS 85-1044/41
 CASTROVIEJO CORNEAL SCISSORS
 RIGHT, STAINLESS-TITANIUM



- 10.5mm blades, stainless
- Blunt tips

Overall length 104 mm, 4.1 inches

- Mat finish
- Flat titanium handle

TWS 85-1045/41
 CASTROVIEJO CORNEAL
 SCISSORS LEFT,
 STAINLESS-TITANIUM



- 10.5 mm blades, stainless
- Blunt tips

Overall length 104 mm, 4.1 inches

- Mat finish
- Flat titanium handle

TWS 85-1032/40
 CASTROVIEJO CORNEAL
 SCISSORS LEFT,
 STAINLESS-TITANIUM



- 7.0mm blades, stainless
- Mat finish/ Blunt tips

Overall length 103 mm, 4 inches

- Flat handle, titanium

TWS 85-1033/40
 CASTROVIEJO CORNEAL
 SCISSORS RIGHT,
 STAINLESS-TITANIUM



- 7.0mm blades, stainless
- Mat finish/ Blunt tips

Overall length 103 mm, 4 inches


- Flat handle, titanium



TW Surgical

CORNEAL SCISSORS

TWS 85-0108/39
KATZIN CORNEAL TRANSPLANT
SCISSORS RIGHT, STAINLESS
MS 017.01 LEFT




Overall length 100 mm, 3.9 inches

- Strongly curved 7 mm blades
- Sharp pointed tips

- Mat finish
- Round handle

CONJUNCTIVAL SCISSORS

TWS 85-1080/42
CONJUNCTIVAL SCISSORS CURVED,
STAINLESS



- 11 mm blades
- Sharp pointed tips

Overall length 108 mm, 4.2 inches

- Mat finish
- Round handle

TWS 85-0011/39
CONJUNCTIVAL SCISSORS CURVED,
STAINLESS
OS014.01 STRAIGHT



- 11.5 mm blades
- Sharp pointed tips

Overall length 100 mm, 3.9 inches

- Mat finish
- Flat handle

TWS 85-0014/55
SCISSORS FOR PLASTIC SURGERY
CURVED, STAINLESS- TITANIUM



- 20mm serrated blades, stainless
- Blunt tips

Overall length 140 mm, 5.5 inches

- Mat finish
- Flat handle, titanium

TWS 85-1089/42
SCISSORS STRAIGHT. STAINLESS MS
062.01 CURVED



- 9.7 mm blades
- Sharp pointed tips

Overall length 108 mm, 4.2 inches

- Mat finish
- Flat handle

TWS 85-0121/44
SCISSORS ANGLED, STAINLESS-
TITANIUM



- 11.5 mm blades, stainless
- Blunt tips

Overall length 112 mm, 4.4 inches

- Mat finish
- Round handle, titanium

TWS 85-5700/22
BARRAQUER IRIS SCISSORS
ANGLED, STAINLESS (BUTTERFLY)



- 7 mm blades
- Sharp tips

Overall length 57 mm, 2.2 inches

- Mat finish
- Flat handle

ANTERIOR CAPSULOTOMY SCISSORS

TWS 85-9000/35
ANTERIOR CAPSULOTOMY
SCISSORS ANGLED, STAINLESS



- Overall length 90 mm, 3.5 inches**
- 10 mm blades
 - Sharp pointed tips
 - Mat finish
 - Flat handle

TWS 85-9001/35
ANTERIOR CAPSULOTOMY
SCISSORS, STAINLESS



- Overall length 90 mm, 3.5 inches**
- 10 mm blades
 - Sharp pointed tips
 - Mat finish
 - Flat handle

TWS 85-8209/32
ANTERIOR CAPSULOTOMY
SCISSORS ANGLED, STAINLESS



- Overall length 82 mm, 3.2 inches**
- 10 mm blades
 - Sharp pointed tips
 - Mat finish
 - Round handle

TWS 85-8208/32
ANTERIOR CAPSULOTOMY
SCISSORS ANGLED, STAINLESS-
TITANIUM



- Overall length 82 mm, 3.2 inches**
- 10 mm blades, stainless
 - Sharp pointed tips
 - Mat finish
 - Round handle, titanium

TWS 85-8010/42
MUSCLE SCISSORS ANGLED RIGHT,
STAINLESS-TITANIUM



- Overall length 108 mm, 4.2 inches**
- 11.5 mm blades, stainless
 - Blunt tips
 - Mat finish
 - Flat handle, titanium


TWS 85-8011/42
MUSCLE SCISSORS ANGLED LEFT,
STAINLESS-TITANIUM



- Overall length 108 mm, 4.2 inches**
- 11.5 mm blades, stainless
 - Blunt tips
 - Mat finish
 - Flat handle, titanium

LENS CUTTERS


TWS 85-2118/44
LENS CUTTER FORCEPS,
STAINLESS-TITANIUM



Overall length 112 mm, 4.4 inches

- 8.5 mm blades, stainless
- Flat handle with two hinge joints, titanium
- Designed to cut and remove all types of foldable intraocular lens through a small incision


TWS 85-2119/44
LENS CUTTER FORCEPS, STAINLESS-TITANIUM



Overall length 112 mm, 4.4 inches

- 9.5 mm blades, stainless
- Designed to cut and remove all types of foldable intraocular lens through a small incision
- Flat handle with two hinge joints, titanium
- Upper blade of the scissors has five bevels for effective cutting and grasping of the IOL


TWS 85-2117/44
LENS CUTTER FORCEPS



Overall length 112mm, 4.4 inches

- 7 mm blades, stainless
- Upper blade has 5 bevels for effective cutting and grasping of IOL
- Flat titanium handle with two hinge joints


TWS 85-1018/43
LENS CUTTER FORCEPS,
STAINLESS-TITANIUM



Overall length 110 mm, 4.3 inches

- 7.8 mm blades, stainless
- Flat handle with two hinge joints, titanium


TWS 85-1017/43
LENS REMOVING FORCEPS



Overall length 110 mm, 4.3 inches

- 11 mm blades, stainless
- Flat handle with two hinge joints, titanium
- Designed to grasp and remove all types of foldable intraocular lens through a small incision

TWS 85-1145/48
LENS FORCEPS



Overall length 114mm, 4.48 inches

- 7 mm blades, stainless
- Upper blade has 1 tooth
- Lower blade has 2 teeth
- Flat titanium handle with two hinge joints

EYE SPECULUMS



TW Surgical

TWS 85-0080/31

EYE SPECULUM TEMPORAL WITH MOVABLE SUPPORTS



Overall length 80 mm, 3.1 inches

- Closed wire blades 15 mm wide with movable supports, stainless
- Mechanism, titanium

TWS 85-0078/30

EYE SPECULUM, TITANIUM




Overall length 78 mm, 3.0 inches

- Closed wire blades 15 mm wide

TWS 85-0070/28

MURDOCK EYE SPECULUM, RIGHT, TITANIUM




Overall length 70 mm, 2.8 inches

- Closed wire blades 18 mm wide
- Locking mechanism

TWS 85-0071/28

MURDOCK EYE SPECULUM, LEFT, TITANIUM



Overall length 70 mm, 2.8 inches

- Closed wire blades 18 mm wide
- Locking mechanism

TWS 85-5400/21

SPRING EYE SPECULUM, TITANIUM



Overall length 54 mm, 2.1 inches

- Open wire blades 15mm wide

TWS 85-5401/21

WRIGHT SPRING EYE SPECULUM, TITANIUM



Overall length 54 mm, 2.1 inches

- Solid blades 15mm wide


EYE SPECULUMS



TW Surgical

TWS 85-0343/17

ARRAQUER EYE SPECULUM ADULT, STAINLESS




Overall length 43 mm, 1.7 inches

- Closed wire blades 15 mm wide
- Dull finish

TWS 85-0404/16

BARRAQUER EYE SPECULUM PEDIATRIC, STAINLESS




Overall length 40 mm, 1.6 inches

- Closed wire blades 10 mm wide
- Dull finish

TWS 85-0383/15

KRATZ- BARRAQUER EYE SPECULUM PEDIATRIC, STAINLESS



Overall length 38 mm, 1.5 inches

- Open wire blades 10 mm wide
- Dull finish

TWS 85-0404/17

KRATZ- BARRAQUER EYE SPECULUM ADULT, STAINLESS




Overall length 44 mm, 1.7 inches

- Open wire blades 14 mm wide
- Dull finish

TWS 85-0440/17

BARRAQUER EYE SPECULUM STRONG ADULT, STAINLESS




Overall length 44 mm, 1.7 inches

- Closed wire blades 15 mm wide
- Dull finish

TWS 85-0353/14

WIRE EYE SPECULUM NEWBORN, STAINLESS



Overall length 35 mm, 1.4 inches

- Closed wire blades 5 mm wide
- Dull finish

TWS 85-0141/16

KRATZ- BECHERT EYE SPECULUM ADULT, STAINLESS



Overall length 41 mm, 1.6 inches

- Open wire blades 14.5 mm wide
- Dull finish

TWS 85-0411/16

BARRAQUER EYE SPECULUM PEDIATRIC, STAINLESS



Overall length 41 mm, 1.6 inches

- Solid blades 11 mm wide
- Dull finish

TWS 85-0242/17

BARRAQUER EYE SPECULUM ADULT, STAINLESS



Overall length 42 mm, 1.7 inches

- Solid blades 15 mm wide
- Dull finish

TWS 85-0474/19

BARRAQUER EYE SPECULUM ADULT, STAINLESS




Overall length 47 mm, 1.9 inches

- Solid blades 20 mm wide
- Dull finish

EYE SPECULUMS

TWS 85-0434/17

BARRAQUER EYE SPECULUM ADULT, STAINLESS




Overall length 43 mm, 1.7 inches

- Closed wire blades 15 mm wide
- Dull finish

TWS 85-0409/16

BARRAQUER EYE SPECULUM PEDIATRIC, STAINLESS




Overall length 40 mm, 1.6 inches

- Closed wire blades 10 mm wide
- Dull finish

TWS 85-8339/15

KRATZ- BARRAQUER EYE SPECULUM PEDIATRIC, STAINLESS




Overall length 38 mm, 1.5 inches

- Open wire blades 10 mm wide
- Dull finish

TWS 85-0944/17

KRATZ- BARRAQUER EYE SPECULUM ADULT, STAINLESS



Overall length 44 mm, 1.7 inches

- Open wire blades 14 mm wide
- Dull finish

TWS 85-0844/17

BARRAQUER EYE SPECULUM STRONG ADULT, STAINLESS




Overall length 44 mm, 1.7 inches

- Closed wire blades 15 mm wide
- Dull finish

TWS 85-0935/14

WIRE EYE SPECULUM NEWBORN, STAINLESS



Overall length 35 mm, 1.4 inches

- Closed wire blades 5 mm wide
- Dull finish

TWS 85-0941/16

KRATZ- BECHERT EYE SPECULUM ADULT, STAINLESS



Overall length 41 mm, 1.6 inches

- Open wire blades 14.5 mm wide
- Dull finish

TWS 85-0841/16

BARRAQUER EYE SPECULUM



PEDIATRIC, STAINLESS

Overall length 41 mm, 1.6 inches

- Solid blades 11 mm wide
- Dull finish

TWS 85-0742/16

BARRAQUER EYE SPECULUM ADULT, STAINLESS



Overall length 42 mm, 1.7 inches

- Solid blades 15 mm wide
- Dull finish

TWS 85-0747/19

BARRAQUER EYE SPECULUM ADULT, STAINLESS



Overall length 47 mm, 1.9 inches

- Solid blades 20 mm wide
- Dull finish

EYE SPECULUMS

TWS 85-0988/35
LIEBERMAN ADJUSTABLE EYE SPECULUM, ADULT, TITANIUM



Overall length 88 mm, 3.5 inches

- Open wire blades 16 mm wide

TWS 85-0979/31
LIEBERMAN ADJUSTABLE EYE SPECULUM, PEDIATRIC, TITANIUM



Overall length 79 mm, 3.1 inches

- Open wire blades 12 mm wide

TWS 85-0982/32
ADJUSTABLE EYE SPECULUM, ADULT, TITANIUM



Overall length 82 mm, 3.2 inches

- Solid blades 15 mm wide

TWS 85-0780/31
ADJUSTABLE EYE SPECULUM, PEDIATRIC, TITANIUM



Overall length 80 mm, 3.1 inches

- Solid blades 11 mm wide

TWS 85-0782/32
KERSHNER REVERSIBLE EYE SPECULUM, ADULT, TITANIUM



Overall length 82 mm, 3.2 inches

- Solid blades 15 mm wide
- Can be used for both nasal & temporal approach

TWS 85-0680/31
LIEBERMAN ADJUSTABLE EYE SPECULUM, TITANIUM OR 013.02 STAINLESS



Overall length 80 mm, 3.1 inches

- Open wire blades 15 mm wide

TWS 85-0480/31
LIEBERMAN ADJUSTABLE EYE SPECULUM, STAINLESS OR 013.011 TITANIUM



Overall length 80 mm, 3.1 inches

- Closed wire blades 15 mm wide

TWS 85-0380/31
LIEBERMAN ADJUSTABLE EYE SPECULUM, REVERSIBLE, TITANIUM



Overall length 80 mm, 3.1 inches

- Solid blades 15 mm wide

EYE SPECULUMS

TWS 85-0171/28
 UNIVERSAL EYE SPECULUM WITH
 REMOVABLE BLADES,
 STAINLESS- TITANIUM
 MS 015.01 TITANIUM



Overall length 71 mm, 2.8 inches

- Removable blades, stainless
- Adult wire blades 16.5 mm
- Pediatric wire blades 13 mm
- Adult solid blades 15 mm
- Pediatric solid blades 12 mm

TWS 85-0172/28
 ASPIRATION ADJUSTABLE EYE
 SPECULUM, LIEBERMAN
 STYLE, STAINLESS-TITANIUM



Overall length 72 mm, 2.8 inches

- Open wire blades 16 mm wide, stainless
- Maximum opening 50 mm

TWS 85-0272/28
 ASPIRATION ADJUSTABLE EYE
 SPECULUM, LIEBERMAN
 STYLE, STAINLESS-TITANIUM



Overall length 72 mm, 2.8 inches

- Open wire blades 16 mm wide, stainless
- Maximum opening 50 mm
- Designed for LASIK

TWS 85-0178/31
 ASPIRATION ADJUSTABLE EYE
 SPECULUM, LIEBERMAN
 STYLE, STAINLESS-TITANIUM



Overall length 78 mm, 3.1 inches

- Open wire blades 16 mm wide, stainless
- Maximum opening 50 mm

TWS 85-0375/29
 ASPIRATION ADJUSTABLE
 EYE SPECULUM, PEDIATRIC,
 STAINLESS-TITANIUM



Overall length 75 mm, 2.9 inches

- Open wire blades 12 mm wide, stainless
- Maximum opening 50 mm

TWS 85-758/23
 ASPIRATION ADJUSTABLE EYE
 SPECULUM, STAINLESS



Overall length 58 mm, 2.3 inches

- Open wire blades 18 mm wide
- Maximum opening 38 mm

EYE SPECULUMS

TWS 85-0688/35
SHAPIRA EYE SPECULUM, TITANIUM



- Overall length 88 mm, 3.5 inches**
- Open wire blades 14 mm wide

TWS 85-0169/27
EYE SPECULUM LATH, STAINLESS




- Overall length 69 mm, 2.7 inches**
- Closed wire blades 15 mm wide

TWS 85-0290/35
ADJUSTABLE EYE SPECULUM, SOLID, TITANIUM



- Overall length 90 mm, 3.5 inches**
- Solid blades 15 mm wide
 - Maximum opening 34 mm

TWS 85-0348/18
ADJUSTABLE EYE SPECULUM, SOLID, TITANIUM



- Overall length 48 mm, 1.8 inches**
- TWS 85-0311/05 Solid blades 11 mm and 5

TWS 85-0648/18
EYE SPECULUM LATH PREMATURE, TITANIUM



- Overall length 48 mm, 1.8 inches**
- TWS 85-0311/05 Solid blades 11 mm and 5
 - TWS 85-0305/03 Solid blades 5 mm and 3 mm wide
 - TWS 85-0306/04 Solid blades 6 mm and 4 mm wide
 - TWS 85-0307/05 Solid blades 7 mm and 5 mm wide

TWS 85-0273/29
OLIVIERI-CANTERA ADJUSTABLE FEMTO EYE SPECULUM,
STAINLESS
OR 038 TITANIUM



- Overall length 73 mm, 2.9 inches**
- Solid blades 14 mm, strong action with ideal curvature for accommodating pivoting microkeratomes
 - Maximum opening 56 mm

TWS 85-0310/39
ADJUSTABLE EYE SPECULUM ADULT, TITANIUM
STAINLESS
OR 038 TITANIUM



- Overall length 100 mm, 3.9 inches**
- Closed wire blades 15 mm wide accommodating pivoting microkeratomes

TWS 85-0311/39
ADJUSTABLE EYE SPECULUM ADULT, TITANIUM



- Overall length 100 mm, 3.9 inches**
- Closed wire blades 15 mm wide

RETRACTORS

TWS 85-1920/47
DESMARRES LID RETRACTOR -SIZE 0



Overall length 120 mm, 4.7 inches

- Round handle, titanium
- Working part, stainless
- **OR 006** - SIZE 0, 10 mm wide
- **OR 006.01** - SIZE 1, 12 mm wide
- **OR 006.02** - SIZE 2, 14 mm wide
- **OR 006.03** - SIZE 3, 16 mm wide
- **OR 006.04** - SIZE 00, 8 mm wide
- **OR 006.05** - SIZE 0, 10 mm wide, titanium
- **OR 025** - SIZE 00, 7.5 mm wide

TWS 85-0915/59
LID RETRACTOR, DOUBLE-END



Overall length 150 mm, 5.9 inches

- Curved working parts, stainless
- 7 mm and 9 mm wide
- Round handle, titanium

TWS 85-0519/58
LID RETRACTOR, DOUBLE-END



Overall length 149 mm, 5.8 inches

- Working parts, stainless
- 10 mm and 12 mm wide
- Round handle, titanium

TWS 85-0381/73
IRIS RETRACTOR, DOUBLE- ENDE



Overall length 183 mm, 7.3 inches

- Working parts, stainless
- Retractor 3.5 mm wide
- Mushroom-shape tip 2.2 mm dia
- Round handle, titanium

TWS 85-0715/58
SCEPHENS ORBITAL RETRACTOR



Overall length 150 mm, 5.8 inches

- Working part with a slit 4 mm wide
- Tube on the top stainless
- Flat handle, titanium

TWS 85-0669/27
STEVENSON LACRIMAL OCULAR
RETRACTOR, TITANIUM



Overall length 69 mm, 2.7 inches

- 4 x 4 blunt prongs

RETRACTORS

TWS 85-0991/47
TISSUE RETRACTOR



Overall length 119 mm, 4.7 inches

- Curved working part 7 × 0.4 mm, stainless
- Round handle, titanium
- To retract conjunctival tissue and maximize exposure of the sclera for making and placing the sutures

TWS 85-6671/69
DOUBLE RETRACTOR-ELEVATED ENDED
ORTMS, BITAL STAIEGLOBNLESS



Overall length 176 mm, 6.9 inches

- 10 mm wide and 15.1 mm wide ends
- Thin ribbon retractor with a gentle "S" curve

TWS 85-6672/69
DOUBLERETRACTOR-ELE ENDED
VAORTMS, BITAL ESTAIGLOBNLESS



Overall length 176 mm, 6.9 inches

- 12 mm wide and 17.2 mm wide ends
- Thin ribbon retractor with a gentle "S" curve

TWS 85-6673/69
DOUBLE RETRACTOR-ELEVATED ENDED
ORTMS, BITAL STAIEGLOBNLESS



Overall length 176 mm, 6.9 inches

- 8 mm wide and 13.6 mm wide ends
- Thin ribbon retractor with a gentle "S" curve

TWS 85-6631/54
KNAPP LACHRYMAL SAC
RETRACTOR, STRAIGHT



- Working part with 8.1 mm wide, stainless
- 4 blunt prongs
- Round handle, titanium

TWS 85-6613/52
KNAPP LACHRYMAL SAC
RETRACTOR, STRAIGHT



Overall length 133 mm, 5.2 inches

- Working part 7.6 mm wide, stainless
- 4 blunt prongs
- Round handle, titanium

TWS 85-6621/64
DOUBLE END RETRAC TITANIUM




Overall length 162 mm, 6.4 inches

- One tip with 3 blunt prongs 9 mm wide
- Other tip 7.6 mm wide
- Round handle

VARIOUS INSTRUMENTS

TWS 85-0312/47

CASTROVIEJO BLADE HOLDER
STRAIGHT, TITANIUM




Overall length 120 mm, 4.7 inches

- TC conical jaws 13 mm long
- Round handle

TWS 85-0495/37

CASTROVIEJO BLADE HOLDER
STRAIGHT, TITANIUM




Overall length 95 mm, 3.74 inches

- TC conical jaws 10 mm long
- Round handle

TWS 85-0324/55

SCHOCKET SCLERAL DOUBLE
ENDED DEPRESSOR



Overall length 142 mm, 5.5 inches

- One side - sphere shaped tip 2.9 mm dia, stain essight round handle titanium
- Other side - 6.5 mm wide stainless

TWS 85-0551/20

BEAVER KNIFE HANDLE, SMALL,
TITANIUM



Overall length 51 mm, 2 inches

- Round handle

TWS 85-0676/30

BEAVER KNIFE HANDLE, MEDIUM,
TITANIUM



Overall length 76 mm, 3.0 inches

- Round handle

TWS 85-0102/40

BEAVER KNIFE HANDLE, LARGE,
TITANIUM




Overall length 102 mm, 4.0 inches

- Round handle


VARIOUS INSTRUMENTS

TWS 85-0515/45
LASEK BOW DISSECTOR, TITANIUM




- Working part 10 mm wide
- Overall length 15 mm, 4.5 inches
- Round handle

TWS 85-0531/51
BEEHLER PUPIL DILATOR




- Slightly curved tube 1.5 mm (15Ga) dia
- 3 working prongs 0.7 mm wide, stainless
- Overall length 130 mm, 5.1 inches
- Round handle, titanium

TWS 85-0516/45
CORNEAL DISSECTOR, ANGLED




- 11.5 mm long and 3 mm wide working part, stainless
- Overall length 16 mm, 4.5 inches
- Round handle, titanium

TWS 85-0518/46
CORNEAL DISSECTOR, CURVED



- 12 mm long and 3 mm wide working part, stainless
- Overall length 18 mm, 4.6 inches
- Round handle, titanium

TWS 85-2101/47
BLUMENTHAL CONJUNCTIVAL DISSECTOR




- Working tip olive type stainless
- 2 x 1.4 mm tip
- Overall length 121 mm, 4.7 inches
- Round handle, titanium

HOOKS



TW Surgical


TWS 85-0611/47
GRAETHER COLLAR BUTTON,
STRAIGHT



Overall length 120 mm, 4.7 inches

- Tip for polishing capsule and retracting iris
- Round handle, titanium
- stainless


TWS 85-0612/47
SINSKEY IRIS & IOL HOOK, ANGLED



Overall length 120 mm, 4.7 inches

- Angled 0.2 mm hook, stainless
- Round handle, titanium


TWS 85-0613/47
KUGLEN IRIS & LENS MANIPULATOR,
ANGLED



Overall length 120 mm, 4.7 inches

- Angled 0.65 mm hook, stainless
- Round handle, titanium


TWS 85-0614/47
BONN IRIS HOOK, ANGLED



Overall length 120 mm, 4.7 inches

- Angled 0.2 mm hook, stainless
- Round handle, titanium


TWS 85-0615/47
RENTSCH LENS MANIPULATOR,
STRAIGHT



Overall length 120 mm, 4.7 inches

- Multifunction tip to retract iris/ capsule, reposition IOL, stainless
- Round handle, titanium


TWS 85-0616/47
OSHER "Y" MANIPULATOR, ANGLED



Overall length 120 mm, 4.7 inches

- Angled "Y" shaped tip, stainless
- Round handle, titanium

TWS 85-0617/50
OSHER "Y" MANIPULATOR, ANGLED



Overall length 127 mm, 5.0 inches


- Vertical Y-hook 0.55mm wide, stainless
- Round handle, titanium

■ 14mm from tip to handle
 TW Medical Wright (U.K.) Ltd.
 128 City Road, London, EC1V 2NX, UK
 info@twrightmedical.com

HOOKS

TWS 85-0810/47

**BECHERT NUCLEUS MANIPULATOR
ANGLED**




Overall length 120 mm, 4,7 inches

- Angled "Y" shaped tip, stainless
- 12 mm from tip to bend
- Round handle, titanium

TWS 85-0811/47

SINSKEY IRIS & IOL HOOK STRAIGHT




Overall length 120 mm, 4,7 inches

- 0.2 mm hook, stainless
- Round handle, titanium

TWS 85-0812/47

**SINSKEY IRIS & IOL HOOK
STRAIGHT, TITANIUM**




Overall length 120 mm, 4,7 inches

- Straight 0.2 mm hook
- Round handle

TWS 85-0812/47

**SINSKEY IRIS & IOL HOOK ANGLED,
TITANIUM**




Overall length 120 mm, 4,7 inches

- Angled 0.2 mm hook
- Round handle

TWS 85-0813/47

4 IN 1 COLLAR BUTTON, ANGLED




Overall length 120 mm, 4,7 inches

- Tip to polishing capsule and retracting iris, stainless
- Round handle, titanium

TWS 85-0814/47

IRIS HOOK, ANGLED




Overall length 120 mm, 4,7 inches

- Angled 0,2 mm hook, stainless
- Round handle, titanium

TWS 85-0815/47

**MICROHOOK, ANGLED RIGHT
MS 002.01 LEFT**




Overall length 120 mm, 4,7 inches

- Curved to the right/left hook stainless
- Round handle, titanium

HOOKS

TWS 85-0410/47

MAUMENE IRIS HOOK SINGLE




■ Blunt tip 0.6 mm, stainless

Overall length 121 mm, 4.7 inches

■ Round handle, titanium

TWS 85-0411/47

IOL MANIPULATION HOOK.
ANGLED OE 005.01 STRAIGHT




■ Mushroom tip, stainless

Overall length 120 mm, 4.7 inches

■ Round handle, titanium

TWS 85-0412/47

FENZL INSERTION HOOK, STRAIGHT




■ Hook, stainless

Overall length 120 mm, 4.7 inches

■ Round handle, titanium

TWS 85-0413/47

FENZL INSERTION HOOK, ANGLED
MS 006.0412 mm from tip to bend




■ 10 mm from tip to bend, stainless

Overall length 118 mm, 4.6 inches

■ Round handle, titanium

TWS 85-0418/47

TWIST HOOK SMALL, RIGHT




■ Anticlockwise tips, stainless

Overall length 120 mm, 4.7 inches

■ Round handle, titanium

TWS 85-0419/47

TWIST HOOK SMALL, LEFT




■ Clockwise tips, stainless

Overall length 120 mm, 4.7 inches

■ Round handle, titanium

TWS 85-0417/47

MICROFINGER HOOK, ANGLED




■ 1.5 mm tip, stainless

Overall length 120 mm, 4.7 inches

■ Round handle, titanium

HOOKS


TWS 85-7410/45
NUCLEUS SUSTAINER, ANGLED



Overall length 14 mm, 4,5 inches

- Tip with ball 0.55mm dia, stainless
- Round handle, titanium
- Instrument for second hand, to be used with prechoppers (for hard cataracts)


TWS 85-5410/45
MANIPULATION HOOK, CURVED



Overall length 14 mm, 4,5 inches

- Tip with ball 0.55mm dia, stainless
- Round handle, titanium
- Instrument for second hand, to be used with prechoppers (for hard cataracts)


TWS 85-5411/47
NUCLEUS SUSTAINER, ANGLED



Overall length 120 mm, 4,7 inches

- 1.3mm tip, stainless
- Round handle, titanium
- Instrument for second hand, to be used with prechoppers (for hard cataracts)


TWS 85-5412/47
NUCLEUS ROTATMS, CURVED



Overall length 120 mm, 4,7 inches

- Curvedtip with notch 0.23mm, stainless
- Round handle, titanium


TWS 85-5421/44
TRIDENT SPLITTER, RIGHT
MS 022.02 LEFT



Overall length 12 mm, 4,4 inches

- Tip 0,85 × 0,6 mm stainless
- Round handle, titanium


TWS 85-5422/44
TRIDENT SPLITTER RIGHT, TITANIUM
MS 022.03 LEFT



Overall length 12 mm, 4,4 inches

- Tip 0.85 × 0.6 mm
- Round handle

TWS 85-5427/46
SINSKEY HOOK, TITANIUM




Overall length 17 mm, 4,6 inches

- Straight hook 0.2 mm
- Round handle

HOOKS


TWS 85-9418/46
ELEVATMS, TITANIUM



Overall length 18 mm, 4,6 inches

- Angled working part 8,4 mm long
- Tip 0.05 mm dia
- Round handle


TWS 85-1209
KUGLEN HOOK, STRAIGHT



Overall length 120 mm, 4,7 inches

- Tip 0,5 mm wide stainless
- Round handle, titanium


TWS 85-1218/46
KUGLEN HOOK, ANGLED



Overall length 18 mm, 4,6 inches

- Tip 0,5 mm wide stainless
- Round handle, titanium


TWS 85-1210/43
AKAHOSHI NUCLEUS
RINGSUSTAINER



Overall length 10 mm, 4,3 inches

- Ring tip 1.0 mm outer dia stainless
- Used to manipulate the hard nucleus
- Provides maximum protection of the posterior capsule
- Round handle, titanium


TWS 85-1137/54
HOOK STRAIGHT, TITANIUM



Overall length 137 mm, 5,4 inches

- 5,6mm working part
- Round handle


TWS 85-1214/45
LESTER LENS
MANIPULATOR, STRAIGHT



Overall length 14 mm, 4,5 inches

- Tip 0.2 × 0.5 mm stainless
- Round handle, titanium

TWS 85-2512/44
LESTER LENS
MANIPULATOR, ANGLED




Overall length 12 mm, 4,4 inches

- Tip 0.2 × 0.5 mm, stainless
- Round handle, titanium

HOOKS


TWS 85-9512/47
IRIS HOOK WITH A BALL TIP



Overall length 121 mm, 4.7 inches

- Tip with a ball 0.25 mm diameter, stainless
- Round handle, titanium


TWS 85-9125/49
IRIS HOOK WITH A BALL TIP



Overall length 125 mm, 4.9 inches

- Tip with a ball 0.25 mm diameter, stainless
- Round handle, titanium


TWS 85-2512/47
LAMELLAR DONOR INSERTION
INSTRUMENT, ANGLED



Overall length 120 mm, 4.7 inches

- 9.5 × 4.2 mm working tip, stainless
- Round handle, titanium
- For scoring the recipient bed and to position the donor lamella


TWS 85-9515
SINSKEY REVERSE HOOK
(DSAEK), TITANIUM



Overall length 15 mm, 4.5 inches

- Angled 0.2 mm hook
- Round handle
- For scoring the recipient bed and to position the donor lamella


TWS 85-9111/43
MANIPULATOR FOR
DLEK PROCEDURE



Overall length 111 mm, 4.3 inches

- Z-hook with blunt tip 0.15 mm diameter, stainless
- Round handle, titanium
- Designed to tuck the edges of the donor lamella


TWS 85-9113/44
MANIPULATOR FOR ICL



Overall length 113 mm, 4.4 inches

- Curved working part 12mm long, stainless
- Round handle, titanium
- 0.7 × 1.15 mm TC tip

TWS 85-9102/47
DOUBLE FIXATION HOOK, LARGE




Overall length 120 mm, 4.7 inches

- 2.0 mm spread between sharp tips, stainless
- Round handle, titanium

HOOKS

TWS 85-9612/47

DOUBLE FIXATION HOOK, SMALL




Overall length 120 mm, 4.7 inches

- 1.5 mm spread between sharp prongs, stainless
- Round handle titanium

TWS 85-1539/60

DOUBLE FIXATION HOOK STRAIGHT, SPACED




Overall length 153 mm, 6.0 inches

- 2.5 mm spread between sharp prongs, stainless
- Round handle titanium


TWS 85-1659/65

DOUBLE FIXATION HOOK STRAIGHT, (FREER SKIN HOOK)



TWS 85-1669/46

GASS RETINAL DETACHMENT HOOK




Overall length 16 mm, 4.6 inches

- 10 mm long tip with hole, stainless
- Round handle titanium

TWS 85-1496/56

GASS RETINAL DETACHMENT HOOK




Overall length 143 mm, 5.6 inches

- 13 mm long tip with hole, stainless
- Round handle titanium

TWS 85-1319/52

GASS RETINAL DETACHMENT HOOK




Overall length 131 mm, 5.2 inches

- 12 mm long with a hole, stainless
- Round handle titanium

TWS 85-1316/52

JAMESON MUSCLE HOOK




Overall length 131 mm, 5.2 inches

- Bulbous tip 3 mm, stainless
- Round handle titanium

HOOKS


TWS 85-1427/56
GREEN STRABISMUS HOOK



Overall length 142 mm, 5,6 inches

- Working part 10 mm long stainless
- Round handle titanium


TWS 85-7127/05
GREEN STRABISMUS HOOK



Overall length 127 mm, 5 inches

- Flat working part 9 mm long, stainless
- Round handle, titanium


TWS 85-7142/56
STRABISMUS HOOK



Overall length 142 mm, 5,6 inches

- Working part 10 mm long stainless
- Round handle titanium


TWS 85-7140/55
GRAEFE MUSCLE HOOK



Overall length 140 mm, 5,5 inches

- Working part 12.3 mm long stainless
- Round handle titanium


TWS 85-7141/55
WRIGHT GROOVED MUSCLE HOOK, RIGHT, TITANIUM



Overall length 140 mm, 5,5 inches

- Tip 14 mm long with groove
- Light round handle


TWS 85-7142/55
WRIGHT MUSCLE DOUBLE HOOK, TITANIUM



Overall length 141 mm, 5,5 inches

- Working part 11.5 mm dia
- Light round handle

TWS 85-7144/55
WRIGHT GROOVED MUSCLE HOOK, LEFT, TITANIUM




Overall length 140 mm, 5,5 inches

- Tip 14 mm long with groove
- Light round handle


HOOKS

TWS 85-1419/55
STEVENS CURVED TENOTOMY HOOK



- Tip 6.2 mm long stainless
- Overall length 141 mm, 5.5 inches
- Round handle titanium

TWS 85-1188/46
STEVENS CURVED TENOTOMY HOOK




- Working part 5 mm long, stainless
- Overall length 118 mm, 4.6 inches
- Tip 1.0 mm dia
- Round handle, titanium


HOOKS




TW Surgical




- Overall length 150 mm, 5.9 inches**
- Angled Y-shaped tip, stainless
 - Round handle, titanium
 - Straight tip with button for polishing capsule and retracting iris, stainless




- Overall length 16 mm, 4.5 inches**
- Angled working parts, stainless
 - Curved tip 0.3 × 0.1 mm
 - Angled tip, 0.3 × 0.1 mm
 - Round handle, titanium



- Overall length 152 mm, 6.0 inches**
- Graefe strabismus hook 2.8mm, stainless
 - Stevens tenotomy hook 6.2mm, stainless
 - Round handle, titanium



- Overall length 152 mm, 6.0 inches**
- Graefe strabismus hook 4mm, stainless
 - Strabismus hook 4.5mm, stainless
 - Round handle, titanium




- Overall length 153 mm, 6.0 inches**
- Jameson muscle hook with bulbous tip 3.3mm, stainless
 - Green strabismus hook 0.2mm, stainless
 - Round handle, titanium

CHOPPERS



TW Surgical

TWS 85-9720/47
SATO KNIFE




- 1 mm blade stainless

Overall length 120 mm, 4,7 inches

- Round handle, titanium



TWS 85-9712/47
LENS CRUSHING KNIFE
OK 002.01 Tip 0.3mm wide




- Blade 1.5 mm stainless
- Tip 0,9 mm wide

Overall length 120 mm, 4,7 inches

- Round handle, titanium



TWS 85-9721/47
CYSTOTOME ANGLED



- 0.2mm blade, stainless

Overall length 120 mm, 4,7 inches

- Round handle, titanium



TWS 85-2079/47
CYSTOTOME CURVED
OK 003.02 11.7mm from tips to bend




- 0.2mm blade, stainless
- 11,2 mm from tip to bend

Overall length 120 mm, 4,7 inches

- Round handle, titanium



TWS 85-9710/47
DRYSDALE LENS CRUSHING KNIFE,
ANGLED




- 1.8mm blade, stainless

Overall length 120 mm, 4,7 inches

- Round handle, titanium



TWS 85-1290/47
DRYSDALE LENS CRUSHING KNIFE,
CURVED




- 1,8mm blade stainless

Overall length 120 mm, 4,7 inches

- Round handle, titanium



TWS 85-1720/47
DRYSDALE LENS CRUSHING KNIFE,
CURVED



- 1.8mm blade stainless
- TC tip

Overall length 120 mm, 4,7 inches

- Round handle, titanium




CHOPPERS



TW Surgical

TWS 85-1020/47
PHACO CHOPPER ANGLED LARGE,
TITANIUM




■ 1.5 mm blade

Overall length 120 mm, 4,7 inches

■ Round handle



TWS85-1021/47
PHACO CHOPPER, ANGLED MEDIUM
TITANIUM




■ 1.3 mm blade

Overall length 120 mm, 4,7 inches

■ Round handle



TWS 85-1022/47
PHACO CHOPPER, ANGLED SMALL
TITANIUM




■ 1,0 mm blade

Overall length 120 mm, 4,7 inches

■ Round handle



TWS 85-1016/45
GOLF CLUB FOREIGN BODY SPUD




■ 1.4mm blade, stainless

Overall length 16 mm, 4,5 inches

■ Round handle, titanium



TWS 85-1023/47
FOREIGN BODY SPUD




■ 2.2mm gently curved blade, stainless

Overall length 120 mm, 4,7 inches

■ Round handle, titanium




TWS 85-1024/47
DAVIS FOREIGN BODY SPUD




■ 4.2mm blade, stainless

Overall length 120 mm, 4,7 inches

■ Round handle, titanium



TWS 85-1025/43
PHACO CHOPPER ANGLED



■ 0.9 mm blade, stainless


Overall length 110 mm, 4,3 inches

■ Round handle, titanium




CHOPPES

TWS 85-1077/43
PHACO CHOPPER ANGLED




- 0.9 mm blade, stainless
 - Round handle, titanium
- Overall length 10 mm, 4.3 inches**

TWS 85-1078/43
PHACO CHOPPER ANGLED





- 0.9 mm blade, stainless
 - Round handle, titanium
- Overall length 10 mm, 4.3 inches**

TWS 85-1079/43
DANIEL - WALKOW PHACO CHOPPER
ANGLED





- 1.3 mm blade with ball 0.4mm dia, stainless
 - Round handle, titanium
- Overall length 10mm, 4.3 inches**

TWS 85-1279/47
HOOK AND PHACO CHOPPER





- 0.43 mm blade, stainless
 - Round handle, titanium
- Overall length 120 mm, 4.7 inches**

TWS 85-1266/47
ROSEN PHACO CHOPPER
UNIVERSAL



- 1.0 mm blade stainless
 - Round handle, titanium
- Overall length 120 mm, 4.7 inches**

TWS 85-1159/45
NICHAMIN CHOPPER, TITANIUM



- 1.1 mm blade with ball 0.5 mm dia
 - Light round handle
- Overall length 115 mm, 4.5 inches**


TWS 85-1159/45
SEIBEL CHOPPER



- 1.7mm blade with ball 0.5 mm dia, stainless
 - Round handle, titanium
- Overall length 110 mm, 4.3 inches**

CHOPPERS


TWS 85-3310/43
SEIBEL CHOPPER



Overall length 10 mm, 4.3 inches

- 1.5mm blade with ball 0.45 mm dia, stainless
- Round handle, titanium


TWS 85-3312/44
PHACO SPLITTER



Overall length 12 mm, 4.4 inches

- 0.85 × 0.7 mm tip, stainless
- Round handle, titanium


TWS 85-1269/47
DANIEL-WALKOW PHACO
CHOPPER ANGLED



Overall length 120 mm, 4.7 inches

- 2.5 × 6.0 mm blade stainless
- Round handle, titanium


TWS 85-1267/47
HOOK AND PHACO CHOPPER



Overall length 120 mm, 4.7 inches

- 2.5 × 8.5 mm semi-sharp blade, stainless
- Round handle, titanium


TWS 85-1265/47
KNIFE HOCKEY TYPE



Overall length 120 mm, 4.7 inches

- 1.4 × 6.0 mm semi-sharp blade, stainless
- Round handle, titanium


TWS 85-1185/43
CHANG CHOPPER RIGHT OK
034.03 LEFT



Overall length 110 mm, 4.3 inches

- Curved microfinger tip, stainless
- Round handle, titanium

TWS 85-1036/45
QUICK CHOPPER




Overall length 16 mm, 4.5 inches

- 1.3mm blade, stainless
- 1.75mm tip
- Round handle, titanium

CHOPPERS

TWS 85-4007/42
QUICK CHOPPER
 OK035.03 1.2mm blade, 1.5mm tip




- 0.8mm blade, stainless
- 1.3mm tip

Overall length 107 mm, 4.2 inches

- Round handle, titanium



TWS 85-4006/41
DRYSDALE LENS CRUSHING KNIFE,
ANGLED




- 1.3mm blade, stainless

Overall length 106 mm, 4.1 inches

- Round handle, titanium



TWS 85-1411/43
QUICK CHOPPER, RIGHT OK
044.01 LEFT




- 1.3 mm blade angled aside, stainless

Overall length 11 mm, 4.3 inches

- Round handle, titanium



TWS 85-1511/43
QUICK CHOPPER (IN LINE WITH
PHACO INCISION), RIGHT
OK 045.01 LEFT




- 1.2mm blade stainless

Overall length 11 mm, 4.3 inches

- Round handle, titanium



TWS 85-0417/42
KANSAS- ALFONSO NUCLEUS
BISECTOR, ANGLED




- 0.8 x 0.5 mm tip stainless

Overall length 107 mm, 4.2 inches

- Round handle, titanium



TWS 85-0515/41
KANSAS-ALFONSO NUCLEUS VECTIS,
RIGHT
OK 043.01 LEFT




- 2 mm wide and 10 mm long tip stainless

Overall length 105 mm, 4.1 inches

- Round handle, titanium



TWS 85-0512/44
PHACONIT QUICK CHOPPER,
ANGLED



- 0.35mm blade, stainless


Overall length 12 mm, 4.4 inches

- Round handle, titanium



CHOPPERS


TWS 85-5106/41
KANSAS-ALFONSO NUCLEUS TRISECTOR,
RIGHT OK 046.01 LEFT



Overall length 106 mm, 4.1 inches

- 1.5 mm wide and 10 mm long tip stainless
- Round handle, titanium


TWS 85-5021/47
LASEK TREPHINE



Overall length 121 mm, 4.7 inches

- To create epithelial flap
- For 8 mm optical zone
- Depth of incision 70 mkm, stainless
- Round handle, titanium


TWS 85-1278/47
LASEK TREPHINE



Overall length 121 mm, 4.7 inches

- To create epithelial flap
- For 9 mm optical zone
- Depth of incision 70 mkm, stainless
- Round handle, titanium


TWS 85-5014/45
LASEK EPITHELIAL MICRO HOE



Overall length 14 mm, 4.5 inches

- Can be used to detach and lift edge of the epithelium lap, stainless
- Round handle, titanium


TWS 85-5015/45
CHOPPER, CURVED



Overall length 15 mm, 4.5 inches

- Tip 2.3 mm long, stainless
- Round handle, titanium


TWS 85-5019/46
NAGAHARA PHACO CHOPPER ANGLED
RIGHT. TITANIUM
OK 052.01 LEFT



Overall length 19 mm, 4.6 inches

- 1.25 mm blade
- Round handle

TWS 85-5011/43
NAGAHARA CHOPPER ANGLED



Overall length 11 mm, 4.3 inches


- 1.0mm sharp blade stainless
- Round handle, titanium

CHOPPERS



TW Surgical


TWS 85-6111/43
NAGAHARA CHOPPER ANGLED



Overall length 11 mm, 4.3 inches

- 1.0 mm sharp blade, stainless
- Round handle, titanium

TWS 85-6112/44
CHOPPER, RIGHT



Overall length 112mm, 4.4 inches

- 1.1mm cutting blade with ball 0.3mm dia, stainless
- Round handle, titanium


TWS 85-6211/44
CHOPPER, LEFT



Overall length 112mm, 4.4 inches

- 1.1mm cutting blade with ball 0.3mm dia, stainless
- Round handle, titanium


TWS 85-1010/43
CONNOR WAND ANGLED



Overall length 10 mm, 4.3 inches

- Tip with ball 0.5 mm dia stainless
- Round handle, titanium

TWS 85-1260/44
CONNOR WAND STRAIGHT




Overall length 12 mm, 4.4 inches

- Tip with ball 0.5 mm dia stainless
- Round handle, titanium

CHOPPERS

TWS 85-7134/53

SEIBEL NUCLEUS CHOPPER/TRIDENT SPLITTER




Overall length 134 mm, 5.3 inches

- 1.7mm blade with ball 0.5 mm dia, stainless
- Round handle, titanium
- 0.85 × 0.6mm splitter, stainless

TWS 85-7129/51

CHANG CHOPPER.
RIGHT OK034.01 LEFT




Overall length 129 mm, 5.1 inches

- Curved micronger tip, stainless
- Round handle, titanium
- Quick chopper tip, stainless

TWS 85-7133/52

SEIBEL VERTICAL&HORIZONTAL CHOPPER
DOUBLE-ENDED
OK038.01LEFT




Overall length 133 mm, 5.2 inches

- Seibel horizontal and vertical, stainless
- Round handle, titanium

TWS 85-7143/56

BOONE CHOPPER DOUBLE-ENDED,
ANGLED




Overall length 143 mm, 5.6 inches

- Vertical and Drisdale, stainless
- Round handle, titanium

TWS 85-1337/52

CHOPPER DOUBLE-ENDED




Overall length 133 mm, 5.2 inches

- Two curved choppers, stainless
- Round handle, titanium

TWS 85-1347/53

DOUBLE-ENDED CHOPPER, RIGH (CHANG
& QUICK)




Overall length 134 mm, 5.3 inches

- Working parts, stainless
- Round handle, titanium

TWS 85-1349/53

DOUBLE-ENDED CHOPPER, LEF
(CHANG & QUICK)



Overall length 134 mm, 5.3 inches

- Working parts, stainless
- Round handle, titanium


CHOPPERS



TW Surgical

TWS 85-8132/52


ROSEN PHACO CHOPPER+ BECHERT NUCLEUS ROTATOR



- Overall length 132mm, 5.2 inches**
- Angled Y-shaped tip, stainless
 - Chopper 1.0 mm long, stainless
 - Round handle, titanium

TWS 85-1328/52


ROSEN PHACO CHOPPER/CONNOR WAND



- Overall length 132 mm, 5.2 inches**
- Tip with ball 0.5mm dia, stainless
 - 1.0mm chopper blade, stainless
 - Round handle, titanium

TWS 85-8142/56


NICHAMIN CHOPPER/ SINSKEY HOOK



- Overall length 142 mm, 5.6 inches**
- 1 mm chopper blade, stainless
 - 0.2mm dia hook, stainless
 - Round handle, titanium

TWS 85-8151/59


CHOPPER/ SPATULA



- Overall length 151 mm, 5.9 inches**
- 1.2mm chopper blade, stainless
 - Spatula tip 0.25 x 0.2 mm, stainless
 - Round handle, titanium

TWS 85-8159/62


CHOPPER/ SPATULA



- Overall length 159 mm, 6.2 inches**
- 1.2mm chopper blade, stainless
 - Spatula tip 0.4 x 0.4 mm stainless
 - Round handle, titanium

TWS 85-1428/56


CHOPPER/ SPATULA



- Overall length 142 mm, 5.6 inches**
- Spatula tip 0.25 x 0.2 mm, stainless
 - 1.7mm blade with ball 0.5mm dia, stainless
 - Round handle, titanium

TWS 85-1498/58

SEIBEL CHOPPER/ SPATULA




- Overall length 149 mm, 5.8 inches**
- Spatula tip 0.4 x 0.4 mm stainless
 - 1.7mm blade with ball 0.5mm dia, stainless
 - Round handle, titanium

CHOPPERS

TWS 85-9157/62

HOCKEY EPITHELIUM REMOVAL KNIFE FOR PRK/
SPATULA





Overall length 157 mm, 6.2 inches

- Spatula tip 0.5 x 0.5 mm, stainless
- 6.0 x 2.5 mm semi sharp blade, stainless
- Round handle, titanium

TWS 85-9138/54

SPLITTER / SPATULA




Overall length 138 mm, 5.4 inches

- 0.6 mm dia flat disk, stainless
- Chopper 1.0 mm blade, stainless
- Round handle, titanium

TWS 85-9108/43

SEIBEL VERTICAL CHOPPER & DONNEFELD
FEMTO SPATULA, RIGHT




Overall length 108 mm, 4.3 inches

- Chopper 0.2 x 0.5 mm, stainless
- Spatula 0.15 x 0.3 mm, stainless
- Round handle, titanium

TWS 85-1089/43

SEIBEL VERTICAL CHOPPER & DONNEFELD
FEMTO SPATULA, RIGHT




Overall length 108 mm, 4.3 inches

- Chopper 0.2 x 0.5 mm, stainless
- Spatula 0.15 x 0.3 mm, stainless
- Round handle, titanium

TWS 85-9012/43

SEIBEL CHOPPER & SPATULA




Chopper 0.25 x 0.5 mm, stainless
Spatula 0.3 mm dia, stainless

- Round handle, titanium

TWS 85-9133/52

STEINERT CLAW CHOPPER & KOCH SPATULA




Overall length 133 mm, 5.2 inches

- Curved chopper 1.6 mm, stainless
- Curved spatula 0.5 mm, stainless
- Round handle, titanium

TWS 85-9154/61

LASEK TREPHINE / FUNNEL 8.0&8.5 MM



Overall length 154 mm, 6.1 inches



- To create epithelial flap
- For 8 mm optical zone
- Depth of incision 70 micrometers
- To apply alcohol solution on epithelium
- Working parts, stainless
- Round handle, titanium

TREPHINES



TW Surgical

TWS 85-1544/61
LASEK TREPHINE / FUNNEL
9.0&9.5 MM



- To create epithelial flap
- For 9 mm optical zone
- Depth of incision 70 mkm

TWS 85-1555/61
EYE TREPHINES CONIC



Overall length 154 mm, 6.1 inches


- To apply alcohol solution on epithelium
- Working parts, stainless
- Round handle, titanium

SPATULAS



TW Surgical


TWS 85-5220/47
EYE SPATULA



Overall length 120 mm, 4.7 inches

- Angled blade 5mm long, stainless
- 0.5 × 0.3 mm wide
- Round handle, titanium


TWS 85-5221/47
EYE SPATULA



Overall length 120 mm, 4.7 inches

- Angled 15 mm long stainless
- 1.0 × 0.3 mm wide
- Round handle, titanium


TWS 85-5222/47
EYE SPATULA



Overall length 120 mm, 4.7 inches

- Angled blade 15 mm long stainless
- 1.5 × 0.3 mm wide
- Round handle, titanium



TWS 85-5223/47
BARRAQUER CYCLODIALYSIS EYE SPATULA



Overall length 120 mm, 4.7 inches

- Curved blade 15 mm long stainless
- 0.5 mm dia
- Round handle, titanium


TWS 85-5224/47
DEMENTIEV MICRO SPATULA, STAINLESS



Overall length 121 mm, 4.7 inches

- Angled blade 0 mm long
- 1.8 × 0.4 mm wide
- Round handle


TWS 85-5216/45
EYE SPATULA



Overall length 16mm, 4.5 inches

- Angled blade 10 mm long stainless
- 0.5 × 0.3 mm wide
- Round handle, titanium

TWS 85-1251/47
IOL SPATULA NOTCHED




Overall length 120 mm, 4.7 inches

- Angled blade 13 mm long stainless
- 0.8 × 0.3 mm wide
- Round handle, titanium


SPLITTERS

TWS 85-2150/59
KIMURA PLATINUM SPATULA




- 4 mm wide blade stainless
- Overall length 150 mm, 5.9 inches
- Round handle, titanium

TWS 85-5712/47
CASTROVIEJO CYCLODIALYSIS SPATULA ANGLED




- 14 mm long blade stainless
- Overall length 120 mm, 4.7 inches
- Tip 0.25 x 0.2 mm wide
- Round handle titanium

TWS 85-1247/47
CASTROVIEJO CYCLODIALYSIS SPATULA ANGLED




- Blade 14 mm long stainless
- Overall length 120 mm, 4.7 inches
- Tip 0.5 x 0.2 mm wide
- Round handle titanium

TWS 85-1268/47
CASTROVIEJO CYCLODIALYSIS SPATULA ANGLED




- Blade 14 mm long stainless
- Overall length 120 mm, 4.7 inches
- Tip 0.75 x 0.2 mm wide
- Round handle titanium

TWS 85-1262/47
CASTROVIEJO CYCLODIALYSIS SPATULA ANGLED




- Blade 14 mm long stainless
- Overall length 120 mm, 4.7 inches
- Tip 1.0 x 0.2 mm wide
- Round handle titanium

TWS 85-1266/47
CASTROVIEJO CYCLODIALYSIS SPATULA ANGLED TITANIUM



- Blade 14 mm long
- Overall length 120 mm, 4.7 inches
- Tip 0.75 x 0.2 mm wide
- Round handle


TWS 85-1263/47
CASTROVIEJO CYCLODIALYSIS SPATULA ANGLED TITANIUM



- Blade 6 mm long
- Overall length 120 mm, 4.7 inches
- Tip 0.3 x 0.25 mm wide
- Round handle

SPATULAS


TWS 85-1150/43
KOCH NUCLEUS SPATULA
 OT 009.01 Tip 0.6 mm wide



Overall length 11 mm, 4.3 inches

- Tip 0.5 mm wide stainless
- Round handle, titanium


TWS 85-1550/45
EYE SPATULA



Overall length 15 mm, 4.5 inches

- Angled and curved 12 mm long blade, stainless
- Round handle, titanium
- Tip 0.25 mm dia


TWS 85-3141/55



Overall length 141 mm, 5.5 inches

- Curved blade 5 mm wide stainless
- 2 teeth 0.5 mm height
- Round handle, titanium


TWS 85-1194/43
PHACO SPATULA
 TITANIUM OT 017.01 STAINLESS



Overall length 110 mm, 4.3 inches

- Blade 11.5 mm long
- Round handle
- Tip 0.8 x 0.2 mm wide


TWS 85-1125/45
BAHAMOND LASIK SPATULA



Overall length 16 mm, 4.5 inches

- Angled blade 8.5 mm long, stainless
- Round handle, titanium
- Tip 2.6 x 0.8 mm wide


TWS 85-1626/45
BAHAMOND LASIK SPATULA



Overall length 16 mm, 4.5 inches

- Angled and curved blade 8.5 mm long stainless
- Round handle, titanium
- Tip 1.9 mm wide

TWS 85-1144/44
LASEK EPITHELIUM
DETACHMENT SPATULA




Overall length 11 mm, 4.4 inches

- Can be used to detach edge of the epithelium flap
- Round handle, titanium

SPATULAS


TWS 85-1288/50
KANSAS NUCLEUS SUSTAINER,
RIGHT



Overall length 128 mm, 5.0 inches

- Blade 18.5 mm long stainless
- Tip 1.2 × 0.4 mm wide
- Round handle, titanium


TWS 85-1258/50
CASTROVIEJO CYCLODIALYSIS EYE
SPATULA



Overall length 128 mm, 5.0 inches

- Angled and curved blade 15 mm long stainless
- Tip 1.9 × 0.25 mm wide
- Round handle, titanium


TWS 85-1064/42
LEHMANN INTRALASEK SPATULA



Overall length 106 mm, 4.2 inches

- Angled and curved 10 mm long blade stainless
- Tip 0.05 × 0.4 mm wide
- Round handle, titanium

TWS 85-1714/46
FEMTOLASIK FLAP SPATULA



Overall length 17 mm, 4.6 inches

- Angled and curved blade 10 mm long, stainless
- Tip 0.8 × 1.0 mm wide
- Round handle, titanium


TWS 85-1713/44
ZALDIVAR FEMTOLASIK SPATULA



Overall length 13 mm, 4.4 inches

- Curved blade 12 mm long stainless
- Tip 0.5 mm wide
- Round handle, titanium


TWS 85-1224/48
LASEK EPITHELIAL FLAP
REPOSITIONING SPATULA



Overall length 122 mm, 4.8 inches

- Blade 1.9 × 0.3 mm wide stainless
- Round handle, titanium

TWS 85-1287/47
MELLES DSAEK PLK SCRAPER
TITANIUM




Overall length 122 mm, 4.7 inches

- Angled 90 degree tip 0.2 × 1.0 mm wide
- Round handle

SPATULAS


TWS 85-5821/47
MELLES DSAEK PLK
SCRAPER TITANIUM



Overall length 122mm, 4.7 inches

- Angled 45 degree tip 0.2 × 1.0 mm wide
- Round handle


TWS 85-5112/47
CASTROVIEJO CYCLODIALYSIS
SPATULA



Overall length 120 mm, 4.7 inches

- Blade 14 mm long stainless
- Tip 1 mm with oval hole 2.0 × 0.5 mm
- Round handle, titanium


TWS 85-1851/46
EYE SPATULA



Overall length 18 mm, 4.6 inches

- Blade 13 mm long stainless
- Tip 0.3 mm wide
- Round handle, titanium


TWS 85-1292/47
EYE SPATULA



Overall length 120 mm, 4.7 inches

- Blade 11 mm long stainless
- Tip 1.7 × 0.2 mm wide
- Round handle, titanium


TWS 85-1293/47
CAPSULE POLISHER ANGLED



Overall length 120 mm, 4.7 inches

- Blade 11 mm long, stainless
- TC sphere shaped tip 1.5 mm dia
- Round handle, titanium


TWS 85-1895/46
PHACOEMULSIFICATION
SPATULA



Overall length 18 mm, 4.6 inches

- Blade 11 mm long stainless
- Round tip 1 mm dia and 0.2 mm thick
- Used to rotate and crack the nucleus
- Round handle, titanium

TWS 85-1128/50
PATON SPOON




Overall length 128 mm, 5.0 inches

- Curved spoon 7.4 × 0.45 mm with 2° angle stainless
- Flat handle, titanium

SPATULAS


TWS 85-1827/46
CAPSULE POLISHER



Overall length 18 mm, 4.6 inches

- Tip 1.6 external dia × 1.3 internal dia, stainless
- Round handle, titanium


TWS 85-1827/46
CAPSULE POLISHER



Overall length 18 mm, 4.6 inches

- Tip 1.6 external dia × 1.3 internal dia, stainless
- Round handle, titanium


TWS 85-1926/46
SPATULA



Overall length 19 mm, 4.6 inches

- Angled and curved blade 12 mm long, stainless
- Tip 3 mm wide
- Round handle, titanium


TWS 85-1097/43
SLADE FEMTO SPATULA



Overall length 10 mm, 4.3 inches

- Blade 95 mm long, stainless
- Flattened tip 1.4 mm long, 0.45 mm wide
- Round handle, titanium
- Designed to open femtosecond created cataract incision and aids in identifying and entering the interface and opening the wound with a sweeping motion

TWS 85-1145/43
DANNENFELD FEMTO SPATULA



Overall length 11 mm, 4.3 inches

- Blade 9.8 mm long stainless
- Flattened tip 1.25 mm long, 0.3 mm wide
- Round handle, titanium
- Designed to open femtosecond created cataract incision and aids in identifying and entering the interface and opening the wound with a sweeping motion

TWS 85-1877/46
BUSIN SPATULA-GUIDE FOR
CORNEAL ENDOTHELIUM
IMPLANTATION, TITANIUM




Overall length 18mm, 4.6 inches

- Suitable for incisions as small as 3.2 mm
- Round handle


SPATULAS

TWS 85-1673/64
DOUBLE ENDED EYE SPATULA




- Overall length 163 mm, 6.4 inches**
- Ang ed blades 12 mm ong stainless
 - Tips 0.5 × 0.3 mm and 1.0 × 0.3 mm wide
 - Round handle, titan'um

TWS 85-1543/06
CASTROVIEJO CYCLODIALYSIS SPATULA ANGLED




- Overall length 153 mm, 6 inches**
- 0.25 × 0.2 mm and 0.5 × 0.2 mm wide
 - 14 mm long blades stainless
 - Round handle, titan'um

TWS 85-1573/06
CASTROVIEJO CYCLODIALYSIS SPATULA ANGLED



- Overall length 153 mm, 6 inches**
- Tips 0.5 × 0.2 mm and 0.75 × 0.2 mm wide
 - 14 mm long blades, stainless
 - Round handle, titan'um

TWS 85-1563/06
CASTROVIEJO CYCLODIALYSIS SPATULA ANGLED




- Overall length 153 mm, 6 inches**
- Tips 0.75 × 0.2 mm and 1.0 × 0.2 mm wide
 - 14 mm long blades stainless
 - Round handle, titan'um

TWS 85-1675/65
CASTROVIEJO DOUBLE-ENDED CYCLODIALYSIS SPATULA OT 007.01 TITANIUM




- Overall length 165 mm, 6.5 inches**
- Blade 15 mm ong and 0.4 × 0.4 mm wide, stainless
 - Blade 11 mm ong and 0.5 × 0.5 mm wide, stainless
 - Round handle, titan'um

TWS 85-1466/55
LINDSTROM LASIK/ PRK SPATULA



- Overall length 140 mm, 5.5 inches**
- Semi - sharp tip is de 'gned to remove epithel'um, stainless
 - Spatula tip 0.5 × 0.5 mm wide, stainless
 - Round handle, titan'um


TWS 85-1472/56
KOCH NUCLES SPATULA/ SINSKEY HOOK



- Overall length 142 mm, 5.6 inches**
- Koch spatula blade 0.5 mm wide, stainless
 - Hook 0.2 mm dia, stainless
 - Round handle, titan'um


SPATULAS

TWS 85-1925/75
ELEVATRS PERIOSTEAL FREER'S




- Overall length 192 mm, 7.5 inches
- Blades 4.0 × 0.4 mm and 4.5 × 0.4 mm, stainless steel
- Round handle, titanium

TWS 85-1433/56
EYE SPATULA/ HOOK TITANIUM




- Overall length 143 mm, 5.6 inches
- Spatula tip 9mm long and 0.6 × 0.25 mm wide, stainless steel
- Hook tip tip 0.2 mm dia, stainless steel
- Round handle, titanium

TWS 85-1216/46
SPATULA/ NUCLEUS ROTATOR




- Overall length 126 mm, 4.6 inches
- Spatula tip 0.75 × 0.2 mm wide, stainless steel
- Nucleus rotator 0.8 mm wide, stainless steel
- Round handle, titanium

TWS 85-1597/60
PATON DOUBLE-ENDED SPATULA AND SPOON




- Overall length 159 mm, 6.0 inches
- Angled spatula tip 6.0 mm long and 0.3 mm wide with 3 lobes, stainless steel
- Spoon 6.0 × 0.45 mm wide, stainless steel
- Round handle, titanium

TWS 85-8145/57
BAHAMOND FLAP LASIK SPATULA




- Overall length 144 mm, 5.7 inches
- Angled blade 8.5 mm long, 1.9 mm wide, stainless steel
- Curved blade 8.5 mm long, 1.9 mm wide, stainless steel
- Round handle, titanium

TWS 85-1336/53
DOUBLE-ENDED LASIK SPATULA



- Overall length 136 mm, 5.3 inches
- Angled blade 11.3 mm long and 1.5 × 0.8 mm wide with a groove 0.5 mm, stainless steel
- Curved blade 0.95 mm wide, stainless steel
- Round handle, titanium

TWS 85-1382/52
FEMTOSECOND THORLAKSON LASIK FLAP LIFTER / SEPARATOR



- Overall length 132 mm, 5.2 inches
- Hook 0.25 mm dia, stainless steel
- Round tip 0.8 mm dia, stainless steel
- Round handle, titanium

Treasure Wright (U.K.) Ltd.

128 City Road, London, EC1V 2NX, UK

info@twrightmedical.com

SPATULAS

TWS 85-5240/55

THORLAKSON LASIK FLAP SEPARATOR/
LIFTER



Overall length 140mm, 5.5 inches

- Spatula and Sinskey Hook, stainless
- Round handle, titanium

TWS 85-5426/49

THORLAKSON LASIK FLAP SEPARATOR/
LIFTER



Overall length 126mm, 4.9 inches

- Spatula 2.2mm with beveled notches, stainless
- Round handle, titanium
- Sinskey hook, stainless

TWS 85-5436/46

DOUBLE-ENDED LASIK SPATULA AND FLAP
RETREATMENT INSTRUMENT



Overall length 136 mm, 4.6 inches

- Angled blade 10 mm long and 0.5 x 0.6mm wide, stainless
- Angled tip 1.2 mm long and 0.8 x 0.2mm wide, stainless
- Round handle, titanium

TWS 85-5431/51

FEMTOLASIK FLAP SPATULA,
TITANIUM



Overall length 131 mm, 5.1 inches

- Curved spatula 12.5 mm long; tip 0.15 mm dia
- Instrument to open and manipulate flap
- Delicate pick 1 mm long
- Round handle

TWS 85-5438/53

FEMTOLASIK FLAP SPATULA



Overall length 135mm, 5.3 inches

- Curved spatula 12.3 mm long, stainless
- Instrument to open and manipulate flap
- Delicate pick 1.5 mm long, stainless
- Round handle, titanium

TWS 85-7740/55

SPATULA



Overall length 140 mm, 5.5 inches

- Angled and curved blade 9.7 mm long, 0.3 mm dia, stainless
- Angled blade 11.2 mm long, 0.3 mm dia, stainless
- Round handle, titanium

TWS 85-7732/52

BOONE DOUBLE-ENDED
FEMTOSPATULA




Overall length 132 mm, 5.2 inches

- Spatula blade 9 mm long, stainless
- Round handle, titanium
- Delicate pick 1 mm long, titanium

SPATULAS

TWS 85-2023/48

DOUBLE-ENDED CURETTE CAPSULE POLISHER




Overall length 123 mm, 4.8 inches

- Sharp working tips 1.6 mm dia, stainless
- Round handle, titanium

TWS 85-2015/59

PHACO SPATULA/MICRO FINGER




Overall length 151 mm, 5.9 inches

- Micro finger, stainless
- Round handle, titanium
- Spatula tip 14 mm long, 0.25 × 0.2 mm wide, stainless

TWS 85-2036/53

CONNOR WAND/SPATULA




Overall length 136 mm, 5.3 inches

- Connor wand, stainless
- Round handle, titanium
- Curved IntraLASIK spatula 1 mm wide, stainless

TWS 85-2016/46

EMTO DOUBLE-ENDED SPATULA / HOOK (CHANSUE)




Overall length 116 mm, 4.6 inches

- Lenticule spatula – blunt spoon, stainless
- Round handle, titanium
- Angled blunt hook, stainless

TWS 85-2031

SEKUNDO PIK SIMLE DISSECT FOR 2.5MM INSICION




Overall length 123mm, 4.8 inches

- Pak tip 1.8mm wide, stainless
- Round handle, titanium
- Sinsky hook 1mm, stainless


FIXATION RINGS

TWS 85-2135/53
THORNTON REVOLVING GLOBE
FIXATION RING, TITANIUM




- Overall length 135 mm, 5,3 inches
- 12 teeth to stabilize and rotate eye globe
 - 16 mm internal dia
 - Round handle

TWS 85-2113/53
THORNTON REVOLVING 3/4 GLOBE
FIXATION RING, TITANIUM




- Overall length 135 mm, 5,3 inches
- 9 teeth to stabilize and rotate eye globe
 - 16 mm internal dia
 - Can be used for AK, RK and clear corneal incision
 - Round handle

TWS 85-1235/53
THORNTON REVOLVING GLOBE
FIXATION RING, TITANIUM




- Overall length 135 mm, 5,3 inches
- 12 teeth to stabilize and rotate eye globe
 - 16 mm internal dia
 - Round handle

TWS 85-1135/53
THORNTON REVOLVING 3/4 GLOBE
FIXATION RING, TITANIUM




- Overall length 135 mm, 5,3 inches
- 9 teeth to stabilize and rotate eye globe
 - 16 mm internal dia
 - Can be used for AK, RK and clear corneal incision
 - Round handle

TWS 85-2129/51
THORNTON GLOBE FIXATION RING,
TITANIUM




- Overall length 129 mm, 5,1 inches
- 12 teeth to stabilize and rotate eye globe
 - 16 mm internal dia
 - Round handle

TWS 85-2128/05
THORNTON 3/4 GLOBE FIXATION
RING, TITANIUM



- Overall length 128 mm, 5 inches
- 9 teeth to stabilize and rotate eye globe
 - 16 mm internal dia
 - Can be used for AK, RK, and clear corneal incision
 - Round handle


TWS 85-2112/48
THORNTON 1/2 GLOBE
FIXATION RING, TITANIUM



- Overall length 123 mm, 4,8 inches
- 7 teeth to stabilize and rotate eye globe
 - 16 mm internal dia
 - Can be used for AK, RK, and clear corneal incision
 - Round handle

FIXATION RINGS


TWS 85-2235/53
THORNTON REVOLVING 3/4 GLOBE
FIXATION RING, TITANIUM



Overall length 135 mm, 5.3 inches

- 9 teeth to stabilize and rotate eye globe
- 14 mm internal dia
- Can be used for PK, RK and clear corneal incision
- Round handle


TWS 85-2229/51
REVOLVING GLOBE FIXATION RING,
TITANIUM



Overall length 129 mm, 5.1 inches

- 9 teeth to stabilize and rotate eye globe
- 13 mm internal dia
- Can be used for PK, RK and clear corneal incision
- Round handle


TWS 85-1722/46
3/4 GLOBE FIXATION RING,
TITANIUM



Overall length 17 mm, 4.6 inches

- Concentric atraumatic rings to stabilize and rotate eye globe
- 15.3 mm internal dia
- Round handle


TWS 85-2213/51
REVOLVING 3/4 GLOBE FIXATION
RING, TITANIUM



Overall length 130 mm, 5.1 inches

- 11.5 mm internal dia
- Concentric atraumatic rings to stabilize and rotate eye globe
- Can be used for PK, RK and clear corneal incision and trabeculectomy
- Round handle


TWS 85-2230/51
HOFMANN-THORNTON REVOLVING
GLOBE FIXATION RING, TITANIUM



Overall length 130 mm 5.1 inches

- 9 teeth to stabilize and rotate eye globe
- 14 mm internal dia
- Round handle


TWS 85-2227/05
THORNTON REVOLVING 3/4 GLOBE
FIXATION RING, TITANIUM



Overall length 127 mm, 5 inches

- 9 teeth to stabilize and rotate eye globe
- 12 mm internal dia
- Round handle

TWS 85-2228/05
THORNTON REVOLVING GLOBE
FIXATION RING, TITANIUM




Overall length 128 mm, 5 inches

- 9 teeth to stabilize and rotate eye globe
- 13 mm internal dia
- Round handle

FIXATION RINGS

TWS 85-2312/44

3/4 GLOBE FIXATION RING, TITANIUM



Overall length 12 mm, 4.4 inches

- 15 mm external dia
- Other side with trapezoidal teeth
- Round handle

TWS 85-2326/49

ATRAUMATIC GLOBE FIXATION RING FOR LASER AND REFRACTIVE SURGERT, TITANIUM




Overall length 12 mm, 4.4 inches

- OZ 015(13 mm inner dia)
- OZ 015.01(14 mm inner dia)
- OZ 015.0'(15 mm inner dia)
- Round handle

TWS 85-2327/49


FLIERINGA FIXATION RING, TITANIUM



- OZ 019 13.2 mm external dia
- OZ 019.01 15.2 mm external dia
- OZ 019.02 17.2 mm external dia
- OZ 019.03 19.2 mm external dia

TWS 85-23132/53

MENDEZ DEGREE GAUGE FOR MARKING AXIS MS ASTIGMATISM, TITANIUM




Overall length 135 mm, 5.3 inches

- 12.2 mm internal dia & 20 mm external dia
- Used to determine incision angles
- Calibrated every 10 degrees from 0° to 180°
- Round handle

TWS 85--2330/51

MENDEZ DEGREE GAUGE FOR MARKING AXIS MS ASTIGMATISM, TITANIUM




Overall length 130 mm, 5.1 inches

- 12.2 mm internal dia & 20 mm external dia
- Used to determine incision angles
- Calibrated every 10 degrees from 0° to 180°
- Round handle

TWS 85-2332/52

MENDEZ DEGREE GAUGE FOR MARKING AXIS MS ASTIGMATISM, TITANIUM




Overall length 132 mm, 5.2 inches

- 12.1 mm internal dia & 16.1 mm external dia
- Used to determine incision angles
- Calibrated every 10 degrees from 0° to 180°
- Round handle

TWS 85-2316/49

MENDEZ DEGREE GAUGE FOR MARKING AXIS MS ASTIGMATISM, TITANIUM




Overall length 126 mm 4.9 inches

- 11 mm internal dia & 14 mm external dia
- Used to determine incision angles
- Calibrated every 10 degrees from 0° to 180°
- Round handle

FIXATION RINGS

TWS 85-2413/53

MENDEZ DEGREE GAUGE FOR MARKING
AXIS MS ASTIGMATISM, TITANIUM




- 18.1 mm external dia
- 12.1 mm internal dia

Overall length 134 mm, 5.3 inches

- Round handle

TWS 85-2427/50

DEGREE FIXATION RING, TITANIUM




- 12 teeth to stabilize and rotate eye globe
- 17 mm external dia
- 14 mm internal dia
- Round handle

Overall length 127 mm, 5.0 inches

TWS 852417/46


KOCH-MENDEZ DEGREE GAUGE



- Degree ring with atraumatic teeth, stainless
- 16 mm external dia, 12 mm inner dia
- Round handle, titanium

Overall length 117 mm, 4.6 inches


TWS 85-2412/48
LASEK FUNNEL, TITANIUM 8.5 MM



- To apply alcohol solution on epithelium
- For 8 mm optical zone
- Round handle

Overall length 122 mm, 4.8 inches

TWS 85-2422/48
LASEK FUNNEL, TITANIUM 9.5 MM




- To apply alcohol solution on epithelium
- For 9 mm optical zone
- Round handle

Overall length 122 mm, 4.8 inches

CALIPERS

TWS 85-8625/34

CASTROVIEJO CALIPER, TITANIUM




Overall length 86 mm, 3.4 inches

- Measures up to 20 mm on both sides

TWS 85-7625/03


CASTROVIEJO CALIPER, TITANIUM



Overall length 76 mm, 3 inches

TWS 85-2586/34

CASTROVIEJO CALIPER CURVED, TITANIUM




Overall length 86 mm, 3.4 inches

- Measures up to 20 mm on both sides

TWS 85-2512/48

CALIPER



Overall length 121 mm, 4.8 inches

- Working tip 0.15 × 0.75 × 3.0 mm with a (mark) on the distance 0.5 mm stainless
- Round handle, titanium

TWS 85-2526/49

PHE HOLDER, TITANIUM




Overall length 126 mm, 4.9 inches

- Tip with slot for blades 0.45 mm wide

MARKERS

TWS 85-2612/47
THORNTON RADIAL MARKER LOW PROFILE




- **OM 002** 3 blades
- **OM 002.01** 4 blades
- **OM 002.02** 5 blades

Overall length 120mm, 4.7 inches

- Working part 12 mm external dia, stainless
- Round handle, titanium

TWS 85-2620/47
THORNTON RADIAL MARKER LOW PROFILE




- **OM003** 4 blades
- **OM003.01** 6 blades
- **OM003.02** 8 blades
- **OM003.03** 10 blades
- **OM 003.04** 12 blades
- **OM003.05** 16 blades

Overall length 120mm, 4.7 inches


- Working part 14 mm external dia, stainless
- Round handle, titanium

TWS 85-2018/47
HÖFFER OPTIC ZONE MARKER WITH CENTERING PIN LOW PROFILE




- Overall length 120-128 mm, 4.7-5.0 inches**
- | | | | |
|----------------------------|----------------------------|----------------------------|----------------------------|
| ■ MS 004 2.00 mm | ■ OZ 004.07 3.75 mm | ■ MS 004.14 5.50 mm | ■ OZ 004.21 8.50 mm |
| ■ MS 004.01 2.25 mm | ■ OZ 004.08 4.00 mm | ■ OZ 004.15 5.75 mm | ■ OZ 004.22 9.00 mm |
| ■ OZ 004.02 2.50 mm | ■ MS 004.09 4.25 mm | ■ OZ 004.16 6.00 mm | ■ OZ 9.50 mm |
| ■ OZ 004.03 2.75 mm | ■ MS 004.10 4.50 mm | ■ OZ 004.17 6.50 mm | ■ Round handle, titanium |
| ■ MS 004.04 3.00 mm | ■ MS 004.11 4.75 mm | ■ OZ 004.18 7.00 mm | ■ Heal, stainless |
| ■ OZ 004.05 3.25 mm | ■ OZ 004.12 5.00 mm | ■ MS 004.19 7.50 mm | |
| ■ OZ 004.06 3.50 mm | ■ OZ 004.13 5.25 mm | ■ OZ 004.20 8.00 mm | |

TWS 85-2602/47
BORES AXIS MARKER




- Overall length 120mm, 4.7 inches**
- Used with Mendez degree gauge to mark the axis
 - Working tip 12 mm long, stainless
 - Round handle, titanium

TWS 85-2622/48
NUIJTS- SOLOMON TORIC AXIS MARKER WITH 5MM



- Overall length 122 mm, 4.8 inches**
- Used with Mendez Degree Gauge
 - Center Opening, stainless
 - Circular opening for maximum visibility
 - Round handle, titanium

TWS 85-2622/48
NUIJTS- SOLOMON TORIC AXIS MARKER WITH 4.5MM




- Overall length 122 mm, 4.8 inches**
- Used with Mendez Degree Gauge
 - Center Opening, stainless
 - Circular opening for maximum visibility
 - Round handle, titanium

MARKERS

TWS 85-2723/48

TORIC AXIS MARKER, TITANIUM




Overall length 123 mm, 4.8 inches

- Used with Mendez Degree Gauge
- 2 blade radial marks
- 12 mm external dia, 6 mm inner dia
- Round handle

TWS 85-2709/43

AXIS MARKER




Overall length 109 mm, 4.3 inches

- 11.4 mm spread between 2 blades, stainless
- 2.5 mm inner dia of ring
- 3.2 mm external dia of ring
- Round handle, titanium

TWS 85-2728/49

OUTSIDE OPTIC ZONE MARKER
LOW PROFILY




Overall length 123- 128 mm, 4.9 – 5.1 inches

- | | | | |
|---------------------|---------------------|----------------------|--------------------------|
| ■ OM 005.01 2.50 mm | ■ OM 005.11 5.25 mm | ■ OM 005.21 9.25 mm | ■ OM 005.31 11.75 mm |
| ■ OM 005.02 2.75 mm | ■ OM 005.12 5.50 mm | ■ OM 005.22 9.50 mm | ■ OM 005.32 12.00 mm |
| ■ OM 005.03 3.00 mm | ■ OM 005.13 5.75 mm | ■ OM 005.23 9.75 mm | ■ OM 005.33 6.25 mm |
| ■ OM 005.04 3.25 mm | ■ OM 005.14 6.00 mm | ■ OM 005.24 10.00 mm | ■ OM 005.34 6.75 mm |
| ■ OM 005.05 3.50 mm | ■ OM 005.15 6.50 mm | ■ OM 005.25 10.25 mm | ■ OM 005.35 7.25 mm |
| ■ OM 005.06 3.75 mm | ■ OM 005.16 7.00 mm | ■ OM 005.26 10.50 mm | ■ OM 005.36 7.75 mm |
| ■ OM 005.07 4.00 mm | ■ OM 005.17 7.50 mm | ■ OM 005.27 10.75 mm | ■ OM 005.37 8.25 mm |
| ■ OM 005.08 4.25 mm | ■ OM 005.18 8.00 mm | ■ OM 005.28 11.00 mm | ■ OM 005.38 8.75 mm |
| ■ OM 005.09 4.50 mm | ■ OM 005.19 8.50 mm | ■ OM 005.29 11.25 mm | ■ Head, stainless |
| ■ OM 005.10 5.00 mm | ■ OM 005.20 9.00 mm | ■ OM 005.30 11.50 mm | ■ Round handle, titanium |

TWS 85-2725/51

OUTSIDE OPTIC ZONE MARKER
LOW PROFILY




Overall length 123- 128 mm, 4.9 – 5.1 inches

- | | | | |
|---------------------|---------------------|---------------------|--------------------------|
| ■ OM 006.01 2.50 mm | ■ OM 006.06 4.00 mm | ■ OM 006.12 5.50 mm | ■ OM 006.18 8.00 mm |
| ■ OM 006.02 2.75 mm | ■ OM 006.07 4.25 mm | ■ OM 006.13 5.75 mm | ■ OM 006.19 8.50 mm |
| ■ OM 006.03 3.00 mm | ■ OM 006.08 4.50 mm | ■ OM 006.14 6.00 mm | ■ OM 006.20 9.00 mm |
| ■ OM 006.04 3.25 mm | ■ OM 006.09 4.75 mm | ■ OM 006.15 6.50 mm | ■ OM 006.21 10.00 mm |
| ■ OM 006.05 3.50 mm | ■ OM 006.10 5.00 mm | ■ OM 006.16 7.00 mm | ■ Head, stainless |
| ■ OM 006.06 3.75 mm | ■ OM 006.11 5.25 mm | ■ OM 006.17 7.50 mm | ■ Round handle, titanium |

TWS 85-2724/47

THORNTON OPTIC ZONE
MARKER WITH CENTERING PIN
LOW PROFILE




Overall length 120-124 mm, 4.7- 4.9 inches

- | | | | |
|---------------------|---------------------|---------------------|--------------------------|
| ■ OM 007.01 2.50 mm | ■ OM 007.06 4.00 mm | ■ MS 007.12 5.50 mm | ■ MS 007.18 8.00 mm |
| ■ OM 007.02 2.75 mm | ■ OM 007.07 4.25 mm | ■ MS 007.13 5.75 mm | ■ MS 007.19 8.50 mm |
| ■ OM 007.03 3.00 mm | ■ OM 007.08 4.50 mm | ■ MS 007.14 6.00 mm | ■ Head, stainless |
| ■ OM 007.04 3.25 mm | ■ OM 007.09 4.75 mm | ■ OM 007.15 6.50 mm | ■ Round handle, titanium |
| ■ OM 007.05 3.50 mm | ■ OM 007.10 5.00 mm | ■ OM 007.16 7.00 mm | |
| ■ OM 007.06 3.75 mm | ■ OM 007.11 5.25 mm | ■ OM 007.17 7.50 mm | |

TWS 85-2735/13

MENDEZ MARKER "COIN",
TITANIUM OM009.01 4.0 x
5.0mm




Overall length 35 mm, 1.3 inches

- Flat working tips 3.5 x 4.0 mm

MARKERS

TWS 85-2812/47
OPTIC ZONE MARKER



Overall length 121 mm, 4.7 inches

- **OM 012** working part 5.2 × 2.5 mm, stainless
- **OM 012.0** working part 5.2 × 3.0 mm stainless
- **OM 012.0** working part 5.7 × 2.5 mm, stainless
- **OM 012.0** working part 5.7 × 3.0 mm stainless
- Round handle, titanium


TWS 85-2850/59
THORNTON DOUBLE-ENDED MARKER



Overall length 150 mm, 5.9 inches

- 8 blades, 14 mm external dia, stainless
- 6 blades, 14 mm external dia, stainless
- Round handle, titanium


TWS 85-2857/22
HARMS TRABECULOTOME,
LEFT OM 018.01 RIGHT



Overall length 57 mm, 2.2 inches

- Tip 9.2 mm long stainless
- Distance between tips 3.5 mm
- Round handle, titanium


TWS 85-2824/49
GULANI LASIK MARKER



Overall length 124 mm, 4.9 inches

- Two rings 2.5 mm and 4 mm inner dia, stainless
- Round handle, titanium


TWS 85-2820/47
OSHER-NEUMANN CORNEAL LOW
PROFILE MARKER



Overall length 120 mm, 4.7 inches

- 2 blades, stainless
- External ring 14 mm dia
- Round handle, titanium


TWS 85-2821/48
LRI MARKER, 6.0&a8.0 mm CHORD LENGTHS,
WITH DEGREE GAUGES, TITANIUM



Overall length 121 mm, 4.8 inches

- Automatically creating marks for 6 and 8 mm chord length of Limbal Relaxing Incisions
- Designed for Limbal Relaxing Incisions and IOL implantation
- Round handle

TWS 85-2812/48
LRI MARKER, 40-60-80, WITH
DEGREE GAUGES, TITANIUM



Overall length 121 mm, 4.8 inches

- Automatically creating marks at 40-60-80 degrees of Limbal Relaxing Incisions
- Designed for Limbal Relaxing Incisions and IOL implantation
- Round handle

MARKERS

TWS 85-2928/50

INSIDE OPTIC ZONE MARKER LOW PROFILE, TITANIUM



Overall length 128 mm, 5.0 inches

- Angled 45 degree working part internal dia 7.00 mm
- Round handle

TWS 85-2921/50

INSIDE OPTIC ZONE MARKER LOW PROFILE, TITANIUM



Overall length 128mm, 5.0 inches

- Angled 45 degree working part, internal dia 8.00 mm
- Round handle

TWS 85-2912/50

INSIDE OPTIC ZONE MARKER LOW PROFILE DOUBLE-ENDED, TITANIUM



Overall length 128mm, 5.0 inches

- Internal dia 7.00 mm
- Internal dia 8.00 mm
- Round handle

TWS 85-2982/50

TORIC IOL MARKER, TITANIUM




Overall length 18 mm, 4.6 inches

- Internal dia 7 mm
- 4 marks
- Round handle

TWS 85-2918/46

TORIC IOL MARKER, TITANIUM



Overall length 18 mm, 4.6 inches

- Internal dia 8 mm
- 4 marks
- Round handle

TWS 85-2920/47

OUTSIDE OPTIC ZONE MARKER LOW PROFILE, DOUBLE-ENDED



Overall length 120 mm, 4.7 inches

- Heads stainless
- External dia 7.75 mm
- External dia 8.00 mm
- Round handle, titanium

TWS 85-2911/47

OUTSIDE OPTIC ZONE MARKER LOW PROFILE DOUBLE-ENDED




Overall length 120 mm, 4.7 inches

- Heads, stainless
- External dia 8.25 mm
- External dia 8.50 mm
- Round handle, titanium

VARIOUS INSTRUMENTS


TWS 85-3012/47
WEBER LENS LOOP



Overall length 120 mm, 4.7 inches

- Loop 6 × 9 mm with cross, stainless
- Round handle, titanium


TWS 85-3020/47
WILDER LENS LOOP OL 002.013.5
mm wide



Overall length 120 mm, 4.7 inches

- 4 × 15 mm curved loop, stainless
- Round handle, titanium


TWS 85-3013/50
LOOP



Overall length 130 mm, 5.0 inches

- OL 003 – 6 × 11.5 mm, stainless
- OL 003.03 – 4 × 11.5 mm, stainless
- OL 003.04 – 3 × 11.5 mm, stainless
- OL 003.05 – 2.7 × 6.5 mm, stainless
- OL 003.06 – 4.5 × 11.5 mm, stainless


TWS 85-3014/57
EYE THERMOCAUTER



Overall length 146 mm, 5.7 inches

- Olive 5 mm dia
- Working tip 0.5 mm dia, stainless
- Round handle, titanium

TWS 85-3021/47
CORNEOSCLERAL PUNCH



Overall length 120 mm, 4.7 inches

- Set of tips: 0.50 mm, 0.75 mm, 1.00 mm, 1.50 mm
- Flat handle, titanium


TWS 85-3082/50
KELLY PUNCH
OP 0050.75mm dia



Overall length 128 mm, 5.0 inches

- Punch 1.0 mm dia, stainless
- Flat handle, titanium

TWS 85-3061/49
MEMBRANE PEELER




Overall length 126 mm, 4.9 inches

- Working tip 0.9 mm dia, stainless
- Spherical 105 degree tip
- Round handle, titanium

PROBES

TWS 85-3114/55

CASTROVIEJO LACHRYMAL DILATOR
DOUBLE-ENDED, TITANIUM




Overall length 140 mm, 5.5 inches

- Needle - shaped tip
- Middle sharp tip (like candle)
- Round handle

TWS 85-1310/39

WILDER LACHRYMAL DILATOR,
MEDIUM TAPER, TITANIUM




Overall length 100 mm, 3.9 inches

- OPT 002 needle shape tip 18 mm long
- OPT 002.01 needle shape tip 22 mm long
- OPT 002.02 needle shape tip 30 mm long
- OPT 002.03 needle shape tip 16 mm long
- OPT 002.04 needle shape tip 20 mm long
- OPT 002.05 needle shape tip 24 mm long
- Round handle

TWS 85-3170/27

RUEDERMANN LACHRYMAL
DILATOR, STAINLESS




Overall length 70 mm, 2.7 inches

- 0.2 mm shoot needle type tip 3.5 mm dia
- Round smooth handle

TWS 85-3138/54

CASTROVIEJO LACHRYMAL
DILATOR DOUBLE-ENDED,
STAINLESS




Overall length 138 mm, 5.4 inches

- Sharp and blunt tips
- Round handle

TWS 85-3108/42

WILDER LACHRYMAL DILATOR,
TITANIUM




Overall length 108 mm, 4.2 inches

- OPT 011 Tip 0.35 mm dia
- OPT 011.01 Tip 0.4 mm dia
- Round handle

TWS 85-3118/54

CASTROVIEJO LACHRYMAL DILATOR
DOUBLE-ENDED, TITANIUM




Overall length 138 mm, 5.4 inches

- Tip 21.5 mm, 0.3 mm dia
- Tip 24 mm long, 0.3 mm dia
- Round handle

PROBES/DILATORS

TWS 85-3230/51

BOWMAN LACHRYMAL PROBE
DOUBLE-ENDED,
STAINLESS




Overall length 130 mm, 5.1 inches

- **OPT 006** tips 0.5 mm and 0.6 mm dia (Size 0000-000)
- **OPT 006.01** tips 0.7 mm and 0.8 mm dia (Size 00-0)
- **OPT 006.02** tips 0.9 mm and 1.1 mm dia (Size 1-2)
- **OPT 006.03** tips 1.3 mm and 1.4 mm dia (Size 3-4)
- **OPT 006.04** tips 1.5 mm and 1.6 mm dia (Size 5-6)
- **OPT 006.05** tips 1.8 mm and 1.9 mm dia (Size 7-8)

TWS 85-1032/43

QUICKERT LACHRYMAL
INTUBATION PROBE, STAINLESS




Overall length 10 mm, 4.3 inches

- **OPT 007** tip 0.4 mm, 0.7 mm dia
- **OPT 007.01** tip 0.8 mm, 0.8 mm dia
- **OPT 007.02** tip 0.9 mm, 0.9 mm dia
- **OPT 007.03** tip 1.1 mm, 1.1 mm dia
- **OPT 007.04** tip 1.3 mm, 1.3 mm dia
- **OPT 007.05** tip 1.4 mm, 1.4 mm dia

TWS 85-3215/59

PROBE, STAINLESS




Overall length 150 mm, 5.9 inches

- Working tips 0.4 mm and 0.5 mm dia (Size 0000-0000)

TWS 85-3214/55

PROBE, STAINLESS




Overall length 140 mm, 5.5 inches

- **OPT 013** Tips 0.6 mm and 0.7 mm dia (Size 000-00)
- **OPT 013.01** Tips 0.7 mm and 0.8 mm dia (Size 00-0)

TWS 85-3214/55

PROBE, STAINLESS




Overall length 130 mm, 5.1 inches

- Working tips 0.5 mm and 0.6 mm dia (Size 0000-000)

TWS 85-3246/57

WORST PIGTAIL PROBE




Overall length 146 mm, 5.7 inches

- Working parts with holes 0.6 mm dia, stainless
- Used in plastic surgery
- Round handle, titanium

PROBES

TWS 85-3310/43

HOOK DOUBLE-ENDED




Overall length 10 mm, 4.3 inches

- Curved tips 0.8 mm with oval holes, stainless
- Round handle titanium

TWS 85-3311/46


IL LACRIMAL PROBE



Overall length 117 mm, 4.6 inches

- Curved tips with holes 0.35 mm dia, stainless
- Can be also used in plastic surgery
- Round handle, titanium

CANNULAS




TWS 85-3448/18
POSTERIOR
CHAMBER CANNULA, STAINLESS



Overall length 48 mm, 1.8 inches

- TC tip 3 mm
- Tube 23Ga




TWS 85-3449/19
CANNULA CURVED WITH LATERAL PORT,
STAINLESS



Overall length 49 mm, 1.9 inches

- Angled 45 degree 13 mm long
 - Dead rounded tip
- Tube 19Ga




TWS 85-3447/18
ASPIRATING CANNULA FOR SUBRETINAL
CONTENTS, STAINLESS



Overall length 47 mm, 1.8 inches

- Flattened tip 0.46 mm wide and 3 mm long
- Tube 21Ga




TWS 85-3452/20
ASPIRATING CANNULA FOR SUBRETINAL
CONTENTS, STAINLESS



Overall length 52 mm, 2.0 inches

- Round tip 0.45 mm dia and 2mm long
- Tube 21Ga




TWS 85-3451/20
ASPIRATING CANNULA FOR SUBRETINAL
CONTENTS, STAINLESS



Overall length 51 mm, 2.0 inches

- Flattened tip 0.46 x 0.3 mm wide
- Angled part of tip 3 mm long
 - Tube 27 Ga




TWS 85-3443/16
ANTERIOR CHAMBER ASPIRATING
WASHOUT CANNULA, STAINLESS



Overall length 43 mm, 1.6 inches

- Angled tube 30 degree, 23Ga
- Spatulated tip 6 mm long




TWS 85-3439/15
BURATTO LASIK IRRIGATING CANNULA,
STAINLESS




Overall length 39mm, 1.5 inches

- Front port 0.2 mm
- Two side ports 0.2 mm
- Tube lattened at the end up to 0.4 mm
 - Tube 23Ga




TWS 85-3444/17
IRRIGATING
CANNULA, STAINLESS



Overall length 44 mm, 1.7 inches

- Front port 0.5 mm dia
- Tube 21 Ga




TWS 85-4434/17
IRRIGATING
CANNULA, STAINLESS



Overall length 44 mm, 1.7 inches

- **MS 021.01** - 0.5mm dia
- **OC 021.06** - 0.6mm dia
- **S 021.0** - 0.5mm dia (ov.length 60mm, 2.3 inches)
- **MS 021.0** - 0.8mm dia




TWS 85-4334/17
IRRIGATING CANNULA,
STAINLESS



Overall length 44 mm, 1.7 inches

- **MS 021.03** - 0.5mm dia
- **MS 021.04** - 0.7mm dia




TWS 85-3834/15
MCINTYRE ANTERIOR CHAMBER MAINTAINER,
STAINLESS



Overall length 38 mm, 1.5 inches

- Curved tip 8.5 mm long with notch 0.75 mm long
- Flattened tip




TWS 85-4134/16
LASIK IRRIGATION CANNULA, STAINLESS



Overall length 41 mm, 1.6 inches

- Curved tip 15 mm long
- Side port 0.3 mm dia
- Tube 21Ga


CANNULAS



TWS 85-4135/16
LASIK IRRIGATION CANNULA, STAINLESS




- Overall length 41 mm, 1.6 inches
- Curved tip 15 mm long
- Front port 0.3 mm dia
- Tube 21Ga



TWS 85-3541/16
LASIK IRRIGATION CANNULA, STAINLESS




- Overall length 41 mm, 1.6 inches
- Curved tip 15 mm long
- Two asymmetric side ports 0.3 mm dia
- Tube 21Ga



TWS 85-4641/16
LASIK IRRIGATION CANNULA, STAINLESS




- Overall length 41 mm, 1.6 inches
- Curved part 10 mm long
- 6 side ports 0.35 mm dia
- Tube 21G



TWS 85-4741/16
RAININ NUCLEUS HYDRODISSECTOR
ANGLED, STAINLESS




- Overall length 41 mm, 1.6 inches
- Angled part 8 mm long
- Oval front port 0.58 x 0.4 mm dia
- Tube flattened at the end
- Tube 25G



TWS 85-4841/16
RAININ NUCLEUS HYDRODISSECTOR,
CURVED, STAINLESS




- Overall length 41 mm, 1.6 inches
- Oval front port 0.58 x 0.4 mm dia
- Tube flattened at the end
- Tube 25G



TWS 85-3551/02
NUCLEUS HYDRODISSECTOR,
STAINLESS




- Overall length 51 mm, 2 inches
- Front port 0.9 mm dia
- Curved tube 33 mm long
- Tube 18Ga



TWS 85-3546
NUCLEUS HYDRODISSECTOR,
STAINLESS



- Overall length 46 mm, 1.8 inches
- Front port 0.7 mm dia
- Curved tube 33 mm long
- Tube 19G



TWS 85-3533/13
PEARCE J-SHAPED MICRO
HYDRODISSECTOR, STAINLESS




- Overall length 33 mm, 1.3 inches
- Front port 0.3 mm dia
- Hook 1.5 x 2 mm
- Tube 25Ga

TWS 85-3555/18
ADAPTER FROM "LUER" TO "RECORD",
TITANIUM

TWS 85-3534/13
ADAPTER FROM "RECORD" TO "LUER",
TITANIUM


CANNULAS



TWS 85-3638/15
SIMCOE IRRIGATING/ASPIRATING
CANNULA, STAINLESS


Overall length 38 mm, 1.5 inches

- Aspirating port 0.2 mm dia
- Silicone tube 160 mm long
 - Tube 23Ga




TWS 85-1806
HESLIN ANTERIOR CHAMBER MAINTAINER,
STAINLESS

Overall length 180 mm, 7.08 inches
Beveled TC blunt tip 6 mm long
Tip 23Ga




TWS 85-1836/78
HESLIN ANTERIOR CHAMBER MAINTAINER,
STAINLESS

Overall length 180 mm, 7.08 inches
Beveled polished blunt tip 6 mm long
Tip 23Ga



TWS 85-3618/78
ANTERIOR CHAMBER MAINTAINER,
STAINLESS


Overall length 180 mm, 7.08 inches
Working tip 1.4 mm long
Beveled tip 20G



TWS 85-1809/78
LEWICKY ANTERIOR CHAMBER
MAINTAINER, STAINLESS

Overall length 180 mm, 7.08 inches

- Working tip 6 mm long
- Screw type tip 20Ga



TWS 85-1801/78
ANTERIOR CHAMBER
MAINTAINER, STAINLESS


Overall length 180 mm, 7.08 inches

- Straight blunt tip with 6 mm long working end
 - Suturing flange
 - Tip 20Ga



CANNULAS


TWS 85-3710/41
 ASPIRATION HANDPIECE FOR BIMANUAL
 TECHNIQUE
 MS 003.01 Tube 0.7 mm (22Ga)



- One top port 0.35 mm
- Tube 0.8 mm (21Ga), stainless
- Smooth tip
- Round handle titanium
- Can be produced with TC tip

Overall length 105 mm, 4.1 inches


TWS 85-3715/41
 IRRIGATION HANDPIECE FOR
 BIMANUAL TECHNIQUE
 ON 004.01 Tube 0.7 mm (22Ga)



- Two side ports 0.35 mm
- Tube 0.8 mm (21Ga), stainless
- Smooth tip
- Round handle titanium
- Can be produced with TC tip

Overall length 105 mm, 4.1 inches


TWS 85-3751/41
 IRRIGATION HANDPIECE FOR
 BIMANUAL TECHNIQUE
 ON 025 Tube 1.2 mm (19Ga)



- One central port
- Tube 0.8 mm (21Ga), stainless
- Round handle titanium

Overall length 105 mm, 4.1 inches


TWS 85-3721/40
 ASPIRATION HANDPIECE FOR
 BIMANUAL TECHNIQUE



- One top port 0.3 mm
- Tube 0.6 mm (23Ga), stainless
- Smooth tip
- Round handle titanium
- Can be produced with TC tip

Overall length 102 mm, 4.0 inches


TWS 85-3712/40
 IRRIGATION HANDPIECE FOR
 BIMANUAL TECHNIQUE



- Two side ports 0.3 mm
- Tube 0.6 mm (23Ga), stainless
- Smooth tip
- Round handle titanium
- Can be produced with TC tip

Overall length 102 mm, 4.0 inches

TWS 85-3750/41
 IRRIGATION ENDOTHELIAL STRIPPER




- Tube 0.9 mm (20Ga), stainless
- Designed for Descemet's stripping
- Round handle titanium

Overall length 105 mm, 4.1 inches

IRRIGATION CHOPPERS


TWS 85-3810/42
NAGAHARA IRRIGATION CHOPPER



Overall length 108 mm, 4.2 inches

- Nagahara tip 1.35 mm long with rounded platform
- Two side ports 0.6 mm
- Tube 0.9 mm (20G), stainless
- Round handle, titanium


TWS 85-3818/42
IRRIGATION CHOPPER



Overall length 108 mm, 4.2 inches

- Tip 1.8 mm long with sharp pointed platform
- Two side ports 0.6 mm
- Tube 0.9 mm (20G), stainless
- Round handle, titanium


TWS 85-3881/42
IRRIGATION CHOPPER



Overall length 108 mm, 4.2 inches

- Angled tip, stainless
- Round handle, titanium


TWS 85-3808/42
IRRIGATION CHOPPER



Overall length 108 mm, 4.2 inches

- Tip with a ball 0.5 mm dia
- Two side ports 0.6 mm
- Tube 0.9 mm (20G), stainless
- Round handle, titanium


TWS 85-3815/41
IRRIGATION CHOPPER



Overall length 105 mm, 4.1 inches

- Triangle shape knife on the tip
- One central port with beveled tip
- Tube 0.8 mm (21G), stainless
- Round handle, titanium

TWS 85-3851/41
CHANG IRRIGATION CHOPPER,
RIGHT




Overall length 105 mm, 4.1 inches

- Tip, stainless
- For right-handed surgeon
- Round handle, titanium

IRRIGATION CHOPPERS

TWS 85-3915/41
CHANG IRRIGATION CHOPPER,
LEFT




- Tip, stainless
- For left-handed surgeon

Overall length 105 mm, 4.1 inches

- Round handle, titanium

TWS 85-3951/41
NAGAHARA IRRIGATION
CHOPPER, RIGHT




- One central port
- Sharp blade 1 mm long

Overall length 105 mm, 4.1 inches

- Tube 0.9 mm (20G), stainless
- Round handle, titanium

TWS 85-3905/41
NAGAHARA IRRIGATION
CHOPPER, LEFT




- One central port
- Sharp blade 1 mm long

Overall length 105 mm, 4.1 inches

- Tube 0.9 mm (20G), stainless
- Round handle, titanium

TWS 85-3950/41
UNIVERSAL IRRIGATION CHOPPER




- One central port
- Beveled 30 degree tip
- Sharp blade 1 mm long

Overall length 105 mm, 4.1 inches

- Tube 0.9 mm (20G), stainless
- Round handle, titanium

TWS 85-3995/41
ROSEN IRRIGATION CHOPPER




- One central port
- Blade 1 mm long

Overall length 105 mm, 4.1 inches

- Tube 0.9 mm (20G), stainless
- Round handle, titanium

TWS 85-3959/41
IRRIGATION CHOPPER




- One central port
- Blade 1.1 mm long

Overall length 105 mm, 4.1 inches

- Tube 0.9 mm (20G), stainless
- Round handle, titanium

CURETTES


TWS 85-9127/50
CHALAZION CURETTE



Overall length 127 mm, 5,0 inches

- Working part 1.8 mm stainless steel
- Round handle, titanium


TWS 85-9129/50
CHALAZION CURETTE



Overall length 129 mm, 5,0 inches

- Working part 2.2 mm stainless steel
- Round handle, titanium


TWS 85-6620/47
MEYERHOEFER CHALAZION
CURETTE



Overall length 120 mm, 4,7 inches

- Working part 1,5 mm stainless steel
- Round handle, titanium


TWS 85-6720/47
MEYERHOEFER CHALAZION
CURETTE



Overall length 120 mm, 4,7 inches

- Working part 1,75 mm stainless steel
- Round handle, titanium


TWS 85-6820/47
MEYERHOEFER CHALAZION
CURETTE



Overall length 120 mm, 4,7 inches

- Working part 2,0 mm stainless steel
- Round handle, titanium


TWS 85-6920/47
MEYERHOEFER CHALAZION
CURETTE



Overall length 120 mm, 4,7 inches

- Working part 2,5 mm stainless steel
- Round handle, titanium

TWS 85-6520/47
MEYERHOEFER CHALAZION
CURETTE



Overall length 120 mm, 4,7 inches

- Working part 3,5 mm stainless steel
- Round handle, titanium



TW
Surgical

Treasure Wright (U.K.) Ltd.

**Address: 128-160 , City Road,
London EC1V 2NX. United Kingdom**

www.twrightmedical.com

info@twrightmedical.com



TW Surgical



**UK
CA CE**