




TW
Surgical

MEASURING DEVICES



**UK
CA** **CE**



TW Surgical

Treasure Wright (U.K.) Ltd.



www.twrightmedical.com



info@twrightmedical.com




TW Surgical

DRILL BITS



TOWNLEY CALIPER

TWS 84-2039/48
4" (102mm)
takes inside and outside measurements up to 4" (10cm)



BONE COMPASS

TWS 84-2039/49
8" (203mm)

Used to determine screw length required during internal fracture fixation (mm graduated). Gauges snap apart for easy cleaning.



MINI SCREW DEPTH GAUGE

TWS 84-2039/50
4 3/4" (121mm)
up to 30mm



SMALL SCREW DEPTH GAUGE

TWS 84-2039/51
6 1/4" (159mm) up to 50mm



LARGE SCREW DEPTH GAUGE

TWS 84-2039/52
8 3/4" (222mm) up to 100mm




X-LARGE SCREW DEPTH GAUGE

TWS 84-2039/53
11" (279mm)
up to 120mm

Treasure Wright (U.K.) Ltd.


128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com

MEASURING DEVICES




CASTROVIEJO CALIPER

TWS 84-2039/54
3 1/2" (89mm)
graduated to 20mm
for hand surgery




CASTROVIEJO CALIPER

TWS 84-2039/55
6 1/2" (165mm)
graduated to 40mm




POLK FINGER GONIOMETER

TWS 84-2039/56
6" (152mm)
cm and inch graduated



ALUMINUM CALIPER

TWS 84-2039/57
9 1/8" (232mm)
inch and mm graduated
up to 6" (180mm)
inside and outside measuring
faces and depth bar



Inside measuring faced
extra long (60mm) for
neurosurgery.

NEURO-CALIPER

TWS 84-2039/58
9" (229mm)
stainless steel
inch and mm graduated
up to 5 3/4" (145mm)
inside and outside measuring
faces and depth bar



TW Surgical

MEASURING DEVICES



- TWS 84-2039/59 6" (152mm)
- TWS 84-2039/60 8" (203mm)
- TWS 84-2039/61 12" (305mm)
- TWS 84-2039/62" (508mm)

RULER

flexible
inch and mm graduated



Ruler for x-rays.


X-RAY SCALE

TWS 84-2039/63
50cm length
mm graduated



MEASURING TAPE

TWS 84-2039/65
linen, inch and cm
graduated up to 150cm



Front

Back

Reverse side has handy gauge to measure K-Wires, Steinmann Pins, twist drills, and pins and rods.

Note:

Gauge Measures
K-Wires .028"-.062"
.7mm-1.6mm

Pins 5/64"-3/16"
2.0mm-4.76mm

Rods .197"-1/4"
5.0mm-6.35mm

K-Wire RULER & PIN GAUGE

TWS 84-2039/64
6" (152mm)
inch and mm graduated

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



EASY CLEAN DEPTH GAUGE

	Screw Size	Hook Probe	Measuring Range
TWS 84-2039/66	1.2-1.5mm	0.8mm	0-30mm
TWS 84-2039/67	2.0-2.4mm	1.2mm	0-30mm
TWS 84-2039/68	2.0-2.4mm	1.2mm	0-50mm
TWS 84-2039/69	2.7mm	1.7mm	0-50mm
TWS 84-2039/70	3.5-4.0mm	2.2mm	0-50mm
TWS 84-2039/71	3.5-4.0mm	2.2mm	0-80mm



TW
Surgical

Treasure Wright (U.K.) Ltd.

**Address: 128-160 , City Road,
London EC1V 2NX. United Kingdom**

www.twrightmedical.com

info@twrightmedical.com



TW Surgical



**UK
CA CE**