

TW Surgical Z Clips



TW Surgical Z clips are vital surgical instruments employed to secure and handle sterile items. They feature a spring-loaded mechanism and serrated teeth jaws allowing for a firm grip and secure closure, ridged gripping surfaces designed to provide holding and minimizing slippage.

TW Surgical Z clips are designed for smooth operation and a secure grip, making them ideal for use with our drape holders and versatile enough for any clipping needs.

TW Surgical Z clips aid in maintaining a sterile field by securely holding items in place, preventing contamination during surgical procedures.

Surgeons heavily rely on these **TW Surgical Z clips** to organize and manage surgical equipment, optimizing the surgical environment and promoting effective infection control.


High Quality Plastic Copolymer is highly suitable for medical sterilization, offering excellent dimensional stability and chemical resistance to hot steam

- **Steam Autoclaving:** Suitable for repeated sterilization cycles
- **Ethylene Oxide (Eto):** Excellent compatibility, allowing for low-temperature sterilization.
- **Plasma/Hydrogen Peroxide:** Suitable, allowing for low-temperature sterilization
- **Grade:** Medical-grade Plastic Copolymer is recommended for healthcare applications.
- **Biocompatibility:** Resistant to bodily fluids.

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom.
info@twrightmedical.com

TW Surgical Z Clips



TW Surgical Z clips can hold up to one or 2 tubing lines to organize all your patient lines and keep tubing off the floor to avoid accidental pulls and falls.

- Secure and handle sterile items in surgical field
- Spring-loaded mechanism and serrated teeth jaws allowing for a firm grip
- Provide holding and minimizing slippage.
- Manufactured from high-quality stainless steel and Plastic Copolymer
- Offer durability and rust resistance
- Durable and high impact,
- Reusable with autoclavable material
- Hold up to 1 or 2 tubing lines to organize all your patient lines

TW Surgical Z clips organize and manage surgical equipment, optimizing the surgical environment and promoting effective infection control.

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom.
info@twrightmedical.com

TW Surgical Z Clips




TWS 00-5326/5006	TW Surgical Z Mini clips, Sharp Teeth Secure and handle sterile items in surgical field 53mm Long, 26mm Width, 50mm Height Teeth Width 35mm Circular Hole 10mm, Open Ended Tubing Hole 10mm x 14mm Height Spring-loaded mechanism secure holding power Serrated 6 nos. of teeth jaws for a firm grip Closing device for bags Material: Stainless Steel and Plastic Copolymer Cleanable, Reusable and Autoclavable
TWS 00-5728/5509	TW Surgical Z Mini clips, Smooth Teeth Secure and handle sterile items in surgical field 57mm Long, 28mm Width, 55mm Height Teeth Width 28mm Oval Hole 10mm x 14mm, Open Ended Tubing Hole 10mm x 16mm Height Spring-loaded mechanism secure holding power Serrated 8 nos. of teeth jaws for a firm grip Closing device for bags Material: Stainless Steel and Plastic Copolymer Cleanable, Reusable and Autoclavable

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom.
info@twrightmedical.com

TW Surgical Z Clips



TWS 00-8660/8509	TW Surgical Z Large Clips, Sharp Teeth Secure and handle sterile items in surgical field 86mm Long, 60mm Width, 85mm Height Teeth Width 40mm Circular Hole 10mm, Open ended Tubing Hole 15mm x 17mm Height Spring-loaded mechanism secure holding power Serrated 9 nos. of teeth jaws for a firm grip Closing device for bags Material: Stainless Steel and Plastic Copolymer Cleanable, Reusable and Autoclavable
TWS 00-8660/8512	TW Surgical Z Large Clips, Smooth Teeth Secure and handle sterile items in surgical field 86mm Long, 60mm Width, 85mm Height Teeth Width 40mm Circular Hole 10mm, Open ended Tubing Hole 15mm x 17mm Height Spring-loaded mechanism secure holding power Serrated 12 nos. of teeth jaws for a firm grip Closing device for bags Material: Stainless Steel and Plastic Copolymer Cleanable, Reusable and Autoclavable

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom.
info@twrightmedical.com