


# Nerve - Retractors - Nerve Hooks

## Graham Nerve Hooks




**Graham Nerve Hook** is a specialized surgical tool to use for elevating and broadening the edges of the dissected sites

- Round Blunt Hook Promoting Atraumatic Retraction.
- Straight Tapering Shaft for Deep Surgical Access
- Ergonomic Solid Handle Ensuring Superior Control.

TWS 25-5221/08	Graham Nerve Hook, Straight, Blunt, 8 mm round width 16.5cm (6 1/2") Length Material: Stainless Steel
TWS 25-5221/09	Graham Nerve Hook, Straight, Blunt, 9 mm round width 16.5cm (6 1/2") Length Material: Stainless Steel
TWS 25-5221/10	Graham Nerve Hook, Straight, Blunt, 10 mm round width 16.5cm (6 1/2") Length Material: Stainless Steel
TWS 25-5221/11	Graham Nerve Hook, Straight, Blunt 11 mm round width 16.5cm (6 1/2") Length Material: Stainless Steel

# Nerve - Retractors - Nerve Hooks

## Flat Nerve Hooks



**Flat Nerve Hook** is a specialized surgical tool to use for elevating and broadening the edges of the dissected sites

- Flat Blunt Hook Promoting Atraumatic Retraction.
- Straight Tapering Shaft for Deep Surgical Access
- Ergonomic Solid Handle Ensuring Superior Control.

TWS 25-5223/08	Flat Nerve Hook, Straight, Blunt, 8 mm L Ended Reach width 16.5cm (6 1/2") Length Material: Stainless Steel
TWS 25-5223/09	Flat Nerve Hook, Straight, Blunt, 9 mm L Ended Reach width 16.5cm (6 1/2") Length Material: Stainless Steel
TWS 25-5223/10	Flat Nerve Hook, Straight, Blunt, 10 mm L Ended Reach width 16.5cm (6 1/2") Length Material: Stainless Steel
TWS 25-5223/11	Flat Nerve Hook, Straight, Blunt 11 mm L Ended Reach width 16.5cm (6 1/2") Length Material: Stainless Steel