

Orthopaedic - Lead Hand

is used to maintain hand and finger position in hand surgery. Manufactured from malleable lead for immobilization. Oxidation of lead surfaces is a normal process. Avoid using chemical sterilization as it can cause severe corrosion.



TWS 10-1080/3625	Lead Hands, Adult 36mm x 25mm Material: Lead (Autoclavable)
TWS 10-1080/2516	Lead Hands, Child 25mm x 16mm Material: Lead (Autoclavable)
TWS 10-1080/3625B	Elastic Silicone Bands (Adult) for fixation to Lead Hands (Autoclavable) (Pack of 10 Nos.)
TWS 10-1080/2516B	Elastic Silicone Bands (Adult) for fixation to Lead Hands (Autoclavable) (Pack of 10 Nos.)

Contains lead: HANDLE WITH CARE!