


Arm Supports



*Arm supports for
any procedure...
Choose your style!*



A. Kross Arm Support

CAT# TWS-7500

Perfect for arm chair, prone, lateral and neuro procedure positioning
SS arm cradle measures 6" x 20"
Disposable Kross Kloth sleeve makes this a sanitary choice, one roll included
Requires one clark socket for mounting
(Kross Cloth available, See pg 21)

B. Multi-Position Armboard

CAT# TWS -6000

6"W x 20"L Phenolic operating surface
Permits complete range of positioning
Stainless steel support bar
Comes w/ 2" Pressure Management pad
Provides arm support for Lateral, Prone, Neuro, Chair and Park Bench positioning
Requires clark socket or super clamp for mounting to siderail of surgical table
For contoured armboard pad - see pg 8

C. Shoulder Mount Armboard

CAT# TWS-4970

100% Radiolucent Carbon Fiber
6" x 26" X-ray translucent surface
Slides under table top cushion
No mounting device, no metal parts
Non-locking swivel gear allows for quick and easy 180° positioning
Comes w/ 2" Pressure Management pad

D. Surgical Table Armguard

CAT# TWS-8000

X-Ray translucent
8" High x 16" Long
Protects patient extremities
Easily slips under table cushion
Prevents nerve damage to patients
Armguard Pads available - see pg 7, 21

E. Standard Armboards

Phenolic operating surface
Rotates 180°, one-finger operation
Standard 24" or 26" length
Easy to attach to all North American surgical tables
Pressure Management and standard pads available in 1", 2", 3", or 4"
(See pg 7)

CAT# TWS -2001

6" x 24" Standard Armboard w/o pad

CAT# TWS-4001

6" x 26" Standard Armboard w/o pad

Standard Armboard and Pad Sets

CAT# TWS-2000

24" Standard Armboard with 2" pad

CAT# TWS-4000

26" Standard Armboard with 2" pad

CAT# TWS-4000-1

26" Standard Armboard with 1" pad

CAT# TWS-4000-3

26" Standard Armboard with 3" pad

CAT# TWS-4000-4

26" Standard Armboard with 4" pad

F. Armboard

CAT# TWS -4501

NEW TWS operating surface
Standard 26" length
Rotates 180°; 21 locking increments
Easy to attach to all surgical tables

Board and pad combination available
Pressure Management and standard pads available in 1", 2", 3", or 4"
(See pg 7)

G. Quick Release Armboard

CAT# TWS-4950

100% Radiolucent Carbon Fiber
6" x 26" X-ray translucent surface
Quick release mounting device
Non-locking swivel gear allows for quick and easy 180° positioning
Comes w/ 2" Pressure Management pad
Mounts to any surgical table siderail

H. Wide Armboard

CAT# TWS-4110

11" x 26" phenolic operating surface
Easy one-finger positioning, in 21 locking increments
Wide surface provides for ultimate patient comfort
Comes w/ 2" Pressure Management pad
Mounts to any surgical table siderail

I. Clamping Armboard

CAT# TWS-4051

Allows for secure positioning on siderail during procedures requiring the table back to be raised w/ a spring loaded rail-lock clamping mechanism
6" x 26" Phenolic operating surface
Rotates 180°; 21 locking increments
Pads available in 1", 2", 3", or 4"
(See pg 7)

J. Pediatric Armboard

CAT# TWS-1200

6" x 20" Phenolic operating surface
Provides close efficient access for staff when working with small patients
Comes with 2" thick pad
Rotates 180°; 21 locking increments
Effortlessly mounts to siderail of any surgical table




B



E

J



I



F



G

H



Surgical Accessories

Arm & Hand Support



TWS-1600
(Shown with Pad)

Arm & Hand Tables

A. Operating Surface

Provides excellent X-ray clarity and C-Arm visualization
3 sizes available (see below)
Supplied with 1" standard pad
Requires 2 clark sockets or round bar clamps for mounting to siderails of all tables

CAT# TWS-1200 [not shown]
12"x30" Table w/ Pad

CAT# TWS-1600
16"x30" Table w/ Pad

CAT# TWS-1800 [not shown]
18"x30" Duraclear™ Table w/ Pad

CAT# TWS-2600
16"x30" Table w/ Pad
Plus Siderail Attachment for additional accessories needed in minor hand surgical procedures
Plus an adjustable/removable telescoping support leg for added strength and stability.

B. New Lightweight Tables – only 7LB's

Provides excellent X-ray clarity and C-Arm visualization
2 sizes available (see below)
Supplied with 1" standard pad
Requires 2 straight bar or round bar clamps for mounting to siderails of all surgical tables

CAT# TWS-3200
16" x 30" Lightweight Table w/ Pad – Round Bar Mount

CAT# TWS-3400
16" x 30" Lightweight Table w/ Pad – Flat Bar Mount

C. PHENOLIC Operating Surface

Constructed from sturdy radiolucent material with built-in clamps for easy attachment and removal on all surgical tables
Supplied with a 2" pressure management pad

CAT# TWS-5200

'T' or Winged shape Arm & Hand Table w/ 2" pressure management pad
Legless design provides unimpeded C-arm access and easy access for surgeon comfort

CAT# TWS-6600

Rectangle Arm & Hand Table w/ 2" pressure management pad
Table comes equipped with stainless steel adjustable/retractable support leg for superior stability

CAT# TWS-6900

Hourglass Arm & Hand Table w/ 2" pressure management pad
Hourglass shape provides ample support and comfort for surgeon
Table comes equipped with stainless steel adjustable/retractable support leg for superior stability

D. Wall Mounting Brackets

For Arm & Hand Table Storage

CAT# TWS-5800

Stainless Steel wall mounting brackets keep Arm & Hand tables up out of harms way

E. Arm & Hand Table

Siderail Attachment

Stainless Steel siderail attachment provides mounting space for additional accessories needed in special surgical procedures.
Request TWS-6050 to be added to table

CAT# TWS-6050-18 [not shown]
18" Stainless Steel siderail

CAT# TWS-6050-24
24" Stainless Steel siderail

F. PHENOLIC Under Mount Table

CAT# TWS-8200

Constructed from sturdy 100% radiolucent phenolic
Easily slides under table cushion
Provides 14" x 25" operating surface of unimpeded imaging
Comes with 2" Pressure Management pad

*Arm & Hand
Minor Surgery Tables*


A

TWS-1600
shown w/o pad



C

TWS-5200



TWS-2600
shown w/o pad

TWS Supplies
Replacement pads for
all makes and models
of minor surgery tables
(See page 7)



TWS-6600
shown w/o pad



TWS-6900



B

TWS-3200
shown w/o pad



D



TWS-5800



E



TWS-6050-24



Round Bar
TWS-3200



Flat Bar
TWS-3400

F



TWS-8200
shown w/o pad

Surgical Accesso

Arm Support Pads

Multi-Purpose Custom Made Arm board Pads

Custom replacement pads are available for all styles of arm Supports. Specify foam style, fabric, thickness and attachment device requirements when ordering

Flame retardant

Anti-microbial and fluid resistant

CAT# TWS 2215-03HPD

22" Length x 15" Width x 3" Thickness
Multi-purpose/ Multi-task Arm board Pad
Hard Pressure management pads – 3"
Flame retardant
Anti-microbial and fluid resistant



CAT# TWS 2316-03HPD

23" Length x 16" Width x 3" Thickness
Multi-purpose/ Multi-task Arm board Pad
Hard Pressure management pads – 3"
Flame retardant
Anti-microbial and fluid resistant



CAT# TWS 0616-03HPD

6" Length x 16" Width x 3" Thickness
Multi-purpose/ Multi-task Arm board Pad
Hard Pressure management pads – 3"
Flame retardant
Anti-microbial and fluid resistant

