


Hand Instrument

TENDON STRIPPERS




TW SURGICAL



- TWS 53-1303/11**
Tendon Stripper 3mm inner dia. X 11 inch length
- TWS 53-1303/14**
Tendon Stripper 3mm inner dia. X 14 inch length
- TWS 53-1303/24**
Tendon Stripper 3mm inner dia. X 24 inch length



- TWS 53-1304/14**
Tendon Stripper 4mm inner dia. X 14 inch length
- TWS 53-1304/24**
Tendon Stripper 4mm inner dia. X 24 inch length



- TWS 53-1305/11**
Tendon Stripper 5mm inner dia. X 11 inch length
- TWS 53-1305/14**
Tendon Stripper 5mm inner dia. X 14 inch length
- TWS 53-1305/20**
Tendon Stripper 5mm inner dia. X 20 inch length
- TWS 53-1305/24**
Tendon Stripper 5mm inner dia. X 24 inch length



- TWS 53-1306/14**
Tendon Stripper 6mm inner dia. X 14 inch length
- TWS 53-1306/24**
Tendon Stripper 6mm inner dia. X 24 inch length

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, UK
info@twrightmedical.com

Hand Instrument



TW SURGICAL



TENDON WEAVING FORCEPS

Provide valuable assistance during tendon weaving.
The conically shaped fine point on the jaw separates tendon without causing long splits.
The jaw can grasp and hold a small tendon stump securely.




TWS 53-5302/05
Tendon Weaving Forcep 4-3/4 inch (12cm) length
Curved,
Delicate sharp point
Single pin tooth

TWS 53-5302/06
Tendon Weaving Forcep 6 inch (15cm) length
Curved,
Delicate sharp point
Single pin tooth



TWS 53-5303/05
Tendon Weaving Forcep 4-3/4 inch (12cm) length
Straight,
Delicate sharp point
Single pin tooth

TWS 53-5303/06
Tendon Weaving Forcep 6 inch (15cm) length
Straight,
Delicate sharp point
Single pin tooth



Treasure Wright (U.K.) Ltd.



128-160 City Road, London, EC1V 2NX, UK
info@twrightmedical.com

Hand Instrument




TW SURGICAL


TENDON GRASPING FORCEPS



TWS 53-5401/06
TENDON GRASPING FORCEPS
WITH NEEDLE TEETH
6 inch (15cm) Length
two pairs of pin teeth at jaws
Closure is adjustable by the thumb-screw
Ratchet lock holds jaws firmly



TWS 53-5401/07
TENDON GRASPING FORCEPS
WITH 4x5 ATRAUMATIC TEETH
5-3/4 inch (14.5cm) Length



TWS 53-5401/08
TENDON GRASPING FORCEPS
WITH 4x5 ATRAUMATIC NEEDLE TEETH
AND O SHAPE JAW
6 1/4 inch (16cm) Length

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, UK
info@twrightmedical.com

Hand Instrument


TENDON GRASPING FORCEPS



TW SURGICAL



**TWS 53-5401/09
TENDON GRASPING FORCEPS**
reduce permanent tendon damage
when applied to tendon ends.
6 inch (15cm) Length,
2.5mm wide jaws



**TWS 53-5355/02
DIETRICH TENDON FIXATION NEEDLE**
Temporarily immobilizes the tendon,
restricting muscular action during the braiding process.,

Treasure Wright (U.K.) Ltd.


128-160 City Road, London, EC1V 2NX, UK
info@twrightmedical.com

Hand Instrument


TENDON GRASPING FORCEPS



TW SURGICAL



TWS 53-5401/10
BRAND TENDON TUNNELING/PASSING FORCEPS
9-7/8inch (25cm) Length,



TWS 53-5368/02
TENDON STRIPPER 2.0mm dia.
(25cm) Length,

TWS 53-5368/03
TENDON STRIPPER 3.0mm dia.
(25cm) Length,

TWS 53-5368/04
TENDON STRIPPER 4.0mm dia.
(25cm) Length,

TWS 53-5368/05
TENDON STRIPPER 5.0mm dia.
(25cm) Length,

TWS 53-5368/06
TENDON STRIPPER 6.0mm dia.
(25cm) Length,

TWS 53-5368/07
TENDON STRIPPER 7.0mm dia.
(25cm) Length,

Treasure Wright (U.K.) Ltd.




128-160 City Road, London, EC1V 2NX, UK
info@twrightmedical.com

Hand Instrument

TENDON GRASPING FORCEPS






TW SURGICAL



TWS 53-6600/10
Tendon Retrievers,
Rigid Shaft, Small Jaw
Tip dia: 1mm
1 x 2 Teeth. Jaw length: 8mm

TWS 53-6600/15
Tendon Retrievers,
Rigid Shaft, Large Jaw
Tip dia: 1.5mm
1 x 2 Teeth. Jaw length: 8mm



TWS 53-6610/10
Tendon Retrievers,
Flexible Shaft, Small Jaw
Tip dia: 1mm
1 x 2 Teeth. Jaw length: 8mm


Treasure Wright (U.K.) Ltd.
128-160 City Road, London, EC1V 2NX, UK
info@twrightmedical.com

Hand Instrument

TENDON GRASPING FORCEPS



TW SURGICAL



TWS 53-5401/15

Carroll Tendon Retriever

7 x 2 mm jaws

Latitudinal serrations

Curve of shaft approximates arc of the fingers and palm at rest

4-1/4 inch (11cm) Length,



TWS 53-5355/02

TENDON FIXATION NEEDLE

**Temporarily immobilizes the tendon,
restricting muscular action during the braiding process.,**

Treasure Wright (U.K.) Ltd.


128-160 City Road, London, EC1V 2NX, UK
info@twrightmedical.com

Hand Instrument

TENDON GRASPING FORCEPS



TW SURGICAL




TWS 53-5421/05
HALSTED-S TENDON GRASPING FORCEPS

Light curve

5 inch (12.5cm) Length,

Simplify the task of maneuvering a tendon around a phalanx




TWS 53-5421/06
HALSTED-S TENDON GRASPING FORCEPS

Medium curve

5 inch (12.5cm) Length,

Simplify the task of maneuvering a tendon around a phalanx




TWS 53-5421/07
HALSTED-S TENDON GRASPING FORCEPS

Strong curve

5 inch (12.5cm) Length,

Simplify the task of maneuvering a tendon around a phalanx



Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, UK
info@twrightmedical.com

Orthopaedic Instrument – Tendon Strippers

The tendon stripper is used to harvest the semitendinosus tendon and deliver it through the incision.



TWS 53-1340/02M	Tendon Stripper 2mm inner dia. x 400mm length, Open End, Length Marking Material: Stainless Steel and Phenolic
TWS 53-1340/03M	Tendon Stripper 3mm inner dia. x 400mm length, Open End, Length Marking Material: Stainless Steel and Phenolic
TWS 53-1340/04M	Tendon Stripper 4mm inner dia. x 400mm length, Open End, Length Marking Material: Stainless Steel and Phenolic
TWS 53-1340/05M	Tendon Stripper 5mm inner dia. x 400mm length, Open End, Length Marking Material: Stainless Steel and Phenolic
TWS 53-1340/06M	Tendon Stripper 6mm inner dia. x 400mm length, Open End, Length Marking Material: Stainless Steel and Phenolic
TWS 53-1340/07M	Tendon Stripper 7mm inner dia. x 400mm length, Open End, Length Marking Material: Stainless Steel and Phenolic

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW SURGICAL

Tendon Strippers

Brand Tendon Strippers

Wire loop handle, shafts from 285 through 600 mm, circular cutting blade for harvesting grafts



Whole Cutting Circle, measured from the inner diameter



3 mm Ø



4 mm Ø



5 mm Ø



6 mm Ø

TWS 58 - 8069/285	285 mm	TWS 58 - 8070/285	285 mm	TWS 58 - 8071/285	285 mm	TWS 58 - 8072/285	285 mm
TWS 58 - 8069/350	350 mm	TWS 58 - 8070/350	350 mm	TWS 58 - 8071/350	350 mm	TWS 58 - 8072/350	350 mm
TWS 58 - 8069/500	500 mm	TWS 58 - 8070/500	500 mm	TWS 58 - 8071/500	500 mm	TWS 58 - 8072/500	500 mm
TWS 58 - 8069/600	600 mm	TWS 58 - 8070/600	600 mm	TWS 58 - 8071/600	600 mm	TWS 58 - 8072/600	600 mm

Slotted Cutting Circle, measured from the inner diameter



3 mm Ø



4 mm Ø



5 mm Ø



6 mm Ø

TWS 58 - 8073/285	285 mm	TWS 58 - 8074/285	285 mm	TWS 58 - 8075/285	285 mm	TWS 58 - 8076/285	285 mm
TWS 58 - 8073/350	350 mm	TWS 58 - 8074/350	350 mm	TWS 58 - 8075/350	350 mm	TWS 58 - 8076/350	350 mm
TWS 58 - 8073/500	500 mm	TWS 58 - 8074/500	500 mm	TWS 58 - 8075/500	500 mm	TWS 58 - 8076/500	500 mm
TWS 58 - 8073/600	600 mm	TWS 58 - 8074/600	600 mm	TWS 58 - 8075/600	600 mm	TWS 58 - 8076/600	600 mm

Bunnell Tendon Strippers

150 mm, knurled handle, slotted to encircle tendon without severing



TWS 58 - 8077/25	2.5 mm Ø
TWS 58 - 8077/30	3.0 mm Ø
TWS 58 - 8077/40	4.0 mm Ø
TWS 58 - 8077/50	5.0 mm Ø
TWS 58 - 8077/60	6.0 mm Ø
TWS 58 - 8077/70	7.0 mm Ø

Treasure Wright (U.K.) Ltd.

160 City Road, London, EC1V 2NX, UK

info@twrightmedical.com