




TW
Surgical

BONE FILES AND RASPS



**UK
CA** **CE**



TW Surgical

Treasure Wright (U.K.) Ltd.



www.twrightmedical.com



info@twrightmedical.com

BONE FILES AND RASPS




TWS 84-2015/05
6" (152mm)
single ended
2mm width angled
plain serrations
delicate



TWS 84-2015/06
6" (152mm)
single ended
2mm width angled
plain serrations
delicate


BONE FILE #92B

BONE FILE #92A



TWS 84-2015/07
6 3/4" (172mm)
double ended
2mm width angled


NAIL RASP #93



TWS 84-2015/08
7" (178mm)
double ended
3mm and 4mm width

BONE RASP #64

BONE FILES AND RASPS




TWS 84-2015/09
7" (178mm)
double ended 3mm and 4mm width


TWS 84-2015/10
8 1/2" (216mm)
double ended 6mm width

POLOKOFF RASP

Plain serrations



TWS 84-2015/13
7" (178mm)
double ended
5mm width
plain serrations



TWS 84-2015/14
7" (178mm)
double ended
7mm and 5mm width
plain serrations


BONE FILE #10

BONE FILE #12




TW Surgical

BONE FILES AND RASPS




TWS 84-2015/15
7" (178mm)
double ended
5mm width
plain and cross serrations

BONE FILE #12A




TWS 84-2015/16
7" (178mm)
double ended
5mm and 7mm width
cross serrations

BONE FILE #12CA



TWS 84-2015/17
7" (178mm)
double ended
6mm and 5mm width
plain serrations


BONE FILE #33



TWS 84-2015/18
7 1/4" (178mm)
double ended
6mm and 5mm width
plain serrations

BONE FILE #45


BONE FILES AND RASPS



TWS 84-2015/19
 8 1/2" (216mm)
 double ended
 13mm width
 fine and coarse
 cross serrations

For rapid shaping of bone.

BONE RASP



TWS 84-2015/20
 6 1/4" (159mm)
 single ended
 8mm at widest point
 fine cross serrations


JOSEPH RASP



TWS 84-2015/23
 7 1/2" (191mm)
 fine
 single ended
 8mm width

TWS 84-2015/24
 7 1/2" (191mm)
 coarse
 single ended
 8mm width

LEWIS RASP
 Plain serrations



TWS 84-2015/25
 8" (203mm)
 single ended
 backward cutting
 7.5mm width
 plain serrations

COTTLE RASP



TW Surgical

BONE FILES AND RASPS



TWS 84-2015/26
7" (178mm) single ended
backward cutting
10mm width

MALTZ RASP



TWS 84-2015/27
8 1/4" (210mm)
double ended
8mm width
forward and
backward cutting

MALTZ RASP IMPROVED



TWS 84-2015/29
8 1/4" (210mm)
single ended
forward cutting
9mm width

TWS 84-2015/30
8 1/4" (210mm)
single ended
backward cutting
9mm width


AUFRICHT RASP



TWS 84-2015/31
8 1/2" (216mm)
double ended
fine and coarse teeth
4-sided, convex
and flat blades
8mm width

FOMON RASP

BONE FILES AND RASPS




TWS 84-2015/32
 8 1/2" (216mm)
 double ended
 fine and coarse teeth
 4-sided, convex and
 flat blades
 8mm width


KLEINERT KUTZ RASP



BONE FILE




TWS 84-2015/33
 9 1/2" (241mm)
 single ended
 bayonet shank
 half round blade
 fine teeth
 15mm width



PUTTI BONE RASP

TWS 84-2015/34
 10 1/2" (267mm)
 double ended
 flat blades
 18mm width tapers to 4mm
 one end curved up
 one end curved sideways




PUTTI BONE RASP

TWS 84-2015/35
 12" (305mm)
 double ended
 one round blade
 10mm width tapers to 4mm
 and one half round blade
 15mm width tapers to 5mm



TW Surgical

BONE FILES AND RASPS



TWS 84-2015/36
11" (280mm)
single ended
one round blade
2-sided, convex and flat blade
23mm width tapers to 6mm

BONE FILE



TWS 84-2015/37
9.5" (241mm)
3.5mm width
angled up
plain serrations

RASP SALEHI LUMBAR FORAMINAL



Patented system

TWS 84-2015/38 Complete set
TWS 84-2015/39 Handle single ended
TWS 84-2015/40 Handle double ended
TWS 84-2015/41 Key

Individual Rasps

TWS 84-2015/42	0.8	extra fine
TWS 84-2015/43	1.0	very fine
TWS 84-2015/44	1.2	fine
TWS 84-2015/45	1.4	fine
TWS 84-2015/46	1.6	medium fine
TWS 84-2015/47	1.8	medium
TWS 84-2015/48	2.0	medium
TWS 84-2015/49	2.2	medium course
TWS 84-2015/50	2.4	course/rough
TWS 84-2015/51	2.6	very course
TWS 84-2015/52	2.6	diamond coated

RASP SET w/INSERTS

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TWS 84-2015/53
11" (280mm) Curved

CANAL FINDER RASP



TW
Surgical

Treasure Wright (U.K.) Ltd.

**Address: 128-160 , City Road,
London EC1V 2NX. United Kingdom**

www.twrightmedical.com

info@twrightmedical.com



TW Surgical



**UK
CA CE**