

All Made of Stainless Steel



TWS 2832-19 Flat Nose Plier
18.5 cm / 7 1/4" Parallel Grip SS



TWS 2832-20 Flat Nose Plier / Cutter
18.5 cm / 7 1/4" Parallel Grip
Max 1.6mm Cutting Capacity



TWS 2832-21 Long Nose Plier 22 cm / 8 1/4"



TWS 2832-22
Vise Grip 20.5 cm / 8"
Stainless Steel



TWS 2833-00 SLAP HAMMER w/ Threaded Shaft Head



TWS 2834-00 SLAP HAMMER w/ Hexa Shaft Head



TWS 2835-20 Regular Jaw Locking Plier
A: 18 cm / 7" B: 20 cm / 8"
C: 23 cm / 9"



TWS 2833-41
Vise Grip 23 cm / 9"
Locking Plier



TWS 2833-42
Vise Grip 20.5 cm / 8"
Locking Plier



TWS 2833-43
Vise Grip 18 cm / 7"
Locking Plier



TWS 2834-31
Needle Nose Locking Plier
30 cm / 12"



TWS 2834-32
Needle Nose Locking Plier
25.5 cm / 10"



TWS 2834-33
Needle Nose
Locking Plier
22 cm / 8.5"



TWS 2835-31
Needle Nose Locking Plier
A: 22 cm / 8.5"
B: 25.5 cm / 10"
C: 30 cm / 12"

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com




TWS 2845-15
Combi-Flat Nose Pliers 15 cm / 6"
for Bending & Cutting



TWS 2846-18
Round / Flat Nose Pliers
18 cm / 7"



TWS 2847-17
Flat Nose Pliers
with crosswise groove
17 cm / 6 3/4"



TWS 2848-18
Flat Nose Pliers
18 cm / 7"
Parallel Grip



TWS 2849-14 Flat Nose Pliers
13.5 cm / 5 1/4"




with longitudinal groove
TWS 2849-15 Hackstock Flat
Nose Pliers 15cm / 6 1/4"



TWS 2847-18
Flat Nose Pliers 17.5 cm / 7"



TWS 2847-19
Flat Nose Pliers Pointed 19.5 cm / 7 3/4"



TWS 2850-16
Platypus Nail
Pulling Forceps
15.5 cm / 6 1/4"

TWS 2850-17
Platypus Nail
Pulling Forceps
15.5 cm / 6 1/4"



TWS 2850-18
Round Nose Pliers
15.5 cm / 6 1/4"




TWS 2850-19
Needle Nose Pliers
15.5 cm / 6 1/4"



TWS 2850-20
Flat Nose Pliers
15 cm / 6"



TWS 2851-14
Needle Nose Pliers
13.5 cm / 5 1/4"



TWS 2851-15
Needle Nose Pliers
Taper to 2mm Tip
14 cm / 5 1/2"
A: Fine Cross
Serration
B: Smooth



TWS 2851-16
18.5 cm / 7 1/4"
Wire Cutting Pliers
Soft wires up to 1.6mm Ø
Hard wires up to 1.0mm Ø



TWS 2851-18
Martin "Diamant"
16 cm / 6 1/4"

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical



TWS 2860-14
14 cm / 5 1/2"
Cirdlage Pliers / Cutter



TWS 2860-15
20 cm / 8"
Wire Bending Pliers



TWS 2860-16
15 cm / 6"
Wire Bending Pliers




TWS 2860-17
20 cm / 8"
Wire Bending Pliers
TUNGSTEN CARBIDE



TWS 2860-18
15 cm / 6"
Needle Nose Pliers



TWS 2860-19
14 cm / 5 1/2"
Flat Nose Pliers



TWS 2860-20
K-Wire Pliers
A: 13 cm / 5"
B: 17.5 cm / 7"



TWS 2860-21
19 cm / 7 1/2"
Screw/Pin
Removing Pliers



TWS 2860-22
14 cm / 5 1/2"
Pin Puller
Grips Small Pins




TWS 2860-23
16 cm / 6 1/2"
Flat Nose Pliers
w/ Cutter
for Cerclage Wire



TWS 2860-24
19 cm / 7 1/2"
Narrow
Nose Pliers
Heavy Jaws
for Pulling
Wire / Pins




TWS 2860-25
17.5 cm / 7"
Long Jaw Pliers
Zinch Long Jaw
Delicate




TWS 2860-26A TWS 2860-26B TWS 2860-26C
UNIVERSAL Pliers
Serrated w/ Heavyduty Spring and Cutting Edge
Max Cutting Cap: 1.1mm

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



ALL TUNGSTEN CARBIDE



- A: 2mm Delicate Jaw TC
- B: 3mm Jaw TC
- C: 2mm Long Jaw TC
- D: 6mm Jaw TC

TWS 2865-22A TWS 2865-22B
TWS 2865-22C TWS 2865-22D
18 cm / 7"
Double Action Wire Extraction Plier

