



**TW**  
Surgical

# MIS INSTRUMENTS



**UK  
CA** **CE**



TW Surgical

Treasure Wright (U.K.) Ltd.



[www.twrightmedical.com](http://www.twrightmedical.com)



[info@twrightmedical.com](mailto:info@twrightmedical.com)



**POSTERIOR COLLATERAL RETRACTOR**

TWS 84-2039/72 8" (203mm) small  
 TWS 84-2039/73 8 1/2" (216mm) medium  
 TWS 84-2039/74 9" (229mm) large



**COLLATERAL RETRACTOR**

TWS 84-2039/75  
 8" (203mm) length



**TIBIAL RETRACTOR**

TWS 84-2039/76  
 7 3/4" (197mm) length



TW Surgical

## MIS INSTRUMENTS



### DOANE KNEE RETRACTOR

6" (152mm) a.k.a. "Z" Retractor

TWS 84-2039/77

57mm blades, 69mm deep, 11mm wide

TWS 84-2039/78

63mm blades, 76mm deep, 11mm wide

TWS 84-2039/79

Overall length of retractor is 8.5" (216mm)

Overall length of each blade is 2.934" (75mm)

Small blade width is .375" (10mm)

Standard blade width is .5" (13mm)



TWS 84-2039/80

Overall length of retractor is 8.5" (216mm)

Overall length of each blade is 2.5" (64mm)

Large blade length is 3.2" and width is .75" (19mm)

Standard blade length is 2.934" and width is .5" (13mm)



### EXTENDED DUAL WIDTH "Z" RETRACTOR

TWS 84-2039/81 small

TWS 84-2039/82 standard

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom  
info@twrightmedical.com



TW Surgical

## MIS INSTRUMENTS



### PCL RETRACTOR

TWS 84-2039/83  
10 1/2" (267mm) length



### COLLATERAL SOFT TISSUE RETRACTOR

TWS 84-2039/84  
9 1/8" (232mm) length



### TIBIAL RETRACTOR

TWS 84-2039/85  
9" (229mm) 3 hole



### RESECTED PATELLA PROTECTOR

TWS 84-2039/86  
3 1/2" (83mm) length

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom  
info@twrightmedical.com



TW Surgical

## MIS INSTRUMENTS



Designed to straddle the cruciate ligament.

### PCL RETRACTOR

TWS 84-2039/87  
10" (254mm)



### TIBIAL RETRACTOR

TWS 84-2039/88  
10 1/4" (260mm)



### BENT HOHMANN

TWS 84-2039/89  
11" (280mm) right angle  
20mm blade



### RIGHT ANGLE HOHMANN RETRACTOR LONG

TWS 84-2039/90  
12 1/2" (318mm) right angle  
18mm blade width

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom  
info@twrightmedical.com



**TW**  
Surgical

## Treasure Wright (U.K.) Ltd.

---

**Address: 128-160 , City Road,  
London EC1V 2NX. United Kingdom**

---

[www.twrightmedical.com](http://www.twrightmedical.com)

[info@twrightmedical.com](mailto:info@twrightmedical.com)



TW Surgical



**UK  
CA CE**