



# Wire Device

Twist cerclage wire quickly, smoothly and evenly and wrap wire around bone and bring wire ends together.



TWS 82-4790	<p>Wire Tightener            11" Long, 10mm serrated jaw            Twist cerclage wire quickly, smoothly and evenly            Engage ratchet and clamp wire firmly.            Pull back on t-handle until desired tension is reached.            Release jaw by disengaging ratchet</p>
-------------	--



TWS 82-4780	<p>Wire Tightener            6 1/4" Long, 3mm serrated jaw, Tungsten carbide inserts            Twist cerclage wire quickly, smoothly and evenly            Engage ratchet and clamp wire firmly.</p>
-------------	---