


Wire Device

Twist cerclage wire quickly, smoothly and evenly and wrap wire around bone and bring wire ends together.



TWS 82-4790	Wire Tightener 11" Long, 10mm serrated jaw Twist cerclage wire quickly, smoothly and evenly Engage ratchet and clamp wire firmly. Pull back on t-handle until desired tension is reached. Release jaw by disengaging ratchet
-------------	--



TWS 82-4780	Wire Tightener 6 1/4" Long, 3mm serrated jaw, Tungsten carbide inserts Twist cerclage wire quickly, smoothly and evenly Engage ratchet and clamp wire firmly.
-------------	---

TWS Open and Close Wire Passer



Applications:

- Fracture fixation
- Cerclage wiring
- Guiding K-wires around bones or implants

Key Features:

- **Hollow channel** for smooth wire passage
- **Durable and autoclavable** stainless-steel design
- **Ergonomic and precise** control for minimal invasiveness

TWS 57-9040/2735	TWS Open and Close Wire Passer Guide Wire Hole 4mm in both sides 27cm Long, 35mm Wide Material: Stainless Steel,
TWS 57-9040/2735	TWS Open and Close Wire Passer Guide Wire Hole 4mm in both sides 29cm Long 45 Wide Material: Stainless Steel,
TWS 57-9040/2735	TWS Open and Close Wire Passer Guide Wire Hole 4mm in both sides 36cm Long 60 Wide Material: Stainless Steel,

Close Loop Wire Passer



TWS 57-3040/03C	Close Loop Wire Passer, C Shape Guide Wire Hole 4mm in both sides 30cm Long, 3cm Wide for Small Bones Material: Stainless Steel and Phenolic Handle
TWS 57-3040/04C	Close Loop Wire Passer, C Shape Guide Wire Hole 4mm in both sides 30cm Long, 4cm Wide for Small Bones and Medium Bones Material: Stainless Steel and Phenolic Handle
TWS 57-3040/05C	Close Loop Wire Passer, C Shape Guide Wire Hole 4mm in both sides 30cm Long, 5cm Wide for Medium Bones Material: Stainless Steel and Phenolic Handle
TWS 57-3040/06C	Close Loop Wire Passer, C Shape Guide Wire Hole 4mm in both sides 30cm Long, 6cm Wide for Medium Bones and Large Bones Material: Stainless Steel and Phenolic Handle
TWS 57-3040/07C	Close Loop Wire Passer, C Shape Guide Wire Hole 4mm in both sides 30cm Long, 7cm Wide for Large Bones Material: Stainless Steel and Phenolic Handle

Close Loop Wire Passer



TWS 57-3040/03CL	Close Loop Wire Passer, C Shape, Curved Left Guide Wire Hole 4mm in both sides 30cm Long, 3cm Wide for Small Bones Material: Stainless Steel and Phenolic Handle
TWS 57-3040/04CL	Close Loop Wire Passer, C Shape, Curved Left Guide Wire Hole 4mm in both sides 30cm Long, 4cm Wide for Small Bones and Medium Bones Material: Stainless Steel and Phenolic Handle
TWS 57-3040/05CL	Close Loop Wire Passer, C Shape, Curved Left Guide Wire Hole 4mm in both sides 30cm Long, 5cm Wide for Medium Bones Material: Stainless Steel and Phenolic Handle
TWS 57-3040/06CL	Close Loop Wire Passer, C Shape, Curved Left Guide Wire Hole 4mm in both sides 30cm Long, 6cm Wide for Medium Bones and Large Bones Material: Stainless Steel and Phenolic Handle
TWS 57-3040/07CL	Close Loop Wire Passer, C Shape, Curved Left Guide Wire Hole 4mm in both sides 30cm Long, 7cm Wide for Large Bones Material: Stainless Steel and Phenolic Handle

Close Loop Wire Passer



TWS 57-3040/03CL	Close Loop Wire Passer, C Shape, Curved Left Guide Wire Hole 4mm in both sides 30cm Long, 3cm Wide for Small Bones Material: Stainless Steel and Phenolic Handle
TWS 57-3040/04CL	Close Loop Wire Passer, C Shape, Curved Left Guide Wire Hole 4mm in both sides 30cm Long, 4cm Wide for Small Bones and Medium Bones Material: Stainless Steel and Phenolic Handle
TWS 57-3040/05CL	Close Loop Wire Passer, C Shape, Curved Left Guide Wire Hole 4mm in both sides 30cm Long, 5cm Wide for Medium Bones Material: Stainless Steel and Phenolic Handle
TWS 57-3040/06CL	Close Loop Wire Passer, C Shape, Curved Left Guide Wire Hole 4mm in both sides 30cm Long, 6cm Wide for Medium Bones and Large Bones Material: Stainless Steel and Phenolic Handle
TWS 57-3040/07CL	Close Loop Wire Passer, C Shape, Curved Left Guide Wire Hole 4mm in both sides 30cm Long, 7cm Wide for Large Bones Material: Stainless Steel and Phenolic Handle

Close Loop Wire Passer



TWS 57-3040/03CUD	Close Loop Wire Passer, C Shape, Curved Up and Down Guide Wire Hole 4mm in both sides 30cm Long, 3cm Wide for Small Bones Material: Stainless Steel and Phenolic Handle
TWS 57-3040/04CUD	Close Loop Wire Passer, C Shape, Curved Up and Down Guide Wire Hole 4mm in both sides 30cm Long, 4cm Wide for Small Bones and Medium Bones Material: Stainless Steel and Phenolic Handle
TWS 57-3040/05CUD	Close Loop Wire Passer, C Shape, Curved Up and Down Guide Wire Hole 4mm in both sides 30cm Long, 5cm Wide for Medium Bones Material: Stainless Steel and Phenolic Handle
TWS 57-3040/06CUD	Close Loop Wire Passer, C Shape, Curved Up and Down Guide Wire Hole 4mm in both sides 30cm Long, 6cm Wide for Medium Bones and Large Bones Material: Stainless Steel and Phenolic Handle
TWS 57-3040/07CUD	Close Loop Wire Passer, C Shape, Curved Up and Down Guide Wire Hole 4mm in both sides 30cm Long, 7cm Wide for Large Bones Material: Stainless Steel and Phenolic Handle

Wire Guide Passer



TWS 57-2020/20S	Wire Guide Passer Guide Wire Hole 2mm in both sides 20cm Long, 2cm Wide for Small Bones Material: Stainless Steel
TWS 57-2020/25S	Wire Guide Passer Guide Wire Hole 2mm in both sides 20cm Long, 2.5cm Wide for Small Bones Material: Stainless Steel
TWS 57-2020/30S	Wire Guide Passer Guide Wire Hole 3mm in both sides 20cm Long, 3cm Wide for Small Bones Material: Stainless Steel
TWS 57-2020/35S	Wire Guide Passer Guide Wire Hole 3mm in both sides 20cm Long, 3.5cm Wide for Small and Medium Bones Material: Stainless Steel
TWS 57-2020/40S	Wire Guide Passer Guide Wire Hole 4mm in both sides 20cm Long, 4cm Wide for Medium Bones Material: Stainless Steel
TWS 57-2020/45S	Wire Guide Passer Guide Wire Hole 4mm in both sides 20cm Long, 4.5cm Wide for Medium and Large Bones Material: Stainless Steel
TWS 57-2020/50S	Guide Wire Hole 4mm in both sides 20cm Long, 5cm Wide for Large Bones Material: Stainless Steel
TWS 57-2020/55S	Guide Wire Hole 4mm in both sides 20cm Long, 5.5cm Wide for Large Bones Material: Stainless Steel