



TW
Surgical

MAGNETIC NEEDLE FINDER



TW Surgical

Treasure Wright (U.K.) Ltd.



www.twrightmedical.com



info@twrightmedical.com

UK
CA

CE



MAGNETIC NEEDLE FINDER

A Rolling Magnetic Needle Finder guaranteed to pick up the smallest needles – even plastic needles with tiny metal tips. Increase or safety with the easy to use, rolling magnetic needle finder. Needle finder has a strong magnetic force that picks up both needles and sharps with ease. Aluminum Housing makes it easier to see what has been picked up.

- Extremely fine suction control due to tear drop shaped suction interrupter
- Individual adaption of the suction tube
- Conical suction tube prevents plugging
- Atraumatic distal tip
- Conical suction tube
- Picks up sutures, needles, screws, clips, and metal waste



<p>TWS 4000-3/34</p>	<p>Magnetic Needle Finder Features 40 in. pole for easy use and portability Magnetic pulling force: 4.5 lbs Telescopic handle adjusts from 23 in. to 40 in. 2-Polyvinyl chloride wheels: 3 in. diameter x 3/4 in. wide Rectangle Magnetic Tube</p>
----------------------	--



TW
Surgical

Treasure Wright (U.K.) Ltd.

**Address: 128-160 , City Road,
London EC1V 2NX. United Kingdom**

www.twrightmedical.com

info@twrightmedical.com