

## TW Surgical positioning pads and cushions for operating room in memory foam with polyurethane coated fabric cover

**TW Surgical positioning pads and cushions** is our **memory foam / high-density foam** pads and positioning cushions with polyurethane coated fabric cover.

The 55 kg/m<sup>3</sup> shape memory foam reduces significantly contact pressure. For large items, thanks to the dual layered memory foam and high-density foam core, our firm cover

**TW Surgical positioning pads and cushions** yield a high rate of pressure relief and comfort. Easy to handle, ideal for bedsores prevention and radiolucent.

**TW Surgical positioning pads and cushions** are covered with a polyurethane coated fabric with waterproof seams.

### Protection cushions for surgery tables


|                                                                                     | Cat. No.      | Designation<br>Sizes (mm)                                  |
|-------------------------------------------------------------------------------------|---------------|------------------------------------------------------------|
|  | TWS 2351-03RE | Protection cushions for surgery tables<br>230 x 510 x 30mm |

### Head protection

|                                                                                     | Cat. No.      | Designation<br>Sizes (mm)        |
|-------------------------------------------------------------------------------------|---------------|----------------------------------|
|  | TWS 5030-15AP | Anatomical pillow<br>300 x 500mm |

## TW Surgical positioning pads and cushions for operating room in memory foam with polyurethane coated fabric cover


### Arm protection

|                                                                                     | Cat. No.        | Designation<br>Sizes (mm)                        |
|-------------------------------------------------------------------------------------|-----------------|--------------------------------------------------|
|    |                 | Round flat armboard pad                          |
|                                                                                     | TWS 5615/06RARP | 560 x 160mm- Height 60mm                         |
|                                                                                     | TWS 6015/06RARP | 600 x 160mm - Height 60mm                        |
|                                                                                     | TWS 6015/06RARP | 660 x 160mm - Height 60mm                        |
|  |                 | Rectangular flat arm board<br>pad                |
|                                                                                     | TWS 5116/06RARP | 510 x 160mm - Height 60mm                        |
|                                                                                     | TWS 5616/06RARP | 560 x 160mm- Height 60mm                         |
|                                                                                     | TWS 6016/06RARP | 600 x 160mm - Height 60mm                        |
|                                                                                     | TWS 6616/06RARP | 660 x 160mm- Height 60mm                         |
|  |                 | Curved arm board pad                             |
|                                                                                     | TWS 4016/06CAB  | 400 x 160mm - Height 60mm                        |
|                                                                                     | TWS 5016/06CAB  | 500 x 160mm - Height 60mm                        |
|                                                                                     | TWS 6616/06CAB  | 660 x 160mm - Height 60mm                        |
|  |                 | curved armboard pad                              |
|                                                                                     | TWS 6016/03RCAB | 600 x 160mm - Height 30mm                        |
|  | TWS 6016/06CAB  | curved armboard pad<br>600 x 160mm - Height 60mm |

## TW Surgical positioning pads and cushions for operating room in memory foam with polyurethane coated fabric cover

### Positioning pads

Cat. No.                      Designation  
                                         Sizes (mm)




#### Cylinder pad

|                |                    |
|----------------|--------------------|
| TWS 3085/85PCP | 85 mm Dia. x 300mm |
| TWS 4085/85PCP | 85 mm Dia. x 400mm |
| TWS 3010/10PCP | 100mm Dia. x 300mm |
| TWS 4010/10PCP | 100mm Dia. x 400mm |
| TWS 5010/10PCP | 100mm Dia. x 500mm |
| TWS 3015/15PCP | 150mm Dia. x 300mm |
| TWS 4015/15PCP | 150mm Dia. x 400mm |
| TWS 5015/15PCP | 150mm Dia. x 500mm |
| TWS 6015/15PCP | 150mm Dia. x 600mm |

#### Half-cylinder pad

|                 |                  |
|-----------------|------------------|
| TWS 3070/35PHCP | 70 x 35 x 300mm  |
| TWS 4070/35PHCP | 70 x 35 x 400mm  |
| TWS 5070/35PHCP | 70 x 35 x 500mm  |
| TWS 3010/50PHCP | 100 x 50 x 300mm |
| TWS 4010/50PHCP | 100 x 50 x 400mm |
| TWS 5010/50PHCP | 100 x 50 x 500mm |
| TWS 6010/50PHCP | 100 x 50 x 600mm |
| TWS 3012/60PHCP | 120 x 60 x 300mm |
| TWS 4012/60PHCP | 120 x 60 x 400mm |
| TWS 5012/60PHCP | 120 x 60 x 500mm |
| TWS 3015/75PHCP | 150 x 75 x 300mm |
| TWS 4015/75PHCP | 150 x 75 x 400mm |
| TWS 5015/75PHCP | 150 x 75 x 500mm |
| TWS 6015/75PHCP | 150 x 75 x 600mm |



Cat. No.

Rectangle pad



|               |                   |
|---------------|-------------------|
| TWS 1515/30RP | 150 x 150 x 30mm  |
| TWS 1515/40RP | 150 x 150 x 40mm  |
| TWS 2020/40RP | 200 x 200 x 40mm  |
| TWS 3015/40RP | 300 x 150 x 40mm  |
| TWS 3020/50RP | 300 x 200 x 50mm  |
| TWS 3030/50RP | 300 x 300 x 50mm  |
| TWS 4020/50RP | 400 x 200 x 50mm  |
| TWS 1520/10RP | 150 x 200 x 100mm |
| TWS 2020/10RP | 200 x 200 x 100mm |
| TWS 2030/10RP | 200 x 300 x 100mm |

Triangle pad



|               |                   |
|---------------|-------------------|
| TWS 2214/30TP | 220 x 140 x 300mm |
| TWS 2214/40RP | 220 x 140 x 400mm |
| TWS 2222/30RP | 220 x 220 x 300mm |
| TWS 2222/40RP | 220 x 220 x 400mm |
| TWS 3220/30RP | 320 x 200 x 300mm |
| TWS 3220/40RP | 320 x 200 x 400mm |
| TWS 3220/50RP | 320 x 200 x 500mm |
| TWS 4226/40RP | 420 x 260 x 400mm |
| TWS 4226/50RP | 420 x 260 x 500mm |
| TWS 5232/40RP | 520 x 320 x 400mm |
| TWS 5232/50RP | 520 x 320 x 500mm |
| TWS 2222/30RP | 220 x 220 x 300mm |
| TWS 2222/40RP | 220 x 220 x 400mm |

Thoracic pad



|               |                   |
|---------------|-------------------|
| TWS 1520/30TP | 150 x 200 x 300mm |
| TWS 1520/40TP | 150 x 200 x 400mm |
| TWS 1520/50TP | 150 x 200 x 500mm |