



**TW**  
Surgical

# HOLLOWWARE



**UK  
CA** **CE**



TW Surgical

Treasure Wright (U.K.) Ltd.



[www.twrightmedical.com](http://www.twrightmedical.com)



[info@twrightmedical.com](mailto:info@twrightmedical.com)



## Round Sterilizer Drum

### ROUND STERILIZER DRUM

Manufactured from 18/10 stainless steel. Domed top and bottom with sliding removable band and drop down clasps. Hinged lid.

## Round Sterilizer Drum "Filter"

### ROUND STERILIZER DRUMS "FILTER"

Manufactured from 18/10 stainless steel. Perforated hinged lid and bottom with interchangeable filters. Designed for vertical steam flow.

#### MODEL SIZE / DIMENSIONS

TWS 84-1000/03	mm. Ø 125 x 80 H	TWS 84-1000/18	mm. Ø 240 x 145 H
TWS 84-1000/04	mm. Ø 125 x 100 H	TWS 84-1000/19	mm. Ø 240 x 165 H
TWS 84-1000/05	mm. Ø 125 x 125 H	TWS 84-1000/20	mm. Ø 240 x 190 H
TWS 84-1000/06	mm. Ø 150 x 80 H	TWS 84-1000/21	mm. Ø 240 x 240 H
TWS 84-1000/07	mm. Ø 150 x 100 H	TWS 84-1000/22	mm. Ø 240 x 290 H
TWS 84-1000/08	mm. Ø 150 x 130 H	TWS 84-1000/23	mm. Ø 240 x 340 H
TWS 84-1000/09	mm. Ø 150 x 150 H	TWS 84-1000/24	mm. Ø 240 x 390 H
TWS 84-1000/10	mm. Ø 150 x 190 H	TWS 84-1000/25	mm. Ø 290 x 145 H
TWS 84-1000/11	mm. Ø 190 x 125 H	TWS 84-1000/26	mm. Ø 290 x 165 H
TWS 84-1000/12	mm. Ø 190 x 145 H	TWS 84-1000/27	mm. Ø 290 x 190 H
TWS 84-1000/13	mm. Ø 190 x 165 H	TWS 84-1000/28	mm. Ø 290 x 240 H
TWS 84-1000/14	mm. Ø 190 x 190 H	TWS 84-1000/29	mm. Ø 290 x 290 H
TWS 84-1000/15	mm. Ø 190 x 240 H	TWS 84-1000/30	mm. Ø 290 x 340 H
TWS 84-1000/16	mm. Ø 190 x 290 H	TWS 84-1000/31	mm. Ø 290 x 390 H
TWS 84-1000/17	mm. Ø 240 x 120 H	TWS 84-1000/32	mm. Ø 340 x 145 H

#### MODEL SIZE / DIMENSIONS

TWS 84-1000/33	mm. Ø 340 x 165H	TWS 84-1000/47	mm. Ø 125 x 80 H	TWS 84-1000/57	mm. Ø 240 x 290 H	TWS 84-1000/72	mm. Ø 390 x 140 H
TWS 84-1000/34	mm. Ø 340 x 190 H	TWS 84-1000/48	mm. Ø 125 x 100 H	TWS 84-1000/58	mm. Ø 240 x 340 H	TWS 84-1000/73	mm. Ø 390 x 165 H
TWS 84-1000/35	mm. Ø 340 x 240 H	TWS 84-1000/49	mm. Ø 125 x 125 H	TWS 84-1000/59	mm. Ø 240 x 390 H	TWS 84-1000/74	mm. Ø 390 x 190 H
TWS 84-1000/36	mm. Ø 340 x 290 H	TWS 84-1000/45	mm. Ø 150 x 100 H	TWS 84-1000/60	mm. Ø 290 x 145 H	TWS 84-1000/75	mm. Ø 390 x 240 H
TWS 84-1000/37	mm. Ø 340 x 340 H	TWS 84-1000/46	mm. Ø 150 x 130 H	TWS 84-1000/61	mm. Ø 290 x 165 H	TWS 84-1000/76	mm. Ø 390 x 290 H
TWS 84-1000/38	mm. Ø 340 x 390 H	TWS 84-1000/47	mm. Ø 150 x 150 H	TWS 84-1000/62	mm. Ø 290 x 190 H	TWS 84-1000/77	mm. Ø 390 x 340 H
TWS 84-1000/39	mm. Ø 390 x 140 H	TWS 84-1000/48	mm. Ø 150 x 190 H	TWS 84-1000/63	mm. Ø 290 x 240 H	TWS 84-1000/78	mm. Ø 390 x 490 H
TWS 84-1000/40	mm. Ø 390 x 165 H	TWS 84-1000/49	mm. Ø 190 x 145 H	TWS 84-1000/64	mm. Ø 290 x 290 H		
TWS 84-1000/41	mm. Ø 390 x 190 H	TWS 84-1000/50	mm. Ø 190 x 165 H	TWS 84-1000/65	mm. Ø 290 x 340 H		
TWS 84-1000/42	mm. Ø 390 x 240 H	TWS 84-1000/51	mm. Ø 190 x 190 H	TWS 84-1000/66	mm. Ø 290 x 390 H		
TWS 84-1000/43	mm. Ø 390 x 290 H	TWS 84-1000/52	mm. Ø 190 x 240 H	TWS 84-1000/67	mm. Ø 340 x 145 H		
TWS 84-1000/44	mm. Ø 390 x 340 H	TWS 84-1000/53	mm. Ø 240 x 145 H	TWS 84-1000/68	mm. Ø 340 x 165 H		
TWS 84-1000/45	mm. Ø 390 x 390 H	TWS 84-1000/54	mm. Ø 240 x 165 H	TWS 84-1000/69	mm. Ø 340 x 190 H		
TWS 84-1000/46	mm. Ø 390 x 490 H	TWS 84-1000/55	mm. Ø 240 x 190 H	TWS 84-1000/70	mm. Ø 340 x 240 H		
		TWS 84-1000/56	mm. Ø 240 x 240 H	TWS 84-1000/71	mm. Ø 340 x 290 H		



**Instrument Boxes**

### INSTRUMENT BOXES

Manufactured from 18/10 stainless steel. Body and lid pressed without welding.

MODEL	SIZE / DIMENSIONS
TWS 84-1000/79	mm. 180 x 80 x 40 H
TWS 84-1000/80	mm. 200 x 100 x 60 H
TWS 84-1000/81	mm. 220 x 120 x 60 H
TWS 84-1000/82	mm. 250 x 120 x 60 H
TWS 84-1000/83	mm. 300 x 150 x 60 H
TWS 84-1000/84	mm. 320 x 150 x 60 H
TWS 84-1000/85	mm. 400 x 160 x 75 H
TWS 84-1000/86	mm. 500 x 200 x 120 H




**Instrument Boilers**

### INSTRUMENT BOILERS

Manufactured from 18/10 stainless steel. Body, lid and perforated internal tray pressed without welding. Hooks to lift the internal tray.

MODEL	SIZE / DIMENSIONS
TWS 84-1000/87	mm. 180 x 80 x 40 H
TWS 84-1000/88	mm. 200 x 100 x 60 H
TWS 84-1000/89	mm. 220 x 120 x 60 H
TWS 84-1000/90	mm. 250 x 120 x 60 H
TWS 84-1000/90	mm. 300 x 150 x 60 H
TWS 84-1000/91	mm. 320 x 150 x 60 H
TWS 84-1000/92	mm. 400 x 160 x 75 H
TWS 84-1000/93	mm. 500 x 200 x 120 H



**Wire Mesh Basket**

### WIRE MESH BASKETS

Manufactured from 18/10 stainless steel

MODEL	L	I	H
TWS 84-1000/96	mm. 210	145	50
TWS 84-1000/97	mm. 240	240	50
TWS 84-1000/98	mm. 240	240	70
TWS 84-1000/99	mm. 240	240	90
TWS 84-1001/00	mm. 240	320	50
MODEL	L	I	H
TWS 84-1001/01	mm. 480	240	50
TWS 84-1001/02	mm. 540	240	50
TWS 84-1001/03	mm. 540	240	70
TWS 84-1001/04	mm. 540	240	90



## INSTRUMENT TRAYS WITH COVER

Manufactured from 18/10 stainless steel. Pressed without welding.

MODEL	SIZE / DIMENSIONS		MODEL	SIZE / DIMENSIONS
TWS 84-1001/06	mm. 200 x 150 x 25		TWS 84-1001/10	mm. 220 x 170 x 50
TWS 84-1001/07	mm. 260 x 180 x 30		TWS 84-1001/11	mm. 280 x 200 x 50
TWS 84-1001/08	mm. 300 x 220 x 30		TWS 84-1001/12	mm. 320 x 240 x 50
TWS 84-1001/09	mm. 350 x 250 x 40		TWS 84-1001/13	mm. 370 x 270 x 50
			TWS 84-1001/14	mm. 460 x 360 x 55

Instrument Tray with cover



## INSTRUMENT TRAYS WITH OUT COVER

Manufactured from 18/10 stainless steel. Pressed without welding.

MODEL	SIZE / DIMENSIONS		MODEL	SIZE / DIMENSIONS
TWS 84-1001/15	mm. 200 x 150 x 25		TWS 84-1001/19	mm. 220 x 170 x 50
TWS 84-1001/16	mm. 260 x 180 x 30		TWS 84-1001/20	mm. 280 x 200 x 50
TWS 84-1001/17	mm. 300 x 220 x 30		TWS 84-1001/21	mm. 320 x 240 x 50
TWS 84-1001/18	mm. 350 x 250 x 40		TWS 84-1001/22	mm. 370 x 270 x 50
			TWS 84-1001/23	mm. 460 x 360 x 55

Instrument Tray without cover



**MODEL      SIZE / DIMENSIONS**

TWS 84-1001/25	mm. 210 x 110 x 50 W/internal tray
TWS 84-1001/26	mm. 240 x 160 x 50 W/internal tray
TWS 84-1001/27	mm. 300 x 180 x 50 W/internal tray
TWS 84-1001/28	mm. 360 x 210 x 50 W/internal tray
TWS 84-1001/29	mm. 400 x 230 x 55 W/internal tray

## Instrument Trays with & without Cover

**MODEL      SIZE / DIMENSIONS**

TWS 84-1001/30	mm. 210 x 110 x 50
TWS 84-1001/31	mm. 240 x 160 x 50
TWS 84-1001/32	mm. 300 x 180 x 50
TWS 84-1001/33	mm. 360 x 210 x 50
TWS 84-1001/34	mm. 400 x 230 x 55

### INSTRUMENT TRAYS WITH COVER

Manufactured from 18/10 stainless steel. Pressed without welding. Two side handles

### INSTRUMENT TRAYS WITH OUT COVER

Manufactured from 18/10 stainless steel. Pressed without welding. Two side handles



Gallipot

### GALLIPOT

Manufactured from 18/10 stainless steel. Pressed without welding.

**MODEL      SIZE / DIMENSIONS**

TWS 84-1001/36	Ø 40 mm ml.30	TWS 84-1001/41	Ø 150 mm ml.750
TWS 84-1001/37	Ø 60 mm ml.50	TWS 84-1001/42	Ø 180 mm ml.1500
TWS 84-1001/38	Ø 80 mm ml.150	TWS 84-1001/43	Ø 200 mm ml.2100
TWS 84-1001/39	Ø 100 mm ml.250	TWS 84-1001/44	Ø 250 mm ml.3500
TWS 84-1001/40	Ø 120 mm ml.450		



## Round Bowls With & Without Cover

### ROUND BOWLS WITH & WITH OUT COVER

Manufactured from 18/10 stainless steel. Pressed without welding.

#### MODEL SIZE / DIMENSIONS

##### WITH COVER

TWS 84-1001/47	Ø 40 mm ml.30	TWS 84-1001/52	Ø 150 mm ml.750
TWS 84-1001/48	Ø 60 mm ml.50	TWS 84-1001/53	Ø 180 mm ml.1500
TWS 84-1001/49	Ø 80 mm ml.150	TWS 84-1001/54	Ø 200 mm ml.2100
TWS 84-1001/50	Ø 100 mm ml.250	TWS 84-1001/55	Ø 250 mm ml.3500
TWS 84-1001/51	Ø 120 mm ml.450		

##### WITH OUT COVER

TWS 84-1001/56	Ø 40 mm ml.30	TWS 84-1001/61	Ø 150 mm ml.750
TWS 84-1001/57	Ø 60 mm ml.50	TWS 84-1001/62	Ø 180 mm ml.1500
TWS 84-1001/58	Ø 80 mm ml.150	TWS 84-1001/63	Ø 200 mm ml.2100
TWS 84-1001/59	Ø 100 mm ml.250	TWS 84-1001/64	Ø 250 mm ml.3500
TWS 84-1001/60	Ø 120 mm ml.450		



TWS 84-1001/46  
**Sputum Cup**

### SPUTUM CUP

Manufactured from 18/10 stainless steel. Pressed without welding, complete with handle and lid.




**Kidney bowls with Cover**

### KIDNEY BOWLS WITH COVER

Manufactured from 18/10 stainless steel. Pressed without welding.

#### MODEL SIZE / DIMENSIONS

TWS 84-1001/67	160 mm	TWS 84-1001/70	250 mm
TWS 84-1001/68	200 mm	TWS 84-1001/71	280 mm
TWS 84-1001/69	300 mm		



**Kidney bowls**

### KIDNEY BOWLS

Manufactured from 18/10 stainless steel. Pressed without welding.

#### MODEL DIMENSIONI / SIZE / DIMENSIONS / MEDIDAS / ROZMERY

TWS 84-1001/72	160 mm	TWS 84-1001/75	250 mm
TWS 84-1001/73	200 mm	TWS 84-1001/76	280 mm
TWS 84-1001/74	300 mm		



## Irrigator

### IRRIGATOR

Manufactured from 18/10 stainless steel, without welding. with handle and holder

**MODEL**      **DIMENSIONI / SIZE / DIMENSIONS /MEDIDAS /ROZMERY**

TWS 84-1001/78	1000 ml
TWS 84-1001/79	10 1500 ml
TWS 84-1001/80	20 2000 ml



## Forceps Jar

### FORCEP JARS

Manufactured from 18/10 stainless steel.

**MODEL**      **DIMENSIONI / SIZE / DIMENSIONS /MEDIDAS /ROZMERY**

TWS 84-1001/86	Ø 30 x 90 mm
TWS 84-1001/87	Ø 50 x 130 mm
TWS 84-1001/88	Ø 50 x 180 mm



## Measuring Jar

### GRADUATED MEASURING JAR

Manufactured from 18/10 stainless steel, without welding. With handle.

**MODEL**      **DIMENSIONI / SIZE / DIMENSIONS /MEDIDAS /ROZMERY**

TWS 84-1001/84	500 ml
TWS 84-1001/85	2000 ml



## Dressing Jar

### DRESSING JAR

Manufactured from 18/10 stainless steel, with cover.

**MODEL**      **DIMENSIONI / SIZE / DIMENSIONS /MEDIDAS /ROZMERY**

TWS 84-1001/89	mm. 75 x 75	TWS 84-1001/93	mm. 180 x 180
TWS 84-1001/90	mm. 100 x 100	TWS 84-1001/94	mm. 200 x 200
TWS 84-1001/91	mm. 125 x 125	TWS 84-1001/95	mm. 250 x 250
TWS 84-1001/92	mm. 150 x 150	TWS 84-1001/96	mm. 300 x 300



Bed Pan

### BED PAN SLEEPER

Manufactured from 18/10 stainless steel, pressed without welding. With handle.

MODEL    DIMENSIONI / SIZE / DIMENSIONS /MEDIDAS /ROZMERY

TWS 84-1001/97  
Adult

TWS 84-1001/98  
Child



Bed Pan

### BED PAN ROUND

Manufactured from 18/10 stainless steel, pressed without welding. With handle.

MODEL    DIMENSIONI / SIZE / DIMENSIONS /MEDIDAS /ROZMERY

TWS 84-1002/01  
Without cover

TWS 84-1002/02  
With cover

## Bed Pan "Fracture" Type

BED PAN "FRACTURE" TYPE  
Manufactured from 18/10 stainless steel.



MODEL    DIMENSIONI / SIZE / DIMENSIONS / MEDIDAS / ROZMERY

TWS 84-1002/04  
Adult  
TWS 84-1002/05  
Child



TWS 84-1002/06  
**Female Urinal**

FEMALE URINAL  
Manufactured from 18/10 stainless steel. With handle.



TWS 84-1002/07  
**Male Urinal**

MALE URINAL  
Manufactured from 18/10 stainless steel. With handle.



**TW**  
Surgical

## Treasure Wright (U.K.) Ltd.

---

**Address: 128-160 , City Road,  
London EC1V 2NX. United Kingdom**

---

[www.twrightmedical.com](http://www.twrightmedical.com)

[info@twrightmedical.com](mailto:info@twrightmedical.com)



TW Surgical



**UK  
CA CE**