

# Gillies Skin Hook

is a fine, precise surgical tool predominantly used in **plastic surgery, dermatology, neurosurgery, and general wound care**. It features a slender handle with a single, sharp or semi-sharp curved hook at the working end, specifically designed to handle delicate tissues under high precision.



## 1. Skin and Soft Tissue Retraction

- **Exposing the surgical site:** It is used to hook and gently pull back skin edges, providing the surgeon with an unobstructed view of the underlying structures.
- **Minimizing tissue trauma:** Unlike heavy surgical forceps or wide retractors that can crush tissue, the fine single prong inflicts minimal point-puncture. This preserves local blood supply and reduces post-operative scarring.

## 2. Elevating and Manipulating Skin Flaps

- **Flap dissection:** During reconstructive procedures, skin grafts, or cosmetic surgeries, the hook safely lifts and stabilizes skin flaps. This allows the surgeon to perform clean cutting and precise cauterization underneath.

## 3. Minor Surgery and Biopsies

- **Lesion removal:** It is heavily utilized in outpatient procedures for retracting small incisions during the removal of moles, cysts, skin tumors, or when harvesting tissue samples for biopsies.


## 4. Wound Alignment and Closure

- **Precise suturing:** When closing long or facial lacerations, the hook gently stretches the skin edges. This ensures the wound margins align perfectly without inversion or overlapping, achieving optimal cosmetic results.

## 5. Specialized Neurosurgical Use


- **Dural retraction:** Certain variations of the Gillies hook are utilized in craniotomies to safely hook and retract the dura mater (the tough outer membrane covering the brain and spinal cord) to expose deeper neurological structures.

# Gillies Hook



|                   |                                                                                      |
|-------------------|--------------------------------------------------------------------------------------|
| TWS 101-0165/65B  | Gillies Hook, 165cm Length, 6.5mm Inner Diameter, Blunt<br>Material: Stainless Steel |
| TWS 101-0165/75B  | Gillies Hook, 165cm Length, 7.5mm Inner Diameter, Blunt<br>Material: Stainless Steel |
| TWS 101-0165/85B  | Gillies Hook, 165cm Length, 8.5mm Inner Diameter, Blunt<br>Material: Stainless Steel |
| TWS 101-0165/110B | Gillies Hook, 165cm Length, 11mm Inner Diameter, Blunt<br>Material: Stainless Steel  |
| TWS 101-0180/65B  | Gillies Hook, 180cm Length, 6.5mm Inner Diameter, Blunt<br>Material: Stainless Steel |
| TWS 101-0180/75B  | Gillies Hook, 180cm Length, 7.5mm Inner Diameter, Blunt<br>Material: Stainless Steel |
| TWS 101-0180/85B  | Gillies Hook, 180cm Length, 8.5mm Inner Diameter, Blunt<br>Material: Stainless Steel |
| TWS 101-0180/110B | Gillies Hook, 180cm Length, 11mm Inner Diameter, Blunt<br>Material: Stainless Steel  |

## Gillies Hook



|                   |                                                                                      |
|-------------------|--------------------------------------------------------------------------------------|
| TWS 101-0165/65S  | Gillies Hook, 165cm Length, 6.5mm Inner Diameter, Sharp<br>Material: Stainless Steel |
| TWS 101-0165/75S  | Gillies Hook, 165cm Length, 7.5mm Inner Diameter, Sharp<br>Material: Stainless Steel |
| TWS 101-0165/85S  | Gillies Hook, 165cm Length, 8.5mm Inner Diameter, Sharp<br>Material: Stainless Steel |
| TWS 101-0165/110S | Gillies Hook, 165cm Length, 11mm Inner Diameter, Sharp<br>Material: Stainless Steel  |
| TWS 101-0180/65S  | Gillies Hook, 180cm Length, 6.5mm Inner Diameter, Sharp<br>Material: Stainless Steel |
| TWS 101-0180/75S  | Gillies Hook, 180cm Length, 7.5mm Inner Diameter, Sharp<br>Material: Stainless Steel |
| TWS 101-0180/85S  | Gillies Hook, 180cm Length, 8.5mm Inner Diameter, Sharp<br>Material: Stainless Steel |
| TWS 101-0180/110S | Gillies Hook, 180cm Length, 11mm Inner Diameter, Sharp<br>Material: Stainless Steel  |



**TW**  
Surgical

## **Treasure Wright (U.K.) Ltd.**

---

**Address: 128-160 , City Road,  
London EC1V 2NX. United Kingdom**

---

[www.twrightmedical.com](http://www.twrightmedical.com)

[info@twrightmedical.com](mailto:info@twrightmedical.com)