



TW
Surgical

LIPOSUCTION



**UK
CA** **CE**



TW Surgical

Treasure Wright (U.K.) Ltd.



www.twrightmedical.com



info@twrightmedical.com



TW Surgical

INFILTRATION CANNULAS



Tumescent Infiltrator

This cannula is used specifically for the infusion of tumescent solution. The configuration of the ports allows even distribution of the fluid.



Tommey Hub

- TWS 84-3037/53 Ø 2.5 mm X 15 cm
- TWS 84-3037/54 Ø 2.5 mm X 20 cm
- TWS 84-3037/55 Ø 3.0 mm X 20 cm
- TWS 84-3037/56 Ø 3.0 mm X 25 cm



Treaded Hub

- TWS 84-3037/58 Ø 2.0 mm X 15 cm
- TWS 84-3037/59 Ø 3.0 mm X 15 cm
- TWS 84-3037/60 Ø 3.0 mm X 20 cm
- TWS 84-3037/61 Ø 3.0 mm X 30 cm
- TWS 84-3037/62 Ø 4.0 mm X 20 cm
- TWS 84-3037/63 Ø 4.0 mm X 30 cm



Ultra Luer Lock

- TWS 84-3037/78 Ø 1.0 mm X 5 cm
- TWS 84-3037/79 Ø 1.0 mm X 7 cm
- TWS 84-3037/80 Ø 1.5 mm X 7 cm
- TWS 84-3037/81 Ø 1.5 mm X 9 cm



Luer Lock

- TWS 84-3037/82 Ø 1.0 mm X 4 cm
- TWS 84-3037/83 Ø 1.0 mm X 5 cm
- TWS 84-3037/84 Ø 1.5 mm X 7 cm
- TWS 84-3037/85 Ø 1.5 mm X 9 cm
- TWS 84-3037/86 Ø 1.5 mm X 10 cm
- TWS 84-3037/87 Ø 1.5 mm X 12 cm
- TWS 84-3037/88 Ø 2.0 mm X 10 cm
- TWS 84-3037/89 Ø 2.0 mm X 12 cm
- TWS 84-3037/90 Ø 2.0 mm X 15 cm
- TWS 84-3037/91 Ø 2.0 mm X 20 cm
- TWS 84-3037/92 Ø 2.5 mm X 15 cm
- TWS 84-3037/93 Ø 2.5 mm X 20 cm
- TWS 84-3037/94 Ø 3.0 mm X 15 cm
- TWS 84-3037/95 Ø 3.0 mm X 20 cm



Luer lock (Stainless Steel)

- TWS 84-3037/57 Ø 1.0 mm X 4 cm
- TWS 84-3037/64 Ø 1.0 mm X 5 cm
- TWS 84-3037/65 Ø 1.0 mm X 7 cm
- TWS 84-3037/66 Ø 1.5 mm X 7 cm
- TWS 84-3037/67 Ø 1.5 mm X 9 cm
- TWS 84-3037/68 Ø 1.5 mm X 10 cm
- TWS 84-3037/69 Ø 1.5 mm X 12 cm
- TWS 84-3037/70 Ø 2.0 mm X 10 cm
- TWS 84-3037/71 Ø 2.0 mm X 12 cm
- TWS 84-3037/72 Ø 2.0 mm X 15 cm
- TWS 84-3037/73 Ø 2.0 mm X 20 cm
- TWS 84-3037/74 Ø 2.5 mm X 15 cm
- TWS 84-3037/75 Ø 2.5 mm X 20 cm
- TWS 84-3037/76 Ø 3.0 mm X 15 cm
- TWS 84-3037/77 Ø 3.0 mm X 20 cm



Power Cannula

- TWS 84-3037/96 Ø 2.0 mm X 15 cm
- TWS 84-3037/97 Ø 2.0 mm X 20 cm
- TWS 84-3037/98 Ø 2.0 mm X 25 cm
- TWS 84-3037/99 Ø 2.0 mm X 30 cm
- TWS 84-3038/00 Ø 2.0 mm X 35 cm
- TWS 84-3038/01 Ø 3.0 mm X 15 cm
- TWS 84-3038/02 Ø 3.0 mm X 20 cm
- TWS 84-3038/03 Ø 3.0 mm X 25 cm
- TWS 84-3038/04 Ø 3.0 mm X 30 cm
- TWS 84-3038/05 Ø 3.0 mm X 35 cm
- TWS 84-3038/06 Ø 4.0 mm X 15 cm
- TWS 84-3038/07 Ø 4.0 mm X 20 cm
- TWS 84-3038/08 Ø 4.0 mm X 25 cm
- TWS 84-3038/09 Ø 4.0 mm X 30 cm
- TWS 84-3038/10 Ø 4.0 mm X 35 cm
- TWS 84-3038/11 Ø 5.0 mm X 15 cm
- TWS 84-3038/12 Ø 5.0 mm X 20 cm
- TWS 84-3038/13 Ø 5.0 mm X 25 cm
- TWS 84-3038/14 Ø 5.0 mm X 30 cm
- TWS 84-3038/15 Ø 5.0 mm X 35 cm
- TWS 84-3038/16 Ø 6.0 mm X 15 cm
- TWS 84-3038/17 Ø 6.0 mm X 20 cm
- TWS 84-3038/18 Ø 6.0 mm X 25 cm
- TWS 84-3038/19 Ø 6.0 mm X 30 cm
- TWS 84-3038/20 Ø 6.0 mm X 35 cm

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical

INFILTRATION CANNULAS



Tumescent Bayonet Infiltrator

The bayonet shaped handles allow an easier manipulation



Tommey Hub

- TWS 84-3038/21 Ø 2.0 mm X 15 cm
- TWS 84-3038/22 Ø 2.5 mm X 20 cm
- TWS 84-3038/23 Ø 3.0 mm X 20 cm
- TWS 84-3038/24 Ø 3.0 mm X 25 cm



Treaded Hub

- TWS 84-3038/30 Ø 2.0mm X 15 cm
- TWS 84-3038/25 Ø 3.0mm X 15 cm
- TWS 84-3038/26 Ø 3.0mm X 20 cm
- TWS 84-3038/27 Ø 3.0mm X 30 cm
- TWS 84-3038/28 Ø 4.0mm X 20 cm
- TWS 84-3038/29 Ø 4.0mm X 30 cm



Ultra Luer Lock

- TWS 84-3038/31 Ø 1.0 mm X 5 cm
- TWS 84-3038/32 Ø 1.0 mm X 7 cm
- TWS 84-3038/33 Ø 1.5 mm X 7 cm
- TWS 84-3038/34 Ø 1.5 mm X 9 cm



Luer Lock

- TWS 84-3038/35 Ø 1.0 mm X 4 cm
- TWS 84-3038/36 Ø 1.0 mm X 5 cm
- TWS 84-3038/37 Ø 1.0 mm X 7 cm
- TWS 84-3038/38 Ø 1.5 mm X 7 cm
- TWS 84-3038/39 Ø 1.5 mm X 9 cm
- TWS 84-3038/40 Ø 1.5 mm X 10 cm
- TWS 84-3038/41 Ø 1.5 mm X 12 cm
- TWS 84-3038/42 Ø 2.0 mm X 10 cm
- TWS 84-3038/43 Ø 2.0 mm X 12 cm
- TWS 84-3038/44 Ø 2.0 mm X 15 cm
- TWS 84-3038/45 Ø 2.0 mm X 20 cm
- TWS 84-3038/46 Ø 2.5 mm X 15 cm
- TWS 84-3038/47 Ø 2.5 mm X 20 cm
- TWS 84-3038/48 Ø 3.0 mm X 15 cm
- TWS 84-3038/49 Ø 3.0 mm X 20 cm



Luer lock (Stainless Steel)

- TWS 84-3038/64 Ø 1.0 mm X 4 cm
- TWS 84-3038/50 Ø 1.0 mm X 5 cm
- TWS 84-3038/51 Ø 1.0 mm X 7 cm
- TWS 84-3038/52 Ø 1.5 mm X 7 cm
- TWS 84-3038/53 Ø 1.5 mm X 9 cm
- TWS 84-3038/54 Ø 1.5 mm X 10 cm
- TWS 84-3038/55 Ø 1.5 mm X 12 cm
- TWS 84-3038/56 Ø 2.0 mm X 10 cm
- TWS 84-3038/57 Ø 2.0 mm X 12 cm
- TWS 84-3038/58 Ø 2.0 mm X 15 cm
- TWS 84-3038/59 Ø 2.0 mm X 20 cm
- TWS 84-3038/60 Ø 2.5 mm X 15 cm
- TWS 84-3038/61 Ø 2.5 mm X 20 cm
- TWS 84-3038/62 Ø 3.0 mm X 15 cm
- TWS 84-3038/63 Ø 3.0 mm X 20 cm



Power Cannula1.5

- TWS 84-3038/66 Ø 2.0 mm X 15 cm
- TWS 84-3038/67 Ø 2.0 mm X 20 cm
- TWS 84-3038/65 Ø 2.0 mm X 25 cm
- TWS 84-3038/67 Ø 2.0 mm X 30 cm
- TWS 84-3038/68 Ø 2.0 mm X 35 cm
- TWS 84-3038/69 Ø 3.0 mm X 15 cm
- TWS 84-3038/70 Ø 3.0 mm X 20 cm
- TWS 84-3038/71 Ø 3.0 mm X 25 cm
- TWS 84-3038/72 Ø 3.0 mm X 30 cm
- TWS 84-3038/73 Ø 3.0 mm X 35 cm
- TWS 84-3038/74 Ø 4.0 mm X 15 cm
- TWS 84-3038/75 Ø 4.0 mm X 20 cm
- TWS 84-3038/76 Ø 4.0 mm X 25 cm
- TWS 84-3038/77 Ø 4.0 mm X 30 cm
- TWS 84-3038/78 Ø 4.0 mm X 35 cm
- TWS 84-3038/79 Ø 5.0 mm X 15 cm
- TWS 84-3038/80 Ø 5.0 mm X 20 cm
- TWS 84-3038/81 Ø 5.0 mm X 25 cm
- TWS 84-3038/82 Ø 5.0 mm X 30 cm
- TWS 84-3038/83 Ø 5.0 mm X 35 cm
- TWS 84-3038/84 Ø 6.0 mm X 15 cm
- TWS 84-3038/85 Ø 6.0 mm X 20 cm
- TWS 84-3038/86 Ø 6.0 mm X 25 cm
- TWS 84-3038/87 Ø 6.0 mm X 30 cm
- TWS 84-3038/88 Ø 6.0 mm X 35 cm

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical

INFILTRATION CANNULAS



Tumescent Infiltrator Angled

This cannula is used specifically for the infusion of tumescent solution. The configuration of the ports allows even distribution of the fluid.



Tommey Hub

- TWS 84-3038/89 Ø 2.0 mm X 15 cm
- TWS 84-3038/90 Ø 2.5 mm X 20 cm
- TWS 84-3038/91 Ø 3.0 mm X 20 cm
- TWS 84-3038/92 Ø 3.0 mm X 25 cm



Treaded Hub

- TWS 84-3038/93 Ø 2.0mm X 15 cm
- TWS 84-3038/94 Ø 3.0mm X 15 cm
- TWS 84-3038/95 Ø 3.0mm X 20 cm
- TWS 84-3038/96 Ø 3.0mm X 30 cm
- TWS 84-3038/97 Ø 4.0mm X 20 cm
- TWS 84-3038/98 Ø 4.0mm X 30 cm



Ultra Luer Lock

- TWS 84-3038/99 Ø 1.0 mm X 5 cm
- TWS 84-3039/00 Ø 1.0 mm X 7 cm
- TWS 84-3039/01 Ø 1.5 mm X 7 cm
- TWS 84-3039/02 Ø 1.5 mm X 9 cm



Luer Lock

- TWS 84-3039/03 Ø 1.0 mm X 4 cm
- TWS 84-3039/04 Ø 1.0 mm X 5 cm
- TWS 84-3039/05 Ø 1.0 mm X 7 cm
- TWS 84-3039/06 Ø 1.5 mm X 7 cm
- TWS 84-3039/07 Ø 1.5 mm X 9 cm
- TWS 84-3039/08 Ø 1.5 mm X 10 cm
- TWS 84-3039/09 Ø 1.5 mm X 12 cm
- TWS 84-3039/10 Ø 2.0 mm X 10 cm
- TWS 84-3039/11 Ø 2.0 mm X 12 cm
- TWS 84-3039/12 Ø 2.0 mm X 15 cm
- TWS 84-3039/13 Ø 2.0 mm X 20 cm
- TWS 84-3039/14 Ø 2.5 mm X 15 cm
- TWS 84-3039/15 Ø 2.5 mm X 20 cm
- TWS 84-3039/16 Ø 3.0 mm X 15 cm
- TWS 84-3039/17 Ø 3.0 mm X 20 cm



Luer lock (Stainless Steel)

- TWS 84-3039/18 Ø 1.0 mm X 4 cm
- TWS 84-3039/19 Ø 1.0 mm X 5 cm
- TWS 84-3039/20 Ø 1.0 mm X 7 cm
- TWS 84-3039/21 Ø 1.5 mm X 7 cm
- TWS 84-3039/22 Ø 1.5 mm X 9 cm
- TWS 84-3039/23 Ø 1.5 mm X 10 cm
- TWS 84-3039/24 Ø 1.5 mm X 12 cm
- TWS 84-3039/25 Ø 2.0 mm X 10 cm
- TWS 84-3039/26 Ø 2.0 mm X 12 cm
- TWS 84-3039/27 Ø 2.0 mm X 15 cm
- TWS 84-3039/28 Ø 2.0 mm X 20 cm
- TWS 84-3039/29 Ø 2.5 mm X 15 cm
- TWS 84-3039/30 Ø 2.5 mm X 20 cm
- TWS 84-3039/31 Ø 3.0 mm X 15 cm
- TWS 84-3039/32 Ø 3.0 mm X 20 cm



Power Cannula

- TWS 84-3039/34 Ø 2.0 mm X 15 cm
- TWS 84-3039/35 Ø 2.0 mm X 20 cm
- TWS 84-3039/33 Ø 2.0 mm X 25 cm
- TWS 84-3039/36 Ø 2.0 mm X 30 cm
- TWS 84-3039/37 Ø 2.0 mm X 35 cm
- TWS 84-3039/38 Ø 3.0 mm X 15 cm
- TWS 84-3039/39 Ø 3.0 mm X 20 cm
- TWS 84-3039/40 Ø 3.0 mm X 25 cm
- TWS 84-3039/41 Ø 3.0 mm X 30 cm
- TWS 84-3039/42 Ø 3.0 mm X 35 cm
- TWS 84-3039/43 Ø 4.0 mm X 15 cm
- TWS 84-3039/44 Ø 4.0 mm X 20 cm
- TWS 84-3039/45 Ø 4.0 mm X 25 cm
- TWS 84-3039/46 Ø 4.0 mm X 30 cm
- TWS 84-3039/47 Ø 4.0 mm X 35 cm
- TWS 84-3039/48 Ø 5.0 mm X 15 cm
- TWS 84-3039/49 Ø 5.0 mm X 20 cm
- TWS 84-3039/50 Ø 5.0 mm X 25 cm
- TWS 84-3039/51 Ø 5.0 mm X 30 cm
- TWS 84-3039/52 Ø 5.0 mm X 35 cm
- TWS 84-3039/53 Ø 6.0 mm X 15 cm
- TWS 84-3039/54 Ø 6.0 mm X 20 cm
- TWS 84-3039/55 Ø 6.0 mm X 25 cm
- TWS 84-3039/56 Ø 6.0 mm X 30 cm
- TWS 84-3039/57 Ø 6.0 mm X 35 cm

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical

INFILTRATION CANNULAS



Facelift Infiltration Cannula

This cannula is used specifically for the infusion of tumescent solution. The configuration of the ports allows even distribution of the fluid.



Tomme Hub

- TWS 84-3039/58 Ø 2.0 mm X 15 cm
- TWS 84-3039/59 Ø 2.5 mm X 20 cm
- TWS 84-3039/60 Ø 3.0 mm X 20 cm
- TWS 84-3039/61 Ø 3.0 mm X 25 cm



Treaded Hub

- TWS 84-3039/61 Ø 2.0mm X 15 cm
- TWS 84-3039/62 Ø 3.0mm X 15 cm
- TWS 84-3039/63 Ø 3.0mm X 20 cm
- TWS 84-3039/64 Ø 3.0mm X 30 cm
- TWS 84-3039/65 Ø 4.0mm X 20 cm
- TWS 84-3039/66 Ø 4.0mm X 30 cm



Ultra Luer Lock

- TWS 84-3039/67 Ø 1.0 mm X 5 cm
- TWS 84-3039/68 Ø 1.0 mm X 7 cm
- TWS 84-3039/69 Ø 1.5 mm X 7 cm
- TWS 84-3039/70 Ø 1.5 mm X 9 cm



Luer Lock

- TWS 84-3039/71 Ø 1.0 mm X 4 cm
- TWS 84-3039/72 Ø 1.0 mm X 5 cm
- TWS 84-3039/73 Ø 1.0 mm X 7 cm
- TWS 84-3039/74 Ø 1.5 mm X 7 cm
- TWS 84-3039/75 Ø 1.5 mm X 9 cm
- TWS 84-3039/76 Ø 1.5 mm X 10 cm
- TWS 84-3039/77 Ø 1.5 mm X 12 cm
- TWS 84-3039/78 Ø 2.0 mm X 10 cm
- TWS 84-3039/79 Ø 2.0 mm X 12 cm
- TWS 84-3039/80 Ø 2.0 mm X 15 cm
- TWS 84-3039/81 Ø 2.0 mm X 20 cm
- TWS 84-3039/82 Ø 2.5 mm X 15 cm
- TWS 84-3039/83 Ø 2.5 mm X 20 cm
- TWS 84-3039/84 Ø 3.0 mm X 15 cm
- TWS 84-3039/85 Ø 3.0 mm X 20 cm



Luer lock (Stainless Steel)

- TWS 84-3039/86 Ø 1.0 mm X 4 cm
- TWS 84-3039/87 Ø 1.0 mm X 5 cm
- TWS 84-3039/88 Ø 1.0 mm X 7 cm
- TWS 84-3039/89 Ø 1.5 mm X 7 cm
- TWS 84-3039/90 Ø 1.5 mm X 9 cm
- TWS 84-3039/91 Ø 1.5 mm X 10 cm
- TWS 84-3039/92 Ø 1.5 mm X 12 cm
- TWS 84-3039/93 Ø 2.0 mm X 10 cm
- TWS 84-3039/94 Ø 2.0 mm X 12 cm
- TWS 84-3039/95 Ø 2.0 mm X 15 cm
- TWS 84-3039/96 Ø 2.0 mm X 20 cm
- TWS 84-3039/97 Ø 2.5 mm X 15 cm
- TWS 84-3039/98 Ø 2.5 mm X 20 cm
- TWS 84-3039/99 Ø 3.0 mm X 15 cm
- TWS 84-3040/00 Ø 3.0 mm X 20 cm



Power Cannula

- TWS 84-3040/01 Ø 2.0 mm X 15 cm
- TWS 84-3040/02 Ø 2.0 mm X 20 cm
- TWS 84-3040/03 Ø 2.0 mm X 25 cm
- TWS 84-3040/04 Ø 2.0 mm X 30 cm
- TWS 84-3040/05 Ø 2.0 mm X 35 cm
- TWS 84-3040/06 Ø 3.0 mm X 15 cm
- TWS 84-3040/07 Ø 3.0 mm X 20 cm
- TWS 84-3040/08 Ø 3.0 mm X 25 cm
- TWS 84-3040/09 Ø 3.0 mm X 30 cm
- TWS 84-3040/10 Ø 3.0 mm X 35 cm
- TWS 84-3040/11 Ø 4.0 mm X 15 cm
- TWS 84-3040/12 Ø 4.0 mm X 20 cm
- TWS 84-3040/13 Ø 4.0 mm X 25 cm
- TWS 84-3040/14 Ø 4.0 mm X 30 cm
- TWS 84-3040/15 Ø 4.0 mm X 35 cm
- TWS 84-3040/16 Ø 5.0 mm X 15 cm
- TWS 84-3040/17 Ø 5.0 mm X 20 cm
- TWS 84-3040/18 Ø 5.0 mm X 25 cm
- TWS 84-3040/19 Ø 5.0 mm X 30 cm
- TWS 84-3040/20 Ø 5.0 mm X 35 cm
- TWS 84-3040/21 Ø 6.0 mm X 15 cm
- TWS 84-3040/22 Ø 6.0 mm X 20 cm
- TWS 84-3040/23 Ø 6.0 mm X 25 cm
- TWS 84-3040/24 Ø 6.0 mm X 30 cm
- TWS 84-3040/25 Ø 6.0 mm X 35 cm

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com

GENERAL SUCTION CANNULAS



StevensStyleCannula



Tommey Hub

- TWS 84-3040/26 Ø 2.0 mm X 15 cm
- TWS 84-3040/27 Ø 2.5 mm X 20 cm
- TWS 84-3040/28 Ø 3.0 mm X 20 cm
- TWS 84-3040/29 Ø 3.0 mm X 25 cm



Treaded Hub

- TWS 84-3040/30 Ø 3.0mm X 20 cm
- TWS 84-3040/31 Ø 3.0mm X 30 cm
- TWS 84-3040/32 Ø 4.0mm X 30 cm
- TWS 84-3040/33 Ø 5.0mm X 30 cm
- TWS 84-3040/34 Ø 6.0mm X 30 cm



Ultra Luer Lock

- TWS 84-3040/35 Ø 1.0 mm X 5 cm
- TWS 84-3040/36 Ø 1.0 mm X 7 cm
- TWS 84-3040/37 Ø 1.5 mm X 7 cm
- TWS 84-3040/38 Ø 1.5 mm X 9 cm



Luer Lock

- TWS 84-3040/39 Ø 1.0 mm X 4 cm
- TWS 84-3040/40 Ø 1.0 mm X 5 cm
- TWS 84-3040/41 Ø 1.0 mm X 7 cm
- TWS 84-3040/42 Ø 1.5 mm X 7 cm
- TWS 84-3040/43 Ø 1.5 mm X 9 cm
- TWS 84-3040/44 Ø 1.5 mm X 10 cm
- TWS 84-3040/45 Ø 1.5 mm X 12 cm
- TWS 84-3040/46 Ø 2.0 mm X 10 cm
- TWS 84-3040/47 Ø 2.0 mm X 12 cm
- TWS 84-3040/48 Ø 2.0 mm X 15 cm
- TWS 84-3040/49 Ø 2.0 mm X 20 cm
- TWS 84-3040/50 Ø 2.5 mm X 15 cm
- TWS 84-3040/51 Ø 2.5 mm X 20 cm
- TWS 84-3040/52 Ø 3.0 mm X 15 cm
- TWS 84-3040/53 Ø 3.0 mm X 20 cm



Luer lock (Stainless Steel)

- TWS 84-3040/54 Ø 1.0 mm X 4 cm
- TWS 84-3040/55 Ø 1.0 mm X 5 cm
- TWS 84-3040/56 Ø 1.0 mm X 7 cm
- TWS 84-3040/57 Ø 1.5 mm X 7 cm
- TWS 84-3040/58 Ø 1.5 mm X 9 cm
- TWS 84-3040/59 Ø 1.5 mm X 10 cm
- TWS 84-3040/60 Ø 1.5 mm X 12 cm
- TWS 84-3040/61 Ø 2.0 mm X 10 cm
- TWS 84-3040/62 Ø 2.0 mm X 12 cm
- TWS 84-3040/63 Ø 2.0 mm X 15 cm
- TWS 84-3040/64 Ø 2.0 mm X 20 cm
- TWS 84-3040/65 Ø 2.5 mm X 15 cm
- TWS 84-3040/66 Ø 2.5 mm X 20 cm
- TWS 84-3040/67 Ø 3.0 mm X 15 cm
- TWS 84-3040/68 Ø 3.0 mm X 20 cm



Power Cannula

- TWS 84-3040/69 Ø 2.0 mm X 15 cm
- TWS 84-3040/70 Ø 2.0 mm X 20 cm
- TWS 84-3040/71 Ø 2.0 mm X 25 cm
- TWS 84-3040/72 Ø 2.0 mm X 30 cm
- TWS 84-3040/73 Ø 2.0 mm X 35 cm
- TWS 84-3040/74 Ø 3.0 mm X 15 cm
- TWS 84-3040/75 Ø 3.0 mm X 20 cm
- TWS 84-3040/76 Ø 3.0 mm X 25 cm
- TWS 84-3040/77 Ø 3.0 mm X 30 cm
- TWS 84-3040/78 Ø 3.0 mm X 35 cm
- TWS 84-3040/79 Ø 4.0 mm X 15 cm
- TWS 84-3040/80 Ø 4.0 mm X 20 cm
- TWS 84-3040/81 Ø 4.0 mm X 25 cm
- TWS 84-3040/82 Ø 4.0 mm X 30 cm
- TWS 84-3040/83 Ø 4.0 mm X 35 cm
- TWS 84-3040/84 Ø 5.0 mm X 15 cm
- TWS 84-3040/85 Ø 5.0 mm X 20 cm
- TWS 84-3040/86 Ø 5.0 mm X 25 cm
- TWS 84-3040/87 Ø 5.0 mm X 30 cm
- TWS 84-3040/88 Ø 5.0 mm X 35 cm
- TWS 84-3040/89 Ø 6.0 mm X 15 cm
- TWS 84-3040/90 Ø 6.0 mm X 20 cm
- TWS 84-3040/91 Ø 6.0 mm X 25 cm
- TWS 84-3040/92 Ø 6.0 mm X 30 cm
- TWS 84-3040/93 Ø 6.0 mm X 35 cm



TW Surgical

GENERAL SUCTION CANNULAS



Cobra Bibevel Style

More aggressive than the Standard, Pyramid or Mercedes, this cannula has a tapered tip with two ports and works effectively within dense fibrous tissue.



Tomme Hub

- TWS 84-3040/94 Ø 2.5 mm X 15 cm
- TWS 84-3040/95 Ø 2.5 mm X 20 cm
- TWS 84-3040/96 Ø 3.0 mm X 15 cm
- TWS 84-3040/97 Ø 3.0 mm X 20 cm
- TWS 84-3040/98 Ø 3.0 mm X 25 cm
- TWS 84-3040/99 Ø 4.0 mm X 25 cm
- TWS 84-3041/00 Ø 4.0 mm X 30 cm
- TWS 84-3041/01 Ø 4.0 mm X 35 cm
- TWS 84-3041/02 Ø 4.5 mm X 30 cm
- TWS 84-3041/03 Ø 4.5 mm X 35 cm



Treaded Hub

- TWS 84-3041/10 Ø 3.0mm X 15 cm
- TWS 84-3041/11 Ø 3.0mm X 25 cm
- TWS 84-3041/12 Ø 4.0mm X 15 cm
- TWS 84-3041/13 Ø 5.0mm X 25 cm



Ultra Luer Lock

- TWS 84-3041/14 Ø 2.0 mm X 15 cm
- TWS 84-3041/15 Ø 2.5 mm X 15 cm
- TWS 84-3041/16 Ø 2.5 mm X 20 cm
- TWS 84-3041/17 Ø 3.0 mm X 15 cm
- TWS 84-3041/18 Ø 3.0 mm X 20 cm
- TWS 84-3041/19 Ø 3.0 mm X 25 cm



Luer Lock

- TWS 84-3041/20 Ø 2.0 mm X 15 cm
- TWS 84-3041/21 Ø 2.5 mm X 15 cm
- TWS 84-3041/22 Ø 2.5 mm X 20 cm
- TWS 84-3041/23 Ø 3.0 mm X 15 cm
- TWS 84-3041/24 Ø 3.0 mm X 20 cm
- TWS 84-3041/25 Ø 3.0 mm X 25 cm



Luer lock (Stainless Steel)

- TWS 84-3041/04 Ø 2.0 mm X 15 cm
- TWS 84-3041/05 Ø 2.5 mm X 15 cm
- TWS 84-3041/06 Ø 2.5 mm X 20 cm
- TWS 84-3041/07 Ø 3.0 mm X 15 cm
- TWS 84-3041/08 Ø 3.0 mm X 20 cm
- TWS 84-3041/09 Ø 3.0 mm X 25 cm



Power Cannula

- TWS 84-3041/26 Ø 2.0 mm X 15 cm
- TWS 84-3041/27 Ø 2.0 mm X 20 cm
- TWS 84-3041/28 Ø 2.0 mm X 25 cm
- TWS 84-3041/29 Ø 2.0 mm X 30 cm
- TWS 84-3041/30 Ø 2.0 mm X 35 cm
- TWS 84-3041/31 Ø 3.0 mm X 15 cm
- TWS 84-3041/32 Ø 3.0 mm X 20 cm
- TWS 84-3041/33 Ø 3.0 mm X 25 cm
- TWS 84-3041/34 Ø 3.0 mm X 30 cm
- TWS 84-3041/35 Ø 3.0 mm X 35 cm
- TWS 84-3041/36 Ø 4.0 mm X 15 cm
- TWS 84-3041/37 Ø 4.0 mm X 20 cm
- TWS 84-3041/38 Ø 4.0 mm X 25 cm
- TWS 84-3041/39 Ø 4.0 mm X 30 cm
- TWS 84-3041/40 Ø 4.0 mm X 35 cm
- TWS 84-3041/41 Ø 5.0 mm X 15 cm
- TWS 84-3041/42 Ø 5.0 mm X 20 cm
- TWS 84-3041/43 Ø 5.0 mm X 25 cm
- TWS 84-3041/44 Ø 5.0 mm X 30 cm
- TWS 84-3041/45 Ø 5.0 mm X 35 cm
- TWS 84-3041/46 Ø 6.0 mm X 15 cm
- TWS 84-3041/47 Ø 6.0 mm X 20 cm
- TWS 84-3041/48 Ø 6.0 mm X 25 cm
- TWS 84-3041/49 Ø 6.0 mm X 30 cm
- TWS 84-3041/50 Ø 6.0 mm X 35 cm

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical

GENERAL SUCTION CANNULAS



Standard, Single Hole

This is the most common cell friendly cannula design for soft tissue removal. The single port is located 180 degrees from the thumb rest on the hub.



Tommey Hub

TWS 84-3041/51	Ø 2.5 mm X 15 cm
TWS 84-3041/52	Ø 2.5 mm X 20 cm
TWS 84-3041/53	Ø 3.0 mm X 15 cm
TWS 84-3041/54	Ø 3.0 mm X 20 cm
TWS 84-3041/55	Ø 3.0 mm X 25 cm
TWS 84-3041/56	Ø 4.0 mm X 25 cm
TWS 84-3041/57	Ø 4.0 mm X 30 cm
TWS 84-3041/58	Ø 4.0 mm X 35 cm
TWS 84-3041/59	Ø 4.5 mm X 30 cm
TWS 84-3041/60	Ø 4.5 mm X 35 cm



Treaded Hub

SI.2002-31	Ø 3.0mm X 15 cm
SI.2002-32	Ø 3.0mm X 25 cm
SI.2002-33	Ø 4.0mm X 15 cm
SI.2002-34	Ø 5.0mm X 25 cm



Ultra Luer Lock

SI.2003-31	Ø 1.5 mm X 12 cm
SI.2003-32	Ø 2.0 mm X 15 cm
SI.2003-33	Ø 2.0 mm X 20 cm
SI.2003-34	Ø 2.5 mm X 15 cm
SI.2003-35	Ø 2.5 mm X 20 cm
SI.2003-36	Ø 3.0 mm X 15 cm
SI.2003-37	Ø 3.0 mm X 20 cm
SI.2003-38	Ø 3.0 mm X 25 cm



Luer Lock

SI.2004-41	Ø 1.5 mm X 12 cm
SI.2004-42	Ø 2.0 mm X 15 cm
SI.2004-43	Ø 2.0 mm X 20 cm
SI.2004-44	Ø 2.5 mm X 15 cm
SI.2004-45	Ø 2.5 mm X 20 cm
SI.2004-46	Ø 3.0 mm X 15 cm
SI.2004-47	Ø 3.0 mm X 20 cm
	Ø 3.0 mm X 25 cm



Luer lock (Stainless Steel)

TWS 84-3041/61	Ø 1.5 mm X 15 cm
TWS 84-3041/62	Ø 2.0 mm X 15 cm
TWS 84-3041/63	Ø 2.0 mm X 20 cm
TWS 84-3041/64	Ø 2.5 mm X 15 cm
TWS 84-3041/65	Ø 2.5 mm X 20 cm
TWS 84-3041/66	Ø 3.0 mm X 25 cm
TWS 84-3041/67	Ø 3.0 mm X 25 cm
TWS 84-3041/68	Ø 3.0 mm X 25 cm



Power Cannula

TWS 84-3041/69	Ø 2.0 mm X 15 cm
TWS 84-3041/70	Ø 2.0 mm X 20 cm
TWS 84-3041/71	Ø 2.0 mm X 25 cm
TWS 84-3041/72	Ø 2.0 mm X 30 cm
TWS 84-3041/73	Ø 2.0 mm X 35 cm
TWS 84-3041/74	Ø 3.0 mm X 15 cm
TWS 84-3041/75	Ø 3.0 mm X 20 cm
TWS 84-3041/76	Ø 3.0 mm X 25 cm
TWS 84-3041/77	Ø 3.0 mm X 30 cm
TWS 84-3041/78	Ø 3.0 mm X 35 cm
TWS 84-3041/79	Ø 4.0 mm X 15 cm
TWS 84-3041/80	Ø 4.0 mm X 20 cm
TWS 84-3041/81	Ø 4.0 mm X 25 cm
TWS 84-3041/82	Ø 4.0 mm X 30 cm
TWS 84-3041/83	Ø 4.0 mm X 35 cm
TWS 84-3041/84	Ø 5.0 mm X 15 cm
TWS 84-3041/85	Ø 5.0 mm X 20 cm
TWS 84-3041/86	Ø 5.0 mm X 25 cm
TWS 84-3041/87	Ø 5.0 mm X 30 cm
TWS 84-3041/88	Ø 5.0 mm X 35 cm
TWS 84-3041/89	Ø 6.0 mm X 15 cm
TWS 84-3041/90	Ø 6.0 mm X 20 cm
TWS 84-3041/91	Ø 6.0 mm X 25 cm
TWS 84-3041/92	Ø 6.0 mm X 30 cm
TWS 84-3041/93	Ø 6.0 mm X 35 cm

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical

GENERAL SUCTION CANNULAS



Mercedese

This circumferential three port configuration cell friendly liposuction cannula is effective for multi-plane soft tissue harvesting.



Tommey Hub

- TWS 84-3041/94 Ø 2.5 mm X 15 cm
- TWS 84-3041/94 Ø 2.5 mm X 20 cm
- TWS 84-3041/96 Ø 3.0 mm X 15 cm
- TWS 84-3041/97 Ø 3.0 mm X 20 cm
- TWS 84-3041/98 Ø 3.0 mm X 25 cm
- TWS 84-3041/99 Ø 4.0 mm X 25 cm
- TWS 84-3042/00 Ø 4.0 mm X 30 cm
- TWS 84-3042/01 Ø 4.0 mm X 35 cm
- TWS 84-3042/02 Ø 4.5 mm X 30 cm
- TWS 84-3042/03 Ø 4.5 mm X 35 cm



Treaded Hub

- TWS 84-3042/37 Ø 3.0mm X 15 cm
- TWS 84-3042/38 Ø 3.0mm X 25 cm
- TWS 84-3042/39 Ø 4.0mm X 15 cm
- TWS 84-3042/40 Ø 4.0mm X 25 cm



Ultra Luer Lock

- TWS 84-3042/41 Ø 1.5 mm X 12 cm
- TWS 84-3042/42 Ø 2.0 mm X 15 cm
- TWS 84-3042/43 Ø 2.5 mm X 15 cm
- TWS 84-3042/44 Ø 2.5 mm X 20 cm
- TWS 84-3042/45 Ø 3.0 mm X 15 cm
- TWS 84-3042/46 Ø 3.0 mm X 20 cm
- TWS 84-3042/47 Ø 3.0 mm X 25 cm



Luer Lock

- TWS 84-3042/48 Ø 1.5 mm X 12 cm
- TWS 84-3042/49 Ø 2.0 mm X 15 cm
- TWS 84-3042/50 Ø 2.5 mm X 15 cm
- TWS 84-3042/51 Ø 2.5 mm X 20 cm
- TWS 84-3042/52 Ø 3.0 mm X 15 cm
- TWS 84-3042/53 Ø 3.0 mm X 20 cm
- TWS 84-3042/54 Ø 3.0 mm X 25 cm



Luer lock (Stainless Steel)

- TWS 84-3042/04 Ø 1.5 mm X 12 cm
- TWS 84-3042/05 Ø 2.0 mm X 15 cm
- TWS 84-3042/06 Ø 2.0 mm X 20 cm
- TWS 84-3042/07 Ø 2.5 mm X 15 cm
- TWS 84-3042/08 Ø 2.5 mm X 20 cm
- TWS 84-3042/09 Ø 3.0 mm X 15 cm
- TWS 84-3042/10 Ø 3.0 mm X 20 cm
- TWS 84-3042/11 Ø 3.0 mm X 25 cm



Power Cannula

- TWS 84-3042/12 Ø 2.0 mm X 15 cm
- TWS 84-3042/13 Ø 2.0 mm X 20 cm
- TWS 84-3042/14 Ø 2.0 mm X 25 cm
- TWS 84-3042/15 Ø 2.0 mm X 30 cm
- TWS 84-3042/16 Ø 2.0 mm X 35 cm
- TWS 84-3042/17 Ø 3.0 mm X 15 cm
- TWS 84-3042/18 Ø 3.0 mm X 20 cm
- TWS 84-3042/19 Ø 3.0 mm X 25 cm
- TWS 84-3042/20 Ø 3.0 mm X 30 cm
- TWS 84-3042/21 Ø 3.0 mm X 35 cm
- TWS 84-3042/22 Ø 4.0 mm X 15 cm
- TWS 84-3042/23 Ø 4.0 mm X 20 cm
- TWS 84-3042/24 Ø 4.0 mm X 25 cm
- TWS 84-3042/25 Ø 4.0 mm X 30 cm
- TWS 84-3042/26 Ø 4.0 mm X 35 cm
- TWS 84-3042/27 Ø 5.0 mm X 15 cm
- TWS 84-3042/28 Ø 5.0 mm X 20 cm
- TWS 84-3042/29 Ø 5.0 mm X 25 cm
- TWS 84-3042/30 Ø 5.0 mm X 30 cm
- TWS 84-3042/31 Ø 5.0 mm X 35 cm
- TWS 84-3042/32 Ø 6.0 mm X 15 cm
- TWS 84-3042/33 Ø 6.0 mm X 20 cm
- TWS 84-3042/34 Ø 6.0 mm X 25 cm
- TWS 84-3042/35 Ø 6.0 mm X 30 cm
- TWS 84-3042/36 Ø 6.0 mm X 35 cm

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical

GENERAL SUCTION CANNULAS



Pyramid Holes Cannula

This commonly used design has three ports located 180 degrees from the thumb print. It is a slightly more aggressive cannula than the Standard during liposuction procedures



Tommey Hub

- TWS 84-3042/55 Ø 2.5 mm X 15 cm
- TWS 84-3042/56 Ø 2.5 mm X 20 cm
- TWS 84-3042/57 Ø 3.0 mm X 15 cm
- TWS 84-3042/58 Ø 3.0 mm X 20 cm
- TWS 84-3042/59 Ø 3.0 mm X 25 cm
- TWS 84-3042/60 Ø 4.0 mm X 25 cm
- TWS 84-3042/61 Ø 4.0 mm X 30 cm
- TWS 84-3042/62 Ø 4.0 mm X 35 cm
- TWS 84-3042/63 Ø 4.5 mm X 30 cm
- TWS 84-3042/64 Ø 4.5 mm X 35 cm



Treaded Hub

- TWS 84-3042/73 Ø 3.0mm X 15 cm
- TWS 84-3042/74 Ø 3.0mm X 25 cm
- TWS 84-3042/75 Ø 4.0mm X 15 cm
- TWS 84-3042/76 Ø 4.0mm X 25 cm



Ultra Luer Lock

- TWS 84-3042/77 Ø 1.5 mm X 12 cm
- TWS 84-3042/78 Ø 2.0 mm X 15 cm
- TWS 84-3042/79 Ø 2.5 mm X 15 cm
- TWS 84-3042/80 Ø 2.5 mm X 20 cm
- TWS 84-3042/81 Ø 3.0 mm X 15 cm
- TWS 84-3042/82 Ø 3.0 mm X 20 cm
- TWS 84-3042/83 Ø 3.0 mm X 25 cm



Luer Lock

- TWS 84-3042/84 Ø 1.5 mm X 12 cm
- TWS 84-3042/85 Ø 2.0 mm X 15 cm
- TWS 84-3042/86 Ø 2.5 mm X 15 cm
- TWS 84-3042/87 Ø 2.5 mm X 20 cm
- TWS 84-3042/88 Ø 3.0 mm X 15 cm
- TWS 84-3042/89 Ø 3.0 mm X 20 cm
- TWS 84-3042/90 Ø 3.0 mm X 25 cm



Luer lock (Stainless Steel)

- TWS 84-3042/65 Ø 1.5 mm X 12 cm
- TWS 84-3042/66 Ø 2.0 mm X 15 cm
- TWS 84-3042/67 Ø 2.0 mm X 20 cm
- TWS 84-3042/68 Ø 2.5 mm X 15 cm
- TWS 84-3042/69 Ø 2.5 mm X 20 cm
- TWS 84-3042/70 Ø 3.0 mm X 15 cm
- TWS 84-3042/71 Ø 3.0 mm X 20 cm
- TWS 84-3042/72 Ø 3.0 mm X 25 cm



Power Cannula

- TWS 84-3042/91 Ø 2.0 mm X 15 cm
- TWS 84-3042/92 Ø 2.0 mm X 20 cm
- TWS 84-3042/93 Ø 2.0 mm X 25 cm
- TWS 84-3042/94 Ø 2.0 mm X 30 cm
- TWS 84-3042/95 Ø 2.0 mm X 35 cm
- TWS 84-3042/96 Ø 3.0 mm X 15 cm
- TWS 84-3042/97 Ø 3.0 mm X 20 cm
- TWS 84-3042/98 Ø 3.0 mm X 25 cm
- TWS 84-3042/99 Ø 3.0 mm X 30 cm
- TWS 84-3043/00 Ø 3.0 mm X 35 cm
- TWS 84-3043/01 Ø 4.0 mm X 15 cm
- TWS 84-3043/02 Ø 4.0 mm X 20 cm
- TWS 84-3043/03 Ø 4.0 mm X 25 cm
- TWS 84-3043/04 Ø 4.0 mm X 30 cm
- TWS 84-3043/05 Ø 4.0 mm X 35 cm
- TWS 84-3043/06 Ø 5.0 mm X 15 cm
- TWS 84-3043/07 Ø 5.0 mm X 20 cm
- TWS 84-3043/08 Ø 5.0 mm X 25 cm
- TWS 84-3043/09 Ø 5.0 mm X 30 cm
- TWS 84-3043/10 Ø 5.0 mm X 35 cm
- TWS 84-3043/11 Ø 6.0 mm X 15 cm
- TWS 84-3043/12 Ø 6.0 mm X 20 cm
- TWS 84-3043/13 Ø 6.0 mm X 25 cm
- TWS 84-3043/14 Ø 6.0 mm X 30 cm
- TWS 84-3043/15 Ø 6.0 mm X 35 cm

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical

GENERAL SUCTION CANNULAS



3 Spiral Port Cannula (Carraway Style)

Spiral three port design on this Cell Friendly liposuction cannula is effective for multi-plane adipose tissue harvesting.



Tommey Hub

- TWS 84-3043/16 Ø 2.5 mm X 15 cm
- TWS 84-3043/17 Ø 2.5 mm X 20 cm
- TWS 84-3043/18 Ø 3.0 mm X 15 cm
- TWS 84-3043/19 Ø 3.0 mm X 20 cm
- TWS 84-3043/20 Ø 3.0 mm X 25 cm
- TWS 84-3043/21 Ø 4.0 mm X 25 cm
- TWS 84-3043/22 Ø 4.0 mm X 30 cm
- TWS 84-3043/23 Ø 4.0 mm X 35 cm
- TWS 84-3043/24 Ø 4.5 mm X 30 cm
- TWS 84-3043/25 Ø 4.5 mm X 35 cm



Treaded Hub

- SI.2002-61 Ø 3.0mm X 15 cm
- SI.2002-62 Ø 3.0mm X 25 cm
- SI.2002-63 Ø 4.0mm X 15 cm
- SI.2002-64 Ø 4.0mm X 25 cm



Ultra Luer Lock

- SI.2003-61 Ø 1.5 mm X 12 cm
- SI.2003-62 Ø 2.0 mm X 15 cm
- SI.2003-63 Ø 2.5 mm X 15 cm
- SI.2003-64 Ø 2.5 mm X 20 cm
- SI.2003-65 Ø 3.0 mm X 15 cm
- SI.2003-66 Ø 3.0 mm X 20 cm
- SI.2003-67 Ø 3.0 mm X 25 cm



Luer Lock

- TWS 84-3043/59 Ø 1.5 mm X 12 cm
- TWS 84-3043/60 Ø 2.0 mm X 15 cm
- TWS 84-3043/61 Ø 2.5 mm X 15 cm
- TWS 84-3043/62 Ø 2.5 mm X 20 cm
- TWS 84-3043/63 Ø 3.0 mm X 15 cm
- TWS 84-3043/64 Ø 3.0 mm X 20 cm
- TWS 84-3043/65 Ø 3.0 mm X 25 cm



Luer lock (Stainless Steel)

- TWS 84-3043/26 Ø 1.5 mm X 12 cm
- TWS 84-3043/27 Ø 2.0 mm X 15 cm
- TWS 84-3043/28 Ø 2.0 mm X 20 cm
- TWS 84-3043/29 Ø 2.5 mm X 15 cm
- TWS 84-3043/30 Ø 2.5 mm X 20 cm
- TWS 84-3043/31 Ø 3.0 mm X 15 cm
- TWS 84-3043/32 Ø 3.0 mm X 20 cm
- TWS 84-3043/33 Ø 3.0 mm X 25 cm



Power Cannula

- TWS 84-3043/34 Ø 2.0 mm X 15 cm
- TWS 84-3043/35 Ø 2.0 mm X 20 cm
- TWS 84-3043/36 Ø 2.0 mm X 25 cm
- TWS 84-3043/37 Ø 2.0 mm X 30 cm
- TWS 84-3043/38 Ø 2.0 mm X 35 cm
- TWS 84-3043/39 Ø 3.0 mm X 15 cm
- TWS 84-3043/40 Ø 3.0 mm X 20 cm
- TWS 84-3043/41 Ø 3.0 mm X 25 cm
- TWS 84-3043/42 Ø 3.0 mm X 30 cm
- TWS 84-3043/43 Ø 3.0 mm X 35 cm
- TWS 84-3043/44 Ø 4.0 mm X 15 cm
- TWS 84-3043/45 Ø 4.0 mm X 20 cm
- TWS 84-3043/46 Ø 4.0 mm X 25 cm
- TWS 84-3043/47 Ø 4.0 mm X 30 cm
- TWS 84-3043/48 Ø 4.0 mm X 35 cm
- TWS 84-3043/49 Ø 5.0 mm X 15 cm
- TWS 84-3043/50 Ø 5.0 mm X 20 cm
- TWS 84-3043/51 Ø 5.0 mm X 25 cm
- TWS 84-3043/52 Ø 5.0 mm X 30 cm
- TWS 84-3043/53 Ø 5.0 mm X 35 cm
- TWS 84-3043/54 Ø 6.0 mm X 15 cm
- TWS 84-3043/55 Ø 6.0 mm X 20 cm
- TWS 84-3043/56 Ø 6.0 mm X 25 cm
- TWS 84-3043/57 Ø 6.0 mm X 30 cm
- TWS 84-3043/58 Ø 6.0 mm X 35 cm

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical

GENERAL SUCTION CANNULAS



Six Spiral Port (Adjust below 3 Spiral Ports)

Six port spiral configuration on this Cell Friendly liposuction cannula harvests maximum tissue volume in minimal time.



Smart Hub

- TWS 84-3043/66 Ø 2.5 mm X 15 cm
- TWS 84-3043/67 Ø 2.5 mm X 20 cm
- TWS 84-3043/68 Ø 3.0 mm X 15 cm
- TWS 84-3043/69 Ø 3.0 mm X 25 cm
- TWS 84-3043/70 Ø 4.0 mm X 25 cm



Luer lock (Stainless Steel)

- TWS 84-3043/85 Ø 2.5 mm X 15 cm
- TWS 84-3043/86 Ø 3.0 mm X 15 cm
- TWS 84-3043/87 Ø 3.0 mm X 25 cm
- TWS 84-3043/88 Ø 4.0 mm X 15 cm
- TWS 84-3043/89 Ø 4.0 mm X 25 cm



Treaded Hub

- TWS 84-3043/71 Ø 3.0mm X 15 cm
- TWS 84-3043/72 Ø 3.0mm X 25 cm
- TWS 84-3043/73 Ø 4.0mm X 15 cm
- TWS 84-3043/74 Ø 4.0mm X 25 cm



Power Cannula

- TWS 84-3043/90 Ø 2.0 mm X 15 cm
- TWS 84-3043/91 Ø 2.0 mm X 20 cm
- TWS 84-3043/92 Ø 2.0 mm X 25 cm
- TWS 84-3043/93 Ø 2.0 mm X 30 cm
- TWS 84-3043/94 Ø 2.0 mm X 35 cm
- TWS 84-3043/95 Ø 3.0 mm X 15 cm
- TWS 84-3043/96 Ø 3.0 mm X 20 cm
- TWS 84-3043/97 Ø 3.0 mm X 25 cm
- TWS 84-3043/98 Ø 3.0 mm X 30 cm
- TWS 84-3043/99 Ø 3.0 mm X 35 cm
- TWS 84-3044/00 Ø 4.0 mm X 15 cm
- TWS 84-3044/01 Ø 4.0 mm X 20 cm
- TWS 84-3044/02 Ø 4.0 mm X 25 cm
- TWS 84-3044/03 Ø 4.0 mm X 30 cm
- TWS 84-3044/04 Ø 4.0 mm X 35 cm
- TWS 84-3044/05 Ø 5.0 mm X 15 cm
- TWS 84-3044/06 Ø 5.0 mm X 20 cm
- TWS 84-3044/07 Ø 5.0 mm X 25 cm
- TWS 84-3044/08 Ø 5.0 mm X 30 cm
- TWS 84-3044/09 Ø 5.0 mm X 35 cm
- TWS 84-3044/10 Ø 6.0 mm X 15 cm
- TWS 84-3044/11 Ø 6.0 mm X 20 cm
- TWS 84-3044/12 Ø 6.0 mm X 25 cm
- TWS 84-3044/13 Ø 6.0 mm X 30 cm
- TWS 84-3044/14 Ø 6.0 mm X 35 cm



Ultra Luer Lock

- TWS 84-3043/75 Ø 2.5 mm X 15 cm
- TWS 84-3043/76 Ø 3.0 mm X 15 cm
- TWS 84-3043/77 Ø 3.0 mm X 25 cm
- TWS 84-3043/78 Ø 4.0 mm X 15 cm
- TWS 84-3043/79 Ø 4.0 mm X 25 cm



Luer Lock

- TWS 84-3043/80 Ø 2.5 mm X 15 cm
- TWS 84-3043/81 Ø 3.0 mm X 15 cm
- TWS 84-3043/82 Ø 3.0 mm X 25 cm
- TWS 84-3043/83 Ø 4.0 mm X 15 cm
- TWS 84-3043/84 Ø 4.0 mm X 25 cm

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical

GENERAL SUCTION CANNULAS



Three Port Cannula

This body sculpting instrument is excellent for harvesting small aliquots of tissue during liposuction procedures.



Tommey Hub

TWS 84-3044/15	Ø 2.5 mm X 15 cm
TWS 84-3044/16	Ø 2.5 mm X 20 cm
TWS 84-3044/17	Ø 3.0 mm X 15 cm
TWS 84-3044/18	Ø 3.0 mm X 20 cm
TWS 84-3044/19	Ø 3.0 mm X 25 cm
TWS 84-3044/20	Ø 4.0 mm X 25 cm
TWS 84-3044/21	Ø 4.0 mm X 30 cm
TWS 84-3044/22	Ø 4.0 mm X 35 cm
TWS 84-3044/23	Ø 4.5 mm X 30 cm
TWS 84-3044/24	Ø 4.5 mm X 35 cm



Treaded Hub

TWS 84-3044/33	Ø 3.0mm X 15 cm
TWS 84-3044/34	Ø 3.0mm X 25 cm
TWS 84-3044/35	Ø 4.0mm X 15 cm
TWS 84-3044/36	Ø 4.0mm X 25 cm



Ultra Luer Lock

TWS 84-3044/37	Ø 1.5 mm X 12 cm
TWS 84-3044/38	Ø 2.0 mm X 15 cm
TWS 84-3044/39	Ø 2.5 mm X 15 cm
TWS 84-3044/40	Ø 2.5 mm X 20 cm
TWS 84-3044/41	Ø 3.0 mm X 15 cm
TWS 84-3044/42	Ø 3.0 mm X 20 cm
TWS 84-3044/43	Ø 3.0 mm X 25 cm



Luer Lock

TWS 84-3044/44	Ø 1.5 mm X 12 cm
TWS 84-3044/45	Ø 2.0 mm X 15 cm
TWS 84-3044/46	Ø 2.5 mm X 15 cm
TWS 84-3044/47	Ø 2.5 mm X 20 cm
TWS 84-3044/48	Ø 3.0 mm X 15 cm
TWS 84-3044/49	Ø 3.0 mm X 20 cm
TWS 84-3044/50	Ø 3.0 mm X 25 cm



Luer lock (Stainless Steel)

TWS 84-3044/25	Ø 1.5 mm X 12 cm
TWS 84-3044/26	Ø 2.0 mm X 15 cm
TWS 84-3044/27	Ø 2.0 mm X 20 cm
TWS 84-3044/28	Ø 2.5 mm X 15 cm
TWS 84-3044/29	Ø 2.5 mm X 20 cm
TWS 84-3044/30	Ø 3.0 mm X 15 cm
TWS 84-3044/31	Ø 3.0 mm X 20 cm
TWS 84-3044/32	Ø 3.0 mm X 25 cm



Power Cannula

TWS 84-3044/51	Ø 2.0 mm X 15 cm
TWS 84-3044/52	Ø 2.0 mm X 20 cm
TWS 84-3044/53	Ø 2.0 mm X 25 cm
TWS 84-3044/54	Ø 2.0 mm X 30 cm
TWS 84-3044/55	Ø 2.0 mm X 35 cm
TWS 84-3044/56	Ø 3.0 mm X 15 cm
TWS 84-3044/57	Ø 3.0 mm X 20 cm
TWS 84-3044/58	Ø 3.0 mm X 25 cm
TWS 84-3044/59	Ø 3.0 mm X 30 cm
TWS 84-3044/60	Ø 3.0 mm X 35 cm
TWS 84-3044/61	Ø 4.0 mm X 15 cm
TWS 84-3044/62	Ø 4.0 mm X 20 cm
TWS 84-3044/63	Ø 4.0 mm X 25 cm
TWS 84-3044/64	Ø 4.0 mm X 30 cm
TWS 84-3044/65	Ø 4.0 mm X 35 cm
TWS 84-3044/66	Ø 5.0 mm X 15 cm
TWS 84-3044/67	Ø 5.0 mm X 20 cm
TWS 84-3044/68	Ø 5.0 mm X 25 cm
TWS 84-3044/69	Ø 5.0 mm X 30 cm
TWS 84-3044/70	Ø 5.0 mm X 35 cm
TWS 84-3044/71	Ø 6.0 mm X 15 cm
TWS 84-3044/72	Ø 6.0 mm X 20 cm
TWS 84-3044/73	Ø 6.0 mm X 25 cm
TWS 84-3044/74	Ø 6.0 mm X 30 cm
TWS 84-3044/75	Ø 6.0 mm X 35 cm

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical

GENERAL SUCTION CANNULAS



Basket Shaped Cannula

The Cell Friendly Basket Cannula is a versatile liposculpting tool. It can operate as an equalizer without vacuum to treat adipose tissue and skin irregularities in revision liposuction and fat transfer procedures. The Basket Cannula is also adept under vacuum as a fat harvesting cannula. For finishing, this cannula can be used as an injector for equalizing and adding adipose tissue to sculpt smoother contours.



Tommey Hub

- TWS 84-3044/76 Ø 3.0 mm X 15 cm
- TWS 84-3044/77 Ø 3.0 mm X 20 cm
- TWS 84-3044/78 Ø 3.0 mm X 25 cm
- TWS 84-3044/79 Ø 4.0 mm X 20 cm
- TWS 84-3044/80 Ø 4.0 mm X 30 cm
- TWS 84-3044/81 Ø 4.0 mm X 35 cm



Holeontip-Stylell



Power Cannula

- TWS 84-3044/82 Ø 2.0 mm X 15 cm
- TWS 84-3044/83 Ø 2.0 mm X 20 cm
- TWS 84-3044/84 Ø 2.0 mm X 25 cm
- TWS 84-3044/85 Ø 2.0 mm X 30 cm
- TWS 84-3044/86 Ø 2.0 mm X 35 cm
- TWS 84-3044/87 Ø 3.0 mm X 15 cm
- TWS 84-3044/88 Ø 3.0 mm X 20 cm
- TWS 84-3044/89 Ø 3.0 mm X 25 cm
- TWS 84-3044/90 Ø 3.0 mm X 30 cm
- TWS 84-3044/91 Ø 3.0 mm X 35 cm
- TWS 84-3044/92 Ø 4.0 mm X 15 cm
- TWS 84-3044/93 Ø 4.0 mm X 20 cm
- TWS 84-3044/94 Ø 4.0 mm X 25 cm
- TWS 84-3044/95 Ø 4.0 mm X 30 cm
- TWS 84-3044/96 Ø 4.0 mm X 35 cm
- TWS 84-3044/97 Ø 5.0 mm X 15 cm
- TWS 84-3044/98 Ø 5.0 mm X 20 cm
- TWS 84-3044/99 Ø 5.0 mm X 25 cm
- TWS 84-3045/00 Ø 5.0 mm X 30 cm
- TWS 84-3045/01 Ø 5.0 mm X 35 cm
- TWS 84-3045/02 Ø 6.0 mm X 15 cm
- TWS 84-3045/03 Ø 6.0 mm X 20 cm
- TWS 84-3045/04 Ø 6.0 mm X 25 cm
- TWS 84-3045/05 Ø 6.0 mm X 30 cm
- TWS 84-3045/06 Ø 6.0 mm X 35 cm



Hole on Tip-Style I



Power Cannula

- TWS 84-3045/07 Ø 2.0 mm X 15 cm
- TWS 84-3045/08 Ø 2.0 mm X 20 cm
- TWS 84-3045/09 Ø 2.0 mm X 25 cm
- TWS 84-3045/10 Ø 2.0 mm X 30 cm
- TWS 84-3045/11 Ø 2.0 mm X 35 cm
- TWS 84-3045/12 Ø 3.0 mm X 15 cm
- TWS 84-3045/13 Ø 3.0 mm X 20 cm
- TWS 84-3045/14 Ø 3.0 mm X 25 cm
- TWS 84-3045/15 Ø 3.0 mm X 30 cm
- TWS 84-3045/16 Ø 3.0 mm X 35 cm
- TWS 84-3045/17 Ø 4.0 mm X 15 cm
- TWS 84-3045/18 Ø 4.0 mm X 20 cm
- TWS 84-3045/19 Ø 4.0 mm X 25 cm
- TWS 84-3045/20 Ø 4.0 mm X 30 cm
- TWS 84-3045/21 Ø 4.0 mm X 35 cm
- TWS 84-3045/22 Ø 5.0 mm X 15 cm
- TWS 84-3045/23 Ø 5.0 mm X 20 cm
- TWS 84-3045/24 Ø 5.0 mm X 25 cm
- TWS 84-3045/25 Ø 5.0 mm X 30 cm
- TWS 84-3045/26 Ø 5.0 mm X 35 cm
- TWS 84-3045/27 Ø 6.0 mm X 15 cm
- TWS 84-3045/28 Ø 6.0 mm X 20 cm
- TWS 84-3045/29 Ø 6.0 mm X 25 cm
- TWS 84-3045/30 Ø 6.0 mm X 30 cm
- TWS 84-3045/31 Ø 6.0 mm X 35 cm

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical

GENERAL SUCTION CANNULAS



Deep Sharp Cut Dissector



Power Cannula

TWS 84-3045/32 Ø 2.0 mm X 15 cm
 TWS 84-3045/33 Ø 2.0 mm X 20 cm
 TWS 84-3045/34 Ø 2.0 mm X 25 cm
 TWS 84-3045/35 Ø 2.0 mm X 30 cm
 TWS 84-3045/36 Ø 2.0 mm X 35 cm
 TWS 84-3045/37 Ø 3.0 mm X 15 cm
 TWS 84-3045/38 Ø 3.0 mm X 20 cm
 TWS 84-3045/39 Ø 3.0 mm X 25 cm
 TWS 84-3045/40 Ø 3.0 mm X 30 cm
 TWS 84-3045/41 Ø 3.0 mm X 35 cm
 TWS 84-3045/42 Ø 4.0 mm X 15 cm
 TWS 84-3045/43 Ø 4.0 mm X 20 cm
 TWS 84-3045/44 Ø 4.0 mm X 25 cm
 TWS 84-3045/45 Ø 4.0 mm X 30 cm
 TWS 84-3045/46 Ø 4.0 mm X 35 cm
 TWS 84-3045/47 Ø 5.0 mm X 15 cm
 TWS 84-3045/48 Ø 5.0 mm X 20 cm
 TWS 84-3045/49 Ø 5.0 mm X 25 cm
 TWS 84-3045/50 Ø 5.0 mm X 30 cm
 TWS 84-3045/51 Ø 5.0 mm X 35 cm
 TWS 84-3045/52 Ø 6.0 mm X 15 cm
 TWS 84-3045/53 Ø 6.0 mm X 20 cm
 TWS 84-3045/54 Ø 6.0 mm X 25 cm
 TWS 84-3045/55 Ø 6.0 mm X 30 cm
 TWS 84-3045/56 Ø 6.0 mm X 35 cm



Power Cannula

TWS 84-3045/57 Ø 2.0 mm X 15 cm
 TWS 84-3045/58 Ø 2.0 mm X 20 cm
 TWS 84-3045/59 Ø 2.0 mm X 25 cm
 TWS 84-3045/60 Ø 2.0 mm X 30 cm
 TWS 84-3045/61 Ø 2.0 mm X 35 cm
 TWS 84-3045/62 Ø 3.0 mm X 15 cm
 TWS 84-3045/63 Ø 3.0 mm X 20 cm
 TWS 84-3045/64 Ø 3.0 mm X 25 cm
 TWS 84-3045/65 Ø 3.0 mm X 30 cm
 TWS 84-3045/66 Ø 3.0 mm X 35 cm
 TWS 84-3045/67 Ø 4.0 mm X 15 cm
 TWS 84-3045/68 Ø 4.0 mm X 20 cm
 TWS 84-3045/69 Ø 4.0 mm X 25 cm
 TWS 84-3045/70 Ø 4.0 mm X 30 cm
 TWS 84-3045/71 Ø 4.0 mm X 35 cm
 TWS 84-3045/72 Ø 5.0 mm X 15 cm
 TWS 84-3045/73 Ø 5.0 mm X 20 cm
 TWS 84-3045/74 Ø 5.0 mm X 25 cm
 TWS 84-3045/75 Ø 5.0 mm X 30 cm
 TWS 84-3045/76 Ø 5.0 mm X 35 cm
 TWS 84-3045/77 Ø 6.0 mm X 15 cm
 TWS 84-3045/78 Ø 6.0 mm X 20 cm
 TWS 84-3045/79 Ø 6.0 mm X 25 cm
 TWS 84-3045/80 Ø 6.0 mm X 30 cm
 TWS 84-3045/81 Ø 6.0 mm X 35 cm



Hole on Tip-Style I

Luer lock (Stainless Steel)

TWS 84-3045/82 Ø 2.0 mm X 10 cm
 TWS 84-3045/83 Ø 2.5 mm X 10 cm
 TWS 84-3045/84 Ø 3.0 mm X 10 cm
 TWS 84-3045/85 Ø 3.5 mm X 10 cm
 TWS 84-3045/86 Ø 4.0 mm X 10 cm
 TWS 84-3045/87 Ø 4.0 mm X 10 cm
 TWS 84-3045/88 Ø 2.0 mm X 15 cm
 TWS 84-3045/89 Ø 2.5 mm X 15 cm
 TWS 84-3045/90 Ø 3.0 mm X 15 cm
 TWS 84-3045/91 Ø 3.5 mm X 15 cm
 TWS 84-3045/92 Ø 4.0 mm X 15 cm
 TWS 84-3045/93 Ø 2.0 mm X 20 cm
 TWS 84-3045/94 Ø 2.5 mm X 20 cm
 TWS 84-3045/95 Ø 3.0 mm X 20 cm
 TWS 84-3045/96 Ø 3.5 mm X 20 cm
 TWS 84-3045/97 Ø 4.0 mm X 20 cm
 TWS 84-3045/98 Ø 2.0 mm X 25 cm
 TWS 84-3045/99 Ø 2.5 mm X 25 cm
 TWS 84-3046/00 Ø 3.0 mm X 25 cm
 TWS 84-3046/01 Ø 3.5 mm X 25 cm
 TWS 84-3046/02 Ø 4.0 mm X 25 cm

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical

GENERAL SUCTION CANNULAS



Tapered Tip Single Central Hole



Power Cannula

- TWS 84-3046/03 Ø 2.0 mm X 15 cm
- TWS 84-3046/04 Ø 2.0 mm X 20 cm
- TWS 84-3046/05 Ø 2.0 mm X 25 cm
- TWS 84-3046/06 Ø 2.0 mm X 30 cm
- TWS 84-3046/07 Ø 2.0 mm X 35 cm
- TWS 84-3046/08 Ø 3.0 mm X 15 cm
- TWS 84-3046/09 Ø 3.0 mm X 20 cm
- TWS 84-3046/10 Ø 3.0 mm X 25 cm
- TWS 84-3046/11 Ø 3.0 mm X 30 cm
- TWS 84-3046/12 Ø 3.0 mm X 35 cm
- TWS 84-3046/13 Ø 4.0 mm X 15 cm
- TWS 84-3046/14 Ø 4.0 mm X 20 cm
- TWS 84-3046/15 Ø 4.0 mm X 25 cm
- TWS 84-3046/16 Ø 4.0 mm X 30 cm
- TWS 84-3046/17 Ø 4.0 mm X 35 cm
- TWS 84-3046/18 Ø 5.0 mm X 15 cm
- TWS 84-3046/19 Ø 5.0 mm X 20 cm
- TWS 84-3046/20 Ø 5.0 mm X 25 cm
- TWS 84-3046/21 Ø 5.0 mm X 30 cm
- TWS 84-3046/22 Ø 5.0 mm X 35 cm
- TWS 84-3046/23 Ø 6.0 mm X 15 cm
- TWS 84-3046/24 Ø 6.0 mm X 20 cm
- TWS 84-3046/25 Ø 6.0 mm X 25 cm
- TWS 84-3046/26 Ø 6.0 mm X 30 cm
- TWS 84-3046/27 Ø 6.0 mm X 35 cm

Fig. 3F



Power Cannula

- TWS 84-3046/53 Ø 2.0 mm X 15 cm
- TWS 84-3046/54 Ø 2.0 mm X 20 cm
- TWS 84-3046/55 Ø 2.0 mm X 25 cm
- TWS 84-3046/56 Ø 2.0 mm X 30 cm
- TWS 84-3046/57 Ø 2.0 mm X 35 cm
- TWS 84-3046/58 Ø 3.0 mm X 15 cm
- TWS 84-3046/59 Ø 3.0 mm X 20 cm
- TWS 84-3046/60 Ø 3.0 mm X 25 cm
- TWS 84-3046/61 Ø 3.0 mm X 30 cm
- TWS 84-3046/62 Ø 3.0 mm X 35 cm
- TWS 84-3046/63 Ø 4.0 mm X 15 cm
- TWS 84-3046/64 Ø 4.0 mm X 20 cm
- TWS 84-3046/65 Ø 4.0 mm X 25 cm
- TWS 84-3046/66 Ø 4.0 mm X 30 cm
- TWS 84-3046/67 Ø 4.0 mm X 35 cm
- TWS 84-3046/68 Ø 5.0 mm X 15 cm
- TWS 84-3046/69 Ø 5.0 mm X 20 cm
- TWS 84-3046/70 Ø 5.0 mm X 25 cm
- TWS 84-3046/71 Ø 5.0 mm X 30 cm
- TWS 84-3046/72 Ø 5.0 mm X 35 cm
- TWS 84-3046/73 Ø 6.0 mm X 15 cm
- TWS 84-3046/74 Ø 6.0 mm X 20 cm
- TWS 84-3046/75 Ø 6.0 mm X 25 cm
- TWS 84-3046/76 Ø 6.0 mm X 30 cm
- TWS 84-3046/77 Ø 6.0 mm X 35 cm

Tapered Tip Two Spiral Holes



Power Cannula

- TWS 84-3046/28 Ø 2.0 mm X 15 cm
- TWS 84-3046/29 Ø 2.0 mm X 20 cm
- TWS 84-3046/30 Ø 2.0 mm X 25 cm
- TWS 84-3046/31 Ø 2.0 mm X 30 cm
- TWS 84-3046/32 Ø 2.0 mm X 35 cm
- TWS 84-3046/33 Ø 3.0 mm X 15 cm
- TWS 84-3046/34 Ø 3.0 mm X 20 cm
- TWS 84-3046/35 Ø 3.0 mm X 25 cm
- TWS 84-3046/36 Ø 3.0 mm X 30 cm
- TWS 84-3046/37 Ø 3.0 mm X 35 cm
- TWS 84-3046/38 Ø 4.0 mm X 15 cm
- TWS 84-3046/39 Ø 4.0 mm X 20 cm
- TWS 84-3046/40 Ø 4.0 mm X 25 cm
- TWS 84-3046/41 Ø 4.0 mm X 30 cm
- TWS 84-3046/42 Ø 4.0 mm X 35 cm
- TWS 84-3046/43 Ø 5.0 mm X 15 cm
- TWS 84-3046/44 Ø 5.0 mm X 20 cm
- TWS 84-3046/45 Ø 5.0 mm X 25 cm
- TWS 84-3046/46 Ø 5.0 mm X 30 cm
- TWS 84-3046/47 Ø 5.0 mm X 35 cm
- TWS 84-3046/48 Ø 6.0 mm X 15 cm
- TWS 84-3046/49 Ø 6.0 mm X 20 cm
- TWS 84-3046/50 Ø 6.0 mm X 25 cm
- TWS 84-3046/51 Ø 6.0 mm X 30 cm
- TWS 84-3046/52 Ø 6.0 mm X 35 cm



TW Surgical

GENERAL SUCTION CANNULAS



One Central Port-Two Parallel Ports



Power Cannula

- TWS 84-3046/78 Ø 2.0 mm X 15 cm
- TWS 84-3046/79 Ø 2.0 mm X 20 cm
- TWS 84-3046/80 Ø 2.0 mm X 25 cm
- TWS 84-3046/81 Ø 2.0 mm X 30 cm
- TWS 84-3046/82 Ø 2.0 mm X 35 cm
- TWS 84-3046/83 Ø 3.0 mm X 15 cm
- TWS 84-3046/84 Ø 3.0 mm X 20 cm
- TWS 84-3046/85 Ø 3.0 mm X 25 cm
- TWS 84-3046/86 Ø 3.0 mm X 30 cm
- TWS 84-3046/87 Ø 3.0 mm X 35 cm
- TWS 84-3046/88 Ø 4.0 mm X 15 cm
- TWS 84-3046/89 Ø 4.0 mm X 20 cm
- TWS 84-3046/90 Ø 4.0 mm X 25 cm
- TWS 84-3046/91 Ø 4.0 mm X 30 cm
- TWS 84-3046/92 Ø 4.0 mm X 35 cm
- TWS 84-3046/93 Ø 5.0 mm X 15 cm
- TWS 84-3046/94 Ø 5.0 mm X 20 cm
- TWS 84-3046/95 Ø 5.0 mm X 25 cm
- TWS 84-3046/96 Ø 5.0 mm X 30 cm
- TWS 84-3046/97 Ø 5.0 mm X 35 cm
- TWS 84-3046/98 Ø 6.0 mm X 15 cm
- TWS 84-3046/99 Ø 6.0 mm X 20 cm
- TWS 84-3047/00 Ø 6.0 mm X 25 cm
- TWS 84-3047/01 Ø 6.0 mm X 30 cm
- TWS 84-3047/02 Ø 6.0 mm X 35 cm



Two Parallel Port- One Central Port



Power Cannula

- TWS 84-3047/03 Ø 2.0 mm X 15 cm
- TWS 84-3047/04 Ø 2.0 mm X 20 cm
- TWS 84-3047/05 Ø 2.0 mm X 25 cm
- TWS 84-3047/06 Ø 2.0 mm X 30 cm
- TWS 84-3047/07 Ø 2.0 mm X 35 cm
- TWS 84-3047/08 Ø 3.0 mm X 15 cm
- TWS 84-3047/09 Ø 3.0 mm X 20 cm
- TWS 84-3047/10 Ø 3.0 mm X 25 cm
- TWS 84-3047/11 Ø 3.0 mm X 30 cm
- TWS 84-3047/12 Ø 3.0 mm X 35 cm
- TWS 84-3047/13 Ø 4.0 mm X 15 cm
- TWS 84-3047/14 Ø 4.0 mm X 20 cm
- TWS 84-3047/15 Ø 4.0 mm X 25 cm
- TWS 84-3047/16 Ø 4.0 mm X 30 cm
- TWS 84-3047/17 Ø 4.0 mm X 35 cm
- TWS 84-3047/18 Ø 5.0 mm X 15 cm
- TWS 84-3047/19 Ø 5.0 mm X 20 cm
- TWS 84-3047/20 Ø 5.0 mm X 25 cm
- TWS 84-3047/21 Ø 5.0 mm X 30 cm
- TWS 84-3047/22 Ø 5.0 mm X 35 cm
- TWS 84-3047/23 Ø 6.0 mm X 15 cm
- TWS 84-3047/24 Ø 6.0 mm X 20 cm
- TWS 84-3047/25 Ø 6.0 mm X 25 cm
- TWS 84-3047/26 Ø 6.0 mm X 30 cm
- TWS 84-3047/27 Ø 6.0 mm X 35 cm



TW Surgical

GENERAL SUCTION CANNULAS



Fig. 3F

Luer lock (Stainless Steel)



- TWS 84-3047/28 Ø 2.0 mm X 7 cm
- TWS 84-3047/29 Ø 2.5 mm X 7 cm
- TWS 84-3047/30 Ø 3.0 mm X 7 cm
- TWS 84-3047/31 Ø 3.5 mm X 7 cm
- TWS 84-3047/32 Ø 4.0 mm X 7 cm
- TWS 84-3047/33 Ø 2.0 mm X 10 cm
- TWS 84-3047/34 Ø 2.5 mm X 10 cm
- TWS 84-3047/35 Ø 3.0 mm X 10 cm
- TWS 84-3047/36 Ø 3.5 mm X 10 cm
- TWS 84-3047/37 Ø 4.0 mm X 10 cm
- TWS 84-3047/38 Ø 2.0 mm X 15 cm
- TWS 84-3047/39 Ø 2.5 mm X 15 cm
- TWS 84-3047/40 Ø 3.0 mm X 15 cm
- TWS 84-3047/41 Ø 3.5 mm X 15 cm
- TWS 84-3047/42 Ø 4.0 mm X 15 cm
- TWS 84-3047/43 Ø 2.0 mm X 20 cm
- TWS 84-3047/44 Ø 2.5 mm X 20 cm
- TWS 84-3047/45 Ø 3.0 mm X 20 cm
- TWS 84-3047/46 Ø 3.5 mm X 20 cm
- TWS 84-3047/47 Ø 4.0 mm X 20 cm
- TWS 84-3047/48 Ø 2.0 mm X 25 cm
- TWS 84-3047/49 Ø 2.5 mm X 25 cm
- TWS 84-3047/50 Ø 3.0 mm X 25 cm
- TWS 84-3047/51 Ø 3.5 mm X 25 cm
- TWS 84-3047/52 Ø 4.0 mm X 25 cm



Power Cannula

- TWS 84-3047/83 Ø 2.0 mm X 15 cm
- TWS 84-3047/84 Ø 2.0 mm X 20 cm
- TWS 84-3047/85 Ø 2.0 mm X 25 cm
- TWS 84-3047/86 Ø 2.0 mm X 30 cm
- TWS 84-3047/87 Ø 2.0 mm X 35 cm
- TWS 84-3047/88 Ø 3.0 mm X 15 cm
- TWS 84-3047/89 Ø 3.0 mm X 20 cm
- TWS 84-3047/90 Ø 3.0 mm X 25 cm
- TWS 84-3047/91 Ø 3.0 mm X 30 cm
- TWS 84-3047/92 Ø 3.0 mm X 35 cm
- TWS 84-3047/93 Ø 4.0 mm X 15 cm
- TWS 84-3047/94 Ø 4.0 mm X 20 cm
- TWS 84-3047/95 Ø 4.0 mm X 25 cm
- TWS 84-3047/96 Ø 4.0 mm X 30 cm
- TWS 84-3047/97 Ø 4.0 mm X 35 cm
- TWS 84-3047/98 Ø 5.0 mm X 15 cm
- TWS 84-3047/99 Ø 5.0 mm X 20 cm
- TWS 84-3048/00 Ø 5.0 mm X 25 cm
- TWS 84-3048/01 Ø 5.0 mm X 30 cm
- TWS 84-3048/02 Ø 5.0 mm X 35 cm
- TWS 84-3048/03 Ø 6.0 mm X 15 cm
- TWS 84-3048/04 Ø 6.0 mm X 20 cm
- TWS 84-3048/05 Ø 6.0 mm X 25 cm
- TWS 84-3048/06 Ø 6.0 mm X 30 cm
- TWS 84-3048/07 Ø 6.0 mm X 35 cm



Fig. 3Hr



Power Cannula

- TWS 84-3047/58 Ø 2.0 mm X 15 cm
- TWS 84-3047/59 Ø 2.0 mm X 20 cm
- TWS 84-3047/60 Ø 2.0 mm X 25 cm
- TWS 84-3047/61 Ø 2.0 mm X 30 cm
- TWS 84-3047/62 Ø 2.0 mm X 35 cm
- TWS 84-3047/63 Ø 3.0 mm X 15 cm
- TWS 84-3047/64 Ø 3.0 mm X 20 cm
- TWS 84-3047/65 Ø 3.0 mm X 25 cm
- TWS 84-3047/66 Ø 3.0 mm X 30 cm
- TWS 84-3047/67 Ø 3.0 mm X 35 cm
- TWS 84-3047/68 Ø 4.0 mm X 15 cm
- TWS 84-3047/69 Ø 4.0 mm X 20 cm
- TWS 84-3047/70 Ø 4.0 mm X 25 cm
- TWS 84-3047/71 Ø 4.0 mm X 30 cm
- TWS 84-3047/72 Ø 4.0 mm X 35 cm
- TWS 84-3047/73 Ø 5.0 mm X 15 cm
- TWS 84-3047/74 Ø 5.0 mm X 20 cm
- TWS 84-3047/75 Ø 5.0 mm X 25 cm
- TWS 84-3047/76 Ø 5.0 mm X 30 cm
- TWS 84-3047/77 Ø 5.0 mm X 35 cm
- TWS 84-3047/78 Ø 6.0 mm X 15 cm
- TWS 84-3047/79 Ø 6.0 mm X 20 cm
- TWS 84-3047/80 Ø 6.0 mm X 25 cm
- TWS 84-3047/81 Ø 6.0 mm X 30 cm
- TWS 84-3047/82 Ø 6.0 mm X 35 cm



Fig. Riz



Power Cannula

- TWS 84-3048/08 Ø 2.0 mm X 15 cm
- TWS 84-3048/09 Ø 2.0 mm X 20 cm
- TWS 84-3048/10 Ø 2.0 mm X 25 cm
- TWS 84-3048/11 Ø 2.0 mm X 30 cm
- TWS 84-3048/12 Ø 2.0 mm X 35 cm
- TWS 84-3048/13 Ø 3.0 mm X 15 cm
- TWS 84-3048/14 Ø 3.0 mm X 20 cm
- TWS 84-3048/15 Ø 3.0 mm X 25 cm
- TWS 84-3048/16 Ø 3.0 mm X 30 cm
- TWS 84-3048/17 Ø 3.0 mm X 35 cm
- TWS 84-3048/18 Ø 4.0 mm X 15 cm
- TWS 84-3048/19 Ø 4.0 mm X 20 cm
- TWS 84-3048/20 Ø 4.0 mm X 25 cm
- TWS 84-3048/21 Ø 4.0 mm X 30 cm
- TWS 84-3048/22 Ø 4.0 mm X 35 cm
- TWS 84-3048/23 Ø 5.0 mm X 15 cm
- TWS 84-3048/24 Ø 5.0 mm X 20 cm
- TWS 84-3048/25 Ø 5.0 mm X 25 cm
- TWS 84-3048/26 Ø 5.0 mm X 30 cm
- TWS 84-3048/27 Ø 5.0 mm X 35 cm
- TWS 84-3048/28 Ø 6.0 mm X 15 cm
- TWS 84-3048/29 Ø 6.0 mm X 20 cm
- TWS 84-3048/30 Ø 6.0 mm X 25 cm
- TWS 84-3048/31 Ø 6.0 mm X 30 cm
- TWS 84-3048/32 Ø 6.0 mm X 35 cm



TW Surgical

GENERAL SUCTION CANNULAS



Fig. 5Ai



Power Cannula

- TWS 84-3048/33 Ø 2.0 mm X 15 cm
- TWS 84-3048/34 Ø 2.0 mm X 20 cm
- TWS 84-3048/37 Ø 2.0 mm X 25 cm
- TWS 84-3048/39 Ø 2.0 mm X 30 cm
- TWS 84-3048/40 Ø 2.0 mm X 35 cm
- TWS 84-3048/41 Ø 3.0 mm X 15 cm
- TWS 84-3048/42 Ø 3.0 mm X 20 cm
- TWS 84-3048/43 Ø 3.0 mm X 25 cm
- TWS 84-3048/44 Ø 3.0 mm X 30 cm
- TWS 84-3048/45 Ø 3.0 mm X 35 cm
- TWS 84-3048/46 Ø 4.0 mm X 15 cm
- TWS 84-3048/47 Ø 4.0 mm X 20 cm
- TWS 84-3048/48 Ø 4.0 mm X 25 cm
- TWS 84-3048/49 Ø 4.0 mm X 30 cm
- TWS 84-3048/50 Ø 4.0 mm X 35 cm
- TWS 84-3048/51 Ø 5.0 mm X 15 cm
- TWS 84-3048/52 Ø 5.0 mm X 20 cm
- TWS 84-3048/53 Ø 5.0 mm X 25 cm
- TWS 84-3048/54 Ø 5.0 mm X 30 cm
- TWS 84-3048/55 Ø 5.0 mm X 35 cm
- TWS 84-3048/56 Ø 6.0 mm X 15 cm
- TWS 84-3048/57 Ø 6.0 mm X 20 cm
- TWS 84-3048/58 Ø 6.0 mm X 25 cm
- TWS 84-3048/59 Ø 6.0 mm X 30 cm
- TWS 84-3048/60 Ø 6.0 mm X 35 cm



Fig. Mr



Power Cannula

- TWS 84-3048/35 Ø 2.0 mm X 15 cm
- TWS 84-3048/36 Ø 2.0 mm X 20 cm
- TWS 84-3048/38 Ø 2.0 mm X 25 cm
- TWS 84-3048/61 Ø 2.0 mm X 30 cm
- TWS 84-3048/62 Ø 2.0 mm X 35 cm
- TWS 84-3048/63 Ø 3.0 mm X 15 cm
- TWS 84-3048/64 Ø 3.0 mm X 20 cm
- TWS 84-3048/65 Ø 3.0 mm X 25 cm
- TWS 84-3048/66 Ø 3.0 mm X 30 cm
- TWS 84-3048/67 Ø 3.0 mm X 35 cm
- TWS 84-3048/68 Ø 4.0 mm X 15 cm
- TWS 84-3048/69 Ø 4.0 mm X 20 cm
- TWS 84-3048/70 Ø 4.0 mm X 25 cm
- TWS 84-3048/71 Ø 4.0 mm X 30 cm
- TWS 84-3048/72 Ø 4.0 mm X 35 cm
- TWS 84-3048/73 Ø 5.0 mm X 15 cm
- TWS 84-3048/74 Ø 5.0 mm X 20 cm
- TWS 84-3048/75 Ø 5.0 mm X 25 cm
- TWS 84-3048/76 Ø 5.0 mm X 30 cm
- TWS 84-3048/77 Ø 5.0 mm X 35 cm
- TWS 84-3048/78 Ø 6.0 mm X 15 cm
- TWS 84-3048/79 Ø 6.0 mm X 20 cm
- TWS 84-3048/80 Ø 6.0 mm X 25 cm
- TWS 84-3048/81 Ø 6.0 mm X 30 cm
- TWS 84-3048/82 Ø 6.0 mm X 35 cm



Fig. P3G



Power Cannula

- TWS 84-3048/83 Ø 2.0 mm X 15 cm
- TWS 84-3048/84 Ø 2.0 mm X 20 cm
- TWS 84-3048/83 Ø 2.0 mm X 25 cm
- TWS 84-3048/85 Ø 2.0 mm X 30 cm
- TWS 84-3048/86 Ø 2.0 mm X 35 cm
- TWS 84-3048/87 Ø 3.0 mm X 15 cm
- TWS 84-3048/88 Ø 3.0 mm X 20 cm
- TWS 84-3048/89 Ø 3.0 mm X 25 cm
- TWS 84-3048/90 Ø 3.0 mm X 30 cm
- TWS 84-3048/91 Ø 3.0 mm X 35 cm
- TWS 84-3048/92 Ø 4.0 mm X 15 cm
- TWS 84-3048/93 Ø 4.0 mm X 20 cm
- TWS 84-3048/94 Ø 4.0 mm X 25 cm
- TWS 84-3048/95 Ø 4.0 mm X 30 cm
- TWS 84-3048/96 Ø 4.0 mm X 35 cm
- TWS 84-3048/97 Ø 5.0 mm X 15 cm
- TWS 84-3048/98 Ø 5.0 mm X 20 cm
- TWS 84-3048/99 Ø 5.0 mm X 25 cm
- TWS 84-3049/00 Ø 5.0 mm X 30 cm
- TWS 84-3049/01 Ø 5.0 mm X 35 cm
- TWS 84-3049/02 Ø 6.0 mm X 15 cm
- TWS 84-3049/03 Ø 6.0 mm X 20 cm
- TWS 84-3049/04 Ø 6.0 mm X 25 cm
- TWS 84-3049/05 Ø 6.0 mm X 30 cm
- TWS 84-3049/06 Ø 6.0 mm X 35 cm

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical

GENERAL SUCTION CANNULAS



Fig.Py



Power Cannula

- TWS 84-3049/07 Ø 2.0 mm X 15 cm
- TWS 84-3049/08 Ø 2.0 mm X 20 cm
- TWS 84-3049/09 Ø 2.0 mm X 25 cm
- TWS 84-3049/10 Ø 2.0 mm X 30 cm
- TWS 84-3049/11 Ø 2.0 mm X 35 cm
- TWS 84-3049/12 Ø 3.0 mm X 15 cm
- TWS 84-3049/13 Ø 3.0 mm X 20 cm
- TWS 84-3049/14 Ø 3.0 mm X 25 cm
- TWS 84-3049/15 Ø 3.0 mm X 30 cm
- TWS 84-3049/16 Ø 3.0 mm X 35 cm
- TWS 84-3049/17 Ø 4.0 mm X 15 cm
- TWS 84-3049/18 Ø 4.0 mm X 20 cm
- TWS 84-3049/19 Ø 4.0 mm X 25 cm
- TWS 84-3049/20 Ø 4.0 mm X 30 cm
- TWS 84-3049/21 Ø 4.0 mm X 35 cm
- TWS 84-3049/22 Ø 5.0 mm X 15 cm
- TWS 84-3049/23 Ø 5.0 mm X 20 cm
- TWS 84-3049/24 Ø 5.0 mm X 25 cm
- TWS 84-3049/25 Ø 5.0 mm X 30 cm
- TWS 84-3049/26 Ø 5.0 mm X 35 cm
- TWS 84-3049/27 Ø 6.0 mm X 15 cm
- TWS 84-3049/28 Ø 6.0 mm X 20 cm
- TWS 84-3049/29 Ø 6.0 mm X 25 cm
- TWS 84-3049/30 Ø 6.0 mm X 30 cm
- TWS 84-3049/31 Ø 6.0 mm X 35 cm



Fig.Le



Power Cannula

- TWS 84-3049/32 Ø 2.0 mm X 15 cm
- TWS 84-3049/33 Ø 2.0 mm X 20 cm
- TWS 84-3049/34 Ø 2.0 mm X 25 cm
- TWS 84-3049/35 Ø 2.0 mm X 30 cm
- TWS 84-3049/36 Ø 2.0 mm X 35 cm
- TWS 84-3049/37 Ø 3.0 mm X 15 cm
- TWS 84-3049/38 Ø 3.0 mm X 20 cm
- TWS 84-3049/39 Ø 3.0 mm X 25 cm
- TWS 84-3049/40 Ø 3.0 mm X 30 cm
- TWS 84-3049/41 Ø 3.0 mm X 35 cm
- TWS 84-3049/42 Ø 4.0 mm X 15 cm
- TWS 84-3049/56 Ø 4.0 mm X 20 cm
- TWS 84-3049/43 Ø 4.0 mm X 25 cm
- TWS 84-3049/44 Ø 4.0 mm X 30 cm
- TWS 84-3049/45 Ø 4.0 mm X 35 cm
- TWS 84-3049/46 Ø 5.0 mm X 15 cm
- TWS 84-3049/47 Ø 5.0 mm X 20 cm
- TWS 84-3049/48 Ø 5.0 mm X 25 cm
- TWS 84-3049/49 Ø 5.0 mm X 30 cm
- TWS 84-3049/50 Ø 5.0 mm X 35 cm
- TWS 84-3049/51 Ø 6.0 mm X 15 cm
- TWS 84-3049/52 Ø 6.0 mm X 20 cm
- TWS 84-3049/53 Ø 6.0 mm X 25 cm
- TWS 84-3049/54 Ø 6.0 mm X 30 cm
- TWS 84-3049/55 Ø 6.0 mm X 35 cm



TW Surgical

HARVESTER CANNULAS



Multiport Speed Harvester- Both End Sharp

A multiport harvester with elevated 1mm microports which yield small enough aliquots of fat tissue to allow for easy and smooth injection with the Tulip 0.7, 0.9 and 1.2mm injectors.



Tommy Hub

TWS 84-3049/58	Ø 2.5 mm X 15 cm
TWS 84-3049/59	Ø 2.5 mm X 20 cm
TWS 84-3049/60	Ø 3.0 mm X 15 cm
TWS 84-3049/61	Ø 3.0 mm X 20 cm
TWS 84-3049/62	Ø 3.0 mm X 25 cm
TWS 84-3049/63	Ø 4.0 mm X 25 cm
TWS 84-3049/64	Ø 4.0 mm X 30 cm
TWS 84-3049/65	Ø 4.0 mm X 35 cm
TWS 84-3049/66	Ø 4.5 mm X 30 cm
TWS 84-3049/67	Ø 4.5 mm X 35 cm



Treaded Hub

TWS 84-3049/73	Ø 3.0mm X 15 cm
TWS 84-3049/74	Ø 3.0mm X 25 cm
TWS 84-3049/75	Ø 4.0mm X 15 cm
TWS 84-3049/76	Ø 4.0mm X 25 cm



Ultra Luer Lock

TWS 84-3049/77	Ø 2.0 mm X 15 cm
TWS 84-3049/78	Ø 2.5 mm X 15 cm
TWS 84-3049/79	Ø 2.5 mm X 20 cm
TWS 84-3049/80	Ø 3.0 mm X 15 cm
TWS 84-3049/81	Ø 3.0 mm X 20 cm
TWS 84-3049/82	Ø 3.0 mm X 25 cm



Luer Lock

TWS 84-3049/83	Ø 2.0 mm X 15 cm
TWS 84-3049/84	Ø 2.5 mm X 15 cm
TWS 84-3049/85	Ø 2.5 mm X 20 cm
TWS 84-3049/86	Ø 3.0 mm X 15 cm
TWS 84-3049/87	Ø 3.0 mm X 20 cm
TWS 84-3049/88	Ø 3.0 mm X 25 cm



Luer lock (Stainless Steel)

TWS 84-3049/57	Ø 2.0 mm X 15 cm
TWS 84-3049/68	Ø 2.5 mm X 15 cm
TWS 84-3049/69	Ø 2.5 mm X 20 cm
TWS 84-3049/70	Ø 3.0 mm X 15 cm
TWS 84-3049/71	Ø 3.0 mm X 20 cm
TWS 84-3049/72	Ø 3.0 mm X 25 cm



Power Cannula

TWS 84-3049/89	Ø 2.0 mm X 15 cm
TWS 84-3049/90	Ø 2.0 mm X 20 cm
TWS 84-3049/91	Ø 2.0 mm X 25 cm
TWS 84-3049/92	Ø 2.0 mm X 30 cm
TWS 84-3049/93	Ø 2.0 mm X 35 cm
TWS 84-3049/94	Ø 3.0 mm X 15 cm
TWS 84-3049/95	Ø 3.0 mm X 20 cm
TWS 84-3049/96	Ø 3.0 mm X 25 cm
TWS 84-3049/97	Ø 3.0 mm X 30 cm
TWS 84-3049/98	Ø 3.0 mm X 35 cm
TWS 84-3049/99	Ø 4.0 mm X 15 cm
TWS 84-3050/00	Ø 4.0 mm X 20 cm
TWS 84-3050/01	Ø 4.0 mm X 25 cm
TWS 84-3050/02	Ø 4.0 mm X 30 cm
TWS 84-3050/03	Ø 4.0 mm X 35 cm
TWS 84-3050/04	Ø 5.0 mm X 15 cm
TWS 84-3050/05	Ø 5.0 mm X 20 cm
TWS 84-3050/06	Ø 5.0 mm X 25 cm
TWS 84-3050/07	Ø 5.0 mm X 30 cm
TWS 84-3050/08	Ø 5.0 mm X 35 cm
TWS 84-3050/09	Ø 6.0 mm X 15 cm
TWS 84-3050/10	Ø 6.0 mm X 20 cm
TWS 84-3050/11	Ø 6.0 mm X 25 cm
TWS 84-3050/12	Ø 6.0 mm X 30 cm
TWS 84-3050/13	Ø 6.0 mm X 35 cm

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical

HARVESTER CANNULAS



Multiport Speed Harvester- One End Sharp

A multiport harvester with elevated 1mm microports which yield small enough aliquots of fat tissue to allow for easy and smooth injection with the Tulip 0.7, 0.9 and 1.2mm injectors.



Tommey Hub

TWS 84-3050/14	Ø 2.5 mm X 15 cm
TWS 84-3050/15	Ø 2.5 mm X 20 cm
TWS 84-3050/16	Ø 3.0 mm X 15 cm
TWS 84-3050/17	Ø 3.0 mm X 20 cm
TWS 84-3050/18	Ø 3.0 mm X 25 cm
TWS 84-3050/19	Ø 4.0 mm X 25 cm
TWS 84-3050/20	Ø 4.0 mm X 30 cm
TWS 84-3050/21	Ø 4.0 mm X 35 cm
TWS 84-3050/22	Ø 4.5 mm X 30 cm
TWS 84-3050/23	Ø 4.5 mm X 35 cm



Treaded Hub

TWS 84-3050/30	Ø 3.0mm X 15 cm
TWS 84-3050/31	Ø 3.0mm X 25 cm
TWS 84-3050/32	Ø 4.0mm X 15 cm
TWS 84-3050/33	Ø 4.0mm X 25 cm



Ultra Luer Lock

TWS 84-3050/34	Ø 2.0 mm X 15 cm
TWS 84-3050/35	Ø 2.5 mm X 15 cm
TWS 84-3050/36	Ø 2.5 mm X 20 cm
TWS 84-3050/37	Ø 3.0 mm X 15 cm
TWS 84-3050/38	Ø 3.0 mm X 20 cm
TWS 84-3050/39	Ø 3.0 mm X 25 cm



Luer Lock

TWS 84-3050/40	Ø 2.0 mm X 15 cm
TWS 84-3050/41	Ø 2.5 mm X 15 cm
TWS 84-3050/42	Ø 2.5 mm X 20 cm
TWS 84-3050/43	Ø 3.0 mm X 15 cm
TWS 84-3050/44	Ø 3.0 mm X 20 cm
TWS 84-3050/45	Ø 3.0 mm X 25 cm



Luer lock (Stainless Steel)

TWS 84-3050/24	Ø 2.0 mm X 15 cm
TWS 84-3050/25	Ø 2.5 mm X 15 cm
TWS 84-3050/26	Ø 2.5 mm X 20 cm
TWS 84-3050/27	Ø 3.0 mm X 15 cm
TWS 84-3050/28	Ø 3.0 mm X 20 cm
TWS 84-3050/29	Ø 3.0 mm X 25 cm



Power Cannula

TWS 84-3050/46	Ø2.0 mm X 15 cm
TWS 84-3050/47	Ø2.0 mm X 20 cm
TWS 84-3050/48	Ø 2.0 mm X 25 cm
TWS 84-3050/49	Ø 2.0 mm X 30 cm
TWS 84-3050/50	Ø 2.0 mm X 35 cm
TWS 84-3050/51	Ø 3.0 mm X 15 cm
TWS 84-3050/52	Ø 3.0 mm X 20 cm
TWS 84-3050/53	Ø 3.0 mm X 25 cm
TWS 84-3050/54	Ø 3.0 mm X 30 cm
TWS 84-3050/55	Ø 3.0 mm X 35 cm
TWS 84-3050/56	Ø 4.0 mm X 15 cm
TWS 84-3050/57	Ø 4.0 mm X 20 cm
TWS 84-3050/58	Ø 4.0 mm X 25 cm
TWS 84-3050/59	Ø 4.0 mm X 30 cm
TWS 84-3050/60	Ø 4.0 mm X 35 cm
TWS 84-3050/61	Ø 5.0 mm X 15 cm
TWS 84-3050/62	Ø 5.0 mm X 20 cm
TWS 84-3050/63	Ø 5.0 mm X 25 cm
TWS 84-3050/64	Ø 5.0 mm X 30 cm
TWS 84-3050/65	Ø 5.0 mm X 35 cm
TWS 84-3050/66	Ø 6.0 mm X 15 cm
TWS 84-3050/67	Ø 6.0 mm X 20 cm
TWS 84-3050/68	Ø 6.0 mm X 25 cm
TWS 84-3050/69	Ø 6.0 mm X 30 cm
TWS 84-3050/70	Ø 6.0 mm X 35 cm



TW Surgical

HARVESTER CANNULAS



Multi Port Harvester - Both End Bevel

A sophisticated, precision-designed harvester that is truly Patient Friendly. The ports in the Sorensen Harvester are machined with alternating 60 degree cutting edges. This unique design is less traumatic than other multi-port harvesters and delivers a fine-grained, high fat content/purity yield that can easily pass through 1.2mm, 0.9mm and 0.7mm blunt injectors.



Tommey Hub

TWS 84-3050/71	Ø 2.5 mm X 15 cm
TWS 84-3050/72	Ø 2.5 mm X 20 cm
TWS 84-3050/73	Ø 3.0 mm X 15 cm
TWS 84-3050/74	Ø 3.0 mm X 20 cm
TWS 84-3050/75	Ø 3.0 mm X 25 cm
TWS 84-3050/76	Ø 4.0 mm X 25 cm
TWS 84-3050/77	Ø 4.0 mm X 30 cm
TWS 84-3050/78	Ø 4.0 mm X 35 cm
TWS 84-3050/79	Ø 4.5 mm X 30 cm
TWS 84-3050/80	Ø 4.5 mm X 35 cm



Treaded Hub

TWS 84-3050/87	Ø 3.0mm X 15 cm
TWS 84-3050/88	Ø 3.0mm X 25 cm
TWS 84-3050/89	Ø 4.0mm X 15 cm
TWS 84-3050/90	Ø 4.0mm X 25 cm



Ultra Luer Lock

TWS 84-3050/81	Ø 2.0 mm X 15 cm
TWS 84-3050/82	Ø 2.5 mm X 15 cm
TWS 84-3050/83	Ø 2.5 mm X 20 cm
TWS 84-3050/84	Ø 3.0 mm X 15 cm
TWS 84-3050/85	Ø 3.0 mm X 20 cm
TWS 84-3050/86	Ø 3.0 mm X 25 cm



Luer Lock

TWS 84-3050/91	Ø 2.0 mm X 15 cm
TWS 84-3050/92	Ø 2.5 mm X 15 cm
TWS 84-3050/93	Ø 2.5 mm X 20 cm
TWS 84-3050/94	Ø 3.0 mm X 15 cm
TWS 84-3050/95	Ø 3.0 mm X 20 cm
TWS 84-3050/96	Ø 3.0 mm X 25 cm



Luer lock (Stainless Steel)

TWS 84-3050/97	Ø 2.0 mm X 15 cm
TWS 84-3050/98	Ø 2.5 mm X 15 cm
TWS 84-3050/99	Ø 2.5 mm X 20 cm
TWS 84-3051/00	Ø 3.0 mm X 15 cm
TWS 84-3051/01	Ø 3.0 mm X 20 cm
TWS 84-3051/02	Ø 3.0 mm X 25 cm



Power Cannula

TWS 84-3051/03	Ø 2.0 mm X 15 cm
TWS 84-3051/04	Ø 2.0 mm X 20 cm
TWS 84-3051/05	Ø 2.0 mm X 25 cm
TWS 84-3051/06	Ø 2.0 mm X 30 cm
TWS 84-3051/07	Ø 2.0 mm X 35 cm
TWS 84-3051/08	Ø 3.0 mm X 15 cm
TWS 84-3051/09	Ø 3.0 mm X 20 cm
TWS 84-3051/10	Ø 3.0 mm X 25 cm
TWS 84-3051/11	Ø 3.0 mm X 30 cm
TWS 84-3051/12	Ø 3.0 mm X 35 cm
TWS 84-3051/13	Ø 4.0 mm X 15 cm
TWS 84-3051/14	Ø 4.0 mm X 20 cm
TWS 84-3051/15	Ø 4.0 mm X 25 cm
TWS 84-3051/16	Ø 4.0 mm X 30 cm
TWS 84-3051/17	Ø 4.0 mm X 35 cm
TWS 84-3051/18	Ø 5.0 mm X 15 cm
TWS 84-3051/19	Ø 5.0 mm X 20 cm
TWS 84-3051/20	Ø 5.0 mm X 25 cm
TWS 84-3051/21	Ø 5.0 mm X 30 cm
TWS 84-3051/22	Ø 5.0 mm X 35 cm
TWS 84-3051/23	Ø 6.0 mm X 15 cm
TWS 84-3051/24	Ø 6.0 mm X 20 cm
TWS 84-3051/25	Ø 6.0 mm X 25 cm
TWS 84-3051/26	Ø 6.0 mm X 30 cm
TWS 84-3051/27	Ø 6.0 mm X 35 cm



TW Surgical

MICRO INJECTORS



Sharp Extractor/Injector

This tip works well for needle extraction of soft tissue and as a soft tissue reinjection needle



Ultra Luer Lock

TWS 84-3051/28	Ø 1.0 mm X 5 cm
TWS 84-3051/29	Ø 1.5 mm X 7 cm
TWS 84-3051/30	Ø 2.0 mm X 15 cm
TWS 84-3051/31	Ø 2.5 mm X 15 cm
TWS 84-3051/32	Ø 3.0 mm X 15 cm



Luer Lock

TWS 84-3051/33	Ø 1.0 mm X 5 cm
TWS 84-3051/34	Ø 1.0 mm X 7 cm
TWS 84-3051/35	Ø 2.0 mm X 15 cm
TWS 84-3051/36	Ø 2.5 mm X 15 cm
TWS 84-3051/37	Ø 3.0 mm X 15 cm



Luer lock (Stainless Steel)

TWS 84-3051/60	Ø 1.1 mm X 3 cm
TWS 84-3051/61	Ø 1.4 mm X 3 cm
TWS 84-3051/62	Ø 1.8 mm X 3 cm
TWS 84-3051/63	Ø 2.0 mm X 3 cm
TWS 84-3051/64	Ø 1.1 mm X 5 cm
TWS 84-3051/65	Ø 1.4 mm X 5 cm
TWS 84-3051/66	Ø 1.8 mm X 5 cm
TWS 84-3051/67	Ø 2.0 mm X 5 cm
TWS 84-3051/68	Ø 1.1 mm X 7 cm
TWS 84-3051/69	Ø 1.4 mm X 7 cm
TWS 84-3051/70	Ø 1.8 mm X 7 cm
TWS 84-3051/71	Ø 2.0 mm X 7 cm
TWS 84-3051/72	Ø 1.1 mm X 10 cm
TWS 84-3051/73	Ø 1.4 mm X 10 cm
TWS 84-3051/74	Ø 1.8 mm X 10 cm
TWS 84-3051/75	Ø 2.0 mm X 10 cm
TWS 84-3051/76	Ø 1.1 mm X 12 cm
TWS 84-3051/77	Ø 1.4 mm X 12 cm
TWS 84-3051/78	Ø 1.8 mm X 12 cm
TWS 84-3051/79	Ø 2.0 mm X 12 cm

Blunt Extractor/Injector

This tip works well for needle extraction of soft tissue and as a soft tissue reinjection needle.



Luer lock (Stainless Steel)

TWS 84-3051/38	Ø 1.1 mm X 3 cm
TWS 84-3051/39	Ø 1.4 mm X 3 cm
TWS 84-3051/40	Ø 1.8 mm X 3 cm
TWS 84-3051/41	Ø 2.0 mm X 3 cm
TWS 84-3051/42	Ø 1.1 mm X 5 cm
TWS 84-3051/43	Ø 1.4 mm X 5 cm
TWS 84-3051/44	Ø 1.8 mm X 5 cm
TWS 84-3051/45	Ø 2.0 mm X 5 cm
TWS 84-3051/46	Ø 1.1 mm X 7 cm
TWS 84-3051/47	Ø 1.4 mm X 7 cm
TWS 84-3051/48	Ø 1.8 mm X 7 cm
TWS 84-3051/49	Ø 2.0 mm X 7 cm
TWS 84-3051/50	Ø 1.1 mm X 10 cm
TWS 84-3051/51	Ø 1.4 mm X 10 cm
TWS 84-3051/52	Ø 1.8 mm X 10 cm
TWS 84-3051/53	Ø 2.0 mm X 10 cm



Ultra Luer Lock

TWS 84-3051/54	Ø 2.0 mm X 15 cm
TWS 84-3051/55	Ø 2.5 mm X 15 cm
TWS 84-3051/56	Ø 3.0 mm X 15 cm



Luer Lock

TWS 84-3051/57	Ø 2.0 mm X 15 cm
TWS 84-3051/58	Ø 2.5 mm X 15 cm
TWS 84-3051/59	Ø 3.0 mm X 15 cm



Blunt Tip Injector

Rounded blunt tip design passes through soft tissue easily without the trauma or potential problems of sharp instruments.



Ultra Luer Lock

TWS 84-3051/80	Ø 1.0 mm X 4 cm
TWS 84-3051/81	Ø 1.0 mm X 5 cm
TWS 84-3051/82	Ø 1.0 mm X 7 cm
TWS 84-3051/83	Ø Curved X 12 cm
TWS 84-3051/84	Ø Curved X 15 cm
TWS 84-3051/85	Ø Curved X 20 cm



Luer Lock

TWS 84-3051/86	Ø 1.0 mm X 4 cm
TWS 84-3051/87	Ø 1.0 mm X 5 cm
TWS 84-3051/88	Ø 1.0 mm X 7 cm
TWS 84-3051/89	Ø Curved X 12 cm
TWS 84-3051/90	Ø Curved X 15 cm
TWS 84-3051/91	Ø Curved X 20 cm



Luer lock (Stainless Steel)

TWS 84-3051/92	Ø 1.1 mm X 3 cm
TWS 84-3051/93	Ø 1.4 mm X 3 cm
TWS 84-3051/94	Ø 1.8 mm X 3 cm
TWS 84-3051/95	Ø 2.0 mm X 3 cm
TWS 84-3051/96	Ø 1.1 mm X 5 cm
TWS 84-3051/97	Ø 1.4 mm X 5 cm
TWS 84-3051/98	Ø 1.8 mm X 5 cm
TWS 84-3051/99	Ø 2.0 mm X 5 cm
TWS 84-3052/00	Ø 1.1 mm X 7 cm
TWS 84-3052/01	Ø 1.4 mm X 7 cm
TWS 84-3052/02	Ø 1.8 mm X 7 cm
TWS 84-3052/03	Ø 2.0 mm X 7 cm
TWS 84-3052/04	Ø 1.1 mm X 10 cm
TWS 84-3052/05	Ø 1.4 mm X 10 cm
TWS 84-3052/06	Ø 1.8 mm X 10 cm
TWS 84-3052/07	Ø 2.0 mm X 10 cm
TWS 84-3052/08	Ø 1.1 mm X 12 cm
TWS 84-3052/09	Ø 1.4 mm X 12 cm
TWS 84-3052/10	Ø 1.8 mm X 12 cm
TWS 84-3052/11	Ø 2.0 mm X 12 cm

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical

MICRO INJECTORS



Fig. 2B



Luer lock (Stainless Steel)

TWS 84-3052/11	Ø 1.1 mm X 3 cm
TWS 84-3052/12	Ø 1.4 mm X 3 cm
TWS 84-3052/13	Ø 1.8 mm X 3 cm
TWS 84-3052/14	Ø 2.0 mm X 3 cm
TWS 84-3052/15	Ø 1.1 mm X 5 cm
TWS 84-3052/16	Ø 1.4 mm X 5 cm
TWS 84-3052/17	Ø 1.8 mm X 5 cm
TWS 84-3052/18	Ø 2.0 mm X 5 cm
TWS 84-3052/19	Ø 1.1 mm X 7 cm
TWS 84-3052/20	Ø 1.4 mm X 7 cm
TWS 84-3052/21	Ø 1.8 mm X 7 cm
TWS 84-3052/22	Ø 2.0 mm X 7 cm
TWS 84-3052/23	Ø 1.1 mm X 10 cm
TWS 84-3052/24	Ø 1.4 mm X 10 cm
TWS 84-3052/25	Ø 1.8 mm X 10 cm
TWS 84-3052/26	Ø 2.0 mm X 10 cm
TWS 84-3052/27	Ø 1.1 mm X 12 cm
TWS 84-3052/28	Ø 1.4 mm X 12 cm
TWS 84-3052/29	Ø 1.8 mm X 12 cm
TWS 84-3052/30	Ø 2.0 mm X 12 cm



Fig. 1P



Luer lock (Stainless Steel)

TWS 84-3052/51	Ø 1.1 mm X 3 cm
TWS 84-3052/52	Ø 1.4 mm X 3 cm
TWS 84-3052/53	Ø 1.8 mm X 3 cm
TWS 84-3052/54	Ø 2.0 mm X 3 cm
TWS 84-3052/55	Ø 1.1 mm X 5 cm
TWS 84-3052/56	Ø 1.4 mm X 5 cm
TWS 84-3052/57	Ø 1.8 mm X 5 cm
TWS 84-3052/58	Ø 2.0 mm X 5 cm
TWS 84-3052/59	Ø 1.1 mm X 7 cm
TWS 84-3052/60	Ø 1.4 mm X 7 cm
TWS 84-3052/61	Ø 1.8 mm X 7 cm
TWS 84-3052/62	Ø 2.0 mm X 7 cm
TWS 84-3052/63	Ø 1.1 mm X 10 cm
TWS 84-3052/64	Ø 1.4 mm X 10 cm
TWS 84-3052/65	Ø 1.8 mm X 10 cm
TWS 84-3052/66	Ø 2.0 mm X 10 cm
TWS 84-3052/67	Ø 1.1 mm X 12 cm
TWS 84-3052/68	Ø 1.4 mm X 12 cm
TWS 84-3052/69	Ø 1.8 mm X 12 cm
TWS 84-3052/70	Ø 2.0 mm X 12 cm



Fig. 1C



Luer lock (Stainless Steel)

TWS 84-3052/31	Ø 1.1 mm X 3 cm
TWS 84-3052/32	Ø 1.4 mm X 3 cm
TWS 84-3052/33	Ø 1.8 mm X 3 cm
TWS 84-3052/34	Ø 2.0 mm X 3 cm
TWS 84-3052/35	Ø 1.1 mm X 5 cm
TWS 84-3052/36	Ø 1.4 mm X 5 cm
TWS 84-3052/37	Ø 1.8 mm X 5 cm
TWS 84-3052/38	Ø 2.0 mm X 5 cm
TWS 84-3052/39	Ø 1.1 mm X 7 cm
TWS 84-3052/40	Ø 1.4 mm X 7 cm
TWS 84-3052/41	Ø 1.8 mm X 7 cm
TWS 84-3052/42	Ø 2.0 mm X 7 cm
TWS 84-3052/43	Ø 1.1 mm X 10 cm
TWS 84-3052/44	Ø 1.4 mm X 10 cm
TWS 84-3052/45	Ø 1.8 mm X 10 cm
TWS 84-3052/46	Ø 2.0 mm X 10 cm
TWS 84-3052/47	Ø 1.1 mm X 12 cm
TWS 84-3052/48	Ø 1.4 mm X 12 cm
TWS 84-3052/49	Ø 1.8 mm X 12 cm
TWS 84-3052/50	Ø 2.0 mm X 12 cm



Fig. 2E



Luer lock (Stainless Steel)

TWS 84-3052/71	Ø 1.1 mm X 3 cm
TWS 84-3052/72	Ø 1.4 mm X 3 cm
TWS 84-3052/73	Ø 1.8 mm X 3 cm
TWS 84-3052/74	Ø 2.0 mm X 3 cm
TWS 84-3052/75	Ø 1.1 mm X 5 cm
TWS 84-3052/76	Ø 1.4 mm X 5 cm
TWS 84-3052/77	Ø 1.8 mm X 5 cm
TWS 84-3052/78	Ø 2.0 mm X 5 cm
TWS 84-3052/79	Ø 1.1 mm X 7 cm
TWS 84-3052/80	Ø 1.4 mm X 7 cm
TWS 84-3052/81	Ø 1.8 mm X 7 cm
TWS 84-3052/82	Ø 2.0 mm X 7 cm
TWS 84-3052/83	Ø 1.1 mm X 10 cm
TWS 84-3052/84	Ø 1.4 mm X 10 cm
TWS 84-3052/85	Ø 1.8 mm X 10 cm
TWS 84-3052/86	Ø 2.0 mm X 10 cm
TWS 84-3052/87	Ø 1.1 mm X 12 cm
TWS 84-3052/88	Ø 1.4 mm X 12 cm
TWS 84-3052/89	Ø 1.8 mm X 12 cm
TWS 84-3052/90	Ø 2.0 mm X 12 cm



TW Surgical

MICRO INJECTORS



Coleman Liposuction Cannula Fig. MIC I



Luer lock (Stainless Steel)

TWS 84-3052/91	Ø 1.1 mm X 3 cm
TWS 84-3052/92	Ø 1.4 mm X 3 cm
TWS 84-3052/93	Ø 1.8 mm X 3 cm
TWS 84-3052/94	Ø 2.0 mm X 3 cm
TWS 84-3052/95	Ø 2.5 mm X 3 cm
TWS 84-3052/96	Ø 1.1 mm X 5 cm
TWS 84-3052/97	Ø 1.4 mm X 5 cm
TWS 84-3052/98	Ø 1.8 mm X 5 cm
TWS 84-3052/99	Ø 2.0 mm X 5 cm
TWS 84-3053/00	Ø 2.5 mm X 5 cm
TWS 84-3053/01	Ø 1.1 mm X 7 cm
TWS 84-3053/02	Ø 1.4 mm X 7 cm
TWS 84-3053/03	Ø 1.8 mm X 7 cm
TWS 84-3053/04	Ø 2.0 mm X 7 cm
TWS 84-3053/05	Ø 2.5 mm X 7 cm
TWS 84-3053/06	Ø 1.1 mm X 10 cm
TWS 84-3053/07	Ø 1.4 mm X 10 cm
TWS 84-3053/08	Ø 1.8 mm X 10 cm
TWS 84-3053/09	Ø 2.0 mm X 10 cm
TWS 84-3053/10	Ø 2.5 mm X 10 cm
TWS 84-3053/11	Ø 1.1 mm X 12 cm
TWS 84-3053/12	Ø 1.4 mm X 12 cm
TWS 84-3053/13	Ø 1.8 mm X 12 cm
TWS 84-3053/14	Ø 2.0 mm X 12 cm
TWS 84-3053/15	Ø 2.5 mm X 12 cm
TWS 84-3053/16	Ø 1.1 mm X 15 cm
TWS 84-3053/17	Ø 1.4 mm X 15 cm
TWS 84-3053/18	Ø 1.8 mm X 15 cm
TWS 84-3053/19	Ø 2.0 mm X 15 cm
TWS 84-3053/20	Ø 2.5 mm X 15 cm
TWS 84-3053/21	Ø 1.1 mm X 20 cm
TWS 84-3053/22	Ø 1.4 mm X 20 cm
TWS 84-3053/23	Ø 1.8 mm X 20 cm
TWS 84-3053/24	Ø 2.0 mm X 20 cm
TWS 84-3053/25	Ø 2.5 mm X 20 cm
TWS 84-3053/26	Ø 1.1 mm X 3.5 cm
TWS 84-3053/27	Ø 1.4 mm X 3.5 cm
TWS 84-3053/28	Ø 1.8 mm X 3.5 cm
TWS 84-3053/29	Ø 2.0 mm X 3.5 cm
TWS 84-3053/30	Ø 2.5 mm X 3.5 cm



Coleman Liposuction Cannula Fig. MIC II



Luer lock (Stainless Steel)

TWS 84-3053/31	Ø 1.1 mm X 3 cm
TWS 84-3053/32	Ø 1.4 mm X 3 cm
TWS 84-3053/33	Ø 1.8 mm X 3 cm
TWS 84-3053/34	Ø 2.0 mm X 3 cm
TWS 84-3053/35	Ø 2.5 mm X 3 cm
TWS 84-3053/36	Ø 1.1 mm X 5 cm
TWS 84-3053/37	Ø 1.4 mm X 5 cm
TWS 84-3053/37	Ø 1.8 mm X 5 cm
TWS 84-3053/39	Ø 2.0 mm X 5 cm
TWS 84-3053/40	Ø 2.5 mm X 5 cm
TWS 84-3053/41	Ø 1.1 mm X 7 cm
TWS 84-3053/42	Ø 1.4 mm X 7 cm
TWS 84-3053/43	Ø 1.8 mm X 7 cm
TWS 84-3053/44	Ø 2.0 mm X 7 cm
TWS 84-3053/45	Ø 2.5 mm X 7 cm
TWS 84-3053/46	Ø 1.1 mm X 10 cm
TWS 84-3053/47	Ø 1.4 mm X 10 cm
TWS 84-3053/48	Ø 1.8 mm X 10 cm
TWS 84-3053/49	Ø 2.0 mm X 10 cm
TWS 84-3053/50	Ø 2.5 mm X 10 cm
TWS 84-3053/51	Ø 1.1 mm X 12 cm
TWS 84-3053/52	Ø 1.4 mm X 12 cm
TWS 84-3053/53	Ø 1.8 mm X 12 cm
TWS 84-3053/54	Ø 2.0 mm X 12 cm
TWS 84-3053/55	Ø 2.5 mm X 12 cm
TWS 84-3053/56	Ø 1.1 mm X 15 cm
TWS 84-3053/57	Ø 1.4 mm X 15 cm
TWS 84-3053/58	Ø 1.8 mm X 15 cm
TWS 84-3053/59	Ø 2.0 mm X 15 cm
TWS 84-3053/60	Ø 2.5 mm X 15 cm
TWS 84-3053/61	Ø 1.1 mm X 20 cm
TWS 84-3053/62	Ø 1.4 mm X 20 cm
TWS 84-3053/63	Ø 1.8 mm X 20 cm
TWS 84-3053/64	Ø 2.0 mm X 20 cm
TWS 84-3053/65	Ø 2.5 mm X 20 cm



TW Surgical

MICRO INJECTORS



Coleman Liposuction Cannula Spoon Tip Fig. MIC III

This tip works well for needle extraction of soft tissue and as a soft tissue re injection needle.



Luer lock (Stainless Steel)

TWS 84-3053/66	Ø 1.1 mm X 3 cm
TWS 84-3053/67	Ø 1.4 mm X 3 cm
TWS 84-3053/68	Ø 1.8 mm X 3 cm
TWS 84-3053/69	Ø 2.0 mm X 3 cm
TWS 84-3053/70	Ø 2.5 mm X 3 cm
TWS 84-3053/71	Ø 1.1 mm X 5 cm
TWS 84-3053/73	Ø 1.4 mm X 5 cm
TWS 84-3053/74	Ø 1.8 mm X 5 cm
TWS 84-3053/75	Ø 2.0 mm X 5 cm
TWS 84-3053/76	Ø 2.5 mm X 5 cm
TWS 84-3053/77	Ø 1.1 mm X 7 cm
TWS 84-3053/78	Ø 1.4 mm X 7 cm
TWS 84-3053/79	Ø 1.8 mm X 7 cm
TWS 84-3053/80	Ø 2.0 mm X 7 cm
TWS 84-3053/81	Ø 2.5 mm X 7 cm
TWS 84-3053/82	Ø 1.1 mm X 10 cm
TWS 84-3053/83	Ø 1.4 mm X 10 cm
TWS 84-3053/84	Ø 1.8 mm X 10 cm
TWS 84-3053/85	Ø 2.0 mm X 10 cm
TWS 84-3053/86	Ø 2.5 mm X 10 cm
TWS 84-3053/87	Ø 1.1 mm X 12 cm
TWS 84-3053/88	Ø 1.4 mm X 12 cm
TWS 84-3053/89	Ø 1.8 mm X 12 cm
TWS 84-3053/90	Ø 2.0 mm X 12 cm
TWS 84-3053/91	Ø 2.5 mm X 12 cm
TWS 84-3053/92	Ø 1.1 mm X 15 cm
TWS 84-3053/93	Ø 1.4 mm X 15 cm
TWS 84-3053/94	Ø 1.8 mm X 15 cm
TWS 84-3053/95	Ø 2.0 mm X 15 cm
TWS 84-3053/96	Ø 2.5 mm X 15 cm
TWS 84-3053/97	Ø 1.1 mm X 20 cm
TWS 84-3053/98	Ø 1.4 mm X 20 cm
TWS 84-3053/99	Ø 1.8 mm X 20 cm
TWS 84-3054/00	Ø 2.0 mm X 20 cm
TWS 84-3054/01	Ø 2.5 mm X 20 cm



Little Injector

Open and more pointed blunt tip is very effective in dense and fibrous tissue.



Ultra Luer Lock

TWS 84-3054/02	Ø 1.0 mm X 4 cm
TWS 84-3054/03	Ø 1.0 mm X 5 cm
TWS 84-3054/04	Ø 1.0 mm X 7 cm



Luer Lock

TWS 84-3054/05	Ø 1.0 mm X 4 cm
TWS 84-3054/06	Ø 1.0 mm X 5 cm
TWS 84-3054/07	Ø 1.0 mm X 7 cm



Luer lock (Stainless Steel)

TWS 84-3054/08	Ø 1.1 mm X 3 cm
TWS 84-3054/09	Ø 1.4 mm X 3 cm
TWS 84-3054/11	Ø 1.8 mm X 3 cm
TWS 84-3054/10	Ø 2.0 mm X 3 cm
TWS 84-3054/12	Ø 1.1 mm X 5 cm
TWS 84-3054/13	Ø 1.4 mm X 5 cm
TWS 84-3054/14	Ø 1.8 mm X 5 cm
TWS 84-3054/15	Ø 2.0 mm X 5 cm
TWS 84-3054/16	Ø 1.1 mm X 7 cm
TWS 84-3054/17	Ø 1.4 mm X 7 cm
TWS 84-3054/18	Ø 1.8 mm X 7 cm
TWS 84-3054/19	Ø 2.0 mm X 7 cm
TWS 84-3054/20	Ø 1.1 mm X 10 cm
TWS 84-3054/21	Ø 1.4 mm X 10 cm
TWS 84-3054/22	Ø 1.8 mm X 10 cm
TWS 84-3054/23	Ø 2.0 mm X 10 cm
TWS 84-3054/24	Ø 1.1 mm X 12 cm
TWS 84-3054/25	Ø 1.4 mm X 12 cm
TWS 84-3054/26	Ø 1.8 mm X 12 cm
TWS 84-3054/27	Ø 2.0 mm X 12 cm



TW Surgical

MICRO INJECTORS



Fig. 1B

This tip works well for needle extraction of soft tissue and as a soft tissue re injection needle.



Luer lock (Stainless Steel)

TWS 84-3054/28	Ø 1.1 mm X 3 cm
TWS 84-3054/29	Ø 1.4 mm X 3 cm
TWS 84-3054/30	Ø 1.8 mm X 3 cm
TWS 84-3054/31	Ø 2.0 mm X 3 cm
TWS 84-3054/32	Ø 2.5 mm X 3 cm
TWS 84-3054/33	Ø 1.1 mm X 5 cm
TWS 84-3054/34	Ø 1.4 mm X 5 cm
TWS 84-3054/35	Ø 1.8 mm X 5 cm
TWS 84-3054/36	Ø 2.0 mm X 5 cm
TWS 84-3054/37	Ø 2.5 mm X 5 cm
TWS 84-3054/38	Ø 1.1 mm X 7 cm
TWS 84-3054/39	Ø 1.4 mm X 7 cm
TWS 84-3054/40	Ø 1.8 mm X 7 cm
TWS 84-3054/41	Ø 2.0 mm X 7 cm
TWS 84-3054/42	Ø 2.5 mm X 7 cm
TWS 84-3054/43	Ø 1.1 mm X 10 cm
TWS 84-3054/44	Ø 1.4 mm X 10 cm
TWS 84-3054/45	Ø 1.8 mm X 10 cm
TWS 84-3054/46	Ø 2.0 mm X 10 cm
TWS 84-3054/47	Ø 2.5 mm X 10 cm
TWS 84-3054/48	Ø 1.1 mm X 12 cm
TWS 84-3054/49	Ø 1.4 mm X 12 cm
TWS 84-3054/50	Ø 1.8 mm X 12 cm
TWS 84-3054/51	Ø 2.0 mm X 12 cm
TWS 84-3054/52	Ø 2.5 mm X 12 cm
TWS 84-3054/53	Ø 1.1 mm X 15 cm
TWS 84-3054/54	Ø 1.4 mm X 15 cm
TWS 84-3054/55	Ø 1.8 mm X 15 cm
TWS 84-3054/56	Ø 2.0 mm X 15 cm
TWS 84-3054/57	Ø 2.5 mm X 15 cm
TWS 84-3054/58	Ø 1.1 mm X 20 cm
TWS 84-3054/59	Ø 1.4 mm X 20 cm
TWS 84-3054/60	Ø 1.8 mm X 20 cm
TWS 84-3054/61	Ø 2.0 mm X 20 cm
TWS 84-3054/62	Ø 2.5 mm X 20 cm



TW Surgical

SPECIALITY CANNULAS



V-Dissector

This cannula is specifically designed for use in the superficial compartment of the truncal area. The cutting notch is very useful for undermining in facial procedures and in dense, fibrous tissue.



Tommey Hub

- TWS 84-3054/63 Ø 2.5 mm X 15 cm
- TWS 84-3054/64 Ø 2.5 mm X 20 cm
- TWS 84-3054/65 Ø 3.0 mm X 15 cm
- TWS 84-3054/66 Ø 3.0 mm X 20 cm
- TWS 84-3054/67 Ø 3.0 mm X 25 cm
- TWS 84-3054/68 Ø 4.0 mm X 25 cm
- TWS 84-3054/69 Ø 4.0 mm X 30 cm
- TWS 84-3054/70 Ø 4.0 mm X 35 cm
- TWS 84-3054/71 Ø 4.5 mm X 30 cm
- TWS 84-3054/72 Ø 4.5 mm X 35 cm



Treaded Hub

- TWS 84-3055/04 Ø 3.0mm X 15 cm
- TWS 84-3055/05 Ø 3.0mm X 25 cm
- TWS 84-3055/06 Ø 4.0mm X 15 cm
- TWS 84-3055/07 Ø 4.0mm X 25 cm



Ultra Luer Lock

- TWS 84-3055/08 Ø 2.0 mm X 15 cm
- TWS 84-3055/09 Ø 2.5 mm X 15 cm
- TWS 84-3055/10 Ø 2.5 mm X 20 cm
- TWS 84-3055/11 Ø 3.0 mm X 15 cm
- TWS 84-3055/12 Ø 3.0 mm X 20 cm
- TWS 84-3055/13 Ø 3.0 mm X 25 cm



Luer Lock

- TWS 84-3055/14 Ø 2.0 mm X 15 cm
- TWS 84-3055/15 Ø 2.5 mm X 15 cm
- TWS 84-3055/16 Ø 2.5 mm X 20 cm
- TWS 84-3055/17 Ø 3.0 mm X 15 cm
- TWS 84-3055/18 Ø 3.0 mm X 20 cm
- TWS 84-3055/19 Ø 3.0 mm X 25 cm



Luer lock (Stainless Steel)

- TWS 84-3054/73 Ø 2.0 mm X 15 cm
- TWS 84-3054/74 Ø 2.5 mm X 15 cm
- TWS 84-3054/75 Ø 2.5 mm X 20 cm
- TWS 84-3054/76 Ø 3.0 mm X 15 cm
- TWS 84-3054/77 Ø 3.0 mm X 20 cm
- TWS 84-3054/78 Ø 3.0 mm X 25 cm



Power Cannula

- TWS 84-3054/79 Ø 2.0 mm X 15 cm
- TWS 84-3054/80 Ø 2.0 mm X 20 cm
- TWS 84-3054/81 Ø 2.0 mm X 25 cm
- TWS 84-3054/82 Ø 2.0 mm X 30 cm
- TWS 84-3054/83 Ø 2.0 mm X 35 cm
- TWS 84-3054/84 Ø 3.0 mm X 15 cm
- TWS 84-3054/85 Ø 3.0 mm X 20 cm
- TWS 84-3054/86 Ø 3.0 mm X 25 cm
- TWS 84-3054/87 Ø 3.0 mm X 30 cm
- TWS 84-3054/88 Ø 3.0 mm X 35 cm
- TWS 84-3054/89 Ø 4.0 mm X 15 cm
- TWS 84-3054/90 Ø 4.0 mm X 20 cm
- TWS 84-3054/91 Ø 4.0 mm X 25 cm
- TWS 84-3054/92 Ø 4.0 mm X 30 cm
- TWS 84-3054/93 Ø 4.0 mm X 35 cm
- TWS 84-3054/94 Ø 5.0 mm X 15 cm
- TWS 84-3054/95 Ø 5.0 mm X 20 cm
- TWS 84-3054/96 Ø 5.0 mm X 25 cm
- TWS 84-3054/97 Ø 5.0 mm X 30 cm
- TWS 84-3054/98 Ø 5.0 mm X 35 cm
- TWS 84-3054/99 Ø 6.0 mm X 15 cm
- TWS 84-3055/00 Ø 6.0 mm X 20 cm
- TWS 84-3055/01 Ø 6.0 mm X 25 cm
- TWS 84-3055/02 Ø 6.0 mm X 30 cm
- TWS 84-3055/03 Ø 6.0 mm X 35 cm

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical

SPECIALITY CANNULAS



Flap

The Flap is primarily used as an undermining dissector, and in smaller diameters, for facial tissue removal. The Flap Cannula in a 2.1mm diameter is very useful for undermining the nasolabial fold prior to autologous tissue injections, as well as for tissue removal in this anatomic area.



Tommey Hub

TWS 84-3055/20	Ø 2.5 mm X 15 cm
TWS 84-3055/21	Ø 2.5 mm X 20 cm
TWS 84-3055/22	Ø 3.0 mm X 15 cm
TWS 84-3055/23	Ø 3.0 mm X 20 cm
TWS 84-3055/24	Ø 3.0 mm X 25 cm
TWS 84-3055/25	Ø 4.0 mm X 25 cm
TWS 84-3055/26	Ø 4.0 mm X 30 cm
TWS 84-3055/27	Ø 4.0 mm X 35 cm
TWS 84-3055/28	Ø 4.5 mm X 30 cm
TWS 84-3055/29	Ø 4.5 mm X 35 cm



Treaded Hub

TWS 84-3055/30	Ø 3.0mm X 15 cm
TWS 84-3055/31	Ø 3.0mm X 25 cm
TWS 84-3055/32	Ø 4.0mm X 15 cm
TWS 84-3055/33	Ø 4.0mm X 25 cm



Ultra Luer Lock

TWS 84-3055/34	Ø 2.0 mm X 15 cm
TWS 84-3055/35	Ø 2.5 mm X 15 cm
TWS 84-3055/36	Ø 2.5 mm X 20 cm
TWS 84-3055/37	Ø 3.0 mm X 15 cm
TWS 84-3055/38	Ø 3.0 mm X 20 cm
TWS 84-3055/39	Ø 3.0 mm X 25 cm



Luer Lock

TWS 84-3055/40	Ø 2.0 mm X 15 cm
TWS 84-3055/41	Ø 2.5 mm X 15 cm
TWS 84-3055/42	Ø 2.5 mm X 20 cm
TWS 84-3055/43	Ø 3.0 mm X 15 cm
TWS 84-3055/44	Ø 3.0 mm X 20 cm
TWS 84-3055/45	Ø 3.0 mm X 25 cm



Luer lock (Stainless Steel)

TWS 84-3055/46	Ø 2.0 mm X 15 cm
TWS 84-3055/47	Ø 2.5 mm X 15 cm
TWS 84-3055/48	Ø 2.5 mm X 20 cm
TWS 84-3055/49	Ø 3.0 mm X 15 cm
TWS 84-3055/50	Ø 3.0 mm X 20 cm
TWS 84-3055/51	Ø 3.0 mm X 25 cm



Power Cannula

TWS 84-3055/52	Ø 2.0 mm X 15 cm
TWS 84-3055/53	Ø 2.0 mm X 20 cm
TWS 84-3055/54	Ø 2.0 mm X 25 cm
TWS 84-3055/55	Ø 2.0 mm X 30 cm
TWS 84-3055/56	Ø 2.0 mm X 35 cm
TWS 84-3055/57	Ø 3.0 mm X 15 cm
TWS 84-3055/58	Ø 3.0 mm X 20 cm
TWS 84-3055/59	Ø 3.0 mm X 25 cm
TWS 84-3055/60	Ø 3.0 mm X 30 cm
TWS 84-3055/61	Ø 3.0 mm X 35 cm
TWS 84-3055/62	Ø 4.0 mm X 15 cm
TWS 84-3055/63	Ø 4.0 mm X 20 cm
TWS 84-3055/64	Ø 4.0 mm X 25 cm
TWS 84-3055/65	Ø 4.0 mm X 30 cm
TWS 84-3055/66	Ø 4.0 mm X 35 cm
TWS 84-3055/67	Ø 5.0 mm X 15 cm
TWS 84-3055/68	Ø 5.0 mm X 20 cm
TWS 84-3055/69	Ø 5.0 mm X 25 cm
TWS 84-3055/70	Ø 5.0 mm X 30 cm
TWS 84-3055/71	Ø 5.0 mm X 35 cm
TWS 84-3055/72	Ø 6.0 mm X 15 cm
TWS 84-3055/73	Ø 6.0 mm X 20 cm
TWS 84-3055/74	Ø 6.0 mm X 25 cm
TWS 84-3055/75	Ø 6.0 mm X 30 cm
TWS 84-3055/76	Ø 6.0 mm X 35 cm



TW Surgical

SPECIALITY CANNULAS



Field Flap

With its elevated leading edge, this cannula is designed to open and expose soft tissue for more precise suctioning along the mandible and sub-mental areas.



Tommey Hub

TWS 84-3055/77	Ø 2.5 mm X 15 cm
TWS 84-3055/78	Ø 2.5 mm X 20 cm
TWS 84-3055/79	Ø 3.0 mm X 15 cm
TWS 84-3055/80	Ø 3.0 mm X 20 cm
TWS 84-3055/81	Ø 3.0 mm X 25 cm
TWS 84-3055/82	Ø 4.0 mm X 25 cm
TWS 84-3055/83	Ø 4.0 mm X 30 cm
TWS 84-3055/84	Ø 4.0 mm X 35 cm
TWS 84-3055/85	Ø 4.5 mm X 30 cm
TWS 84-3055/86	Ø 4.5 mm X 35 cm



Treaded Hub

TWS 84-3055/93	Ø 3.0mm X 15 cm
TWS 84-3055/94	Ø 3.0mm X 25 cm
TWS 84-3055/96	Ø 4.0mm X 15 cm
TWS 84-3055/95	Ø 4.0mm X 25 cm



Ultra Luer Lock

TWS 84-3055/96	Ø 2.0 mm X 15 cm
TWS 84-3055/97	Ø 2.5 mm X 15 cm
TWS 84-3055/98	Ø 2.5 mm X 20 cm
TWS 84-3055/99	Ø 3.0 mm X 15 cm
TWS 84-3056/00	Ø 3.0 mm X 20 cm
TWS 84-3056/01	Ø 3.0 mm X 25 cm



Luer Lock

TWS 84-3056/02	Ø 2.0 mm X 15 cm
TWS 84-3056/03	Ø 2.5 mm X 15 cm
TWS 84-3056/04	Ø 2.5 mm X 20 cm
TWS 84-3056/05	Ø 3.0 mm X 15 cm
TWS 84-3056/06	Ø 3.0 mm X 20 cm
TWS 84-3056/07	Ø 3.0 mm X 25 cm



Luer lock (Stainless Steel)

TWS 84-3055/87	Ø 2.0 mm X 15 cm
TWS 84-3055/89	Ø 2.5 mm X 15 cm
TWS 84-3055/90	Ø 2.5 mm X 20 cm
TWS 84-3055/91	Ø 3.0 mm X 15 cm
TWS 84-3055/92	Ø 3.0 mm X 20 cm
TWS 84-3055/88	Ø 3.0 mm X 25 cm



Power Cannula

TWS 84-3056/08	Ø 2.0 mm X 15 cm
TWS 84-3056/09	Ø 2.0 mm X 20 cm
TWS 84-3056/10	Ø 2.0 mm X 25 cm
TWS 84-3056/11	Ø 2.0 mm X 30 cm
TWS 84-3056/12	Ø 2.0 mm X 35 cm
TWS 84-3056/13	Ø 3.0 mm X 15 cm
TWS 84-3056/14	Ø 3.0 mm X 20 cm
TWS 84-3056/15	Ø 3.0 mm X 25 cm
TWS 84-3056/16	Ø 3.0 mm X 30 cm
TWS 84-3056/17	Ø 3.0 mm X 35 cm
TWS 84-3056/18	Ø 4.0 mm X 15 cm
TWS 84-3056/19	Ø 4.0 mm X 20 cm
TWS 84-3056/20	Ø 4.0 mm X 25 cm
TWS 84-3056/21	Ø 4.0 mm X 30 cm
TWS 84-3056/22	Ø 4.0 mm X 35 cm
TWS 84-3056/23	Ø 5.0 mm X 15 cm
TWS 84-3056/24	Ø 5.0 mm X 20 cm
TWS 84-3056/25	Ø 5.0 mm X 25 cm
TWS 84-3056/26	Ø 5.0 mm X 30 cm
TWS 84-3056/27	Ø 5.0 mm X 35 cm
TWS 84-3056/28	Ø 6.0 mm X 15 cm
TWS 84-3056/29	Ø 6.0 mm X 20 cm
TWS 84-3056/30	Ø 6.0 mm X 25 cm
TWS 84-3056/31	Ø 6.0 mm X 30 cm
TWS 84-3056/32	Ø 6.0 mm X 35 cm

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical

SINGLE-USE INSTRUMENTS

Ultra Luer Lock

The Super Luer Lock is on our Single-use instruments.



Tumescent Infiltrator

This cannula is used specifically for the infusion of tumescent solution and for hydraulic dissection. The configuration of the ports allows even distribution of the fluid.

- TWS 84-3056/33 Ø 2.0 mm X 12 cm
- TWS 84-3056/34 Ø 2.0 mm X 20 cm



Tonnard Harvester

A multiport harvester with 1 mm elevated sharp microports. The 20 port configuration allows for harvesting both on the forward and backward strokes.

- TWS 84-3056/35 Ø 3.0 mm X 20 cm
- TWS 84-3056/36 Ø 2.4 mm X 20 cm



Tonnard Harvester

Excellent for harvesting small aliquots of tissues

- TWS 84-3056/37 Ø 2.5 mm X 15 cm
- TWS 84-3056/38 Ø 2.0 mm X 15 cm



Carraway Harvester

Spiral three port design effective for multi-plane tissue harvesting.

- TWS 84-3056/39 Ø 2.5 mm X 15 cm
- TWS 84-3056/40 Ø 2.0 mm X 10 cm
- TWS 84-3056/41 Ø 2.0 mm X 15 cm
- TWS 84-3056/42 Ø 1.5 mm X 8 cm
- TWS 84-3056/43 Ø 1.5 mm X 10 cm



Injector

Round blunt tip design passes through soft tissue easily without the trauma or potential problems of sharp injectors

- TWS 84-3056/44 Ø 1.5 mm X 9 cm
- TWS 84-3056/45 Ø 1.5 mm X 7 cm
- TWS 84-3056/46 Ø 1.0 mm X 5 cm
- TWS 84-3056/47 Ø 1.0 mm X 4 cm



Spoon Tip Injector

Blunt injector designed with a more open port for rapid injection of larger volumes - a torso and buttock workhouse

- TWS 84-3056/48 Ø 3.0 mm X 20 cm
- TWS 84-3056/49 Ø 2.0 mm X 15 cm





Miller Speed Harvester

- TWS 84-3056/50 Ø 2.0 mm X 15 cm

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com

HANDLES FOR CANNULAS



HANDLES FOR CANNULAS

TWS 84-3056/55

with connection thread for cannulas hose connection



TWS 84-3056/56

with connection thread for cannulas hose connection



TWS 84-3056/57

with connection thread for cannulas hose connection



TWS 84-3056/58

with connection thread for cannulas hose connection



TWS 84-3056/59

Pocar

The Pocar is a percutaneous punch that is used to make a circumferential incision prior to infiltration and tissue extraction. It is designed to maintain a tight seal while minimizing scarring.





TW Surgical

LIPOSUCTION ACCESSORIES

NanoTransfer

TWS 84-3056/60

NanoTransfer




TWS 84-3056/61

Mesh for NanoTransfer



Anaerobic Transfers

- TWS 84-3056/62 Luer 1.2mm
- TWS 84-3056/63 Luer 1.4mm
- TWS 84-3056/64 Luer 2.4mm



Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical

LIPOSUCTION ACCESSORIES

TWS 84-3056/67

Luer to Luer



Decant Stands

- TWS 84-3056/65 60cc to Luer
- TWS 84-3056/66 Luer to Luer



SnapLoks

These syringe snap locks fit into the lunger flanges of the syringe. When the plunger is pulled out the lock snaps open on the barrel lip thus holding the syringe under vacuum.

- TWS 84-3056/68 60cc
- TWS 84-3056/69 50cc
- TWS 84-3056/70 35cc
- TWS 84-3056/71 20cc
- TWS 84-3056/72 10cc





TW Surgical

LIPOSUCTION ACCESSORIES



Universal Power Injector

TWS 84-3056/73

The Power Injector allows the surgeon to re-inject soft tissue through large or small diameter injectors with very little effort and maximum control. The Universal Power Injector will hold a 10cc or a 20cc BD syringe or a BD syringe or a 60cc Monoject syringe.



Fat Syringe

TWS 84-3056/74

The Power Injector allows the surgeon to re-inject soft tissue through large or small diameter injectors with very little effort and maximum control. The Universal Power Injector will hold a 10cc or a 20cc BD syringe or a BD syringe or a 60cc Monoject syringe.

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



Sterilization Box

TWS 84-3056/75



Syringe Decanting and Injection

TWS 84-3056/76	60 cc
TWS 84-3056/77	50 cc
TWS 84-3056/78	30 cc
TWS 84-3056/79	20 cc
TWS 84-3056/80	10 cc



TW
Surgical

Treasure Wright (U.K.) Ltd.

**Address: 128-160 , City Road,
London EC1V 2NX. United Kingdom**

www.twrightmedical.com

info@twrightmedical.com



TW Surgical



**UK
CA CE**