



TW
Surgical

IMPACTOR



TW Surgical

Treasure Wright (U.K.) Ltd.



www.twrightmedical.com



info@twrightmedical.com

UK
CA

CE



SPINAL IMPACTOR

Spinal Impactors are used in neurosurgical procedures that require manipulation of bone, specifically lumbar and cervical vertebrae. This instrument is designed for use in reinforcing and impacting bone grafts that are done and is especially well suited for use in the unique shape of lumbar and cervical vertebrae. Different heads making it highly customizable and a good fit for almost any procedure.



<p>TWS 93-1352/1015</p>	<p>Spinal Bone Graft Impactor, Straight 35cm Length, 5mm x 10mm Concave Grained Impactor Head, Material: Phenolic Handle and Stainless Steel</p>
<p>TWS 93-1352/1419</p>	<p>Spinal Bone Graft Impactor, Angled 35cm Length, 14mm x 19mm Rectangle Grained Impactor Head, Material: Phenolic Handle and Stainless Steel</p>

SPINAL IMPACTOR




TWS 93-1203/0916	Spinal Bone Graft Impactor, Straight 20cm Length, 9mm x 16mm Rectangle Grained Impactor Head, Material: Phenolic Handle and Stainless Steel
TWS 93-1203/1523	Spinal Bone Graft Impactor, Straight 20cm Length, 15mm x 23mm Rectangle Grained Impactor Head, Material: Phenolic Handle and Stainless Steel
TWS 93-1204/0916	Spinal Bone Graft Impactor, Bayonet 20cm Length, 9mm x 16mm Rectangle Grained Impactor Head, Material: Phenolic Handle and Stainless Steel
TWS 93-1204/1523	Spinal Bone Graft Impactor, Bayonet 20cm Length, 15mm x 23mm Rectangle Grained Impactor Head, Material: Phenolic Handle and Stainless Steel

SPINAL IMPACTOR



TWS 93-1353/0916	Spinal Bone Graft Impactor, Straight 35cm Length, 9mm x 16mm Rectangle Grained Impactor Head, Material: Phenolic Handle and Stainless Steel
TWS 93-1353/1523	Spinal Bone Graft Impactor, Straight 35cm Length, 15mm x 23mm Rectangle Grained Impactor Head, Material: Phenolic Handle and Stainless Steel
TWS 93-1354/0916	Spinal Bone Graft Impactor, Bayonet 35cm Length, 9mm x 16mm Rectangle Grained Impactor Head, Material: Phenolic Handle and Stainless Steel
TWS 93-1354/1523	Spinal Bone Graft Impactor, Bayonet 35cm Length, 15mm x 23mm Rectangle Grained Impactor Head, Material: Phenolic Handle and Stainless Steel

SPINAL IMPACTOR



TWS 93-1828/0916	A Bone Graft Impactor, 34cm Length, 9mm x 16mm Horizontal Stripes Impactor, Material: Phenolic Handle and Stainless Steel
TWS 93-1828/0923	L Bone Graft Impactor, 34cm Length, 9mm x 23mm Horizontal Stripes Impactor, Material: Phenolic Handle and Stainless Steel



TWS 93-1829/10	Circular Bone Graft Impactor, 43cm Length, 33cm Shaft 1 cm Diameter Diagonal Lattice Impactor, Material: Stainless Steel
----------------	---

BONE TAMPER

Bone Tamper is a specialized orthopedic tool that surgeons can use for compacting and tamping bone structures and bone grafting into fracture spaces, looking towards to grow new tissue and repair the fracture.




<p>TWS 93-1250/0916NH</p>	<p>Spinal Bone Graft Impactor, Bayonet 25cm Length, 9mm x 16mm Rectangle Grained Impactor Head, Material: Stainless Steel</p>
<p>TWS 93-1250/1523NH</p>	<p>Spinal Bone Graft Impactor, Bayonet 25cm Length, 15mm x 23mm Rectangle Grained Impactor Head, Material: Stainless Steel</p>


BONE TAMPER

TWS Bone Tamper is a specialized orthopedic tool that surgeons can use for tamping bone structures and bone grafting into fracture spaces and compacting to grow new tissue and repair the fracture. Surgeons can tamp small and large grafting into the necessary places.

- Round-Shaped Tip ends in a circular flat surface.
- Solid Design for Deep Surgical Access.
- Knurled handle provides a texture for gripping and minimizes fatigue for Superior Control.
- Serrated tip surface Ergonomic bone grafting packing.
- Premium grade stainless steel, granting the highest quality and performance



BONE TAMPER



TWS 93-6150/02	Bone Tamper 2mm Dia., 160mm Length, Material: Stainless Steel
TWS 93-6150/03	Bone Tamper 3mm Dia., 160mm Length, Material: Stainless Steel
TWS 93-6150/04	Bone Tamper 4mm Dia., 160mm Length, Material: Stainless Steel
TWS 93-6150/05	Bone Tamper 5mm Dia., 160mm Length, Material: Stainless Steel
TWS 93-6150/06	Bone Tamper 6mm Dia., 160mm Length, Material: Stainless Steel
TWS 93-6150/07	Bone Tamper 7mm Dia., 160mm Length, Material: Stainless Steel
TWS 93-6150/08	Bone Tamper 8mm Dia., 160mm Length, Material: Stainless Steel
TWS 93-6150/09	Bone Tamper 9mm Dia., 160mm Length, Material: Stainless Steel
TWS 93-6150/10	Bone Tamper 10mm Dia., 160mm Length, Material: Stainless Steel



TW
Surgical

Treasure Wright (U.K.) Ltd.

**Address: 128-160 , City Road,
London EC1V 2NX. United Kingdom**

www.twrightmedical.com

info@twrightmedical.com