









Ivd and Pituitary Rongeurs – Straight Jaws

Product numbers below are for stainless steel finish. Black coating finish is also available.

Pituitary Rongeur is designed to remove tissues, bone and disc in small surgical sites. The instrument Shaft length are in 130mm, 150mm and comes in Straight Jaws Pattern, Jaws Angled Up Pattern and Jaws Angled Down Pattern. The cup-shaped jaws have the dimensions of 1.5mm, 2mm, 3mm, 4mm, 5mm and 6mm width with 10 mm long. It features short and long shanks that help surgeons in getting better access and reach.



					
1.5mm Jaw	2mm Jaw	3mm Jaw	4mm Jaw	5mm Jaw	6mm Jaw

Ivd and Pituitary Rongeurs – Straight Jaws Material: Stainless Steel

IVD AND PITUITARY RONGEURS SHAFT LENGTH						
Shaft Length	13 cm	15 cm	18 cm	20 cm	25 cm	31 cm
1.5 mm Jaw	TWS 56-1313/15	TWS 56-1315/15	TWS 56-1318/15	TWS 56-1320/15	TWS 56-1325/15	TWS 56-1331/15
2 mm Jaw	TWS 56-1313/20	TWS 56-1315/20	TWS 56-1318/20	TWS 56-1320/20	TWS 56-1325/20	TWS 56-1331/20
3 mm Jaw	TWS 56-1313/30	TWS 56-1315/30	TWS 56-1318/30	TWS 56-1320/30	TWS 56-1325/30	TWS 56-1331/30
4 mm Jaw	TWS 56-1313/40	TWS 56-1315/40	TWS 56-1318/40	TWS 56-1320/40	TWS 56-1325/40	TWS 56-1331/40
5 mm Jaw	TWS 56-1313/50	TWS 56-1315/50	TWS 56-1318/50	TWS 56-1320/50	TWS 56-1325/50	TWS 56-1331/50
6 mm Jaw	TWS 56-1313/60	TWS 56-1315/60	TWS 56-1318/60	TWS 56-1320/60	TWS 56-1325/60	TWS 56-1331/60

Treasure Wright (U.K.) Ltd.

128 City Road, London, EC1V 2NX, UK

info@twrightmedical.com



Ivd and Pituitary Rongeurs – Jaws Angled Up

Product numbers below are for stainless steel finish. Black coating finish is also available.

Pituitary Rongeur is designed to remove tissues, bone and disc in small surgical sites. The instrument Shaft length are in 130mm, 150mm and comes in Straight Jaws Pattern, Jaws Angled Up Pattern and Jaws Angled Down Pattern. The cup-shaped jaws have the dimensions of 1.5mm, 2mm, 3mm, 4mm, 5mm and 6mm width with 10 mm long. It features short and long shanks that help surgeons in getting better access and reach.



1.5mm Jaw	2mm Jaw	3mm Jaw	4mm Jaw	5mm Jaw	6mm Jaw

Ivd and Pituitary Rongeurs – Jaws Angled Up Material: Stainless Steel						
	IVD AND PITUITARY RONGEURS SHAFT LENGTH					
Shaft Length	13 cm	15 cm	18 cm	20 cm	25 cm	31 cm
1.5 mm Jaw	TWS 56-1413/15	TWS 56-1415/15	TWS 56-1418/15	TWS 56-1420/15	TWS 56-1425/15	TWS 56-1431/15
2 mm Jaw	TWS 56-1413/20	TWS 56-1415/20	TWS 56-1418/20	TWS 56-1420/20	TWS 56-1425/20	TWS 56-1431/20
3 mm Jaw	TWS 56-1413/30	TWS 56-1415/30	TWS 56-1418/30	TWS 56-1420/30	TWS 56-1425/30	TWS 56-1431/30
4 mm Jaw	TWS 56-1413/40	TWS 56-1415/40	TWS 56-1418/40	TWS 56-1420/40	TWS 56-1425/40	TWS 56-1431/40
5 mm Jaw	TWS 56-1413/50	TWS 56-1415/50	TWS 56-1418/50	TWS 56-1420/50	TWS 56-1425/50	TWS 56-1431/50
6 mm Jaw	TWS 56-1413/60	TWS 56-1415/60	TWS 56-1418/60	TWS 56-1420/60	TWS 56-1425/60	TWS 56-1431/60

Treasure Wright (U.K.) Ltd.

128 City Road, London, EC1V 2NX, UK

info@twrightmedical.com









Ivd and Pituitary Rongeurs – Jaws Angled Down

Product numbers below are for stainless steel finish. Black coating finish is also available.

Pituitary Rongeur is designed to remove tissues, bone and disc in small surgical sites. The instrument Shaft length are in 130mm, 150mm and comes in Straight Jaws Pattern, Jaws Angled Up Pattern and Jaws Angled Down Pattern. The cup-shaped jaws have the dimensions of 1.5mm, 2mm, 3mm, 4mm, 5mm and 6mm width with 10 mm long. It features short and long shanks that help surgeons in getting better access and reach.



					
1.5mm Jaw	2mm Jaw	3mm Jaw	4mm Jaw	5mm Jaw	6mm Jaw

Ivd and Pituitary Rongeurs – Jaws Angled Down Material: Stainless Steel

IVD AND PITUITARY RONGEURS SHAFT LENGTH						
Shaft Length	13 cm	15 cm	18 cm	20 cm	25 cm	31 cm
1.5 mm Jaw	TWS 56-1513/15	TWS 56-1515/15	TWS 56-1518/15	TWS 56-1520/15	TWS 56-1525/15	TWS 56-1531/15
2 mm Jaw	TWS 56-1513/20	TWS 56-1515/20	TWS 56-1518/20	TWS 56-1520/20	TWS 56-1525/20	TWS 56-1531/20
3 mm Jaw	TWS 56-1513/30	TWS 56-1515/30	TWS 56-1518/30	TWS 56-1520/30	TWS 56-1525/30	TWS 56-1531/30
4 mm Jaw	TWS 56-1513/40	TWS 56-1515/40	TWS 56-1518/40	TWS 56-1520/40	TWS 56-155/40	TWS 56-1531/40
5 mm Jaw	TWS 56-1513/50	TWS 56-1515/50	TWS 56-1518/50	TWS 56-1520/50	TWS 56-1525/50	TWS 56-1531/50
6 mm Jaw	TWS 56-1513/60	TWS 56-1515/60	TWS 56-1518/60	TWS 56-1520/60	TWS 56-1525/60	TWS 56-1531/60

Treasure Wright (U.K.) Ltd.

128 City Road, London, EC1V 2NX, UK

info@twrightmedical.com