



TW
Surgical

PLASTIC MALLET



TW Surgical

Treasure Wright (U.K.) Ltd.



www.twrightmedical.com



info@twrightmedical.com


UK
CA

CE



PLASTIC MALLET

A Plastic Mallet is a surgical instrument used in Orthopedic procedures to apply force to osteotomes, chisels, and gouges for manipulating, cutting, or scraping bones. It provides a comfortable and safe grip. It is designed to last through repeated use and sterilization.



<p>TWS 59.7873</p>	<p>Plastic Mallet Size: 8cm x 6cm x 24cm Long Material: Woven Linen Phenolic</p>
--------------------	--



TW
Surgical

Treasure Wright (U.K.) Ltd.

**Address: 128-160 , City Road,
London EC1V 2NX. United Kingdom**

www.twrightmedical.com

info@twrightmedical.com