

## Minimally Invasive Achilles Tendon Repair Device

The MIS Tendon Suture Repair Device is an elegant method to treat acute Achilles tendon ruptures. It is a minimally invasive procedure that allows direct visual control of the repair, as well as percutaneous introduction of the sutures.

The surrounding soft tissues and the tendon itself are always treated with the utmost care to avoid any local trauma.



TWS 10-2640/30	MIS Tendon Suture Repair Device 120mm length, with AO Fitting Handle 1.5mm Hole x 7 holes Material: Stainless Steel and Phenolic Handle
TWS 111-1516/15	Surgical Needle 1.6mm Dia. 15 cm Long