



TW
Surgical

MALLEABLE RETRACTOR



TW Surgical

Treasure Wright (U.K.) Ltd.



www.twrightmedical.com



info@twrightmedical.com



MALLEABLE RETRACTOR

TW Surgical supplies malleable retractor with a malleable blade that can be shaped to the specific contours of a surgical site, allowing surgeons to gently retract tissues and organs for better access.



<p>TWS 150-1220/10</p>	<p>Malleable Retractor with angled Tip 2mm 150mm x 20mm x 1mm Thickness Material: Stainless Steel</p>
------------------------	---



TW
Surgical

Treasure Wright (U.K.) Ltd.

**Address: 128-160 , City Road,
London EC1V 2NX. United Kingdom**

www.twrightmedical.com

info@twrightmedical.com