



TW
Surgical

OSTEOTOMES CHISELS AND GOUGES



**UK
CA CE**



TW Surgical

Treasure Wright (U.K.) Ltd.



www.twrightmedical.com



info@twrightmedical.com



TW Surgical

OSTEOTOMES, CHISELS AND GOUGES




- TWS 84-2027/81 5mm
- TWS 84-2027/82 10mm
- TWS 84-2027/83 16mm
- TWS 84-2027/8 25mm

OSTEOTOME BLADES




TWS 84-2027/85



TWS 84-2027/86


- TWS 84-2027/87 5mm
- TWS 84-2027/88 10mm
- TWS 84-2027/89 16mm
- TWS 84-2027/90 25mm

CHISEL BLADES



KEY
TWS 84-2027/91

Phenolic handle is autoclavable to 290°F (143°C).



- TWS 84-2027/95 1" width, 50mm radius
- TWS 84-2027/96 1" width, 60mm radius
- TWS 84-2027/97 1" width, 70mm radius

Blades should be discarded after use.

GOUGE BLADES

TWS 84-2027/98 complete set with handle, key and all blades
TWS 84-2027/99 7" black phenolic handle handpiece only

INTERCHANGEABLE OSTEOTOME, CHISEL AND GOUGE SET



Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical



OSTEOTOMES CHISELS AND GOUGES

5" Straight

TWS 84-2028/00	2mm	TWS 84-2028/31	14mm
TWS 84-2028/01	4mm	TWS 84-2028/32	15mm
TWS 84-2028/02	6mm	TWS 84-2028/33	18mm
TWS 84-2028/03	8mm	TWS 84-2028/34	19mm
TWS 84-2028/04	10mm	TWS 84-2028/35	20mm
TWS 84-2028/05	12mm		



MINI-LAMBOTTE OSTEOTOME

The unusually long bevel allows the cut bone to bend gradually as it is cut.

7" Straight		7" Curved	
With calibration lines		With calibration lines	
TWS 84-2028/06	10mm	TWS 84-2028/09	5mm
TWS 84-2028/07	20mm	TWS 84-2028/10	10mm
TWS 84-2028/08	25mm		

LONG BEVEL OSTEOTOME


7" Straight		7" Curved	
With calibration lines		With calibration lines	
TWS 84-2028/11	4mm	TWS 84-2028/20	4mm
TWS 84-2028/12	5mm	TWS 84-2028/21	5mm
TWS 84-2028/13	6mm	TWS 84-2028/22	6mm
TWS 84-2028/14	8mm	TWS 84-2028/23	8mm
TWS 84-2028/15	10mm	TWS 84-2028/24	10mm
TWS 84-2028/16	12mm	TWS 84-2028/25	12mm
TWS 84-2028/17	16mm	TWS 84-2028/26	16mm
TWS 84-2028/18	19mm	TWS 84-2028/27	19mm
TWS 84-2028/19	22mm	TWS 84-2028/28	22mm
		TWS 84-2028/29	25mm

LAMBOTTE OSTEOTOME



TW Surgical

OSTEOTOMES, CHISELS AND GOUGES



9" Straight

Without calibration lines


TWS 84-2028/36	1/4"	6mm
TWS 84-2028/37	1/2"	13mm
TWS 84-2028/38	3/4"	19mm
TWS 84-2028/39	1"	25mm
TWS 84-2028/40	1 1/4"	32mm
TWS 84-2128/40	1 1/2"	38mm

9" Curved

Without calibration lines

TWS 84-2028/41	1/4"	6mm
TWS 84-2028/42	1/2"	13mm
TWS 84-2028/43	3/4"	19mm
TWS 84-2028/44	1"	25mm
TWS 84-2028/45	1 1/4"	32mm
TWS 84-2028/46	1 1/2"	38mm

LAMBOTTE OSTEOTOME



Straight

TWS 84-2028/47	1/8"	3mm
TWS 84-2028/48	3/16"	4mm
TWS 84-2028/49	1/4"	6mm
TWS 84-2028/50	5/16"	8mm
TWS 84-2028/51	3/8"	10mm
TWS 84-2028/52	1/2"	13mm
TWS 84-2028/53	5/8"	16mm
TWS 84-2028/54	3/4"	19mm

Curved

TWS 84-2028/55	1/8"	3mm
TWS 84-2028/56	3/16"	4mm
TWS 84-2028/57	1/4"	6mm
TWS 84-2028/58	5/16"	8mm
TWS 84-2028/59	3/8"	10mm
TWS 84-2028/60	1/2"	13mm
TWS 84-2028/61	5/8"	16mm
TWS 84-2028/62	3/4"	19mm

HOKE OSTEOTOME

5 1/4" (140mm) straight



TW Surgical

OSTEOTOMES, CHISELS AND GOUGES



- TWS 84-2028/63 3/16" (4mm)
- TWS 84-2028/64 1/4" (6mm)
- TWS 84-2028/65 3/8" (10mm)
- TWS 84-2028/66 1/2" (13mm)
- TWS 84-2028/67 3/4" (19mm)

HOKE CHISEL


5 1/4" (133mm)
straight



- TWS 84-2028/68 2mm
- TWS 84-2028/69 4mm
- TWS 84-2028/70 6mm
- TWS 84-2028/71 8mm
- TWS 84-2028/72 10mm
- TWS 84-2028/73 12mm
- TWS 84-2028/74 15mm

SHEEHAN OSTEOTOME

6 1/4" (159mm)
straight



- TWS 84-2028/80 straight
- TWS 84-2028/81 curved, guard right
- TWS 84-2028/82 curved, guard left

SILVER OSTEOTOME

7" (178mm)
with single guard



- TWS 84-2028/77 6mm
- TWS 84-2028/78 12mm
- TWS 84-2028/79 18mm

U.S. ARMY PATTERN OSTEOTOME

7" (178mm)




TW Surgical

OSTEOTOMES, CHISELS AND GOUGES




TWS 84-2028/83	6mm
TWS 84-2028/84	12mm
TWS 84-2028/85	14mm

U.S. ARMY PATTERN CHISEL
7" (178mm)




TWS 84-2028/86	6mm
TWS 84-2028/87	12mm
TWS 84-2028/88	14mm

U.S. ARMY PATTERN GOUGE
7" (178mm)



TWS 84-2028/89	4mm
TWS 84-2028/90	6mm
TWS 84-2028/91	8mm
TWS 84-2028/92	10mm
TWS 84-2028/93	12mm
TWS 84-2028/94	14mm

ALEXANDER CHISEL
7" (178mm)



TWS 84-2028/95	4mm
TWS 84-2028/96	6mm
TWS 84-2028/97	8mm
TWS 84-2028/98	10mm
TWS 84-2028/99	12mm
TWS 84-2029/00	14mm


ALEXANDER GOUGE
7" (178mm)

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com




TW Surgical



TWS 84-2029/01	4mm
TWS 84-2029/02	6mm
TWS 84-2029/03	8mm
TWS 84-2029/04	10mm
TWS 84-2029/05	12mm

LEXER MINI CHISEL

7" (178mm)



TWS 84-2029/06	5mm
TWS 84-2029/07	7mm
TWS 84-2029/08	10mm
TWS 84-2029/09	15mm
TWS 84-2029/10	20mm
TWS 84-2029/11	25mm
TWS 84-2029/12	30mm

LEXER CHISEL

8 3/4" (222mm)

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical

OSTEOTOMES, CHISELS AND GOUGES



1/4"



1/2"



1"

8" Straight

TWS 84-2029/13	1/4"	6mm
TWS 84-2029/14	3/8"	10mm
TWS 84-2029/15	1/2"	13mm
TWS 84-2029/16	5/8"	16mm
TWS 84-2029/17	3/4"	19mm
TWS 84-2029/18	1"	25mm
TWS 84-2029/19	1 1/4"	32mm

8" Curved

TWS 84-2029/20	1/4"	6mm
TWS 84-2029/21	3/8"	10mm
TWS 84-2029/22	1/2"	13mm
TWS 84-2029/23	5/8"	16mm
TWS 84-2029/24	3/4"	19mm
TWS 84-2029/25	1"	25mm
TWS 84-2029/26	1 1/4"	32mm



SMITH PETERSON OSTEOTOME



1/4"



1/2"



3/4"

8" Straight

TWS 84-2029/27	1/4"	6mm
TWS 84-2029/28	3/8"	10mm
TWS 84-2029/29	1/2"	13mm
TWS 84-2029/30	5/8"	16mm
TWS 84-2029/31	3/4"	19mm
TWS 84-2029/32	1"	25mm
TWS 84-2029/33	1 1/4"	32mm

8" Curved

TWS 84-2029/34	1/4"	6mm
TWS 84-2029/35	3/8"	10mm
TWS 84-2029/36	1/2"	13mm
TWS 84-2029/37	5/8"	16mm
TWS 84-2029/38	3/4"	19mm
TWS 84-2029/39	1"	25mm
TWS 84-2029/40	1 1/4"	32mm



SMITH PETERSON GOUGE

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



TW Surgical

OSTEOTOMES, CHISELS AND GOUGES



TWS 84-2029/41	10mm
TWS 84-2029/42	14mm
TWS 84-2029/43	18mm
TWS 84-2029/44	22mm

STILLE-TYPE OSTEOTOME
8" (203mm) straight



TWS 84-2029/45	10mm
TWS 84-2029/46	14mm
TWS 84-2029/47	18mm
TWS 84-2029/48	22mm

STILLE-TYPE CHISEL
8" (203mm) straight



TWS 84-2029/49	10mm
TWS 84-2029/50	14mm
TWS 84-2029/51	18mm
TWS 84-2029/52	22mm

STILLE-TYPE GOUGE
8" (203mm) straight



1"



1 1/4"



9" (229mm) Straight

TWS 84-2029/53	1/4" (6mm)
TWS 84-2029/54	3/8" (10mm)
TWS 84-2029/55	1/2" (13mm)
TWS 84-2029/56	5/8" (16mm)
TWS 84-2029/57	3/4" (19mm)
TWS 84-2029/58	7/8" (22mm)
TWS 84-2029/59	1" (25mm)
TWS 84-2029/60	1 1/8" (29mm)
TWS 84-2029/61	1 1/4" (31mm)
TWS 84-2029/62	1 1/2" (38mm)

9" (229mm) Curved


TWS 84-2029/63	1/4" (6mm)
TWS 84-2029/64	3/8" (10mm)
TWS 84-2029/65	1/2" (13mm)
TWS 84-2029/66	5/8" (16mm)
TWS 84-2029/67	3/4" (19mm)
TWS 84-2029/68	7/8" (22mm)
TWS 84-2029/69	1" (25mm)
TWS 84-2029/70	1 1/8" (29mm)
TWS 84-2029/71	1 1/4" (31mm)
TWS 84-2029/72	1 1/2" (38mm)

HIBBS OSTEOTOME



TW Surgical

OSTEOTOMES, CHISELS AND GOUGES




9" (229mm) Straight

TWS 84-2029/73	1/4" (6mm)
TWS 84-2029/74	3/8" (10mm)
TWS 84-2029/75	1/2" (13mm)
TWS 84-2029/76	5/8" (16mm)
TWS 84-2029/77	3/4" (19mm)
TWS 84-2029/78	7/8" (22mm)
TWS 84-2029/79	1" (25mm)
TWS 84-2029/80	1 1/8" (29mm)
TWS 84-2029/81	1 1/4" (31mm)
TWS 84-2029/82	1 1/2" (38mm)

9" (229mm) Curved

TWS 84-2029/83	1/4" (6mm)
TWS 84-2029/84	3/8" (10mm)
TWS 84-2029/85	1/2" (13mm)
TWS 84-2029/86	5/8" (16mm)
TWS 84-2029/87	3/4" (19mm)
TWS 84-2029/88	7/8" (22mm)
TWS 84-2029/89	1" (25mm)
TWS 84-2029/90	1 1/8" (29mm)
TWS 84-2029/91	1 1/4" (31mm)
TWS 84-2029/92	1 1/2" (38mm)

HIBBS CHISEL



9" (229mm) Straight

TWS 84-2029/93	1/4" (6mm)
TWS 84-2029/94	3/8" (10mm)
TWS 84-2029/95	1/2" (13mm)
TWS 84-2029/96	5/8" (16mm)
TWS 84-2029/97	3/4" (19mm)
TWS 84-2029/98	7/8" (22mm)
TWS 84-2029/99	1" (25mm)
TWS 84-2030/00	1 1/8" (29mm)
TWS 84-2030/01	1 1/4" (31mm)
TWS 84-2030/02	1 1/2" (38mm)

9" (229mm) Curved


TWS 84-2030/03	1/4" (6mm)
TWS 84-2030/04	3/8" (10mm)
TWS 84-2030/05	1/2" (13mm)
TWS 84-2030/06	5/8" (16mm)
TWS 84-2030/07	3/4" (19mm)
TWS 84-2030/08	7/8" (22mm)
TWS 84-2030/09	1" (25mm)
TWS 84-2030/10	1 1/8" (29mm)
TWS 84-2030/11	1 1/4" (31mm)
TWS 84-2030/12	1 1/2" (38mm)

HIBBS GOUGE




TW Surgical


OSTEOTOMES, CHISELS AND GOUGES



SHOULDER PENETRATING GOUGE



SWAN NECK GOUGE



11" (280mm) Straight

TWS 84-2030/15	1/4"	(6mm)
TWS 84-2030/16	3/8"	(10mm)
TWS 84-2030/17	1/2"	(13mm)
TWS 84-2030/18	5/8"	(16mm)
TWS 84-2030/19	3/4"	(19mm)
TWS 84-2030/20	1"	(25mm)
TWS 84-2030/21	1 1/4"	(31mm)
TWS 84-2030/22	2"	(51mm)

11" (280mm) Curved

TWS 84-2030/23	1/4"	(6mm)
TWS 84-2030/24	3/8"	(10mm)
TWS 84-2030/25	1/2"	(13mm)
TWS 84-2030/26	5/8"	(16mm)
TWS 84-2030/27	3/4"	(19mm)
TWS 84-2030/28	1"	(25mm)
TWS 84-2030/29	1 1/4"	(31mm)

COBB OSTEOTOME



TW Surgical


OSTEOTOMES, CHISELS AND GOUGES



TWS 84-2030/30
straight
5/32" (7mm) blade width



TWS 84-2030/31
medium curve
5/32" (7mm) blade width



TWS 84-2030/33
strong curve
5/32" (7mm) blade width



TWS 84-2030/34
reverse curve
5/32" (7mm) blade width

COBB-TYPE SPINAL GOUGE

11" (280mm)

All stainless, with knurled handle


Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com




TW Surgical

OSTEOTOMES, CHISELS AND GOUGES



TWS 84-2030/35
18 1/2" (470mm) with stainless steel handle
.28" (7mm) tip width straight
3 1/2" tip length
10 3/4" shaft length
4 1/4" handle length
.45" (11.5mm) shaft diameter

ENDOSCOPIC HIBBS GOUGE



TWS 84-2030/36
18 1/2" (470mm) with stainless steel handle
3/8" (10mm) tip width straight
3 1/2" tip length
10 3/4" shaft length
4 1/4" handle length
.45" (11.5mm) shaft diameter

ENDOSCOPIC HIBBS OSTEOTOME



TW
Surgical

Treasure Wright (U.K.) Ltd.

**Address: 128-160 , City Road,
London EC1V 2NX. United Kingdom**

www.twrightmedical.com

info@twrightmedical.com



TW Surgical



**UK
CA CE**